325mm

HENRY SCHEIN®

900-6438

INFRARED EAR THERMOMETER

Memory • Fever Alarm • Replaceable Battery

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
The Parts of the unit

1. Probe Cover
2. Probe
3. On/Off Button
4. Test Button
5. Battery Cover

Specifications

| | |
|-----------------------|--|
| Point of Measurement | Infrared measurement at the eardrum |
| Measuring range | 32,0 - 43,0 °C (89,6 - 109,4 °F) |
| Operating conditions | 10 - 40 °C (50 - 104 °F), humidity: 15 - 95 % RH |
| Storage conditions | -25 - 55 °C (-13 - 131 °F), humidity: 15 - 95 % RH |
| Resolution | 0.1 °C or 0.1 °F |
| Accuracy (laboratory) | 35.5 °C - 42.0 °C: +/- 0.2 °C (95.9 °F - 107.6 °F: +/- 0.4 °F); other ranges +/- 0.3°C / +/- 0,6°F |
| Memory | Stores up to 10 measurements in the memory |
| Automatic switch-off | 60 seconds |
| Reaction time | 2 seconds |
| Battery | DC 3V lithium battery CR 2032 |
| Battery life | More than 12 months if used for 10 minutes per day |
| Dimensions | 106 mm (l) x 33 mm (b) x 45 mm (h) |
| Weight | 34,5 g (including battery) |
| Quality guarantee: | The unit conforms to the Council Directive 93/42/EEC concerning medical devices and the standard DIN EN 12470-5:2003 Medical Thermometers, Part 5: Standards for Infra-Red Ear Thermometers. A quality management system in compliance with EN ISO 13485 ensures that the required quality criteria are met. |

- Store out of children's reach! Measurement of children's temperatures must be conducted under adult supervision!
- If the unit is dropped, false readings may occur.
- Put on the new probe cover before testing.
- The probe tip should be cleaned before each use. It can be cleaned by wiping it with a cloth or cotton-wool ball soaked in rubbing alcohol.
- If the lens is dirty, false readings may occur.
- Use of the thermometer is no substitute for medical treatment!
- The unit has an off-set function to compensate for the influence of the ambient temperature.
- When it is switched on, the unit conducts an automatic self-testing procedure. Any faults are shown on the display.

WHAT IS BODY TEMPERATURE?

Body temperature refers to the internal temperature of the body (core temperature). "Normal" body temperature fluctuates over the course of 24 hours. It tends to be lowest in the early morning and highest in the late afternoon. Body temperature also varies depending on the part of the body where the temperature is measured and also differs with age. Hence, the body temperature measurement obtained depends on when and where the temperature is taken and on who the patient is.


| Normal temperature ranges - by point of measurement | |
|---|----------------|
| Axillary | 34.7 - 37.3 °C |
| Oral | 35.5 - 37.5 °C |
| Ear | 35.8 - 38.0 °C |
| Rectal | 36.6 - 38.0 °C |

| Normal temperature ranges - by age | |
|------------------------------------|----------------|
| 0 - 2 years | 36.4 - 38.0 °C |
| 3 - 10 years | 36.1 - 37.8 °C |
| 11 - 65 years | 35.9 - 37.6 °C |
| > 65 years | 35.8 - 37.5 °C |

TEMPERATURE MEASUREMENT - EAR

In order to ensure accurate measurement, the patient should first rest for 20 minutes in the same room as the thermometer. Immediately after exercise or a bath, the temperature may be raised.

1. Place a new probe cover on the sensor firmly, making sure the probe and the inner portion of the probe cover fit snugly. To ensure proper hygiene and to avoid contamination, a new probe cover should be applied for each user. If the probe cover is in good condition, it may be used several times by the same individual.
2. Press the On/Off Button or the Test Button. The LCD screen will appear as shown in Diagram 2 (Note: The last recorded temperature reading will be present.) Wait for screen to change to setting shown in Diagram 3. The thermometer is now ready for a new measurement.
3. Place the probe snugly into the ear canal. Make sure the probe tip faces the eardrum as depicted in Diagram 4. Press the Test Button. A long beep will sound approximately 2 seconds after the Test Button is activated. This signals the end of the process as results are displayed on the LCD screen. (See Diagram 5)
4. The thermometer is ready for testing again once an ear mark displays on the LCD. A waiting period of 20 seconds between testing is recommended to ensure complete accuracy.
5. The thermometer will shut off automatically after 60 seconds of inactivity. To prolong battery life, press the On/Off Button to turn the unit off.

It is advisable to take the temperature in the same ear, since the values measured in the left and right ear may be different.
• In the case of patients whose auditory canal is curved, the temperature values measured may deviate owing to anatomical influences.
• The ear should not be obstructed, for example by excessive amounts of ear wax.

Memory

- Press and hold the On/Off button
- The first memory cell number ("0") will then show the stored temperature (the last reading).
- Press the On/Off button again. Now the reading before last will be shown
- The other values stored in the memory can be viewed by pressing the On/Off button again.


Changing the temperature scale

In order to switch over to Fahrenheit or Celsius scale, please turn off the thermometer then press and hold the on/off button until either c or f is present as shown in diagram 9.

Battery Replacement and Insertion

When the symbol appears on the LCD panel, the battery needs to be replaced.

Press down on the part of the cover marked with the arrow in the direction of the arrow shown in the illustration.

Remove the battery.
Attention: Avoid using any sharp metal object in this operation.

Insert a new CR2032 lithium battery, taking care to match the polarities (+) as indicated.

Replace the battery cover by pushing it in the direction indicated by the arrow in the illustration.

If no signal is received after changing the battery, insert the new battery again and press the ON/OFF button.
When disposing of the unit and/or batteries, be sure to adhere to current regulations.


Fault Indication

| FAULT OR FAULT SYMBOL | FAULT DESCRIPTION | CORRECTIVE MEASURE |
|--|---|--|
| No display on the LCD panel. | The battery has run out. Incorrect battery polarity. | Replace the battery. Please note: The (+) side of the battery must face upwards. |
| Measurement not possible (or an abnormal value is displayed) | The thermometer is not ready. | Wait until the Ear symbol is displayed. |
| An abnormal temperature value is displayed. | The probe tip is dirty or damaged. Did you hear the beep after pressing the ON button? | Clean the probe tip or get it repaired. Wait until you hear the beep before removing the thermometer from the ear. |
| LO or HI symbol is displayed | The temperature measured is outside the measuring range (0 °C -100 °C or 32 °F -212 °F) Low: temperature too low High: temperature too high | Check that the probe tip is clean and that the thermometer is properly inserted into the ear. |
| | The battery has run out. | Replace the battery. |
| | The ambient temperature is outside the operating temperature range. | To ensure accurate measurement, let the thermometer rest at operating temperature for 20 minutes prior to use. |

Limited warranty

This thermometer is guaranteed for 2 years from the date of purchase against any manufacturing defect, conditional upon normal household use. If your unit does not work owing to defective parts or assembly, we undertake to repair it without charge. Apart from the battery, all parts of the unit are subject to this warranty. Damage to the unit caused by improper use is not covered by this warranty. We recommend that the accuracy of the unit be checked annually by an authorised laboratory (in Germany this is compulsory). This checking procedure is not a service provided under the warranty. Use of the unit under operating conditions outside the range indicated may result in false readings. The period of manufacture is indicated by the LOT number on the unit. The first two digits indicate the month and the third and fourth digits the year of manufacture.

Description of the parts

1. Probe Cover
2. Probe
3. On/Off Button
4. Test Button
5. Battery Cover

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