



OUR MISSION:

Offer innovative digital solutions in the field of dental shade taking with high-performance, quality products and a service that meets the expectations of international dentistry.

OUR TEAM:

Strong, motivated and skilled team (average year : 33 yo)

GENERAL DATA:

Address : Limoges (France)
Established date : 10th of june 2013
Regulatory : ISO 13485 (2019)



TABLE OF CONTENTS

1) THE CONCEPT

2) THE TARGETS

3) A VALIDATED TECHNOLOGY

- ⇒ Spectrophotometry and the visible spectrum
- ⇒ The CIE L*a*b* / L*C*h color space
- ⇒ Reading values on Vision software
- ⇒ Spectrophotometer vs Colorimeter
- ⇒ Evolution and limiting factors of manual shade guides

4) COBRA

- ⇒ Pack content
- ⇒ A patented measuring head
- ⇒ Data analysis

5) VISION SOFTWARE

- ⇒ Centralization of aesthetic data
- ⇒ Software tools
- ⇒ Shade guide configuration (composites)
- ⇒ Comparison tool
- ⇒ Smartphone photo app

6) THE BOREA CONNECT

- ⇒ What is that ?
- ⇒ Join the community
- ⇒ The advantages of the platform

7) BOREA WEBSITE

8) THE COMPETITORS

- ⇒ Which solutions?
- ⇒ Complementarity of intraoral scanners
- ⇒ VITA- Easyshade V© comparison

9) COBRA MODULE

10) COBRA & MIX CERAM INTUITION

11) WHY CHOOSE BOREA



COBRA
SHADE SCAN

**BEFORE****NOW**

Two type of information are essential to transmit to the laboratory:

- The optimal **dimension** is scanned using your intraoral scanner.
- The **color** is now recorded and then transmitted thanks to Cobra.

For the first time, you benefit from an all-digital workflow and its many advantages.

1900

MANUAL SHADE GUIDE

2016

SILICONE

2022

COBRA

ALGINATE

INTRA-ORAL SCANNER



PRACTITIONERS

Digitally oriented dentists

- Dentists equipped with intraoral scanners, cameras, etc.
- Dentistry in general is going digital (CEREC systems, scanners, 3D printers, etc.)

Specialists: implantologists and prosthodontists.

- Color is crucial and is usually a source of non-compliance.
- The aesthetic expectations of patients lead them to consider the purchase of a spectrophotometer

General practitioners:

- Indirect and composite restorations.
- Bleaching treatments.

**OBJECTIVITY AND COMFORT
WITH AN ACCURATE AND REPEATABLE SOLUTION**

DENTAL TECHNICIANS

- Help to control the quality of the work carried out.
- Shade taking at the laboratory.
- Flexibility (choice of the level of detail of the cartography: 1, 3 or 9 parts, or the detailed mapping)
- Free software (Vision software and Borea connect access)
- Increased visibility for the laboratory
- Standardization of manufacturing processes



Realize

A complete analysis in
a few seconds



Delete


scrap for
color defects



Benefit



from a free software suite
and visibility

**STANDARDIZATION, ELIMINATION OF SCRAP
AND INCREASED VISIBILITY FOR LABORATORIES**



SPECTROPHOTOMETRY

- A spectrophotometer allows the acquisition of spectral data **over the entire tooth**.
- The wavelengths are analyzed and then compared to a database in order to deduce the color of a tooth.
- Each acquisition is **a safe operation** for the patient


 λ_1, t_1
 λ_2, t_2
 λ_3, t_3
 λ_4, t_4
 λ_5, t_5
 λ_6, t_6


TECHNICAL SPECIFICATIONS

- Cobra performs acquisitions under 6 wavelengths.
- These wavelengths are in the visible spectrum (between 380 nm and 780 nm).
- 4 LEDs per wavelength for homogeneous emission on the tooth.
- An acquisition under white light for the tooth image.
- All acquisitions are polarized.

**COBRA ANALYZES THE ENTIRE TOOTH
THANKS TO THE BEST TECHNOLOGY AVAILABLE TODAY : SPECTROPHOTOMETRY**

LCh tooth color space

Complete Lab color space



CIE L*a*b* coordinates are a **colorimetric standard** defined by the Commission Internationale de l'Éclairage (CIE) in 1976. It allows a multitude of functionalities and standardized communication with the laboratory.

- Calculation of colorimetric differences (ΔE).
- Quantification of bleaching treatments.
- Editing help.
- Comparison of natural tooth and prosthesis in the mouth.
- Compatibility with manufacturing recipe applications.

This space expresses the color in three values :

- L* for the brightness indicated between black (0) and white (100)
- a* from green (-) to red (+)
- b* from blue (-) to yellow (+)

Based on a three-dimensional representation, the CIE Lab system allows a more precise description of colors and approaches the properties of the human eye. It is thus easy to quantify color deviations or to dissociate and study the luminosity (L) of a color.

- 1 $L^*a^*b^*$ / L^*c^*h values of the measured tooth
- 2 Choice of reference shade guide
- 4 Closest shade reference according to the selected shade guide
- 3 Delta E – color difference between the two values
- 5 Shade reference area (green dot)

The information transmitted on the software screen is very clear and show the $L^*a^*b^*$ coordinates of the tooth measured in parallel with the values of the selected commercial shade guide.

The screenshot displays the 'SHADE VIEWER' interface for a patient named 'Smith William'. The interface includes a central image of a tooth with a green dot (5) indicating the measurement area. To the right, the 'COLORIMETRIC TOOLS' section shows the selected shade guide as 'Vita 3D Master' (2) and the measured shade as '1M1' (3). The 'Measured value' section displays the L^* (78.4), a^* (0.0), and b^* (10.7) coordinates. The 'Distance (ΔE)' is shown as 1.9 (4), with a color scale from 'Imperceptible' to 'Visible'.

**VISION SOFTWARE OFFERS INTUITIVE
AND CLEAR READING OF THE $L^*a^*b^*$ COORDINATES**

SPECTROPHOTOMETER vs COLORIMETER

Colorimeters and **spectrophotometers** are the two types of color measurement instruments used to capture, analyze, and communicate color.

Confusion exists among many users and health professionals between spectrophotometry and colorimetry.

Spectrophotometer

- Colorimetric measuring instrument that takes into account the light intensity according to the color.
- **Full-spectrum** colorimetric measurement (as opposed to the trichromatic procedure of a colorimeter).
- Generates colorimetric data that is **beyond observation by the human eye**.



Spectrophotometer - Full visible spectrum



Colorimeter - RGB

How does a spectrophotometer work?

- An internal light source emits at different wavelengths to represent the entire color spectrum.
- Each LED representing a specific wavelength of the spectrum emits one after the other so as to diffuse the light on the tooth.
- Under each emission wavelength, the sensor measures the light intensity and the absorbance of the analyzed tooth.
- These data collected are then analyzed by a comparison algorithm which compares the information collected with the spectral signatures of the various previously recorded shade guides.
- The spectrophotometer displays this information in digital form.

Colorimeter

- A trichromatic color measurement tool that provides an objective assessment of color characteristics from light passing through the red, green, and blue primary filters.
- It simulates the way the human eye perceives color.

How does a colorimeter work?

- The tooth is illuminated at a 45° angle by an internal light source.
- The light passes through the trichromatic filters, which indicates the amount of red, green and blue light reflected by the tooth.
- Filter measurements are quantized in RGB values, which simulate the sensitivity of human eyes to light.

Colorimeter = technology for analyzing the color of intraoral scanners



« Intraoral scanners (...) remain **unreliable** because they have **a narrow choice of colors** and yield results that are significantly different from those of spectrophotometers. As a result, **they should not be used for deciding the shade for a final restoration.**

Spectrophotometers are the best-adapted tools for shade matching since they record the L*, a*, and b* values of one color.

«A comparative in vivo study of new shade matching procedures»
The International Journal of Computerized Dentistry

Dr. Hampé-Kautz, Dr. Ali Salehi, Dr. Bernard Senger & Dr. Olivier Etienne.

[Discover the study here](#)

**INTRAORAL SCANNERS USE COLORIMETRY
AND SHOULD NOT BE USED TO DECIDE THE SHADE !**

Evolution and limiting factors of manual shade guides

Ideal situation

- Neutral environment
- Day light
- Without direct sunlight
- End of the morning
- Cloudy sky

Subjectivity

- Dependent on ocular physiology
- Retinal age
- visual education
- Surrounding light quality
- Patient (lipstick, color of clothes,...)




**ANALYSIS WITH A MANUAL SHADE GUIDE
IS TOTALLY SUBJECTIVE!**

Evolution and limiting factors of manual shade guides




VITA classical A1-D4®.

- Created in 1956
- 16 shades
- Practical shade sample sequence based on hues and empirical findings






VITA SYSTEM 3D-MASTER®.

- Created 1998
- 26 Shades
- Shade sample sequence for optimum coverage of the color space of natural teeth
- Classified by L (luminosity), C (chroma) & h (hue)



Protocol of use of the VITA 3D-MASTER shade guide:

VITA Toothguide 3D-MASTER Instructions for Use – Example




Source : VITA ®

- Complex and time consuming protocol
- Subjective analysis
- Difficulty in transmitting true information to the laboratory
- While Cobra performs this analysis **automatically** and **objectively** on each pixel of the tooth !

OCULAR PHYSIOLOGY + ENVIRONNEMENT = SUBJECTIVITY

PACK CONTENT


**COBRA**
SHADE SCAN

- 1 Cobra Spectrophotometer
- 6 Calibration tips (autoclavable)
- Docking station
- Charger
- Dongle
- Mobile app
- Microfiber wipe
- Quick Start Guide
- Disinfection and sterilization procedure

A sterilizable calibration tip and much more...

A multifunction tip combined with a patented measuring head:

- **Precise:** It allows automatic control of all your acquisitions.
- **Objectivity:** It frees from all external light disturbances.
- **Assisted:** It ensures precise positioning of the device on the tooth and control of the working distance for perfect acquisition.
- **Hygienic:** It eliminates all cross-contamination and each tip is designed for 100 sterilization cycles!!



PRODUCT HIGHLIGHT: Up to 100 sterilization cycles per tip!




CMOS sensor + LEDs

**WORK WITH PEACE OF MIND
THANKS TO AUTOMATIC CALIBRATION**


Available shade-guides

- AIDITE 3D Pro Zir
- IVOCLAR A-D Shade Guide
- IVOCLAR Chromascop
- IVOCLAR Vivodent SP E
- KURARAY Noritake
- ZIRKONZAHN Prettau 2 Dispersive
- ZIRKONZAHN Prettau 3 Dispersive
- SHOFU Vintage
- VITA classical A1-D4®
- VITA Classical A1-D4®+ bleach shades
- VITA Classical A1-D4® (2023)
- VITA Classical A1-D4® (2023) + bleach shades
- VITA Toothguide 3D-MASTER® with bleach shades
- VITA Bleachedguide 3D-MASTER®
- VITA VitaBlocs 3D Master


Shade analysis




Polarised picture




Overall shade




3 parts shade



9 parts shade




Detailed mapping



Translucency

**A COMPLETE MAPPING OF THE TOOTH
IN A SINGLE ACQUISITION**



PRODUCT HIGHLIGHTS:

Positioning assistance with a video stream on the device screen and on your computer screen.

A real **didactic and motivational experience** with your patients.

Improve the dental practice **brand image**.




Operation	Time
Creation of a patient in Vision software (Name & surname at least)	8 sec
Shade-taking positionning	6 sec
Acquisition	1 sec
Automatic data transfer by wifi	8 sec
TOTAL	23 sec
If sending the native file to the lab through Borea connect plateform without order form	25
TOTAL	48 sec

**LESS THAN A MINUTE
FOR A COMPLETE DIGITAL WORKFLOW !**



VISION SOFTWARE AND ITS APPLICATIONS



**FREE SOFTWARE FOR THE DENTAL PRACTICE AND ITS LABORATORIES.
NO OPERATING OR UPDATING COSTS !**



Save and centralize the patient's aesthetic data in a dedicated patient folder:


- Shade file (.rpk)
- Bleaching and Shade comparison report (.pdf)
- 3D model (stl, obj, ply)
- Patient pictures

The screenshot displays the VISION software interface for a patient named CAILLAUD DIANE. The main view shows a dental arch with teeth numbered 11 through 48. A legend indicates that blue circles represent 'Teinte présente' (Shade present) and orange circles represent 'Modèle 3D présent' (3D model present). A sidebar on the left contains navigation options: Borea Connect, Patients, Données, Académie, Paramètres, and Notifications. The top right shows the patient's name and a search icon. Below the patient name, there are buttons for 'IMPORTER (.rpk)', 'AJOUTER PHOTO', 'AJOUTER 3D', and 'AJOUTER TEINTE'. A table on the right lists aesthetic data entries under the heading 'ÉLÉMENTS COMMANDES'.

Type	Date / Heure	N° dent
<input type="checkbox"/>	21/10/2021 - 14:57	--
<input type="checkbox"/>	15/10/2021 - 14:36	--
<input type="checkbox"/>	04/10/2021 - 11:17	16 - 26
<input type="checkbox"/>	17/09/2021 - 17:45	--
<input type="checkbox"/>	R 18/08/2021 - 17:32	23
<input type="checkbox"/>	R 19/07/2021 - 09:36	11


SOFTWARE FOR SIMPLY CENTRALIZING AND ARCHIVING ALL THE AESTHETIC DATA OF YOUR PATIENTS

Tooth outline




- Automatically or manually define the tooth contour.
- The contour is essential for a complete analysis of the tooth.

Shade analysis



- Polarized image
- Global analysis
- 3 and 9 parts
- Detailed mapping
- Instant conversion to different shade guides


Translucency analysis



- Easily view a map of the tooth translucency.
- A gray scale indicates the presence of translucency on a scale of 0% to 100%.


**AN INTERFACE FOR ANALYSIS
ADAPTED TO YOUR NEEDS**

Black and white analysis




- Convert color polarized photo to grayscale.
- Change brightness and contrast settings.
- Bring out the distribution of dentin, enamel, characterize the surface condition, highlight the remarkable characteristics of the analyzed tooth.

Measuring tool



- Locate and measure remarkable characteristics of a tooth.


CIE L*a*b* / L*c*h analysis




- Determine pixel-by-pixel color coordinates according to CIE L*a*b* or L*C*h standards.
- Easily calculate to the pixel the colorimetric difference «Delta E or ΔE» between the actual shade of the tooth and the closest shade according to the reference shade guide.

**AN INTERFACE FOR ANALYSIS
ADAPTED TO YOUR NEEDS**

In the case of direct restorations, it is possible to configure the selected shade guide by indicating the shades actually available in the practice.



Mapping with all the shades of the Vita classical



Cartography after removal of shades B2 to D4.




**CONFIGURE YOUR SHADE GUIDE
ACCORDING TO YOUR AVAILABLE COMPOSITE REFERENCES!**

What's the point ?



The Vision software allows the comparison between several acquisitions.

This shade comparator allows a visual and colorimetric comparison ($L^*a^*b^*$) with the referenced shade guides.

It is the ideal tool to evaluate and quantify bleaching treatments but also to help with correction by allowing the lab to adjust the achievements based on objective criteria.



The screenshot shows the 'ELEMENTS' tab of the Vision software interface. At the top, there are tabs for 'ELEMENTS' and 'ORDERS'. Below the tabs, there are group action buttons: 'COMPARE', 'EXPORT', 'ORDER', and a trash icon. An orange arrow points to the 'COMPARE' button. Below the buttons is a table with columns for 'Type', 'Date / Time', and 'Tooth No.'. The table contains four rows of data, with the last two rows selected (checked).

<input type="checkbox"/>	Type	Date / Time	Tooth No.
<input type="checkbox"/>	 	14/01/2022 - 17:07	24 - 14
<input type="checkbox"/>	 	14/01/2022 - 17:07	45 - 31
<input checked="" type="checkbox"/>	R	14/01/2022 - 16:13	22
<input checked="" type="checkbox"/>	R	11/01/2022 - 15:26	21

A COMPARATOR TO EVALUATE ON OBJECTIVE CRITERIA


Shade comparison

- Compare specific areas of the tooth.
- Ability to select the type of analysis area, circular or rectangular.
- Perform analyzes on one or three areas.



Visual comparison

- Allows to compare two or three polarized images.
- Use the mouse to move the cursor to the right or to the left.



Easily generate pdf reports:

- “**shade comparison report**”
- “**Bleaching reports**” based on Vita Bleachguide 3D Master®.



Scan to discover
a clinical case
[Or click here](#)

VISION
SHADE SCAN

PAGE 1/2

MR PIERRE-FRANÇOIS CRENIER
58 avenue de la Gare
87800 La Meyze
pf.crenier@borea-dental.com

RAPPORT D'ÉCLAIRCISSEMENT
21/10/2021

PATIENT
CAILLAUD DIANE, 2, 31 ans

TEINTIER
Vita Bleachedguide 3D Master

DENT DE RÉFÉRENCE
11

VISION
SHADE SCAN

PAGE 2/2

MR PIERRE-FRANÇOIS CRENIER
58 avenue de la Gare
87800 La Meyze
pf.crenier@borea-dental.com

RAPPORT D'ÉCLAIRCISSEMENT
21/10/2021

PATIENT
CAILLAUD DIANE, 2, 31 ans

TEINTIER
Vita Bleachedguide 3D Master

DENT DE RÉFÉRENCE
11


lundi 19 juillet 2021 / vendredi 17 septembre 2021

Teinte	19/07/2021			17/09/2021			Delta		
	1.5M2	1.5M2	-	1M1.5	1M1.5	-	1M1	1M1	-
L	78.8	80.2	1.4	80.2	81.1	0.9	71.3	68.9	-2.4
a	1.1	0.9	-0.2	-1.2	-0.2	0.9	-0.9	-0.1	0.7
b	16.4	18.2	1.8	10.4	11.0	0.6	6.9	6.6	-0.3

ΔE = 1.7 ΔE = 1.3 ΔE = 2.5

**EASILY PRODUCE BLEACHING REPORTS AND IMPROVE
COMMUNICATION WITH PATIENTS!**

Smartphone application



“Pics App solution”

- Quickly add context photos of your patient through the app.
- Automatic transmission in the patient folder.
- Send up to 4 photos.
- Image editing tools (Crop, rotate ...).
- Free download on Play store and Apple store.
- Photos are not saved on the smartphone.

**DESCRIBE THE AESTHETIC CONTEXT THROUGH
A FREE APPLICATION !**




What is Borea Connect?

- A free order management platform (Practice and laboratory access).
- Easy sending of aesthetic data (shades, patient photos, 3D files)
- Total accessibility from a web browser on computer, smartphone, tablet.

**A FREE COMMUNICATION PLATFORM FOR THE DENTAL PRACTICE
AND ITS LABORATORIES PARTNERS !**

The advantages of the Borea Connect community:

- Monitor prosthetic orders in real time.
- Receive source files (.rpk) and generate complete pdf order forms.
- Benefit from a powerful platform, updated regularly and free of charge.
- Bring visibility to its activity on a community platform where everyone can be seen.




Borea Connect – Dashboard

The 'MY ORDERS' section displays a detailed list of orders with columns for Order No., Date, Status, Prosthetic Patient, Prosthetic, Teeth No., Implant, Attached documents, and Messages. The table includes several rows of data, such as orders for 'LABORATORY BOREA' and 'DELAGE Jean'.

Order No.	Date	Status	Prosthetic Patient	Prosthetic	Teeth No.	Implant	Attached documents	Messages
1407	27 October 8, 2022	IN PROGRESS	Laboratory BOREA	Fixed prosthesis	11, 24		0/1	0/1
1408	27 October 9, 2022	COMPLETED	Laboratory BOREA	Fixed prosthesis			0/1	0/1
1409	27 October 9, 2022	COMPLETED	DELAGE Jean	Fixed prosthesis			0/1	0/1
1745	27 October 18, 2022	COMPLETED	Laboratory BOREA	Fixed prosthesis			0/1	0/1
1584	27 October 26, 2022	ACCEPTED	Laboratory BOREA	Fixed prosthesis			0/1	0/1
1586	27 October 26, 2022	ACCEPTED		Fixed prosthesis			0/1	0/1
1433	27 November 1, 2022	ACCEPTED		Fixed prosthesis			0/1	0/1

Borea Connect – Order tab






Borea Connect – Community map

**OPTIMIZE ORDER MANAGEMENT
AND BRING VISIBILITY TO YOUR ACTIVITY!**




Access to digital order form and color information.
Generate a complete pdf order form.






Full description of the solutions offered by Borea.

- Cobra
- Vision software
- Borea Connect



“Academy” page

- Video
- Tutorials
- Webinars
- Studies
- Clinical cases
- Testimonials



A « Resources" page

- Software download
- Brochures
- Guides
- FAQs
- Support page.

**A COMPLETE WEBSITE FOR AN IN-DEPTH UNDERSTANDING
AND TRAINING!**

What solutions exist for shade taking in the dental office?



Technology	Polarised image Smartphones	Polarised image Cameras	Colorimetry Intra oral scanners	Colorimetry Optishade	Spectrocolorimetry Vita - Easysshade	Spectrophotometry Handy	Spectrophotometry Cobra
Digital shade indication	X	X	✓	✓	✓	✓	✓
Polarised image	✓	✓	X	✓	X	✓	✓
Translucency	X	X	X	X	X	✓	✓
CIE L*a*b* Coordinates	X	✓	X	✓	✓	✓	✓
Shade comparison	X	✓	X	✓	X	✓	✓
Bleaching evaluation	X	✓	X	✓	X	✓	✓
Objectivity / Accuracy	X	✓	X	✓	✓	✓	✓
Autocalibration	X	X	X	✓	X	✓	✓
Calibrated system	X	✓	X	✓	✓	✓	✓
Compatibility PC / MAC	X	✓	✓	X	X	✓	X (2024)
Without subscription fee (Operation / update)	X	X	X	X	✓	✓	✓
Compatibility with softwares for recipes	X	✓	X	✓	X	X	✓
Remote video stream	X	X	✓	✓	X	X	✓
Order follow up platform	X	X	✓	X	X	✓	✓
Acquisition cost	525 € Ex. Tax 630 € Incl. Tax (ie: Smile Lite)	Approx. 3300 € Ex. Tax Approx . 4000 € In. Tax	3 ~ 5K € Ex. Tax 4 ~ 6K € In. Tax (ie: Trios, Emerald, CS)	1590 € Ex. Tax 1908 € In. Tax (ie: Optishade)	1600 € Ex. Tax 1920 € In. Tax	2750 € Ex. Tax 3300 € In. Tax	3300 € Ex. Tax 3960 € In. Tax

No longer available



1

Color and
translucency analysis
+ Patient photos
« Pics App Solution »



2

Dental
impression with
intraoral scanner.



3

Vision software for
processing and
sending data via
Borea connect.



4

Manufacture of
dental
protheses in
the laboratory.




5

Integrated crown
in the smile

- Shade taking reports are compatible **with all scanners** and management software (.pdf).
- Cobra native files (.rpk) are compatible with the majority of intraoral scanner platforms (Medit, Planmeca, Carestream, 3DISC...).
- Digital impressions import to Borea Connect for an independent workflow.

**THE SHADE SCANNER IS THE UNIQUE COMPLEMENT TO THE INTRA ORAL
SCANNER FOR A 100% DIGITAL WORKFLOW!**

Comparison with the Vita spectrophotometer - Easyshade V ©

**RAYPLICKER COBRA**

- Complete analysis of the tooth.
- Spectral analysis of about 419,000 Px
- Different level of cartography (from + to – detailed).
- 3 parts mapping represents an average of all pixels present on each third of the tooth and not a single point.
- Multi shade guides.
- Translucency.
- Self-calibration at each acquisition.
- Immediate transfer to the computer.
- Help with video positioning.

VITA EASYSHADE V ©

- Spectral analysis on a restricted area.
- Area of approximately 9px per acquisition
- Result: 3 reference shades that required 3 color shots
- Compatible Vita classic and Vita 3D master only.
- Non-automatic calibration

1 SINGLE ANALYSIS WITH COBRA = HUNDREDS OF EASYSHADE ANALYZES !

upcoming development



MODULE

COBRA
SHADE SCAN



**All Cobra technology,
reliability and accuracy at
your fingertips.**





A device designed to be integrated into the dental unit:

- Shade taking becomes available like any other handpiece.
- Compatibility with all screens available in the dental office.
- The association with intra-oral scanners allows a complete digital workflow.

**THE COLOR SCANNER INTEGRATED IN THE CHAIR
FOR INCREASED ERGONOMICS!**

Everything becomes easier with prosthesis manufacturing recipes for the laboratory

upcoming development



The ideal recipe application for laboratories

- Compatible with most ceramic powder references on the market.
- Personalized recipes according to the user's stock references.
- A high-performance work tool that standardizes manufacturing processes.
- Accessible from an Internet browser. (computer, smartphone, tablet)

Christian CLAVERIE
Dental technician
&
Denis ELKAÏM
Dental technician
Meilleur Ouvrier de France




<https://weapp.mxceram.com>

BOREA
THE SHADE COMPANY

COBRA
INTUITION

COBRA & MX CERAM INTUITION WITH THEIR RECIPE SYSTEM GUIDE YOU TOWARDS THE RIGHT COLOR!



JOIN THE COBRA COMMUNITY

A validated and recognized solution, take advantage of the experience of many users.

AIMING FOR QUALITY WITHOUT ADDITIONAL COST

Our products are guaranteed for 2 years, parts and labour.
Our solution includes the handpiece, the Vision software and the Borea Connect platform.

No hidden costs, no extras. Only consumables: calibration tips
Each tip is validated for 100 autoclave cycles.

ENJOY FRENCH MANUFACTURING & KNOW-HOW

Our products are medical devices made in France.

We favor local and national suppliers.

Our mission: to develop, innovate and design in compliance with quality standards. BOREA is ISO 13485 certified.

MASTER THE TECHNOLOGY

BOREA is made up of a team of experts who works on possible developments to offer you the best of technology. BOREA has also patented all its innovations.

A MULTI PATENTED SOLUTION:

- Dental shade taking device (2015)
- Process for carrying out and checking compliance of a dental prosthesis from the parameters obtained with a shade taking device (2016)
- Lighting device for tooth shade analysis (2016)
- Automatic association of a tooth shade analysis with a digital dental impression file (stl). (2019)



More info at
www.borea-dental.com



CRENIER Pierre-François (Marketing Manager) pf.crenier@borea-dental.com

CHEVALIER Fabrice (Product specialist) f.chevalier@borea-dental.com

Tel: +33 (0)5 44 00 01 70

