



112-7156

Digital Pediatric Table Scale



Contents

Congratulations!	3
Safety	3
Before you start...	3
Unpacking	3
Power supply	3
Operation.....	4
Remove Protective LCD Film	4
Buttons	4
Display Characters	4
User Weighing	4
Scale Components	5
Cleaning	6
Troubleshooting	6
General Problem Resolution	6
Maintenance.....	7
Changing the batteries	7
Technical data	7
Disposal.....	8
Disposing of the device.....	8
Technical Support.....	8
Warranty.....	8

Congratulations!

Stylish, convenient and easy. These are the words you'll think of as soon as you see the new **Henry Schein Digital Pediatric Table Scale**. The gently curved weighing tray keeps the baby safe during every step of the weighing process. If the baby is restless, our Sta-Bilt stabilizing software feature will ensure accurate measurement even if the baby moves. The large, bright back lit 1.5" LCD figures make it easy and convenient to read the result. With the aid of the Auto-Hold feature, the result is displayed even after the baby has been taken off the scale. Powered by 4 AA batteries or AC power adapter. The integrated measuring tape converts the scales into a complete weighing and measuring system.

Safety

Before using the new scale, please take time to read the following safety instructions.

- Follow the instructions for use.
- Do not put your baby unattended on the scale.
- Avoid excessive movement of the scale.
- The Digital Pediatric Table Scale is a high-precision measuring instrument, please avoid using in extreme temperature environments. Use in such environments may shorten the life of the product and or cause unstable readings.
- Avoid sharp objects when touching the buttons.
- Scale platform cannot be imposed by weight or force greater than the maximum weight (45 lbs and 20kg); otherwise it will damage the sensor.
- Avoid the use of corrosive cleaners.

Before you start...

Unpacking


Remove the packaging and place the scale on a firm, level surface. Included with your scale are the following components:

- 4-AA batteries
- Baby Tape Measure
- Operation Manual
- AC adaptor

Power supply

The scale is generally powered by the supplied 4-AA batteries. The included AC adaptor can be used to power the scale. The 4-AA batteries will provide enough power to power the scale for 2 years.

- Install the supplied 4-AA batteries.
- The connection socket for the AC transformer is in the rear of the base of the scale.



Operation

Remove Protective LCD Film

Before use, please remove the LCD protective film.



Buttons

Zero: When the zero button is pressed, the scale will be reset to zero weight.

KG/LB: Press this button to toggle between KG/LB modes.

Reweight: This will allow a weight to be reinitiated without having the infant removed from the platform.

Display Characters

▶0◀

Indicates that the scale has an accurate zero reading.

HOLD

Indicates a stable weight has been achieved and held on the display.

LOW BATT

Indicates when the batteries require replacement.

User Weighing

1. Place the infant on the scale platform.
2. When a stable weight is reached, the display will lock and hold the weight reading.
3. The scale will hold the weight and after 10 seconds shut off automatically.

Scale Components

- ① Weighing Platform
- ② Touch Panel
- ③ LCD
- ④ Baby Tape Measure
- ⑤ Power Adapter Connection
- ⑥ Battery Compartment
- ⑦ Scale Feet



Fig. 1 Front View




Fig. 2 Back View

Cleaning

Clean all external surfaces with a clean damp cloth or tissue. Mild antimicrobial or antibacterial soap and water solution may be used. Dry with a clean soft cloth. Do not immerse the scale into cleaning or other liquid solutions. Do not use abrasive cleaners.

Under no circumstances should abrasive or acid cleaners, white spirit, benzene or the like be used for cleaning. Such substances can damage the high-quality surfaces. Proper care and cleaning is essential to ensure a long life of accurate and effective operation.

Troubleshooting

General Problem Resolution

Problems	What to Do or Check
Weight reading will not repeat or scale does not return to zero when weight is removed.	Person or object other than that being weighed is touching the scale. Verify that there is nothing caught in the feet area or interfering with the feet.
Scale will not come to zero when the ZERO button is pressed.	Before every use, please make sure the scale is stable and not in use. Press the ZERO button and scale will zero.
Scale drifts off zero.	Check for air currents and/or vibration around the scale.
Scale shuts itself off or will not turn on.	Press the On/Off button to turn on the indicator. Check that new batteries are installed properly in the battery tray.

Error Messages

Error Message	What to Do or Check
Er Ad	The A/D communication is not detected. If this problem persists, please contact Henry Schein technical support for assistance.
ov Ld	The scale weight capacity (75 lbs & 20kg) has been exceeded and is in overload mode. The load on the scale exceeds the scale capacity. Remove excess weight from the scale. If this problem persists please contact Henry Schein technical support for assistance.
95 ol	The scale is in an extreme overload condition. The load on the scale has exceeded the scale capacity by more than 150%. Immediately remove excess weight from the scale as damage not covered under the warranty can occur. If this problem persists please contact Henry Schein technical support for assistance.

Maintenance

The product must be carefully set up and regularly maintained. Depending on the frequency of use, it is recommend servicing the scale every 3 to 5 years by a qualified service lab. For any additional questions please contact Henry Schein Customer Service at 1-888-392-8492 for further information.

Changing the batteries

When **LOW BATT** appears in the display, you should change the batteries. You will need 4-AA batteries.

- Carefully turn over the scale so that the base of the scale is accessible.
- Push the battery compartment closure towards “Battery” and open the cover.
- Remove the used batteries and insert new batteries in the holder.

Technical data

Item # 112-7156	
Width:	28"
Height:	4.25"
Depth:	17"
Platform Size	25"L x 13.38" W
Weight	8.2 lbs
Power supply	6V/300mA power adapter (sold separately)
Batteries	4 x 1.5V alkaline AA batteries
LCD	1.5" text character w/orange back light
Capacity	45 lbs x 0.2 oz (20 kg x 5 g)
Accuracy	5g / 0.2oz

Disposal

Disposing of the device

Batteries

Do not dispose of the device in domestic waste. The device must be disposed of properly as electronic waste. Follow the national regulations which apply in your case.


Do not throw used batteries away in domestic waste. Dispose of batteries at collection points in the vicinity. When buying new batteries, select those low in harmful substances and containing no mercury (Hg), cadmium (Cd) or lead (Pb).

Technical Support

For further information, contact our service department at:

1-888-392-8492

Distributed by (in US only):
HENRY SCHEIN INC.
135 DURYE ROAD
MELVILLE, NY 11747 USA



YYWKRR



Warranty

A two year warranty from date of delivery applies to defects attributable to poor materials or workmanship. All moveable parts - batteries, rechargeable batteries etc. - are excluded. Defects which come under warranty will be made good for the customer at no charge on production of the receipt. No further claims can be entertained. The costs of transport in both directions will be borne by the customer should the equipment be located anywhere other than the customer's premises. In the event of transport damage, claims under warranty can be honored only if the complete original packaging was used for any transport and the scale secured and attached in that packaging just as it was when originally packed. All the packaging should therefore be retained.

A claim under warranty will not be honored if the product is opened by persons not authorized by Henry Schein, Inc.

Henry Schein scales' liability under this warranty is limited to the repair or replacement of the defective product and in no event be liable for consequential or indirect damages to equipment or personnel.

REV DATE: 10/14