

MetriZyme™

Dual enzymatic detergent

Low foaming, dual enzymatic detergent designed for easy one-step pre-cleaning by disintegrating blood and protein from delicate instruments and equipment.

KEEP OUT OF REACH OF CHILDREN



WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

INSTRUCTIONS FOR USE

Use fresh MetriZyme for each instrument or set of instruments. Discard diluted MetriZyme solution after each use.

Manual Cleaning: Add 1 oz. of concentrate to one gallon of warm water (68°F-104°F). Soak instruments and equipment immediately after use, until soil is dissolved and removed. Soak for a minimum of 1 minute. Soak longer if matter is dried on, until visibly clean. Immerse instruments completely and aspirate solution into any channels and lumens. Use an instrument brush for light mechanical cleaning if necessary. Following soak period, visually inspect for cleanliness and rinse thoroughly with warm water. Rough dry. After cleaning with MetriZyme, equipment and instruments must be disinfected/sterilized according to appropriate procedures.

Automatic Reprocessors: Add 1/2 oz. of concentrate per gallon of water to the reprocessor reservoir, or attach container to the reprocessor's automatic dispensing unit.

Ultrasonic Units: Add 1 oz. of concentrate per gallon of water. Recommended water temperature is 68°F-104°F. A soak time of 2 minutes or more is recommended for ultrasonic cleaners. Process instruments according to manufacturer's directions.

Washer Sterilizers/Decontaminators: Attach concentrate container to the unit's automatic dispensing unit.

Note: MetriZyme is not a sterilant/disinfectant.

STORAGE AND DISPOSAL

Storage: Store at ambient temperature.

Solution Disposal: For unused solution, flush thoroughly with large quantities of water into sewage disposal system in accordance with federal, state, and local regulations.

For used solution, the waste solution must be characterized by the generator and disposed of in accordance with federal, state and local regulations.

Container Disposal: Do not reuse empty container. Recycle when available, if not, discard in trash.

RISK PHRASE(S):

- May cause eye and skin irritation.
- May cause throat, mouth and stomach irritation if ingested.
- Repeated inhalation may result in sensitization and allergic reaction in hypersensitive individuals.

PRECAUTIONARY MEASURES:

- Avoid contact with strong oxidizing agents and removal of water.
- Avoid spraying or high pressure washing.

PERSONAL PROTECTIVE EQUIPMENT:

- Protective gloves, eye gear, and impervious gown must be worn.

FIRST AID MEASURES:

EYES- Flush immediately with water for 15 minutes. Seek medical attention. Temporary redness of eye may occur.

SKIN- Wash with soap and water.

INHALATION- Remove to fresh air. Consult physician if allergic response exhibited.

INGESTION- Do not induce vomiting. Consult physician.

FOR FURTHER INFORMATION REFER TO

Metrex Research, MetriZyme Safety Data Sheet

For product or technical information, contact Metrex at **1-800-841-1428** or visit our website at **www.metrex.com**

MADE IN USA

Metrex Research, LLC.
1717 W Collins Ave.
Orange, CA 92867