Instruction Manual

570-2091-Pediatric SpO2 Sensor

It is compatible to 900-8275 Henry Schein Handheld Pulse Oximeter.

1. Intended Use

The 570-2091-Pediatric Sensor is intended for spot checking or continuous monitoring of functional arterial oxygen saturation and pulse rate for pediatric patients in non-invasive with oximeter equipment in hospitals, hospital-type facilities, and home environments.

2. Cleaning and Disinfecting

- Please unplug the sensor from the monitor before cleaning or disinfecting.
- Clean or disinfect the sensor before attaching to a new patient to avoid infection.
- Clean the sensor and patient contact surfaces with soft cloth moistened in water or a mild soap solution.
- To disinfect the sensor, wipe the sensor and patient contact surface with disinfecting solution. We recommend 70% isopropyl alcohol as disinfecting solution.

Caution: Don’t sterilize by irradiation, steam or ethylene oxide.

3. Contraindications

The sensor is a fully compatible re-usable replacement sensor. It is contraindicated for use on active patients.

4. Warnings

- This sensor should be used with compatible monitors, instruments or Oximeter modules only. If not, improper performance may be caused.
- Intravascular dyes may affect measurement accuracy of the sensor and monitor.
- Excessive motions may cause inaccurate measurements during measurement.
- Finger with nail polish or fake nail, or poorly perfused finger, too big finger or improper placement of the finger may cause inaccurate measurements.
- Operation may be affected under bright lights. If necessary, please shield the sensor with an opaque material.
- Operation may be affected by electromagnetic sources.
- The sensor is not recommended to use when imaging equipment work, such as Magnetic Resonance Imaging (MRI).
- Do not immerse the sensor in any kind of liquid. Do not sterilize the sensor with ETO.
- Please stop using the sensor immediately if it is damaged in any way.
- Please stop using the sensor immediately if it is wetted, and wipe it with a dry cloth.
- Do not clip or wrap the measuring position too tightly.
- Please check and reposition the sensor to an alternate location every four hours during continuous monitoring.
- Be careful to the cables to reduce the possibility of patient entanglement or strangulation.
- Follow the pulse oximeter or monitor’s operation manual to connect the sensor with them.
- This equipment should not be adjacent to or stacked with other equipment to use.

5. Measurement Principle

The principle of the Sensor relies on the property of spectrum which collects data from hemoglobin and oxyhemoglobin in glow area and area of near- infrared source and applies the formular “Lammert Beer” to establishing the evaluating equation. In fact, the basic calculating technique is to combine electro-photo SpO2 detecting technology with recording mock pulse rate technology to get the result, the process demonstrates as below: Firstly, use glow with wavelength of 660 nm and approximate infrared with wavelength of 940 nm to put light on the finger tip ready to read analogue data.

Secondly, process the data by electric circuit and microprocessor, and then convey the data to LED display for reading.

Diagram of Operation Principle

- Red and Infrared-ray Emission Tube
- Red and Infrared-ray Receipt Tube

6. Instruction for Use

Appearance


Simple Introduction: The pediatric sensor is for children whose age is between 3 and 12 years old.

- Select a suitable site for the sensor. The first finger is the preferred location. Alternative sites recommended are the thumb, little finger or the large toe.
- The patient’s finger must be inserted right to the end of the sensor.
- The sensor should be oriented in such a way that the cable is positioned along the top of the hand.
- Plug the sensor into the monitor and verify proper operation as described in the instrument’s operator manual.
- Hold the connector rather than the cable when connecting or disconnecting the finger

sensor to the oximeter.

Rx only: "Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner".

Note: Do not twist unnecessarily or use excessive force then using, connecting, disconnecting, or storing the sensor.

	Recommended position
	Correct inserting & Cable position
	Connect to monitor

7. Specifications

SpO2: Measuring Range: 70%~100%
Accuracy: 70%≤SpO2≤100%, ±2 digits; <70%, Unspecified
Probe LED Specifications:

	Wavelength	Radiant Power
RED	660±3nm	1.8mW
IR	940±10nm	2.0mW

Pulse Rate: Measurement Range: 30~235bpm
Resolution: 1bpm
Accuracy: 30≤PR<100, ±2bpm; 100≤PR≤235, ±2%

Storage temperature: -20°C~+70°C

Relative humidity: 10%~95%, No condensation

Operation temperature: +5°C~+40°C

Relative humidity: ≤80%, No condensation

If you have any questions regarding any of this information, please contact our Technical Service Department.

Latex Free. No toxicity, which meets biology compatibility requirement.

The functional tester cannot be used to assess the accuracy of the Oximeter.



Index2 that made by Bioteck company is a function tester. The user can use this function tester to measure the Oximeter.

The test methods used to establish the SpO₂ accuracy is clinical testing. The Oximeter used to measure the arterial haemoglobin oxygen saturation levels and these levels are to be compared to the levels determined from arterial blood sampling with a CO-Oximeter.

8. Classification

According to the degree of protection against electric shock: TYPE BF APPLIED PART; According to the degree of protection against ingress of water: IPX1 According to the mode of operation: CONTINUOUS OPERATION

9. Symbol Explanation

Symbol	Explanation	Symbol	Explanation
	Don't throw away		Keep upright
	Attention, consult the accompanying documents		Fragile
	Storage temperature and relative humidity		Protected from moisture
	Date of Manufacture		Manufacturer's information
	Follow instructions for use		Type BF applied part
	Serial number		Protected against dripping water
	Recyclable		

Note: The illustrations used in this manual may differ slightly from the appearance of the actual product.

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ALL RIGHTS RESERVED