

SIEMENS



Datasheet

ACUSON X700 Ultrasound System

Women's Imaging

Release 2.0

www.siemens.com/ultrasound

TABLE OF CONTENTS

GENERAL INFORMATION	4	FREEZE, CINE AND CINE	
System Architecture	4	POST-PROCESSING FUNCTIONS	10
User Interface	5	Cine Review	10
Control Panel	5	Post-Processing Features in Freeze Frame	
System	5	or Cine	10
Monitor	5	STUDY TYPES	10
TRANSDUCER TECHNOLOGY	6	EXAM-SPECIFIC MEASUREMENTS AND	
Specialty Transducers	6	REPORTS	11
Transducer Ports	6	Abdomen	11
Transducer Storage	6	Obstetrics	11
Hard Drive	6	Gynecology	11
OPERATING MODES	6	Cerebrovascular	11
2D	6	Peripheral Vascular	11
Color Doppler	6	Venous	12
Spectral Doppler	6	Thyroid	12
M-mode	6	DIGITAL PATIENT STUDY STORAGE AND	
DISPLAY MODES	7	ARCHIVING	12
2D-mode Image Processing	7	Patient Study Management	12
MultiHertz Multiple Frequency Imaging	7	ADVANCED TECHNOLOGIES	
Tissue Harmonic Imaging (THI)	7	AND FEATURES	13
Focusing	7	Dynamic TCE Technology	13
2D Image Display	7	Tissue Grayscale Optimization (TGO)	13
2D Calipers – Generic Measurements and		SieClear Compounding	13
Calculations	8	Advanced SieClear Spatial Compounding	13
Color Doppler Velocity Imaging	8	DICOM 3.0 Connectivity	13
Power Doppler Imaging/Directional		DICOM Modality Worklist	14
Power Doppler	8	DICOM MPPS – Modality Performed	
Color and Power Doppler Display	9	Procedure Step	14
Pulsed Wave Spectral Doppler	9	Integrated Gel Warmer (Option)	14
Spectral Doppler Display	9	Clarify Vascular Enhancement (VE)	
Spectral Doppler Calipers – Generic Measurements		Technology (Option)	14
and Calculations	9	SieScape Panoramic Imaging	14
M-mode	9	<i>syngo</i> Auto Follicle Measurements (Option)	15
M-mode Image Display	9	eSie Touch Elasticity Imaging (Option)	15
M-mode Calipers – Generic Measurements and		<i>syngo</i> Auto OB Measurements (Option)	15
Calculations	9	<i>syngo</i> Velocity Vector Imaging (VVI)	
QuickSet User-Programmable System		Technology (Option)	15
Parameters	10	Cadence Contrast Agent Imaging* (Option)	16

TABLE OF CONTENTS

Barcode Reader (Option).....	16
DICOM Structured Reporting (Option).....	16
DOCUMENTATION DEVICES.....	17
Optional On-board Video Devices.....	17
SYSTEM INPUT/OUTPUT.....	17
Video and Audio Input/Output (I/O).....	17
System Interface Connections.....	17
Acoustic Output Management.....	17
SYSTEM DIMENSIONS.....	17
ELECTRICAL/ENVIRONMENTAL SPECIFICATIONS.....	17
INTEGRATING THE HEALTHCARE ENTERPRISE (IHE).....	18
STANDARDS COMPLIANCE.....	18
Quality Standards.....	18
Design Standards.....	18
Acoustic Output Standards.....	18
CE DECLARATION.....	19
EU AUTHORIZED REPRESENTATIVE.....	19

ACUSON X700 Ultrasound System

Sometimes, the most valuable ultrasound feature is the one you can't see. With the ACUSON X700™ ultrasound system, we've engineered a system that offers one of the most critical innovations of all: confidence, in both your diagnoses and your investment. With the ACUSON X700 system, we had the luxury of choosing features from some of the most proven, sophisticated imaging systems on the market, and packaging them in a purpose-built platform. We've designed the ACUSON X700 system with built-in features that permit upgrades on your timetable. This means you have the flexible system you need today, and the freedom to grow to meet the demands of tomorrow.

GENERAL INFORMATION

System Architecture

All-digital signal processing and multi-beam formation technology provide best-in-class imaging in all modes and enable seamless integration of features and options to:

- Allows for parallel Quad beam processing of the RF signal data in the time and amplitude domains with new generation all-digital beamformer technology
 - System imaging frequency: 1.0 MHz to 17.0 MHz
 - 2D-mode line density: Up to 512 lines
 - Processing channels: Up to 23,360 channels
 - Total system dynamic range: > 205 dB
- Native™ tissue harmonic imaging, DTO™ Dynamic Tissue Optimization technology, Dynamic TCE™ tissue contrast enhancement technology, SieScape™ panoramic imaging, Advanced SieClear™ spatial compounding,



SieClear™ compounding and Clarify™ vascular enhancement (VE) technology are built into the ACUSON X700 system to provide excellent image quality

- Enhance productivity through application-specific imaging presets, TGO™ tissue grayscale optimization technology, enhanced measurement and report functionality, knowledge-based workflow tools including syngo® Auto OB measurements, syngo® Velocity Vector Imaging™ (VVI) technology, QuikStart rapid boot and standardized imaging protocols
- Streamline connectivity with solutions such as DICOM Print/Store, DICOM Modality Worklist, DICOM MPPS and DICOM structured reporting for OB/GYN, and Vascular
- Increase functionality with fourSight™ 4D transducer technology, eSie Touch™ elastography imaging, syngo® Auto Follicle, Advanced fourSight™ technology and Cadence™ contrast agent imaging technology*

* At the time of publication, the U.S Food and Drug Administration has cleared ultrasound contrast agents only for use in LVO. Check current regulations for the country in which you are using this system for contrast agent clearance.

User Interface

- On-screen text, control panel overlay and operating instructions are available in the following languages:
 - Chinese
 - English
 - French
 - Italian
 - German
 - Spanish
- Russian keyboard and operating instructions supported
- Operating instructions are available in 17 additional languages
- Thumbnail menu provides on-screen thumbnail images and dynamic clips during exams
- On-screen acoustic power indicator (AIUM/NEMA output display standard)

Control Panel

- Simple, intuitive user interface with Home Base design minimizes repetitive hand motions and enables motor-memory learning
- A floating control panel allows for a wide range of adjustments for operator comfort in standing and sitting positions
- Left/Right swivel articulation $\pm 90^\circ$
- Vertical adjustment: 10 cm
- Backlit key illumination
- Variable key brightness indicates active functions
- Multi-level task light: ON, 1, 2, 3, OFF
- Customizable soft keys for easy and immediate viewing of on-screen menus
- Full-size QWERTY keyboard supports text entry, function keys and system programming
- Wrist support to help reduce operator fatigue

System

- Height-adjustable and swivel control panel
- Control panel with lock lever
- Multi-directional articulating monitor arm to help improve ergonomics
 - Arm rotation: -90° to $+90^\circ$
 - FPD rotation: -80° to $+80^\circ$
 - Tilt: 15° to 90°
 - Up: 125 mm
 - Pull: 250 mm
- Integrated high performance audio speakers
- Wheels
 - Castor size: 125 mm
 - Front castor (2 ea): Bi-brake system (direction lock and total lock)
 - Rear castor (2 ea): Total lock

Monitor

- Flat Panel Display (LED/LCD), 20-inch widescreen, color, high resolution, and progressive scan (non-interlaced) with In Plane Switching (IPS) technology
- Resolution: 1600 x 900 pixels
- Total screen area: 1600 x 900 pixels
- Recordable image area clips: 1024 x 768 pixels
- Total screen capture: 1600 x 900 pixels
- Monitor tilt: 15° up, 90° down and swivel of $\pm 80^\circ$
- Digital on-screen display of brightness control
- Energy saving display power management
- (4) Illumination controls including menu, increase, decrease and light
- Angle of view: 178°

TRANSDUCER TECHNOLOGY

Ultra-sensitive wideband transducers, matched with user-selectable MultiHertz™ multiple frequency imaging, improve resolution and penetration. Up to seven 2D and THI frequencies and up to two color Doppler and spectral Doppler frequencies expand the clinical versatility of a single transducer, thereby maximizing transducer investment.

- Innovative ultra low-loss lens materials and microelectronic technologies for efficient performance and increased signal bandwidth
- microCase™ transducer miniaturization technology and SuppleFlex™ transducer cables
- SuppleFlex cables and integrated cable management provide protection during exams and transport
- Independent 2D and color frequencies for optimal resolution and penetration
- Frequency range: 1.3 – 16.0 MHz
- Hanafy lens acoustic technology
- Universal, stainless steel and disposable biopsy guides for specified linear and curved array transducers

Specialty Transducers

The ACUSON X700 system supports a range of transducers and capabilities. For more detailed transducer information, please consult the Transducer Flyer.

Transducer Ports

- (3) Active transducer ports that support high density phased array, curved array and linear array transducers
- Electronic transducer selection (instantaneous switching between transducers)
- Supports 512 micro-pinless (MP) and DL (260) type connectors
- Industrial design provides easy access to the transducer ports

Transducer Storage

- (6) Configurable transducer holders support all transducer designs and provide gel bottle storage
- Special transducer holder provides secure storage and easy access to endocavity transducer
- Transducer holders can be removed for cleaning

Hard Drive

- Internal 500 GB hard drive
- Allows storage of patient studies that include images, clips, reports and measurements
- Image storage capacity up to 300,000 images with compression
- For more detailed information, please consult the Digital Patient Study Storage and Archiving section on page 12

OPERATING MODES

2D

- Fundamental 2D
- Phased THI
- Filtered THI
- Alternating THI

Color Doppler

- Velocity-based color Doppler
- Power Doppler
- Directional power Doppler
- Color M-mode
- Color Doppler tissue imaging

Spectral Doppler

- Pulsed wave
- Spectral Doppler tissue imaging
- Duplex and Triplex modes

M-mode

- Color M-mode
- Anatomical M-mode

DISPLAY MODES

- Selectable split screen display formats in 2D or 2D/color with M-mode and/or spectral Doppler mode: top-bottom or side-by-side in real-time and digital cine replay
- 4B-mode allowing simultaneous display of 4 static B-mode images
- Virtual Format
- Dual from freeze
- Split/Zoom

Flexible combination of imaging modes in side-by-side Dual and Dual Select in real-time, and digital cine replay.

2D-mode Image Processing

- All-digital parallel signal processing with frame rate up to 498 fps (transducer dependent)
- MultiHertz imaging: Up to 7 user-selectable transmit frequencies
- Res/Speed selection: 6 levels
- Persistence: 5 levels
- Edge enhancement: 4 levels
- Dynamic range selection: 30 – 90 dB in 3 or 5 dB increments, application dependent
- Gain: -30 to +30 dB in 1 dB increments
- DTO technology: 3 levels
- Dynamic TCE technology: 3 levels
- Depth/Gain compensation: 8 controls
- User-selectable gain balance: 9 maps
- User-selectable 2D colorization: 16 maps
- Maximum display depth: 30 cm
- Minimum display depth: 1 cm

MultiHertz Multiple Frequency Imaging

Siemens' unique MultiHertz multiple frequency imaging is designed to combine the resolution and penetration of several transducers into one. At the push of a button, the user can independently change frequencies for 2D, THI, color and spectral Doppler to select the optimal combination for image resolution, penetration and sensitivity.

- Transmit frequencies: Up to 7 user-selectable frequencies
 - 2D and M-mode: Up to 4 fundamental frequencies
 - THI: Up to 4 frequencies
 - Color, power, or pulsed wave Doppler modes: Up to 2 frequencies

Tissue Harmonic Imaging (THI)

Selectable harmonic frequencies increase success with difficult-to-image patients, improving diagnostic confidence, and dramatically improving contrast and spatial resolution by reducing noise and clutter in the image.

- MultiHertz imaging capability in THI
- Available on the 4C1, 6C2, EC9-4w, C6F2, C8F3, EV9F3, VF10-5, and VF12-4 transducers

Focusing

- Transmit focal zones: Up to 4 zones
- Digital dynamic receive focusing with dynamic apodization
- Multi-position, user-selectable position
- Can use multiple focal zones simultaneously

2D Image Display

- Full screen, Split, Quad and Dual Select screen formats as well as Dual, Dual seamless, Dual select and Dual from Freeze
- Curved Vector format
- L/R flip and U/D flip for all formats in real-time and digital cine replay

- Split/Zoom
- Image depth: 1 – 30 cm in 1.0 cm increments (transducer dependent)
- Virtual Format Imaging (transducer dependent)
 - Left/right steer
 - Trapezoid Imaging
- Digital read/write Zoom with image pan
 - Available on live and cine replay images
 - At least 2.5x and up to 10x zoom (transducer dependent)
- 4B-mode
- 90° image rotation (all linear transducers)

2D Calipers – Generic Measurements and Calculations

- Multiple cursor sets on frozen, live, dual screen and cine playback images
- Distance measurements: Up to 8 measurements per screen
 - Distance measurement
 - Depth measurement from skin line
 - Angle measurement
 - Area and circumference: ellipse, trace
- Compound Measurements
 - Volume: user-selectable preset by 1 distance, 2 distance, 3 distance, or 1 ellipse and 1 distance
 - Flow volume: 1 velocity and 1 distance, or 1 velocity and 1 ellipse
 - Stenosis: user-selectable preset calculated by 2 ellipse or 2 distance measurements

Color Doppler Velocity Imaging

- Available on all imaging array transducers
- Multi-beam formation technology provides quad signal processing for color Doppler frame rates up to 177 fps (transducer dependent)
- Transmit frequencies: Up to 2 user-selectable frequencies per transducer

- Left/right steer on all linear transducers
- Color Doppler invert
- Advanced processing in color Doppler mode resulting in excellent spatial resolution and superior Flash suppression
- AutoColor flow state optimization with high, medium and low flow settings
- Color Doppler velocity maps: Up to 9 user-selectable maps (7 velocity and 2 velocity/ variance)
- Velocity scale range: ± 0.6 – ± 244.4 cm/sec
- PRF scale range: 100 – 19,500 Hz (transducer dependent)
- Gain: -20 – 20 dB in 1 dB increments
- Color Doppler line density: 6 selections
- Wall filter: 4 selections
- Color smoothing: 4 levels
- Tissue/color priority: 5 selections
- Color Doppler persistence: 5 levels
- Velocity tag
- Peak hold: Off, 1 sec, 2 sec and 3 sec

Power Doppler Imaging/Directional Power Doppler

- Available on all imaging array transducers
- Multi-beam formation technology provides quad signal processing for power Doppler frame rates up to 195 fps (transducer dependent)
- Left/right steer on all linear array transducers
- Transmit frequencies: Up to 2 user-selectable frequencies per transducer
- Power Doppler maps: Up to 16 maps (8 directional and 8 non-directional)
- PRF scale range: 100 – 19,500 Hz (transducer dependent)
- Gain: -20 – 20 dB in 1 dB increments
- Power Doppler line density: 6 selections
- Wall filter: 4 selections
- Power Doppler smoothing: 4 levels
- Tissue/power Doppler priority: 5 selections
- Color persistence: 5 levels

Color and Power Doppler Display

- 2D/C mode, Split 2D-2D/C mode
- Dual real-time 2D/C mode
- 2D/C/D mode (simultaneous triplex), 2D/C/D mode (update)

Pulsed Wave Spectral Doppler

- Available on all imaging array transducers
- Transmit frequencies: Up to 2 user-selectable frequencies per transducer
- DTI capability available on select transducers
- Sweep speed: 8 selections
- Post-processing gray maps: 8 maps
- Doppler colorization maps: 12 user-selectable maps
- Gain: 0 – 90 dB in 1 dB increments
- PRF range: 100 – 19,500 Hz
- Velocity scale range: ± 1.5 – ± 350 cm/sec with 0° angle correction
- Angle correction: 0 – 89° in 1° increments
- Gate size: 0.2 – 20 mm
- Wall filter: 25 – 3906 Hz, 8 steps (transducer dependent)
- Baseline shift: 17 levels
- Spectral invert
- Autotrace function

Spectral Doppler Display

- Full screen Doppler trace, 2D/Doppler mode, triplex or update 2D/C/Doppler
- Imaging display: 4 formats
 - Top-bottom: 1/3-2/3, 1/2-1/2, 2/3-1/3
 - Side-by-side: 40-60

Spectral Doppler Calipers – Generic Measurements and Calculations

- Multiple cursor sets on frozen and cine playback images
- Velocity/Frequency/Pressure Gradient
- Heart rate/Heart cycle/Time

- Autotrace measurements in real time and freeze including calculations for PS, ED, TAMx, TAMn, PI, RI and S/D
- Resistive Index (RI)
- Pulsatility Index (PI), including Peak-to-Peak method
- Time Average Velocity max (TAV)
- Systolic/diastolic ratio (S/D)
- Velocity Time Integral (VTI)
- Acceleration/Deceleration
- Flow volume using combined velocity and distance, or velocity and ellipse measurements
- Doppler angle correction after measurement

M-mode

- Available on all imaging array transducers
- Anatomical M-mode – live, cineloop and DIMAQ image review
- Frequencies: Up to 5 user-selectable frequencies, including fundamental and harmonics
- Edge enhancement: 4 selections
- Display dynamic range: 30 – 70 dB in 5 dB increments
- Gain: -30 – 30 dB in 1 dB increments
- Gray maps: 9 maps
- M-mode colorization maps: 16 maps
- Sweep speed: 8 selections

M-mode Image Display

- Full screen M-mode, 2D/M-mode
- Imaging display: 4 formats
 - Top-bottom: 1/3-2/3, 1/2-1/2, 2/3-1/3
 - Side-by-side: 40-60

M-mode Calipers – Generic Measurements and Calculations

- Multiple cursor sets on frozen and cine playback images
 - Distance
 - Time
 - Slope
 - Heart rate

QuickSet User-Programmable System Parameters

All imaging modes and parameters are customizable and programmable using QuickSet user-programmable system parameters.

- Up to 32 QuickSet parameters supported
- QuickSet parameters combine all preferred imaging mode parameters, annotation, text and measurements into a single user preset

FREEZE, CINE AND CINE POST-PROCESSING FUNCTIONS

Cine Review

The cine feature is standard and offers the ability to review real-time acquired data. All real-time, post-acquisition optimization functions are available in cine review.

- Frame-by-frame cineloop review and continuous cine motion review, including control of playback rate
- Independent cine review in mixed modes (2D/M, 2D/Doppler, 2D/C/Doppler)
- Independent cine review in 2D Dual Select mode with image align playback feature
- Standard cine memory: Up to 2729 frames
- Acoustic clip capture from cine review
- Anatomical M-mode available

Post-Processing Features in Freeze Frame or Cine

- 2D-mode
 - Zoom/pan
 - Gray map
 - Colorization map
 - Measurements/reports/annotations/pictograms
- Color Doppler
 - Zoom/pan
 - Color map
 - Color invert
 - Measurements/reports/annotations/pictograms



- Spectral Doppler
 - Gray map
 - Doppler colorization map
 - Angle correct
 - Measurements/reports/annotations/pictograms
- M-mode
 - Gray map
 - M-mode colorization map
 - Measurements/reports/annotations/pictograms

STUDY TYPES

The ACUSON X700 system is designed to support many women's healthcare imaging applications. Factory-supplied exam and transducer dependent imaging presets have been carefully optimized for each application to provide consistency, reliability and increased productivity. All applications include body markers, text and annotation labels. Selected applications support customized worksheets and reports. The following study types are supported:

- Abdominal
- Renal
- Obstetrics
- Gynecology

- Early Obstetrics
- Vascular (C-Vas, P-Vas, Venous)
- Small Parts (Breast, Testicle, Thyroid)

EXAM-SPECIFIC MEASUREMENTS AND REPORTS

All measurement and report packages are available for use with all exam types.

All exam specific measurement and reports support:

- All general measurements and calculations
- Comprehensive, customizable, patient reports and worksheets
- Customizable anatomy descriptions
- Physician summary utility – supports on-system report generation including customizable letterhead, patient data, results, graphs, images, comments, recommendations and a customizable signature line

The following Measurements and Reports packages are available on the ACUSON X700 system:

Abdomen

Obstetrics

- Early Obstetrics Gestational Age (GA) measurements are MSD, CRL, and Yolk Sac
- Gestational Age parameter labels are MSD, CRL, BPD, OFD, HC, AC, ATD, ASD, FL, HL, UL, TL, FT, FTA and BN
- 20 user-defined measurement labels
- Calculations include: EFW from the selected reference, HC/AC, TCD/AC, LVW/HW, CorBPD, FL/AC, FL/BPD, CI, AFI, AXT
- Comprehensive Fetal Heart measurements and calculations
- Facial Angle
 - Nuchal Translucency and Nuchal Fold measurements

- Calculations for both Gestational Age (GA) and Estimated Date of Confinement (EDC)
- Early OB and Standard OB patient reports include worksheets for viewing the progress of the report and editing during the exam process
- Multiple fetus reporting capabilities
 - Maximum of 4
- Growth Analysis Graphs with exam file linking
- Detailed Fetal Heart report page

Gynecology

- Micturated and residual volume calculation
- Uterus, Right and Left Ovary, Right and Left Follicle, CRL, MSD, GS and Yolk Sac measurements
- Follicle Measurement supports up to 15 follicles
- Follicles measurement methods supported
 - Distance
 - 2Dist + Avg
 - 3Dist + Avg
 - 2Dist Avg
 - 3Dist Avg
 - Area
 - Volume
 - Circumference

Cerebrovascular

- CCA prox, CCA mid, CCA dist, ICA prox, ICA mid, ICA dist, ECA and VA measurements
- Area Percent Stenosis and Diameter Percent Stenosis measurements

Peripheral Vascular

- CIA, EIA, CFA, PFA, SFA prox, SFA mid, SFA dist, POP A, TRUNK ATA, PTA, PER A and DPA measurements
- Right and left extremity measurements

Venous

- Right and left extremity measurements

Thyroid

- Rt. Lobe, Lt. Lobe Isthmus and nodule measurements
- Supports measurement of up to 15 separate nodules
- Measurements can display all 3 measurements plus volume calculation

The ACUSON X700 system supports customizable labeled measurements (20 B-mode, 5 Doppler and 5 M-mode) for the following exam types: Abdomen, Breast, Thyroid, Venous, Renal, and Small Parts. Additionally, all reporting packages support user-defined descriptors (up to 39).

DIGITAL PATIENT STUDY STORAGE AND ARCHIVING

The DIMAQ-IP integrated workstation allows for digital acquisition, storage and review of complete ultrasound studies, including static images, dynamic clips, measurements, calculations and reports.

Studies can be reviewed and quantified on-board, stored on the system hard drive and transferred to the built-in DVD multi-drive (DVD-R/RW & CD-R/RW) or USB Flash drive for cost-effective archival.

Patient Study Management

Playback of digitally stored images in a selectable 1-up, 4-up, 9-up, 16-up or 25-up screen format. The patient study screen allows search, selection and deleting of studies or export to DVD multi-drive (DVD-R/RW and CD-R/RW).

- 300 GB (of the 500 GB internal hard drive) internal hard drive reserved for patient data management
 - Patient database sorting by Name, ID and Study Date
- Hard drive capacity:
 - Approximately 300,000 B/W and color images
- Storage and retrieval of static images
- Storage and retrieval of cine clips
 - Retrospective clip capture during real-time imaging with a selectable duration of 1, 2, or 4 seconds
 - Prospective clip capture during real-time imaging with a selectable duration of 1 to 120 seconds
- Export of patient studies from hard drive
- Storage and retrieval of reports
- Supports measurements and calculations on archived study and on saved and retrieved images
- Acoustic clip capture from cine review
- M-mode still frame scroll and store
- PW spectral Doppler still frame scroll and store
- Export of patient studies from hard drive to DVD-R/RW and CD-R/RW drive. Studies can be individually selected or batched copied
- The system supports the following data export file formats: RTF, PDF, TIFF, AVI, JPG and DICOM. Connectivity to PACS, other off-line storage (such as USB flash drive) or EMR device is achieved via LAN or WLAN connection.
- Compatible with removable 650 MB, 700 MB and 790 MB CD-R and 650 MB or 700 MB CD-RW
- Removable 4.7 GB single layer DVD and 8.5 GB single side double layer DVD
- Supports export to USB Flash Drive

ADVANCED TECHNOLOGIES AND FEATURES

Dynamic TCE Technology

- Dynamic TCE technology is a proprietary, advanced post-processing method for speckle reduction
- Compatible with other advanced imaging modes including Advanced SieClear compounding, THI and TGO technology
- Supports all primary and secondary exam types
- Three levels available: Low, Medium and High

Tissue Grayscale Optimization (TGO)

TGO tissue grayscale optimization technology provides one-button image optimization. It automatically adjusts image brightness to the tissue type being imaged and equalizes the overall image gain. The user-definable threshold accommodates different user preferences for gain settings and various room lighting conditions. TGO technology improves the consistency and quality of the ultrasound images to enhance productivity by removing time-consuming and operator-dependent manual adjustments. TGO technology can be used with every transducer, for every exam type and at every imaging frequency, including THI.

SieClear Compounding

The SieClear compounding option uses multiple lines of sight to increase contrast resolution and improve tissue differentiation of low contrast lesions by reducing image speckle. Tissue boundaries and interfaces appear sharper and more continuous. SieClear compounding is accessible in THI and is compatible with other advanced imaging options including SieScape imaging, TGO technology and Clarify vascular enhancement technology.



Advanced SieClear Spatial Compounding

This feature combines two distinct technologies to create exceptional image quality: Advanced SieClear spatial compounding and SieClear compounding. This combination of technologies provides exceptional improvements in border definition.

- Up to 7 steering angles available on linear transducer, 5 available on curved array transducers
- Supports all primary and secondary exam types

DICOM 3.0 Connectivity

Enables digital data transfer via a DICOM network for both printing and storage. The ACUSON X700 system acts as a DICOM Storage Class User and DICOM Print Class User.

Functionality supported:

- Connectivity to PACS system for storage of all digital images and dynamic clips with patient demographic data
- In-progress store during the exam
- Image printing to DICOM color and grayscale printers
- DICOM storage commitment
- DICOM exchange media export to DVD-R/RW and CD-R/RW

- DICOM region calibration
- DICOM interchange media viewer software SHOWCASE®
- Interchange media database that identifies the CD to which the patient study has been burned

DICOM Modality Worklist

Enables query and direct download of the patient worklist schedule from the Hospital/Radiology Information System (HIS/RIS) to the ACUSON X700 system, automatically populating the “New Patient” screen with patient demographic information.

DICOM MPPS – Modality Performed

Procedure Step

Enables automatic exchange of Modality Performed Procedure Step information with the Hospital/Radiology Information System (HIS/RIS).

Integrated Gel Warmer (Option)

- Temperature control precision: $\pm 1^\circ$
 - Low: 31°C
 - Medium: 34°C
 - High: 37°C
- Easy-to-use power on/off control switch
- LED color for status indicator
 - Standby mode: Off
 - Operating mode: Orange
- Safety protection for electrical overload
- Weight: ~480 g
- Size: 83 mm x 123 mm x 103 mm

QuikStart (Option)

QuikStart reduces time required for power-up and power-down events by allowing the system to enter a specialized standby mode.

- Standby mode power-down in 3 sec; complete boot-up in less than 12 sec

- Standby time: > 40 min
- Full recharge: < 180 min

Clarify Vascular Enhancement (VE) Technology (Option)

Clarify VE technology is a real-time, adaptive, pixel-by-pixel analysis implemented through a simple, time-saving user interface that provides multiple levels of clarification to optimize tissue contrast resolution and definition of both tissue and vessel walls according to user preference. Clarify VE technology is available on all curved and linear transducers.

SieScape Panoramic Imaging

Grayscale panoramic imaging allows acquisition and display of images up to 60 cm in length to a maximum curvature of 360° .

- Available on all curved and linear transducers
- Can be displayed in full field acquisition or zoomed for detail viewing
- Measurements available

3D/4D

- 3-Scape™ real-time 3D imaging
 - Provides a freehand acquisition technique
 - Supported on the 4C1 and 6C2 transducers
- *fourSight*™ 4D Transducer Technology
 - Provides real-time 3D images
 - Utilizes mechanical acquisition
 - Up to 30 vol/sec
 - C6F2, C8F3 and the EV9F3 4D transducers supported
 - Offers an easy-to-use interface for rapid acquisition and volume rendering Curved Top VOI
 - 4D cine
 - MPR Measurements

- Advanced *fourSight* Technology
 - Offers enhanced 3D/4D acquisition, data rendering and post-processing functionality
 - MultiSlice format allows the user to select range, slice spacing and format for viewing each slice. The MultiSlice format supports up to 36 slices at once.
 - Thick Slice Imaging (TSI) enables definition of a view plane and creates a thick slice around the region of interest. TSI delivers improved contrast resolution and provides more information in a single image.
 - Curved MPR enables real-time multiplanar reformatting of images into any linear or curved plane. This permits the user to set points along a curved object in order to bring all objects along this line into the same plane for viewing, such as the fetal spine.

syngo Auto Follicle Measurements (Option)

The *syngo*® Auto Follicle measurement option is an automated measurement technique that enables fast and accurate assessment of multiple follicles. It helps reduce exam time by allowing the user to automatically capture and record measurements.

eSie Touch Elasticity Imaging (Option)

eSieTouch elasticity imaging is a real-time quantitative imaging method that calculates and displays the relative stiffness of the tissue.

eSieTouch elasticity imaging allows the user to generate the elastogram by applying gentle sequential compression cycles during standard B mode imaging. This relative displacement of tissue is displayed as an elastogram in a live dual image display of the grayscale or color image with the standard B mode image.

- Unique mapping options in grayscale and color further enhance the ease of interpreting an elastogram

- Area, Distance and Strain ratio measurement capability allow for quantitative comparison of two images
- A quality factor indication provides feedback on the quality of acquisition and allows more acute selection of most appropriate frame(s) for assessment or measurement

syngo Auto OB Measurements (Option)

The *syngo* Auto OB measurements option contains an innovative algorithm which provides the ability to automatically measure the most common fetal structures required for fetal biometry: BPD, HC, AC, FL and HL. This technology reduces user variance, providing consistency between operators and reduces keystrokes. All measurement results are saved to the OB report.

syngo Velocity Vector Imaging (VVI) Technology (Option)

syngo VVI technology visualizes, measures and displays global and regional myocardial motion and mechanics from a 2D image and sophisticated tracking algorithm. Using individual vectors to display direction and relative velocity (based on length of the vectors) of frame-to-frame tissue movement, *syngo* VVI technology delivers motion measurement at any point in the cardiac cycle. *syngo* VVI technology provides a unique graphical presentation of the mechanics of myocardial function.

- Algorithm processes ultrasound clips obtained in all views of the heart, as well as generic moving tissue (e.g. vessel wall)
- Not limited by Doppler angle dependencies, frame rate or mean velocities
- ECG signal is not required, allowing use with fetal echo
- Tracking algorithm incorporates multiple sources of information, including speckle tracking

- Displays include:
 - Visual assessment of wall motion and vector dynamics
 - Time curves of the selected velocity vector components
 - 3D representations of parametric color M-mode
 - Selected components of tissue velocity (tissue strain and tissue strain rate)
 - Time curves of global and segmental LV volumes and Ejection Fraction
 - Synchrony analysis including time-to-peak
- Compatible with standard acquisition frame rate clips and acoustic clip capture (e.g. isovolumetric contraction and relaxation events)

Cadence Contrast Agent Imaging* (Option)

Provides excellent detection of the responses from contrast agent microbubbles, coupled with a high resolution display. The ability to simultaneously detect both the signals returning from the contrast agent and the signals arising from tissue allows the user to switch between a contrast-only display and a tissue-only display for instantaneous confirmation of contrast presence and position, providing real-time flexibility and enhanced diagnostic confidence.

- Optimized for 4C1 transducer for abdominal applications
- MultiHertz imaging provides improved fine-tuning for low-MI contrast investigations
- Integrated burst/reflow control for destruction, reperfusion investigations
- On-screen stopwatch feature
- Frame rate triggering with extended clip capture times of up to 20 minutes
- Burst (bubble destruction) mode

Barcode Reader (Option)

- Allows fast and accurate patient information data input
- Easy attachment to USB port
- Supports 2D and 1D patient barcode
- Can scan up to 3 individual barcodes: patient, physician and sonographer
- Inputs the following patient identifying data:
 - Patient Name (First and Last)
 - Patient ID
 - Physician ID
 - Sonographer

DICOM Structured Reporting (Option)

DICOM Structured Reporting (SR) provides a standardized report architecture to allow for easy transfer of measurements to offline PCs, workstations and archiving systems. DICOM Structured Reporting will automatically populate measurements to their respective fields in an external software package. This option is available for the following applications:

- OB/GYN
- Vascular

Wireless Data Transfer (Option)

Utilizes USB dongle to enable wireless connectivity between the ultrasound system and the facility's LAN to provide functionality equivalent to a wired network.

The Wireless Option supports connectivity with:

- DICOM services – Modality worklist, print, storage commitment and store
- Siemens Remote Service – Remote update handling for storage distribution and NetViewer for remote application support and remote

* At the time of publication, the U.S Food and Drug Administration has cleared ultrasound contrast agents only for use in LVO. Check current regulations for the country in which you are using this system for contrast agent clearance.

troubleshooting

Technical Specification

- Standards: IEEE 802.11n, 802.11g, 802.11b, 802.11a
- Security features: WEP, WPA, WPA2 personal, WPA and WPA2 Enterprise

ASUS – Supported in all countries

Ultrasound System Security – Virus Protection (Option)

Embedded virus protection solution that protects the system against advanced persistent threats, viruses, malware and other executing software by detecting and preventing any unwanted change to improve IT compliance and security.

DOCUMENTATION DEVICES

Optional On-board Video Devices

- Up to 2 documentation devices (B/W printer and color printer/DVR recorder) can be integrated into the system cart and controlled from the system control panel
- Supported devices:
 - Mitsubishi USB P95DW B/W Printer
 - Sony UP-D897 B/W Printer
 - Mitsubishi USB CP30DW Color Printer
 - Sony UP-D25MD Color Printer
 - TEAC DVR UR-50BDS

SYSTEM INPUT/OUTPUT

Video and Audio Input/Output (I/O)

- VGA out and DVI out
- DVI in/out

- External audio out

System Interface Connections

- Network
 - (1) Ethernet connector, type RJ45 (10/100 BaseT)
- Peripherals
 - (6) USB 2.0 ports (2 USB ports in front, 4 USB ports in the rear)
 - DVI in/out
 - (2) AC Main Outlet

Acoustic Output Management

- On-screen acoustic power indicator (AIUM/NEMA output display standard)

SYSTEM DIMENSIONS

- Height: 138.6 – 166 .3 cm (54.6 – 65.5 in)
- Width: 55.6 cm (21.9 in)
- Depth: 90.5 cm (35.6 in) without rear handle, 97.9 cm (38.5 in) with rear handle
- Weight: 102 kg (225 lbs)
 - 98 kg (216 lbs) without OEM's

ELECTRICAL/ENVIRONMENTAL SPECIFICATIONS

The ACUSON X700 system is available in one mainframe configuration, suitable for use in 100/115V and 230V environments.

- Power connections: 100-120/200-240 VAC, 50/60 Hz
- Built-in AC isolation transformer
- Power consumption: maximum 600 VA with OEM's
- Atmospheric pressure range: 700 – 1060 hPa

(525 – 795 mm Hg) or up to 3000 m (9,483 ft)

- Ambient temperature range (without OEM's): 10 – 40° C (50 – 104° F)
- Humidity: 30 – 80%, non-condensing, during operation
- Maximum heat output: 2150 BTU/hr
- Average noise level output: 46 – 47 dB

INTEGRATING THE HEALTHCARE ENTERPRISE (IHE)

Having all relevant information at one's fingertips is a prerequisite for optimal and efficient patient care. Seamless integration of the hospital's IT and Imaging Systems and their capabilities to exchange information without restriction are key success factors for facilitating daily work. This is why Siemens has been instrumental in launching and advancing the IHE (Integrating the Healthcare Enterprise) Initiative. Our commitment and dedication enable us to provide clinicians with the ACUSON X700 system one of many innovative products embedded with the building blocks necessary in supporting clinicians' need for seamless health information exchange.

For more information on the ACUSON X700 system and the Siemens commitment to the IHE initiative, please visit www.siemens.com/IHE.

STANDARDS COMPLIANCE

The ACUSON X700 system is in compliance with the following standards, including all applicable amendments at the time of product release.

Quality Standards

- FDA QSR 21 CFR Part 820
- ISO 9001
- ISO 13485

Design Standards

- UL 60601-1
- CSA C22.2 No. 601-1
- EN 60601-1 and IEC 60601-1
- EN 60601-1-1 and IEC 60601-1-1
- EN 60601-1-2 and IEC 60601-1-2 (Class B)
- EN 60601-2-18 and IEC 60601-2-18
- EN 60601-2-37 and IEC 60601-2-37
- EN 60601-2-25 and IEC 60601-2-25
- EN 60601-1-4 and IEC 60601-1-4
- EN 60601-1-6 and IEC 60601-1-6
- ISO 14971
- EN 62304 and IEC 62304

Acoustic Output Standards

- IEC 61157 (Declaration of Acoustic Power)
- IEC 62359 ed. 2 (Test Methods for the Determination of TI and MI)
- AIUM/NEMA UD-2, Acoustic Output Measurement Standard for Diagnostic Ultrasound
- AIUM/NEMA UD-3, Standard for Real-time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment

CE DECLARATION

This product is provided with a CE marking in accordance with the regulations stated in Council Directive 93/42/EEC of June 14, 1993 concerning Medical Devices. The CE marking only applies to medical devices that have been put on the market according to the above referenced Council Directive. Unauthorized changes to this product are not covered by the CE marking and the related Declaration of Conformity.



EU AUTHORIZED REPRESENTATIVE

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
phone: +49 9131 84 0

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens Healthcare Sales Representative for the most current information.

The products/features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons their future availability cannot be guaranteed. Please contact your local Siemens organization for further details.

3-Scape, ACUSON X700, AcuNav, Advanced SieClear, Axius, Cadence, Clarify, DTI, DTO, Dynamic TCE, eSie Touch, *fourSight*, microCase, MultiHertz, Native, QuickSet, SieClear, SieScape, SuppleFlex, SwiftLink, TGO, and Velocity Vector Imaging are trademarks of Siemens Medical Solutions USA, Inc., and *syngo* is a trademark of Siemens AG.

SHOWCASE is a registered trademark of Trillium Technology, Inc.

Microsoft and Windows are registered trademarks of the Microsoft Corporation in the United States and other countries

DS 1014 | © 10.2014, Siemens Medical Solutions USA, Inc.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Telephone: +49 9131 84-0
Germany

www.siemens.com/healthcare

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
685 E Middlefield Road
Mountain View, CA 94043
USA
Telephone: 1-888-826-9702
www.siemens.com/ultrasound