

You've Earned the Right to Wear the Best

At the Heart of your 3M™ Littmann® Stethoscope


3M™ Littmann® Stethoscopes Product Comparison Chart




	Electronic Stethoscopes		Cardiology Stethoscopes			Classic Stethoscopes			Select/Lightweight Stethoscopes	
	Electronic Stethoscope Model 3200	Electronic Stethoscope Model 3100	Master Cardiology™ Stethoscope	Cardiology S.T.C. (Soft Touch Chestpiece) Stethoscope	Cardiology III™ Stethoscope	Master Classic II™ Stethoscope	Classic II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscope	Select Stethoscope	Lightweight II S.E. Stethoscope
Performance (1-10 Scale)	10+	10+	10	9	9	8	7	7	6	6
Main Application Areas	Cardiology/ High Performance Bluetooth®	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Blood Pressure, Limited Physical Assessment	Blood Pressure, Limited Physical Assessment
Chestpiece Design and Technology	Chrome-coated metal	Chrome-coated metal	Polished Stainless Steel Single sided	Matte Finish Stainless Steel Single sided with Soft Touch Insert	Machined Stainless Steel Double sided Tunable	Plated and Polished Alloy Single sided	Machined Stainless Steel Double sided	Machined Stainless Steel Double sided	Coated Alloy Single sided	Metal/Resin Composite Double sided
Headset Material	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy
Overall Length	27 inch/69 cm	27 inch/69 cm	22 inch/56 cm 27 inch/69 cm	27 inch/69 cm	22 inch/56 cm 27 inch/69 cm	27 inch/69 cm	28 inch/71 cm	28 inch/71 cm	28 inch/71 cm	28 inch/71 cm
Total Weight	185 grams	185 grams	175, 185 grams	180 grams	175, 180 grams	160 grams	135 grams	105, 95 grams	120 grams	118 grams
Diaphragm Details, Diameter	51.0 mm	51.0 mm	Tunable 1.75 inch/4.4 cm	Tunable 1.88 inch/4.8 cm	Two Tunable Diaphragms 1.7 inch/4.3 cm 1.3 inch/3.3 cm	Tunable 1.75 inch/4.4 cm	Tunable 1.75 inch/4.4 cm 1.25 inch/3.2 cm	1.31 inch/3.3 cm 1 inch/2.5 cm 1.06 inch/2.7 cm 0.75 inch/1.9 cm	Tunable 1.75 inch/4.4 cm	Tunable 2.125 inch/5.4 cm
Warranty	2 Years	2 Years	7 Years	6 Years	5 Years	3 Years	3 Years	3 Years	2 Years	2 Years

Tunable Technology

Littmann tunable diaphragm technology offers enhanced acoustics, versatility and time-saving convenience. Switching between bell and diaphragm modes is as easy as a simple pressure change on the chestpiece. There is no longer a need to remove, index and replace for same-site auscultation.



Bell Mode
(low-frequency)
Light contact with chestpiece allows you to hear low-frequency sounds.



Diaphragm Mode
(high-frequency)
Without turning the chestpiece over, use firm pressure to hear high-frequency sounds.



3M™ Littmann® Stethoscope Ordering Information and Tube Color Guide

		Black	Burgundy	Navy Blue	Hunter Green	Plum	Chocolate	Gray	Pine Green	Ceil Blue	Purple	Raspberry	Caribbean Blue	Orange	Seafoam Green	Red	Lilac	Royal Blue	Ocean Blue	Peach	Pearl Pink	Coral	Dark Olive Green	Sand	Lavender	Turquoise	Lemon-Lime	Emerald	Champagne	
Model 3200	27 in (69 cm)	3200BK27	3200BU	3200NB				3200T																						
Model 3100	27 in (69 cm)	3100BK27	3100BU	3100NB																										
Master Cardiology™	22 in (56 cm)	2159																												
	27 in (69 cm)	2160 2161 ¹ 2175 ² 2176 ⁵	2163	2164	2165 2183 ²	2167						2178											2182 ⁵							
Cardiology S.T.C.	27 in (69 cm)	4471	4472	4473		4475																								
Cardiology III™	22 in (56 cm)	3127																												
	27 in (69 cm)	3128 3131BE ¹ 3128BRS ² 3157SM ⁵ 3152RBW ⁶	3129	3130	3134	3135	3137 3137CPR ³	3136		3146		3148	3138			3140							3149	3166 ⁵	3165 ⁵	3158 ⁷	3159 ⁵	3161 ⁵		
Master Classic II™	27 in (69 cm)	2144L 2141 ¹ 2142G ⁴	2146	2147	2632				2634	2633	2143		2630												2636 ¹					
Classic II S.E.	28 in (71 cm)	2201 2218BE ¹ 2201BRS ² 2827SM ⁵	2211	2205	2208		2820CPR ³	2203	2818	2813	2209	2210 2829RBW ⁶	2206 2823 ⁶	2812				2215	2819	2822	2817			2942 ⁵	2941 ⁵			2836 ⁵	2837 ⁵	2830 ⁵
Classic II Pediatric	28 in (71 cm)	2113		2123					2131			2122	2119 2153 ⁶	2155		2113R		2136		2154	2115									
Classic II Infant	28 in (71 cm)	2114							2132			2125 2157 ⁶	2124	2179		2114R		2156		2158	2120									
Select	28 in (71 cm)	2290	2293						2305	2301	2294	2296	2291					2298	2306	2310	2292									
Lightweight II S.E.	28 in (71 cm)	2450	2451							2454			2452		2455		2453				2456									

Teaching stethoscopes, catalog numbers 2138 and 2139 as well as veterinary stethoscope, catalog number 1392 are also available.



**3M Health Care
Littmann® Brand Products**
3M Center, Building 275-4E-01
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957
www.3m.com/healthcare
www.Littmann.com/us

**3M Health Care
3M Canada**
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1 800 364-3577

Bluetooth® is a registered trademark of Bluetooth SIG. 3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Master Classic II and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2014. All rights reserved. 70-2008-9603-6



Special chestpiece/headset finishes are only available in a limited number of tubing colors.



¹ Black Edition – black finish chestpiece and eartubes.



² Brass Edition – brass finish chestpiece and eartubes.




³ Copper Edition – copper finish chestpiece and eartubes.



⁴ Gold Edition – gold finish chestpiece and eartubes.



⁵ Smoke Edition – smoke finish chestpiece and eartubes.



⁶ Rainbow Edition – rainbow finish chestpiece and sunshine gold eartubes.



⁷ Rainbow Edition – rainbow finish chestpiece and brass finish eartubes.