

RADIOLOGY



Join Our Community
Join Us on Facebook | Follow Us on Twitter
View Us on YouTube

To Order: 1-800-P-SCHEIN (1-800-772-4346) 8am-9pm (et)
To Fax: 1-800-329-9109 24 Hours
www.henryschein.com/medical



GE imagination at work



KONICA MINOLTA



Table of Contents

Pages 3–8	Contrast Media/Barium Supplies
Page 10	Radiopaque Markers
Pages 11–13	Contrast Delivery
Pages 14–17	Equipment
Pages 18–20	Radiation Protection
Page 21	Accessories
Page 22	X-ray Film and Cassettes
Page 23	Processing
Page 24	Illuminators
Page 25	Disposables
Page 27	Equipment
Pages 28–29	Ultrasound
Pages 30–31	Filing and Storage

Henry Schein is committed to being your valued partner for all the Radiologic specialty needs. This guide is intended to serve as an overview of our most popular radiological items.

Don't see what you are looking for? Ask your Sales Consultant!

From X-ray room equipment to filing and storage needs, we have the solution for you!

Call us today!

Magnevist Injection

Magnevist (brand of gadopentetate dimeglumine) injection is the N-methylglucamine salt of the gadolinium complex of diethylenetriamine pentaacetic acid, an injectable contrast medium for magnetic resonance imaging (MRI) to be administered by intravenous injection. Each mL contains 469.01 mg gadopentetate dimeglumine, 0.39 mg meglumine, 0.15 mg diethylenetriamine pentaacetic acid and water for injection. A clear, colorless to slightly yellow solution which contains no antimicrobial preservative.

#1973676, 5-mL Vial (616-0027)	20/case
#1311935, 10-mL Vial (616-0028)	20/case
#1311745, 15-mL Vial (616-0029)	20/case
#3218450, 20-mL Vial (616-0030)	20/case
#1240340, 100-mL Pharmacy Bulk Vial (616-0013)	10/case
#1213321, 10-mL Prefilled Syringe (616-0031)5/pkg
#1213347, 15-mL Prefilled Syringe (616-0032)5/pkg
#1213727, 20-mL Prefilled Syringe (616-0011)5/pkg

Gadavist

Gadavist™ (gadobutrol) injection is a gadolinium-based contrast agent indicated for intravenous use in diagnostic magnetic resonance imaging (MRI) in adults and children (2 years of age and older) to detect and visualize areas with disrupted blood brain barrier (BBB) and/or abnormal vascularity of the central nervous system.

#1658772, 7.5-mL SDV (616-0017)	20/case
#1641851, 10-mL SDV (616-0018)	20/case
#1665637, 15-mL SDV (616-0020)	20/case
#1664754, 30-mL Pharmacy Bulk Pack (616-0021)	10/case
#1684950, 65-mL Pharmacy Bulk Pack (616-0012)	10/case
#1650787, 7.5-mL Prefilled Syringe (616-0014)5/box
#1655026, 10-mL Prefilled Syringe (616-0015)5/box
#1656826, 15-mL Prefilled Syringe (616-0016)5/box

Eovist

Eovist injection is the first and only gadolinium-based MRI contrast agent approved by the FDA for the detection and characterization of known or suspected focal liver lesions in adults.

#3278959, 10-mL SDV (616-0002)5/box
--	--------

Ultravist® Injection

Ultravist Injection is available as a stable, ready-to-use aqueous solution of iopromide in vials or bottles in concentrations of 240, 300, and 370 mg iodine per mL. As a proven X-ray contrast medium, Ultravist has a broad range of indications for intra-arterial procedures and intravenous procedures.

Ultravist 240 #342-10, 100-mL Vial (616-0003)	10/box
#342-21, 200-mL Vial (616-0025)	10/box
Ultravist 300 #344-05, 50-mL Vial (616-0004)	10/box
#344-10, 100-mL Vial (616-0005)	10/box
#344-12, 125-mL Vial (616-0006)	10/box
#344-15, 150-mL Vial (616-0007)	10/box
#344-21, 200-mL Pharmacy Bulk Pack (616-0024)	10/box
#344-58, 500-mL Pharmacy Bulk Pack (616-0023)	8/box
Ultravist 370 #346-05, 50-mL Vial (616-0008)	10/box
#346-10, 100-mL Vial (616-0009)	10/box
#346-15, 150-mL Vial (616-0022)	10/box
#346-20, 200-mL Vial (616-0026)	10/box
#346-58, 500-mL Pharmacy Bulk Pack (616-0010)	10/box

NOW, MORE THAN EVER, RELY ON HENRY SCHEIN BRAND PRODUCTS

SAVE WITHOUT COMPROMISING ON EXCELLENCE, WITH OUR OWN HIGH-QUALITY HENRY SCHEIN BRAND



We offer over 6,000 quality products you can trust to fulfill your practice needs – each bearing the Henry Schein Seal of Excellence – your guarantee of satisfaction.

- Categories for all of your most commonly used Radiological items
- Fulfillment opt-in program for maximizing your supply purchases
- Extra benefits for PRIVILEGES Members
- Consultation on product mix that is right for you

HENRY SCHEIN®
BRAND PRODUCTS



Isovue 200, 41%
 #1314-30, 50-mL Glass Vial
 (113-4783)10/box



Isovue-M 200, 41%
 #1411-11, 10-mL Vials
 (840-6446)10/box
 #1411-25, 20-mL Vials
 (840-8206)10/box



Isovue-250, 51%
 #1317-05, 50-mL Vial
 (115-6902)10/case
 #1317-02, 100-mL Vial
 (115-3175)10/case



Isovue-M 300, 61%
 #1412-15, 15-mL Vial
 15 mL
 (840-6445)10/case



Isovue 300, 61%
 #1315-25, 30-mL Vial
 (110-9291)10/case
 #1315-30, 50-mL Vial
 (840-7106)10/case
 #1315-47, 75-mL Bottle
 (114-7222)10/box
 #1315-35, 100-mL Bottle
 (109-3019)10/box
 #1315-50, 150-mL Bottle
 (112-4806)10/case
 Isovue Multipack 300, 61%
 #1315-41, 200-mL Bottle
 (114-0736)10/case
 #1315-98, 500-mL Bottle
 (117-4869)6/case



Isovue 370, 76%
 #1316-01, 50-mL Bottle
 (116-9187)10/box
 #1316-52, 75-mL Bottle
 (116-9488)10/box
 #1316-35, 100-mL Bottle
 (109-3061)10/box
 #1316-04, 125-mL Bottle
 (116-0092)10/box
 #1316-37, 150-mL Bottle
 (115-3176)10/box
 #1316-41, 200-mL Bottle
 (118-4392)10/box
 Isovue Multipack® 370, 76%
 #1316-98, 500-mL Bottle
 (840-9716)6/case

Cystografin
 #0149-60, 100-mL Bottle
 (113-1376)10/case
 #0149-57, 300-mL Bottle
 (112-4611)10/case
 #1410-30, Dilute, 300-mL Bottle
 (109-2480)10/case



Gastrografen Solution
 A formulation of iopamidol in orange-flavored aqueous solution for gastroenterography by oral administration or, where necessary, by enema. Its non-ionic nature permits a considerable improvement in the radiological exploration of the gastrointestinal tract compared to traditional contrast media.
 #0445-35, 30-mL Bottle
 (111-7943)24/case
 #0445-40, 120-mL Bottle
 (499-5019)ea

Sinografin
 #0523-30, 10-mL SDV
 (115-3786)10/case

Kinevac
 #0556-15, 5mcg Vials
 (634-0011)10/case



ProHance® Contrast Medium
 ProHance® (gadoteridol) is a contrast medium for magnetic resonance imaging (MRI) in the form of a sterile, apyrogenic solution for intravenous injection. In MRI, ProHance® provides contrast enhancement of the brain, spine and surrounding tissues, resulting in improved visualization (compared with unenhanced MRI) of lesions with abnormal vascularity or those thought to cause a disruption of the normal blood-brain barrier. ProHance® can also be used for whole-body MRI, including the head, neck, liver, breast, musculoskeletal system, and soft tissue pathologies.
 #1111-04, 5-mL SDV
 (113-2005)5/box
 #1111-01, 10-mL SDV
 (113-2003)5/box
 #1111-02, 15-mL SDV
 (113-2006)5/box
 #1111-03, 20-mL SDV
 (113-2002)5/box
 #1111-70, 50-mL Multipack
 (113-2036)5/box
 #1111-16, 10-mL Prefilled Syringe
 (113-1998)5/box
 #1111-45, 17-mL Prefilled Syringe
 (111-9432)5/box



MultiHance® Contrast Agent
 MultiHance® (gadobenate dimeglumine) is a paramagnetic MRI contrast agent. The paramagnetic gadolinium ion is complexed with gadobenate, a chelating agent that forms a highly stable coordinating sphere around the gadolinium ion. It is salified with meglumine, forming a 0.5 M solution. MultiHance is endowed with all the physicochemical features typical of gadolinium class compounds, but in addition it has outstanding relaxivity properties.
 #5164-12, 5-mL SDV
 (111-8330)5/box
 #5164-13, 10-mL SDV
 (111-3834)5/box
 #5164-14, 15-mL SDV
 (114-7082)5/box
 #5164-15, 20-mL SDV
 (111-4242)5/box
 MultiHance® MultiPack,
 #5264-16, 50-mL MDV
 (111-4243)5/box
 #5264-17, 100-mL Bottle
 (113-1997)5/box



Oral Contrast Agents

CT Barium

Esopho-Cat Barium Sulfate Suspension, 3%, w/v
 #4501-11, 30-g Jar with Spoon
(728-1241).....24/case

Volumen LHV Thin Liquid Barium Sulfate for Suspension, 40% w/v after Reconstitution
 #4507-01, 450-mL Bottle
(845-7168).....24/case

Redi-Cat® 2 Barium Sulfate Suspension, 2.1% w/v, 2.0% w/v
 #4503-03, Apple, 450 mL
(172-7836).....24/case

#4503-04, Banana, 450 mL
(300-2922).....24/case

#4503-02, Banana, 250 mL
(931-4686).....24/case

Barium Sulfate Suspension, 2.1% w/v
 #4503-05, 450 mL, Vanilla
(904-0151).....24/case

#4503-06, 250-mL Vanilla
(961-1170).....24/case

#4503-01, 450 mL, Berry
(614-2896).....24/case

#4503-07, 450 mL, Mochaccino
(607-0002).....24/case

#4501-04, 450-mL Bottle
(728-6535).....24/case

#4501-06, 1900-mL Jug
(728-5271).....4/case

#4501-10, 900-mL Bottle
(260-1289).....12/case

Barium Sulfate Suspension, 1.3%, w/v, 1.2% w/v
 #4501-3, 900-mL Bottle
(728-9441).....12/case

Barium Sulfate Suspension, 1.3%, w/v
 #4501-05, 1900-mL Jug
(101-7808).....4/case

#4501-09, 450-mL Bottle
(728-2784).....24/case

E-Z Cat® Dry Barium Sulfate for Suspension, 2.0% w/w after mixing
 #4501-07, 23-g Packet
(166-4048).....50/case

EZ Cat® Barium Sulfate for Suspension Concentrate (4.9% w/v, 4.6% w/w)
 #4501-01, 225-mL Bottle mL
(240-3648).....24/case

Capacity Mixing Container with Lid
 #4500-01, 32 oz
(614-2868).....25/case



Varibar®

40% w/v

#900002, Nectar Barium Sulfate Suspension, 240 mL
(695-8941).....24/case

#900005, Honey Barium Sulfate Suspension, 250 mL
(110-1850).....12/case

#900006, Pudding Barium Sulfate Esophageal Paste, 230 mL
(846-3564).....12/case

40% w/v after reconstitution
 #900001, Thin Liquid Barium Sulfate for Suspension, 148g
(630-2266).....24/case



E-Z-Paque® Barium Sulfate for Suspension

#901901, 96% w/w, 176-g Bottle
(728-3052).....24/case
 #902601, 96% w/w, 1200-g Jug
(728-4807).....8/box
 #9030-01, 96% w/v, 10-kg Pail
(102-4420).....ea

Liquid E-Z-Paque® Barium Sulfate Suspension

#902901, 60% w/v, 41% w/w, 1900-mL Jug
(728-2700).....ea
 #902801, 60% w/v, 41% w/w, 355-mL Bottles
(728-4267).....24/case

E-Z-Paste® Barium Sulfate Esophageal Cream

#902101, 60% w/w, 1 lb.
(101-7200).....12/case

Ultra-R® Barium Sulfate for Suspension

#901902, 95% w/w, 6 oz
(386-0665).....24/case



Liquid Polibar® and Liquid Polibar® Plus Barium Sulfate

Contrast agent used to diagnose problems with the gastrointestinal tract.

Polibar® ACB Barium Sulfate (96% w/w)
 #9005-04, Prefilled Disposable Barium Enema Kit Powder for Rectal Suspension, 397-g Bag
(114-8185).....24/case

#9005-06, Polibar ACB Barium in Snap-Cap Enema Bag, 16 oz
(110-3173).....24/case

#9003-02, Empty Enema Kits, Empty Snap Bag
(702-5891).....48/case

#9003-04, Empty Enema Kits, 1-Piece Screw Cap, Exacta Bag
(337-5899).....48/case

#4500-02, CT Enema Kit, Empty Enema Kit with Conventional Rigid Tip
(356-3980).....12/case

E-Z Dose Kit

#9001-01, Unit Dose Container with 650 mL of Liquid Polibar Plus, ½" x 60" Tubing with Miller Air Tip with Silicone Trimline Retention Cuff
(113-9505).....6/case

#9002-02, 64-oz Liquid Polibar®, 100% w/v
(728-0257).....ea

#9002-03, 64-oz Liquid Polibar® Plus, 105% w/v
(670-7082).....4/case

#9013-01, Enema Tip, Rounded Flex-Tip®
(993-3130).....24/case

#9013-03, Enema Tip, Infant Flex-Tip®
(392-4248).....24/case

#9013-05, Enema Tip, 14F Pediatric Flex-Tip®
(110-2372).....24/box

#9012-02, Miller Air Tip
(399-7535).....48/case


#9012-01, Flexi-Cuff with Non-Latex Cuff
(728-7929).....48/case

#9013-06, Cone Colostomy Tip
(108-3925).....12/box

#9013-07, Nipple Colostomy Tip
(333-7672).....12/box

#9004-04, Air Bulb Insufflator
(889-0738).....12/box

#9004-05, One-Shot Cuff Inflator
(889-7114).....12/box



E-Z-HD™
 #9017-02, High Density Barium Sulfate for Suspension 250% w/v, 340-g (12-oz) Bottle
(728-7974)24/case

Digibar® 190 Barium Sulfate for Suspension
 #9017-01, 190% w/v, 12 oz
(158-1192)24/case

Maxibar™ Barium Sulfate Suspension
 #9027-01, 210% w/v, 120 mL
(338-5477)24/case



E-Z-Gas® II Effervecent Granules
 #9020-01, 4 g
(728-1046)50/box

E-Z-Disk™ Barium Sulfate Tablets
 #9021-02, 648g
(302-5261)100/btl

LiquiPrep Magnesium Citrate Oral Solution
 #3902-01
(728-2894)24/case

LoSo Prep™ Bowel Cleansing System
 Contains one packet of magnesium carbonate, citric acid, and potassium citrate for oral solution of 18-g magnesium citrate; four 5-mg bisacodyl tablets, USP; one 10-mg bisacodyl suppository, USP.
 #3902-02
(728-2659)10/case

E-Z-Scan®
 #6010-03, 1 Gal
(325-4766)4/case
 #6010-04, 5 Liter
(672-3032)4/case

EnteroVu™ Barium Sulfate for Suspension
 #9014-01, 81% w/w, 110-g Packets
(687-6031)12/box

Liquid EnteroVu Barium Sulfate Suspension
 #9014-06, 13%, 600 mL
(341-7368)12/box

Liquid EnteroVu Barium Sulfate Suspension
 #9014-07, 24% w/v 600 mL
(102-6811)12/box

Enteroclysis Catheter
 #9015-01, 155cm
(257-5786)ea

#9015-03, 152cm
(362-7717)ea

Entero-H™ Barium Sulfate Suspension
 #9014-08, 80 w/v, 1900 mL
(338-2289)4/case

Cysto-Aid™
 #9010-04, 1000 mL
(457-1313)25/case

Vinyl Retention Ring
 #9011-01
(728-0373)12/box

HAVE YOU JOINED PRIVILEGES?

Enroll today!

**EXCLUSIVE
REWARDS FOR
PREFERRED
CLIENTS**

1.866.MED.VIPS (1.866.633.8477)
www.henryschein.com/privilegesmd

Optiray® 240 Injectable

Nonionic, low-osmolar contrast media agent intended for intravascular administration. Therapeutically and biologically inert when injected into body for use in organ or tissue enhancement in CT, X-ray, and fluoroscopy imaging procedures. Indicated for use in adults for angiography throughout the cardiovascular system. Uses include cerebral, coronary, peripheral, visceral, and renal arteriography, venography, aortography, left ventriculography, CT of head and body, and pediatric excretory urography. Not for intrathecal use. Each mL contains 509 mg of ioversol, 3.6 mg of tromethamine as a buffer, & 0.2 mg of edetate calcium disodium as a stabilizer.

- #1324-06, 50-mL Bottle
(748-0074).....25/box
- #1324-11, 100-mL Bottle
(748-0060).....12/box
- #1324-16, 150-mL Bottle
(748-0069).....12/box
- #1324-21, 200-mL Bottle
(748-0067).....12/box
- #1324-78, 50-mL Handheld Ultraject Prefilled Syringe
(748-0125).....20/box
- #1324-81, 125-mL Ultraject Prefilled Power Injector Syringe
(748-0044).....20/box



Optiray™ 300 Injectable Contrast Agent

Optiray™ 300 contrast agent is a non-ionic, low osmolar, lower viscosity injectable and intended for intravascular administration. It is a prescription drug that is intended to be therapeutically and biologically inert when injected into the body for use in organ or tissue enhancement in CT, X-ray and fluoroscopy imaging procedures for which it is approved. Each mL of Optiray™ 300 contrast agent contains 636 mg of ioversol, 3.6 mg of tromethamine as a buffer and 0.2 mg of edetate calcium disodium as a stabilizer.

- #1332-06, 50-mL Bottle
(748-0054).....25/case
- #1332-11, 100-mL Bottle
(748-0020).....12/case
- #1332-16, 150-mL Bottle
(748-0026).....12/case
- #1332-21, 200-mL Bottle
(748-0077).....12/box
- #1332-78, 50-mL Ultraject™ Prefilled Handheld Syringe
(748-0126).....20/box
- #1332-83, 100-mL Ultraject™ Prefilled Power Injector
(748-0011).....20/case
- #1332-83, 100-mL Ultraject™ Prefilled Power Injector
(748-0011).....20/case
- #1332-00, 100-mL Ultraject™ Prefilled RFID Power Injector
(748-0071).....20/box
- #1332-61, 500-mL Bottle, Multidose Pharmacy Bulk Package
(748-0031).....6/box



Optiray™ 320 Injectables

Nonionic, low-osmolar contrast media agent intended for intravascular administration. Therapeutically and biologically inert when injected into body for use in organ or tissue enhancement in CT, X-ray, and fluoroscopy imaging procedures. Indicated for use in adults for angiography throughout the cardiovascular system. Uses include cerebral, coronary, peripheral, visceral, and renal arteriography, venography, aortography, left ventriculography, CT of head and body, and pediatric excretory urography. Not for intrathecal use. Each mL contains 678 mg of ioversol, 3.6 mg of tromethamine as a buffer, and 0.2 mg of edetate calcium disodium as a stabilizer.

- #1323-06, 50-mL Vial
(748-0029).....25/case
- #1323-75, 50-mL Ultraject™ Prefilled Handheld Syringe
(748-0124).....20/case
- #1323-77, 50-mL Prefilled Syringe
(748-0082).....20/case
- #1323-83, 100-mL Ultraject™ Prefilled Power Injector Syringe
(748-0119).....20/case
- #1323-11, 100-mL Bottle
(748-0012).....12/case
- #1323-81, Ultraject™ Prefilled Power Injector Syringe, 125 mL
(748-0009).....20/case
- #1323-16, 150-mL Bottle
(748-0033).....12/case
- #1323-00, 100-mL Ultraject™ Prefilled Power Injector Syringe
(748-0076).....20/case
- #1323-02, 20-mL Vial
(748-0079).....25/box
- #1323-04, 30-mL Vial
(748-0058).....25/box
- #1323-21, 200-mL Bottle
(748-0038).....12/box
- #1323-27, 125-mL Ultraject™ Prefilled Power Injector Syringe
(748-0045).....20/box
- #1323-41, 75-mL Bottle
(748-0084).....12/box
- #1323-61, 500-mL Bulk Bottle
(748-0030).....6/box
- #1323-85, 75-mL Ultraject™ Prefilled RFID Power Injector Syringe
(748-0064).....20/box
- #1323-91, 75-mL Ultraject™ Prefilled Power Injector Syringe
(748-0025).....20/box

Optiray™ 350 Injectables

Nonionic, low-osmolar contrast media agent intended for intravascular administration. Therapeutically and biologically inert when injected into body for use in organ or tissue enhancement in CT, X-ray, and fluoroscopy imaging procedures. Indicated for use in adults for angiography throughout the cardiovascular system. Uses include cerebral, coronary, peripheral, visceral, and renal arteriography, venography, aortography, left ventriculography, CT of head and body, and pediatric excretory urography. Not for intrathecal use. Each mL contains 741 mg of ioversol, 3.6 mg of tromethamine as a buffer, and 0.2 mg of edetate calcium disodium as a stabilizer.

- #1333-06, 50-mL Bottle
(748-0041).....25/box
- #1333-41, 75-mL Bottle
(748-0050).....12/box
- #1333-11, 100-mL Bottle
(748-0017).....12/box
- #1333-16, 150-mL Bottle
(748-0024).....12/box
- #1333-21, 200-mL Bottle
(748-0035).....12/box
- #1333-51, 250-mL Multidose Pharmacy Bulk Package Bottle
(748-0022).....12/box
- #1333-61, 500-mL Pharmacy Bulk Package Bottle
(748-0016).....6/box
- #1333-78, 50-mL Ultraject Prefilled RFID Power Injector Syringe
(748-0127).....20/box
- #1333-91, 75-mL Ultraject Prefilled Power Injector Syringe
(748-0019).....20/box
- #1333-83, 100-mL Ultraject Prefilled Power Injector Syringe
(748-0010).....20/box
- #1333-81, 125-mL Ultraject Prefilled Power Injector Syringe
(748-0015).....20/box
- #1333-77, 50-mL Syringe
(748-0059).....20/box
- #1333-00, 100-mL Ultraject Prefilled RFID Power Injector Syringe
(748-0042).....20/box
- #1333-85, 7-mL Ultraject Prefilled RFID Power Injector Syringe
(748-0040).....20/box



Omnipaque™ Contrast Media

A low osmolar, nonionic, iodinated contrast agent. Its indications include a broad range of intravascular diagnostic procedures such as a coronary angiography, spinal cord imaging, and body cavity procedures including shoulder and knee joints. Omnipaque is approved for use in adults and children, and is available in a variety of packagings.

- #Y-101, 180 mg, 10-mL Vial (638-9526).....10/box
- #Y-102, 180 mg, 20-mL Vial (393-2273).....10/box
- #Y-203, 240 mg, 10-mL Vial (638-2431).....10/box
- #Y-220, 240 mg, 20-mL Vial (600-4244).....10/box
- #Y-250, 240 mg, 50-mL Vial (106-1173).....10/box
- #Y-526, 240 mg, 200-mL Polymer Bottle (953-1539).....10/box
- #Y-313, 300 mg, 100-mL Bottle (737-6509).....10/box
- #Y-306, 300 mg, 10-mL Vial (931-4685).....10/box
- #Y-308, 300 mg, 50-mL Vial (601-2164).....10/box
- #Y-352, 300 mg, 50-mL Bottle (248-6022).....10/box
- #Y-530, 300 mg, 50 mL (108-5540).....10/box
- #Y-318, 300 mg, 150-mL Bottle (695-8952).....10/box
- #Y-540, 350 mg, 50 mL (108-5551).....10/box
- #Y-413, 350 mg, 100-mL Bottle (605-0199).....10/box
- #Y-524, 240 mg, 150-mL Polymer Bottle (106-2767).....10/box
- #Y-526, 240 mg, 200-mL Polymer Bottle (953-1539).....10/box
- #Y-531, 300 mg, 75-mL Polymer Bottle (114-8309).....10/box
- #Y-532, 300 mg, 100-mL Polymer Bottle (114-9045).....10/box
- #Y-534, 300 mg, 150-mL Polymer Bottle (114-9048).....10/box
- #Y-542, 300 mg, 100-mL Polymer Bottle (112-4076).....10/box
- #Y-415, 350 mg, 150-mL Bottle (113-6186).....10/box

COVIDIEN

positive results for life™

Conray Contrast Media

Ionic, high-osmolar contrast media agent that is therapeutically and biologically inert when injected into body for use in organ or tissue enhancement in CT, X-ray, and fluoroscopy imaging procedures. Indicated for use in excretory urography, cerebral angiography, peripheral arteriography, venography, arthrography, direct cholangiography, endoscopic retrograde cholangiopancreatography (ERCP), CT of head and body, IV-DSA, and IA-DSA.

- Conray™ 43 #318309, 250-mL Bottle (341-7326).....12/case
- Conray™ 30 #095211, 150-mL Vial (748-0078).....12/case
- Conray™ 43 #318315, 50-mL Vial (748-0047).....50/box
- #318309, 250-mL Vial (748-0075).....12/box
- Conray™ 60 #095313, 30-mL Vial (748-0056).....50/case
- #095315, 50-mL Vial (748-0032).....50/case
- #095309, 100-mL Bottle (748-0065).....12/case
- #095311, 150-mL Bottle (748-0062).....12/case
- Cysto-Conway II #086207, 250-mL Bottle (748-0023).....12/case

COVIDIEN

positive results for life™

Optimark

- #1177-04, 10-mL Vial (748-0053).....10/case
- #1177-06, 15-mL Vial (748-0037).....10/case
- #1177-02, 5-mL Vial (748-0080).....10/case
- #1177-08, 20-mL Vial (748-0021).....10/case
- #1177-50, 50-mL Vial (748-0043).....5/case
- #1177-11, 10-mL Prefilled Syringe (748-0120).....10/box
- #1177-16, 15-mL Prefilled Syringe (748-0122).....10/box
- #1177-21, 20-mL Prefilled Syringe (748-0121).....10/box
- #1177-31, 30-mL Prefilled Syringe (748-0123).....10/box

COVIDIEN

positive results for life™

Hexabrix™ 320 Contrast Agent

Each milliliter of this low osmolar ionic dimer contains 393 mg of ioxaglate meglumine, 196 mg of ioxaglate sodium and 0.10 mg of edetate calcium disodium as a stabilizer. The solution contains 3.48 mg (0.15mEq) sodium in each milliliter and provides 32% (320 mg/mL) organically bound iodine.

- #5505-51, 20-mL Vial (748-0073).....10/box
- #5505-06, 50-mL Vial (748-0063).....25/case



GE imagination at work

Visipaque Contrast Media

Of all the iodine concentrations, this is the only iso-osmolar contrast medium available for intravascular use. It is indicated for Contrast Enhanced Computer Tomography (CECT) imaging of the head and body.

- Visipaque 270**
270 mg/mL
#V020, 50-mL Glass Vial (118-3021).....10/box
- #V026, 150-mL Bottle (638-0002).....10/box
- #V552, 100-mL Bottle (106-3773).....10/box
- Visipaque 320**
320 mg/mL
#V041, 50-mL Hanging Bottle (225-6981).....10/box
- #V045, 100-mL Bottle (854-3243).....10/box
- #V047, 200-mL Bottle (116-6190).....10/box
- 150-mL PlusPak (118-2491).....10/box
- #V568B, 500-mL PluPak Bulk (116-6175).....10/box



GE imagination at work

Omniscan (Gadodiamide) Injection

Omniscan Injection is the formulation of the gadolinium complex of diethylenetriamine pentaacetic acid bismethylamide, an injectable, nonionic extracellular enhancing agent for magnetic resonance imaging, administered by intravenous injection. Provided as a sterile, clear to slightly yellow, aqueous solution. Each mL contains 287 mg gadodiamide, 12 mg caldiumide sodium, and water for injection.

- #J-068, 5-mL Vial (115-9513).....10/box
- #J-100, 10-mL Vial (104-8590).....10/box
- #J-120, 15-mL Fill in 20-mL Vial (600-8007).....10/box
- #J-140, 20-mL Vial (381-9006).....10/box
- #J-160, 10mL Prefilled Syringe (638-5345).....10/box
- #J-170, 15-mL Prefilled Syringe (673-9006).....10/box
- #J-180, 20-mL Prefilled Syringe (993-6306).....10/box
- #J-800, 100-mL PlusPak Polymer Bottle (114-0840).....10/box

Our Portfolio of Practice Solutions

Designed to support your business goals and help you deliver the highest quality care possible to your patients.

HENRY SCHEIN DxRx SOLUTIONS

Questions?
Call Our DxRx Solutions Hotline

The Henry Schein DxRx Solutions Hotline provides information on the following topics for diagnostics, pharmaceuticals, and vaccines:

- ACIP recommendations for vaccinations
- CPT codes and J codes†
- Brand-to-generic comparisons
- Indications and usage as provided on a package insert
- Package inserts
- State law requirements for school and day care entry*
- Latest updates on supply issues and much more

†For informational purposes only.

Source: Reimbursement Codes.com and/or www.cms.gov <<http://www.cms.gov>>. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

*For informational purposes only; does not constitute legal advice.

Source: www.immunize.org

DxRxSolutions
DIAGNOSE • TREAT • PREVENT

1.877.523.SHOT 8:30am–5:00pm, et
1.877.523.7468
rxsolutions@henryschein.com

HENRY SCHEIN BRAND PRODUCTS

Rely On Us
For Quality, Value, and Selection

Health care professionals around the world use Henry Schein Brand products to fulfill their practice supply needs. There are over 6,000 quality products to choose from—each bearing the Henry Schein Seal of Excellence—your guarantee of satisfaction.

Henry Schein Brand products:

- Save on the products and supplies you use most
- Selected for high quality and high value
- Categories for disposables, instruments, diagnostics, and more

HENRY SCHEIN
BRAND PRODUCTS

1.800.PSCHEIN
www.henryschein.com/hsibrand

PRIVILEGES REWARDS

Rewarding Loyalty

By using any of our electronic ordering solutions, **PRIVILEGES** members earn **Rewards Points** (one point for each dollar spent) that can be redeemed for name-brand merchandise, travel, event tickets, equipment, and more!

Plus, PRIVILEGES benefits include:

- Value Certificates to reduce your expenses and offer you the opportunity to earn Bonus Points
- **PRIVILEGES** monthly E-promotions
- Exclusive quarterly E-newsletters
- Access to practice marketing services
- Repair services specials

Privileges EXCLUSIVE
REWARDS FOR
PREFERRED
CLIENTS

To Learn More or to Enroll: 1.866.MED.VIPS
www.henryschein.com/privilegesmd
privilegesmd@henryschein.com

PRACTICE MARKETING

Win More Patients With Experienced
Health Care Consultants

Henry Schein Medical and the most experienced health care marketing consultants are here to help you win more patients. Some of the Healthcare Success Strategies' products and services include:

- Education
- Internet marketing
- Expert guidance
- Branding and advertising

As a preferred provider, we are offering Henry Schein Medical customers up to 15% off our products and services.

HEALTHCARE SUCCESS
STRATEGIES
YOUR RESOURCE FOR BUILDING & SUSTAINING PROFITABLE GROWTH

1.800.656.0907 ext. 4
www.healthcaresuccess.com/henryschein

HENRY SCHEIN FINANCIAL SERVICES

Business Solutions
For Health Care Practitioners

Henry Schein Financial Services specializes in strategies for success. We offer financial resources that can be tailored to your practice goals.

Contact us today to review our full array of financial resources developed to move your business forward with the competitive advantage you deserve.

- Equipment financing
- Start-up financing
- Demographic site analysis reports
- Henry Schein credit card
- Credit card acceptance
- Patient collections

HENRY SCHEIN
FINANCIAL SERVICES
1.800.443.2756
hsfs@henryschein.com
www.henryschein.com/hsfs

HENRY SCHEIN RX SAMPLES SERVICE

Our Rx Samples Service
Is Free, Secure, and Easy

Eligible Customers can order pharmaceutical samples online—easily and conveniently.

Sponsored by leading pharmaceutical companies—with no cost to you.

Sampling enables you to:

- Initiate patient treatment faster
- Increase patient compliance
- Access leading brands not currently available to you today
- Increase patient satisfaction; patients appreciate the convenience, ease, and thoughtfulness of samples

HENRY SCHEIN
SAMPLES SERVICE **Rx**

1.800.PSCHEIN
www.henryschein.com/rxsamples

The preferred colonic transit diagnostic test

Reliably and economically diagnose GI conditions such as chronic constipation, colonic inertia, hypomotility, outlet delay, and bowel obstruction.

Two Easy Testing Methods:

The **Simple Sitzmarks Method** helps detect a blockage or colonic inertia in chronically constipated patients by tracking O-rings in the colon.

The **Segmental Sitzmarks Method** uses all rings to help determine the amount of time that it takes food to pass through each section of the colon.



Available in three shapes:



SITZMARKS® RADIOPAQUE MARKERS Colonic Transit Time Test

available through HENRY SCHEIN®

From Konsyl Pharmaceuticals – manufacturer of the #1 doctor-recommended, all-natural psyllium fiber.



Sitzmarks® Radiopaque Markers

A simple, reliable, and affordable method of diagnosis for many gastrointestinal (GI) conditions, including chronic constipation, colonic inertia, hypomotility, outlet delay, and bowel obstruction. This unique colonic transit diagnostic test helps professionals select the best therapeutic option for adult patients with severe constipation who have otherwise negative GI evaluations.

Radiopaque Markers.....10/box
Specify:
#8100, O-Ring.....(330-2596)
#8100-24TC, Tri-Chamber..(111-6709)
#8100-24DD, Double D.....(861-9772)



Our focus is your patient's **COMFORT**

Surshield® Safety Winged Infusion Sets

- Simple, one-handed activation
- Minimal steps for activation
- No change in needle insertion technique required
- Confirmation that needle is securely locked

Surshield Safety Winged Infusion Sets have the soft, flexible, interlocking wings that readily conform to body's contour. The needles are bonded to the wings at just under a 3° downward angle, for less risk of vessel damage or occlusion of the lumen. The memory free tubing reduces recoil during use. These products have no components made of natural rubber latex and are DEHP free.

Surshield® Safety Winged Infusion Sets

Product Code	Item Code	Description
SV*S21BL	315-8301	21G x 3/4" needle, 12" tubing 50/Box
SV*S23BI	315-3535	23G x 3/4" needle, 12" tubing 100/Box
SV*S19BLS	315-6492	19G x 3/4" needle, 3 1/2" tubing Each
SV*S21BLS	315-7074	21G x 3/4" needle, 3 1/2" tubing Each
SV*S23BLS	315-1046	23G x 3/4" needle, 3 1/2" tubing 100/Box
SV*S25BLS	315-2073	25G x 3/4" needle, 3 1/2" tubing 100/Box



® and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURSHIELD is a registered trademark of Terumo Corporation, Tokyo, Japan. © 2012 Accession # 04272012-2293



TERUMO MEDICAL PRODUCTS



SURFLO® Winged Infusion Sets

These soft, flexible, interlocking wings readily conform to the body's contour. The wing is angled for easy positioning of the needle during venipuncture, for less risk of vessel damage or occlusion of the lumen. The precisely honed needle with Ultra-Thin wall design ensures ease of penetration and continued access with a single insertion. DEHP-free.

3 1/2" Tubing

#SV*25BLS, 25 ga x 3/4" Thin Wall Needle (315-2073).....100/box

#SV*21BLS, 21 ga x 3/4" Ultra-Thin Wall Needle (315-7074)ea

3 1/2" Tubing with Intermittent Plug
#SV*23BI, 23 ga x 3/4" Ultra-Thin Wall Needle (315-3535).....100/box

8" Tubing

#SV*27EL, 27 ga x 1/2" Regular Wall Needle (315-5460).....100/box

12" Tubing

#SV-25BLK, 25 ga x 3/4" Thin Wall Needle (100-3454)50/box

#SV-23BLK, 23 ga x 3/4" Ultra-Thin Wall Needle (100-5102).....50/box

#SV-22BLK, 22 ga x 3/4" Ultra-Thin Wall Needle (315-5572)50/box



#SV21-BLK, 21 ga x 3/4" Ultra-Thin Wall Needle (100-2145).....50/box

#SV-19BLK, 19 ga x 3/4", Thin Wall Needle (100-0269)50/box

#SV-18BLK, 18 ga x 3/4" Thin Wall Needle (315-0465)50/box



Surflo® IV Catheters

Consistently safe and simple venipuncture. Eliminates worries about leakage or separation of IV catheter and hub. Flashback chamber in clear, sure-grip hub gives immediate indication of proper placement. Flexible, thin-wall catheter with large inner diameter ensures good blood flow.

Catheters.....Ea.
Specify:

#SR-OX1451CA, 14-ga x 2" Long Catheter, 16-ga Needle, Orange (315-5809)

#SR-OX1632CA, 16-ga x 1 1/4" Long Catheter, 16-ga Needle, Gray (315-5452)

#SR-OX1651CA, 16-ga x 2" Long Catheter, 18-ga Needle, Gray (315-3043)

#SR-OX1851C, 18-ga x 2" Long Catheter, 20-ga Needle, Green (315-8073)

#SR-OX1832CA, 18-ga x 1 1/4" Long Catheter, 20-ga Needle, Green (315-5131)

#SR-OX2025CA, 20-ga x 1" Long Catheter, 22-ga Needle, Pink (315-2314)

#SR-OX2032CA, 20-ga x 1 1/4" Long Catheter, 22-ga Needle, Pink (315-9704)



#SR-OX2051C, 20-ga x 2" Long Catheter, 22-ga Needle, Pink (315-8751)

#SR-OX2225CA, 22-ga x 1" Long Catheter, 24-ga Needle, Blue (315-9675)

#SR-OX2419CA, 24-ga x 3/4" Long Catheter, 27-ga Needle, Yellow (315-8857)



TERUMO

Terumo® Hypodermic Syringes without Needle

Standard, Luer-Lock, Luer-Slip Tips, and Eccentric Luer-Slip Tips

- 3cc**
#SS-03L, 3cc Luer-Lock-Tip Syringe without Needle, Graduations: 0.1cc (3cc Total) (100-5094)100/box
- 5cc**
#SS-05L, 5cc Luer-Lock-Tip Syringe without Needle, Graduations: 0.2cc (6cc Total) (100-7847)100/box
- 10cc**
#SS-10L, 10cc Luer-Lock-Tip Syringe without Needle, Graduations: 0.2cc (12cc Total) (100-6667)100/box
- 20cc**
#SS-20L2, 20cc/mL Luer-Lock Tip, 1cc Graduations (25cc Total), Latex-Free (315-3268)50/box

- 30cc**
#SS-30L, 30cc Luer-Lock-Tip Syringe without Needle, Graduations: 1.0cc (35cc Total) (100-8334)25/box
- 60cc**
#SS-60L, 60cc Luer-Lock-Tip Syringe without Needle, Graduations: 1.0cc (60cc Total) (100-5188)25/box
- 10cc**
#SS-10S, 10cc Slip-Tip Syringe without Needle, Graduations: 0.2cc (12cc Total) (315-1696)100/box
- 20cc**
#SS20S, 20cc Luer-Slip-Tip Syringe without Needle, Graduations: 1.0cc (25cc Total) (371-3810)200/case
- 30cc**
#SS-30S, 30cc Luer-Slip-Tip Syringe without Needle, Graduations: 1.0cc (35cc Total) (315-2891)25/box
- 60cc**
#S-60C, 60cc (2 oz.) Catheter-Tip Syringe without Needle, Graduations: 1.0cc (60cc Total) (315-0045)25/box

B|BRAUN

Macro Drip IV Administration Set with 1 Injection Site-73"

Universal spike, injection site 6" above distal end, roller clamp, slide clamp, and SPIN-LOCK® connector. Latex-free. Priming volume: 16 mL, 15 drops/mL.

- #1416 (507-4189)ea



TERUMO

Disposable Hypodermic Syringes with Needle and Luer-Lock Tip

Offer maximum patient comfort through ultimate sharpness and smoothness. Latex-free. Clear syringe barrel with bold, precise scale markings. Translucent needle hub for quick observation of flashback.

- 3cc**
#SS-03L2325, 23 ga x 1" (100-6223)100/box
- #SS-03L2225, 22 ga x 1" (100-1798)100/box
- #SS-03L2125, 21 ga x 1" (100-1378)100/box
- #SS-03L2025, 20 ga x 1" (100-0707)100/box
- 5cc**
#SS-05L2125, 21 ga x 1" (100-6872)100/box
- #SS-05L2025, 20 ga x 1" (100-4808)100/box
- 10cc**
#SS-10L2125, 21 ga x 1" (100-0358)100/box
- #SS-10L2025, 20 ga x 1" (100-8831)100/box



BD™ Vacutainer® Safety-Lok™ Blood Collection Sets

Safety-engineered winged sets for both blood collection and infusion. Each set features a translucent, integrated, protective shield that allows for one-handed activation immediately after use to help minimize the risk of needlestick injuries and to offer clear visibility of blood flashback.

- #367292, 7" Tube with Luer, 23 ga x 3/4" (987-0529)ea
- #367285, 12" Tube with Luer, 25 ga (840-3447)ea
- #367283, 12" Tube with Luer, 23 ga (987-1301)50/box
- #367298, 12" Tube without Luer, 25 ga (623-4670)50/box
- #367297, 12" Tube without Luer, 23 ga (987-4692)50/box



BD™ Syringe/Needle Combination
BD Luer-Lok™ Tip, Sterile, Single Use

- 3 mL**
#309571, 23 ga x 1" (987-7076)100/box
- #309572, 22 ga x 1" (987-2558)100/box
- #309574, 22 ga x 1 1/2" (987-8154)100/box
- #309575, 21 ga x 1" (987-5994)100/box
- 5 mL**
#309630, 22 ga x 1" (987-0441)100/box
- #309631, 22 ga x 1 1/2" (987-1405)100/box
- #309632, 21 ga x 1" (987-4478)100/box
- #309633, 21 ga x 1 1/2" (987-4275)100/box
- #309634, 20 ga x 1" (987-3628)100/box
- 10 cc**
#309644, 20 ga x 1" (987-1271)100/box



B|BRAUN

Introcane® Safety™ Catheters

Passive design minimizes needlestick injuries without change in technique. Requires minimal in-servicing and a short learning curve. Eliminates risk of inadvertent activation. Compliance virtually ensured since safety mechanism cannot be bypassed.

- 24 ga x 3/4", Straight (507-0341)ea



- 24 ga x 3/4", Wing (507-2343)50/box
- 22 ga x 1", Straight (507-5531)ea
- 22 ga x 1", Wing (507-7254)50/box
- 20 ga x 1", Straight (507-5995)ea
- 20 ga x 1", Straight (507-5328)50/box
- 20 ga x 1", Wing (507-4839)50/box

Polyurethane (PUR)

- 24 ga x 3/4", Straight (507-6928)50/box
- 24 ga x 3/4", Wing (507-0641)ea
- 22 ga x 1", Straight (507-7701)ea
- 20 ga x 1 1/4", Straight (507-7274)50/box
- 18 ga x 1 1/4", Straight (507-3345)50/box



BD Insyte™ Autoguard™ Shielded IV Catheters

- #381423, Blue, 22 ga x 1" (0.9 mm x 25 mm)ea
(987-5015)
- #381523, Blue, 22 ga x 1" (0.9 mm x 25 mm), Notched Needle **(987-5417)**50/box
- #381412, Yellow, 24 ga x 0.75" (0.7 mm x 19 mm)ea
(840-3491)
- #381433, Pink, 20 ga x 1" (1.1 mm x 25 mm)200/case
(987-6926)
- #381434, Pink, 20 ga x 1.16" (1.1 mm x 30 mm)50/case
(987-4421)
- #381444, Green, 18 ga x 1.16" (1.3 mm x 30 mm)ea
(840-3495)



BD™ Saf-T-Intima™ Closed IV Catheter Systems

- Y Adapter and Needle Shield
Combines needlestick safety with a closed-system design. The result is a virtually bloodless procedure and greater safety for practitioner, all from the first and only closed IV catheter system available for vascular access. Made from clinically proven BD Vialon™ biomaterial with wings.
- Systems25/box
- Specify:
- #383312, with Removable PRN Adapter, Yellow, 24 ga x 3/4" (0.7 mm x 19 mm).....**(987-3762)**
 - #383322, with Removable PRN Adapter, Blue, 22 ga x 3/4" (0.9 mm x 19 mm).....**(987-4887)**
 - #383336, with Y-Adapter, Pink, 20 ga x 1" (1.1 mm x 25 mm)**(119-3024)**
 - #383346, with Y-Adapter, Green, 18 ga x 1" (1.3 mm x 25 mm)**(987-9844)**
 - #383313, "Y" Adapter & Needle Shield, 24 ga x 3/4" (0.7 mm x 19 mm)ea
(987-7203)



Baxter

0.9% Sodium Chloride Irrigation, USP

- IV Solutions, in Plastic Bags
- #2B1306, 50 mL **(153-3897)**ea
 - #2B1307, in VIAFLEX Plastic, Single Pack, 100 mL **(153-0110)**ea
 - #2B1321, 150 mL **(153-8937)**ea
 - #2B1322Q, 250 mL **(153-7468)**ea
 - #2B1323Q, 500 mL **(153-7162)**ea
 - #2B1301, in VIAFLEX Plastic, Quad Pack, 50 mL **(153-3826)**4/pkg
 - #2B1308, in VIAFLEX Plastic, 16-Pack, 50 mL **(153-6129)**96/case
Contains: 4 packs of 16.
 - #2B1302, in VIAFLEX Plastic, Quad Pack, 100 mL **(153-7916)**4/pkg
 - #2B1309, in VIAFLEX Plastic, 16-Pack, 100 mL **(153-0731)**16/box



positive results for life™

Contrast Media Administration Sets

Fluid administration set with large bore is appropriate for cystography. Includes Christmas tree-style adapter for connection to multiple types of Foley catheters. Stopcock 3-way I.V. set allows for handheld bolus injection and drip infusion through a single I.V. hookup. Luer locknuts are designed specifically for Ultraject™ prefilled contrast media syringes for power injection.

- #0144-03, Fluid Administration Set **(116-3916)**12/box
- #2657-03, IV Administration Set **(145-0570)**12/box



BD™ Nexiva™ Closed IV Catheter System

- Closed IV catheter system for peripheral access. Reduces life-threatening exposure while enhancing IV therapy success and clinical outcomes.
- Features and benefits:
- No blood exposure or blood mess
 - Built-in stabilization platform to secure the catheter
 - BD Vialon™ biomaterial: proven catheter material to reduce IV therapy complications
 - BD Instaflash™ needle technology designed to improve first-stick success
 - Includes BD Q-syte, proven to reduce catheter blood-stream infections
 - All integrated parts simplify the IV therapy process
 - Reduces contamination at the IV site
 - Larger bore extension tubing for higher flow rates
- #383531, 24 ga x 0.75", Dual Port **(987-5113)**20/box
 - #383532, 22 ga x 1", Dual Port **(987-5112)**20/box
 - #383536, 20 ga x 1" HF Dual Port **(111-6357)**80/case
 - #383537, 20 ga x 1.25" HF Dual Port **(111-6711)**80/case



B. BRAUN

0.9% Sodium Chloride IV Solutions

- Excel® DEHP-Free, PVC-Free, and Latex-Free Containers
- #L8002, 250-mL Bag **(507-8002)**ea
 - #L8001, 500-mL Bag **(100-4737)**ea
 - #L8000, 1000-mL Bag **(100-2808)**ea
 - #S8004-5384, 50-mL Fill/100 mL **(507-4046)**ea
 - #S8004-5264, 100-mL Fill/150 mL **(507-2187)**ea



positive results for life™

Ready-Box™ Contrast Media Warmer

Made of lightweight aluminum. The transparent plastic front door, compartmentalized shelves and partitions allow easy access, visibility and flexible storage of contrast media boxes, vials, and bottles. The key lock limits access to help ensure security of the contents, and the open coil heater aids even heat distribution.

- #1550CW, Ready-Box **(275-1107)**ea



Barium Paper Cups

- 14-ounce paper cups feature printed mixing scales of 0 to 14 oz and 0 to 40cc.
- #15601, 14 oz/40cc100/pkg

Integrity™ 1000 Series

RADIOGRAPHIC SYSTEMS



Call Henry Schein for all of your
Radiographic Room Equipment Needs! We Have it ALL!

A 5-year parts warranty assures long-term value.

Simple...



KONICA MINOLTA

Nothing Makes Your Transition To Digital Easier

ALL-IN-ONE DESIGN

ImagePilot exemplifies the true meaning of an integrated digital radiography system. It combines patient registration, CR acquisition, distributed workflow, image review, and long-term storage in one easy-to-use and maintain system. These advanced features meet private practice workflow needs and eliminate the need for a separate PACS.

Standard ImagePilot Features:

- All-In-One Design—Easiest for the staff to learn. Reduces IT support requirements.
- AutoPilot—Automatic image processing for unsurpassed consistency in image quality.
- Simplified Image Acquisition—Single click workflow increases productivity.
- Same User Interface At Every Workstation—learn only one system.
- Easy-To-Use And Feature Rich Viewing Screen—Image adjustments, annotations, measurements, inputs and outputs, customizable for basic to advanced users.

ImagePilot



ImagePilot Sigma Highlights:

- REGIUS Sigma Tabletop CR Reader
 - 25 or 45 cassettes per hour
 - Flexible plate technology for compact design and cost-effectiveness
- Small 20"W x 24"D x 14"H
- Lightweight—under 70 pounds
- Reliable—very few replacement parts
- Touch Screen capable—easy operation
- Informity ready!



ImagePilot LS Highlights:

- Rapid Output Speed—40 cassettes per hour.
- Redundant Storage—mirrored 500GB hard drives (RAID 1)
- Distribute Workflow with Multiple Users—Two additional licenses (3 total including the server).
- Import Priors or Export Patient Images to CD—Better than film for transportability and viewing.



Viztek Introduces ViZion Straight-Arm DR

- Generator:**
SHF50, 50 kW, 600 mA HF-AP, 3-phase with AEC - Other generator sizes available
- Tube Stand:**
Straight -Arm APR techniques with ViZion 17x17 CsI detector, three field AEC ion chamber, and manual SID display
- Tube:**
Toshiba, E7884FX, 300kHU tube, 150KVP, 0.6 x 1.2mm FS
- Collimator:**
Manual collimator
- Table:**
Standard mobile table
- ViZion U-Arm DR Package:**
- 17x17 Cesium panel
 - Dell 3500 Computer with monitor
 - Opal-Acquire software
 - Modality print, store, and worklist
 - Comprehensive on-site applications

*Weight bearing stand and manual stitching available as options




- Complete Straight-Arm System
- Acquire Software/Computer Hardware + Monitor
- Standard Mobile Table
- On-Site Applications



Sensible Solutions for Refurbished Radiology



OEC 9600 C-Arm




OEC 9800 C-Arm



GE OEC 9800 Super C

**Before you make your selection,
Call your Henry Schein
Sales Consultant
for more information!**



Your Practice. Your Preference.
Your Productivity Solution.

PACS SOLUTIONS FOR EVERY PRACTICE

Feature-Rich PACS– Affordable Solution for Your Practice

Enjoy the benefits of the web-based MXReview PACS paired with the modality of your choice. MXReview PACS solution provides your practice with the technology to optimize workflow at an affordable price.

The scalable software allows users to have the sophisticated tool-sets of an enterprise-level while also providing a sub-set of optional features for a customized fit. The software includes tools for technologists QC as well as a full featured web viewer.

Workflow Your Way

MXReview PACS caters to practices of all sizes. Avoid overspending by customizing your solution's features and functionalities based on your facilities size, needs, and budget.

MXReview PACS solution follows medical standards for DICOM, IHE, and HL7. The PACS seamlessly integrates with any third party vendor RIS, EMR, and Practice Management software.

No matter what the size or type of facility, MXReview can be hand tailored for the most cost effective, workflow driving, business solution.





Easy Wrap Apron

0.5-mm Lead Vinyl
This design is even more comfortable and features 2-sided Velcro® closures so fit can be adjusted more evenly.

Dark Blue
(100-0858)ea



Apronette

0.5-mm Lead Vinyl
Designed to shield the pelvic region. Held by a belt and Velcro® closure and can be worn front or back.

Beige, 18" x 24"
(594-5735)ea



Sensi Flex™ Gloves

The ultimate in sensitivity and flexibility, these gloves fit like cloth gloves. Perfect for Cath/EP labs, Orthopedics, Endoscopy, Urological procedures and a multitude of other radiological procedures. The perfect choice when maximum tactile sensitivity is important.

- Size 7, 0.22 mm (969-4046)1 pair
- Size 7, 0.30 mm (594-9856)1 pair
- Size 7½, 0.22 mm (109-4943)1 pair
- Size 7½, 0.30 mm (594-4253)1 pair
- Size 8, 0.30 mm (594-7711)1 pair
- Size 8½, 0.30 mm (594-4428)1 pair



Thyroid Collar

Protects neck area from radiation during exposure. For neck sizes 9" to 17".

Beige, 0.5 mm
(594-7704)ea



Gonad Shield

0.5-mm Lead Equivalency Protection
Specifically designed to protect the genital area, the gonad shield comes with straps to hold it in place.

(594-6940)ea



X-ray Gloves

0.5-mm Flexible Lead Rubber
Medium

(100-4010)pair



HENRY SCHEIN® BRAND PRODUCTS



X-RAY APRONS

Covers from thyroid to below gonadal area. All-vinyl surface.
0.3-mm lead equivalency.
Approximately 24" x 27".

- Multi-ply construction: pliable and comfortable
- Heavy-duty vinyl fronts: clean easily and are durable
- Non-slip foam backing: added support

Adult with Collar
24" x 27"ea
Specify:
Light Blue(100-2528)
Beige/Tan(100-5725)
Mauve.....(101-4125)
Textured Gray.....(101-5015)

Adult without Collar
24" x 25"ea
Specify:
Light Blue(100-0989)
Beige/Tan(100-3152)
Mauve.....(100-7806)
Textured Gray.....(100-4339)

Child with Collar
21" x 22"ea
Specify:
Light Blue(100-9022)
Beige/Tan(100-6216)
Mauve.....(101-9961)
Textured Gray.....(101-1205)

Panoramic Style
Front: 22" x 27"
Back: 18" x 18"ea
Specify:
Light Blue(100-6321)
Beige/Tan(100-0224)

Technician Style
24" x 37"ea
Specify:
Light Blue(100-7058)
Beige/Tan(100-5444)
Mauve.....(100-5372)
Textured Gray.....(100-3250)

Protective Collar
5" x 20"ea
Specify:
Light Blue(100-0013)
Beige/Tan.....(101-2573)
Mauve.....(101-2591)
Textured Gray(101-3752)



DeBroeck Wall Shoulders™
Durable and attractive X-ray garment storage. Extend the life of your X-ray apron.



Model 2120
For Grommet- or Loop-Style
X-ray Apronsea
Specify:
Linen(365-5954)
Sandstone(365-0111)
Eclipse.....(365-6518)



Model 3120
For Panoramic-Poncho-Style
X-ray Apronsea
Specify:
Linen(365-6730)
Sandstone(365-0112)
Eclipse.....(365-2315)



Model 3130
For Full Front-and-Back Panoramic-Style
X-ray Apronsea
Specify:
Linen.....(365-5146)
Sandstone(365-0113)
Eclipse.....(365-1474)



Apron Hanger
(135-2628)ea



HENRY SCHEIN®
X-RAY APRON HANGER
For adult or child cling shield.
Helps prolong the life of X-ray aprons.
• Heavy-duty, white epoxy-coated steel:
easily mounted to door or wall
• Provides safe storage for aprons
when not in use: helps prolong the life
of X-ray aprons
• Dimensions:
8¼"H x 14⅞"W x 1⅜"D
(100-3530).....ea



PALMERO
HEALTH CARE
X-ray Apron Hanger
White-coated stainless hanger for adult or
child cling shield. 2"H x 12"L.
(134-1752)ea




Wolf
Tri-Rak Apron and Glove Rack
Wall-mounted Tri-Rak holds a protective
apron and a pair of protective gloves.
Dimensions: 16"H x 28"W x 8"D
(594-8562)ea



Wolf
Uni-Hanger
Allows you to hang any protective apron on
a wall hook or in a conventional closet.
(594-3076)ea



BUX
CENTAL
Panoramic Apron Hanger
(135-5413)ea



Xact® Extremity Ring Markers

The 15-mm markers are perfect for marking areas of interest on hands, arms, and feet. They create a clear and easy-to-see circular indicator of known 15-mm diameter. The 10-mm markers are perfect for indicating areas of interest in the fingers and toes. Many facilities use them to mark problematic screws and other hardware prior to surgery, and to identify areas of digital pain.

- Designed for X-rays of bone such as in the hands and feet
- Are easy to use; just peel and stick the self-adhering backing to the patient and take multiple view radiographs.
- Used in podiatry, orthopedics, plastic surgery, and emergency departments
- Lead-free and latex-free
- Thin, flexible, and highly conforming
- Hypoallergenic adhesive—they won't fall off
- Convenient and environmentally friendly dispenser box

Extremity Ring Markers
 10 mm
(452-6918)60/box
 15 mm
(452-0000)60/box



DENSPLY
RINN

Mount Marking Instrument

Unique white pencil that writes easily on most film mounts, including paper, plastic, or vinyl.

(189-2477)3/box



Wolf

“PRO-PEL” Marking Pencils

Wax skin pencils for marking film, glass, and metal.

Black
(102-1330)ea
 Red
(803-6416)ea
 Red Marking Pencil Refills
(803-6417)5/box



Patient Positioning Foam Set

- Foam spinal wedge, 7" x 7" x 21.25"
 - Foam oblique finger, 3.13" x 5.5" x 5"
 - Foam circular disk, 7" diameter
 - Foam wedge, 16", set of two, 2.5" x 10"
- (113-6537)**4/pkg



Cervical Sandbag Set

- Ready to use
- Tear-resistant with sewn seams
- Heat-sealed inner sleeve

#15731
(113-6538)2/pkg
 Contains: 2–10-lb 11" x 11" bags.

Digital Triple Timer

Accurate timing for darkroom, operator, or lab. Displays 3 countdown timers at once in hours, minutes, and seconds. Timer, stopwatch, and clock. Automatic memory.

Model#T189B
(796-2889)ea



Count Up/Count Down Timer

Timer and stopwatch, count down from 99:59, count up from 99:59, displays minutes and seconds, automatic memory, large display, loud alarm, magnet on back. Specifications: Color: white, Digit height: 0.5", Display size: 1.9" x 0.83", Unit size: 2.63" x 2" x 1.13", Power source: 1–AAA battery, Weight: 2.25 oz.

Model #T1170
(796-1927)ea



Electronic Timer

Pocket-sized with battery included. Times from 100 minutes to 1 second; loud tone when time is complete. 1-year warranty.

2⁵/₁₆" x 2⁷/₁₆" x 1/2"
(100-8822)ea



instadose™ Dosimeter

Provides an instant readout when connected to any computer with internet access. Provides unlimited readings for a year at one low price, reassignment of devices online, elimination of the collection and return of badges for processing, confirmation of a safe radiation work environment for your staff and patients, and compliance with state and federal regulations.

Annual Badgeea.
 Contains:
 Specify:
 Black **(116-1424)**
 Blue **(116-1421)**
 Pink **(116-1422)**
 Silver **(116-1423)**
 Green **(114-4896)**



Luxel®+ X-ray Monitoring Badges

Experts agree that there is no safe level of radiation. Even the smallest dose may cause cancer and genetic damage. Any health care facility may be exposing staff to unsafe levels of radiation. That's why accurately measuring and recording radiation exposure over the long term is important for your people and your practice. The Landauer Luxel Plus with state-of-the-art Optically Stimulated Luminescence (OSL) technology, has the highest sensitivity available today—10 times the sensitivity of film badges!

- Features include:
- Fastest reporting in the industry—5-day average turnaround
 - A lifetime exposure record is archived for safekeeping
 - Personalized badges
 - Web-based reporting
 - Automatic renewal program
- Monthly
(683-8079)ea
 Quarterly
(683-2339)ea
 Fetal Badge
 (7-Month)
(683-4113)ea






Adjustable Mobile Shield

Designed to protect the patient against scattered radiation, this shield is ideal for chest, skull, and spine radiography. The screen is made of lead vinyl with a .25-mm lead protection and is attractively covered on both sides by easy-to-maintain vinyl. It is 24" (60cm) wide and is adjustable with a minimum height of 24" (60cm) and a maximum height of 48" (122cm). It's mounted on an all-welded, durable chrome-finished frame and rolls easily on 4 conductive casters.

Specifications: 48"H (Open) [30"H (Closed)] x 26¾"W x 18"D (at Base)
(201-4707)ea




**HENRY SCHEIN®
X-RAY ENVELOPES**

- Printed (104-2559)500/box
- Plain, Unprinted (104-9771)500/box
- 8" x 10" (100-8383)100/box



**HENRY SCHEIN®
FLOATING X-RAY THERMOMETER**

(100-8280)ea



X-ray ID Systems

- Printer (100-2055)ea
- Printer, Automatic with Adjustable Exposure Timer (594-3870)ea
- Printer Cards (100-4790)1000/box



L and R Marker Clips

Fast, convenient method for marking X-ray films. Clips over cassette frame and reveals L or R during patient's exposure.

(100-4294)ea



Swivel Safelight

Perfect for direct spotlighting, this wall-mountable unit swivels to any position. Includes a GBX-type filter, lightbulb, and 6' cord with 3-prong plug. 6-month warranty.

12"H x 6"W
(100-8206)ea



Label/X-ray Tape

¾" Width

Fast, simple way to mark X rays. Pull the length of tape you want, print patient data on tape, remove adhesive backing, place on cassette, and expose.

25' Roll (100-2783)ea

50' Roll (334-0004)ea

100' Roll (100-2069)ea

Precut 3" Strips (110-7137)300/box

MultiPack Precut 3" Strips (268-7358)100/box



Thickness Caliper

Calibrated in inches and centimeters, this aluminum caliper is perfect for measuring areas to be X-rayed. It will provide years of accurate use.

(594-9442)ea



Film Bin

Lightproof film storage bin in a tabletop format. Bin can also be mounted on the side of a desk. Holds up to 150 14" x 17" films. Automatically closes the door after the films are removed.


19"L x 15½"W x 11½"H
(594-0763)ea



HENRY SCHEIN®

X-RAY FILM—BLUE SENSITIVE

A fine-grain, high-contrast, blue-sensitive, all-purpose film, designed for 90-second rapid processing; adaptable to any processing cycle.



- 8" x 10" (900-4282)100/box
- 10" x 12" (900-4283)100/box
- 11" x 14" (900-4291)100/box
- 14" x 17" (900-4284)100/box

HENRY SCHEIN®

X-RAY FILM—BLUE SENSITIVE, HALF SPEED

- 10" x 12" (900-4289)100/box
- 14" x 17" (900-4288)100/box



HENRY SCHEIN®

X-RAY FILM—DUPLICATING

- 10" x 12" (900-4293)100/box

HENRY SCHEIN®

X-RAY FILM—GREEN SENSITIVE

Compatible with all rare-earth systems, this film can be processed in all popular brands of X-ray film developers, both automatic and manual.



- 7" x 17" (900-4292)100/box
- 8" x 10" (900-4285)100/box
- 9⁷/₁₆" x 11³/₁₆" (900-4286)100/box
- 10" x 12" (900-4281)100/box
- 11" x 14" (900-4287)100/box
- 11¹/₁₆" x 13²⁵/₃₂" (900-4290)100/box
- 14" x 17" (900-4280)100/box



Kodak Ektavision® G and L Extraoral Film—Panoramic

Get low-crossover technology that captures crisp radiographic detail with both. Ideal for panoramic and cephalometric radiographs, TMJ cases, implant procedures, and more. Both feature 400-speed film system that minimizes radiation exposure. Ektavision G Extraoral Film: Kodak's sharpest extraoral film to visualize maximum amount of radiographic information for diagnosis. Panoramic radiographs taken with films like this may even assist in detection of calcifications in carotid arteries, a possible predictor of heart attack or stroke. Ektavision L Extraoral Film: low-contrast, wide-latitude film. Ideal tool for enhancing soft tissue visualization for CT and cephalometric radiographs. Process both films automatically or manually, and use current cassettes with Kodak Ektavision or Lanex Regular screens.



- Ektavision G, 5" x 12" (111-1555)50/box
- Ektavision G, 6" x 12" (111-6395)50/box
- Ektavision G, 8" x 10" (111-4508)50/box
- T-Mat G, 5" x 12" (111-2889)50/box
- T-Mat G, 6" x 12" (111-1862)50/box
- T-Mat G, 8" x 10" (111-2257)50/box
- T-Mat H, 5" x 12" (111-2671)50/box
- T-Mat H, 6" x 12" (111-6296)50/box
- T-Mat H, 8" x 10" (111-5758)50/box
- X-Omat DBF, 5" x 12" (111-8466)50/box
- X-Omat DBF, 6" x 12" (111-0278)50/box
- X-Omat DBF, 5" x 7" (111-8135)50/box
- X-Omat DBF, 8" x 10" (111-0020)50/box
- T-Mat L, 8" x 10" (111-5047)50/box

WE SELL
Kodak | Dental Systems

Kodak Ektavision® G and L Extraoral Films—Duplicating

- X-Omat 2
- 1¹/₄" x 1⁵/₈" (111-2968)150/box
- 5" x 12" (111-3147)50/box
- 15cm x 30cm (111-1601)50/box
- 8" x 10" (111-8416)50/box



HENRY SCHEIN®

SURFACE BARRIER—LIGHT HANDLE COVER

Disposable barrier product to aid in the prevention of cross contamination.

- (100-9864) 4"W x 5¹/₂"L500/box
- (100-9860) 2¹/₂"W x 10"L500/box
- (104-7423) 2¹/₂"W x 10"L500/box
- (100-9863) 27¹/₂"W x 24"L225/box



Cover Film—Disposable Sheets with Adhesive Back

- 4" x 6"1200/roll
- Specify:
- Clear (100-9869)
- Blue (101-1838)
- (101-4275) 10¹/₂" x 14"500/box
- (101-6507) 15"W x 26"L250/roll



Radiomat™ SG™ Full Speed Orthochromatic Film

This film provides high contrast and excellent detail throughout the entire diagnostic range and rapid processing. This medium speed, high contrast film is compatible with any green-emitting rare-earth intensifying screen. When used with 400-speed green rare-earth screens, images have high clarity, high contrast, and overall image sharpness. When used with 100-speed green rare-earth screens, the film provides outstanding detail in extremity images. All come in a non-interleaved, reclosable box.

- 8" x 10" (594-1853)100/box
- 10" x 12" (765-0680)100/box
- 11" x 14" (777-5870)100/box



Deluxe Cassettes

Manufactured for the best film contact, die-formed stainless steel frames are hand-welded for strength and durability. Screen not included.



- 8" x 10" (100-3460)ea

Intensifying Screen Cleaner

Used for a safe and convenient way to clean screens. Apply liquid to gauze pad, wipe clean and air dry.

- (111-2615)16 oz



Spray 2000™ Plus

Roller film processor cleaner. Use once a week to remove residue on roller transports. Protects film from roller deposits and prolongs the life of roller transports.

Complete Kit
(698-9587)ea
Contains: 1 liter of Spray 2000™ Plus, sprayer attachment & 4 cleaning webs.



Formula 2000™ Plus

Nontoxic, nonhazardous cleanser. Effectively removes deposits and stains from roller transports, processors, and dip tanks.

(698-5966).....2/pkg



Kodak | Dental Systems

RP X-Omat

#192, Developer Replenisher, 5-gal Unit

(111-7166).....2/pkg

#92, Fixer and Replenisher, 5-gal Unit Automatic processing solution. Provides consistently high image quality of dental film.

(111-9069).....2/pkg

#187, Developer Starter, to Start 8 gal of Developer
(111-8911)24 oz



Kodak | Dental Systems

GBX

Chemical concentrate designed for manual (hand tank) processing of dental X-ray film. Each twin-pack makes 1 gal.

Twin-Pack
(111-3842).....ea

Developer and Replenisher, 28 oz
(111-8408)ea

Developer and Replenisher, 5 gal
(111-6804)ea

Fixer and Replenisher, 28 oz
(111-7864)ea

Fixer and Replenisher, 5 gal
(111-1610)ea



**HENRY SCHEIN®
AUTOMATIC A & B
PROCESSING
SOLUTIONS**

Premixed, ready-to-use automatic processing solutions that are designed for use in all brands of automatic roller-type processors.

- Suitable for all brands of dental X-ray film: cost-effective
- Ready to use: no mixing required; no mess
- Long shelf life: extended use

Developer and Fixer Set
(101-5075)4/case

Contains: 2 gal of ea: developer & fixer.

Developer Only
(101-4405)4/case

Contains: 4 gal.

Fixer Only
(101-7957).....4/case

Contains: 4 gal.



**Air Techniques
Chemistry**

For use with all roller transport processors. For both automatic and manual replenishing.

Developer and Fixer Set
(698-6513)ea

Contains: 2 gal ea. developer & fixer.

Developer Only, 4-gal Case
(698-7446)ea

Fixer Only, 4-gal Case
(698-8081)ea

Starter Solution, 12-oz Bottle
(698-0763)ea



**DENSPLY
RINN**

**RRP Rapid Process
Developer and Fixer**

Long-lasting and quick processing: 30 seconds in developer and 30 seconds in fixer. Provides high-quality results in density, contrast, and archival quality.

Twin-Pack
(189-4910)ea

Contains: 1 qt of ea.: developer & fixer (working-strength solutions).



SDS Southland

**Clear Image
Radiographic
Cleaner**

100% biodegradable and nontoxic. Used weekly, it guarantees sharper images. Safe for all processors. Cleans developer, fixer, oxidation, silver, and algae deposits from rollers, tanks, webs, and transports in 3-5 minutes. Acid-free, chromate-free, and no hazardous shipping charge required.

1 qt with Foamer Spray
(107-4017)ea

1-gal Refill
(107-4189)ea



**HENRY SCHEIN®
FORMULA RTC™**

Use quarterly.
1-liter Bottle
(101-4643).....2/box



Lh Low-hazard Cleaner

One-part concentrate dissolves deposits and stains from roller assemblies and tanks. Environmentally safe. No hazardous duty fee.

Quart
(107-2264)6/case



**Solutek
CORPORATION**

**SDX Concentrate
Developer and Fixer**

These are single-solution concentrates for manual processing of all brands and types of medical X-ray films. The special one-part formulations are easy to mix and eliminate the need for mixing of multiple chemical parts, lessening the possibility of mixing errors.

SDX Developer produces highest speed and contrast, plus low base and fog. Clean-working over a wide range of temperatures, it is environmentally friendly and contains no boron. SDX Fixer contains hardener. It is noncorrosive, and rapidly clears and hardens all X-ray films. Resists formulation of sludges.

1-gal SDX Developer
(237-1702).....ea

1-gal SDX Fixer
(237-9133).....ea



**Solutek
CORPORATION**

**SoluDent Automatic
Developer and Fixer
for Roller-type Processors**

Ready-to-use processing solutions provide consistent results, long shelf life, and exceptional quality. Premixed chemistry prevents solution mixing errors and time-consuming retakes. Assures high quality films with less handling time and fewer problems. Suitable for use in all automatic roller-type dental X-ray film processors.

2-1-gal Bottles of Developer & Fixer
(237-3364)4 gal/case

Contains: 2-1-gal bottles of each: developer & fixer.



**Solutek
CORPORATION**

**SoluDent SDX
Developer and Fixer**

The developer is a one-part liquid concentrate that is packaged in 26-oz bottles. Each bottle makes 1 gallon of working strength developer. The fixer is a one-part liquid concentrate, including hardener, packaged in 26-oz bottles. Each bottle makes 1 gallon of working strength fixer. Both products are used for manual processing of all brands and types of medical and dental X-ray film.

26-oz Bottle SDX Developer
(237-3264)ea

26-oz Bottle SDX Fixer
(237-0771)ea





InSight X-ray Illuminator

This illuminator utilizes LCD technology, providing an X-ray viewbox with a uniquely even light surface in a super-thin cabinet. Just 2" deep, it gives the look of a recessed unit without the hassle of cutting holes in the wall. The panel accommodates one 14" x 17" X-ray. The auto switch activates the viewing area when the films are inserted, and has a manual on/off switch as well.

21"L x 20.25"W x 1.6"D
(106-8228)ea




Econoline Illuminators

18 1/8"H x 3"D

These economical illuminators have a 14" x 17" viewing area per unit, and are 3" thick, making them unobtrusive when hung on walls. They feature 2-bulb illumination (15 W each), 1,713 candela per square meter (bulbs are included with unit), baked-on white polyurethane finish, Magic Grip film retainer, and 8' line cord with 3-prong plug.

Single Unit, 14 1/8"W
(594-2429)ea
2-in-1 Unit, 28 1/4"W
(594-4023)ea
4-in-1 Unit, 56 1/2"W
(594-4457)ea



Trimline Illuminator

Features film-activated "auto switch," which turns illuminator on when X-ray is inserted, as well as standard on/off.

Constructed of steel, finished in baked-on white polyurethane, unit comes standard with 2-lamp (1,713 candela per square meter) illumination (unit comes with bulbs in place), roller grip film retainer, and 8' line cord with hospital-grade plug. Unique to Trimline is Wolf E-"Z"-Bar™, a simple wall-mounting system that lets you install illuminator in minutes.

Viewing area: 14"W x 17"H.
14"W x 22"H x 3 3/4"D
(100-6719)ea



HZ Series Horizontal Illuminator

Perfect for desktop use, this illuminator can also be wall hung. Made of heavy-gauge steel and finished in baked-on white polyurethane, it features two 20-W fluorescent lamps (shipped in place), a Magic Grip film retainer, rocker switch, and 5' line cord with 3-prong plug. 12" x 24" unbreakable Plexiglas® viewing area is perfect for viewing two 10" x 12" or 8" x 10" radiographs together.

13 1/2"H x 22 1/8"W x 4"D
(594-1670)ea



HENRY SCHEIN®

#200 FILM ILLUMINATOR

Stainless steel frame (steel housing)
Magic Grip film retainer. Holds two 18" window mounts or 5" x 7" or 8" x 10" film.
Viewing area: 10" x 13". Overall dimensions: 15 1/4" x 11 1/2" x 2 1/4".
1-year warranty. 110 V.

(313-6330)ea

HENRY SCHEIN®

#TR130 CEPHALOMETRIC TRACING BOX AND VIEWER

Lies flat, stands on desk, or hangs on wall. Includes 2 clips for holding film. Viewing area: 10" x 10". Overall dimensions: 12 1/4" x 11 1/2" x 2 1/4". 1-year warranty. 110 V.

(313-7060)ea

Not available in 220 V.



HENRY SCHEIN®

#100 SLIM LINE FILM ILLUMINATOR

Viewing area: 5" x 14". Overall dimensions: 15 1/4" x 6 1/4" x 2 1/4".
1-year warranty. 110 V.

- Viewing area: 5" x 14"
- Overall dimensions: 15 1/4" x 6 1/4" x 2 1/4"
- 1-year warranty
- 110 V

Illuminatorsea

Specify:

Beige(313-4538)

Light Gray Textured Case... (313-2507)

White(313-8131)

1" x 15" Fluorescent Bulb, White
(101-8301)ea



HENRY SCHEIN®

SLIM LINE PLASTIC FILM ILLUMINATOR

Solid plastic, lightweight construction. For desk or wall mounting. "Instant on" blast: quick on/off without light flicker. Sturdy on/off switch built into the unit: years of trouble-free use. Quick-release feature: easy removal of viewing plate to replace light bulb. Uniform illumination: optimum film viewing; accommodates a full series (up to 18) of intraoral film, as well as panoramic film. White plastic finish: easy to clean and fits into any office decor. 13 1/2" x 5" viewing area. Overall frame: 15" x 6 1/2" x 2". 1-year warranty.

- Solid plastic, lightweight construction
- For desk or wall mounting
- "Instant on" blast: quick on/off without light flicker
- Sturdy on/off switch built into the unit: years of trouble-free use
- Quick-release feature: easy removal of viewing plate to replace light bulb
- Uniform illumination: optimum film viewing; accommodates a full series (up to 18) of intraoral film, as well as panoramic film
- White plastic finish: easy to clean and fits into any office decor
- 13 1/2" x 5" viewing area
- Overall frame: 15" x 6 1/2" x 2"
- 1-year warranty

(100-4669)ea



**Ultimate Fabricel®
X-ray Gown and Cape**

Poly Outside, Blue
Laminated tissue/poly fabric provides superior strength and fluid resistance. Attached ties.

Gown, 30" x 42"50/case
(101-9087)



**Choice Reversible
X-ray Gown**

1-Ply Tissue/Poly
Wear this gown poly side in or out.

Blue, 30" x 42"50/case
(101-5598)



**Ultimate Mammography
Capes**

Dexter material, medium blue, 30" x 21" Front or back opening. Soft, pliable, leakproof and durable. Crepe exterior absorbs and deters runoff.

(308-7256)100/case



Ultra Mammography Cape, SMS nonwoven, dark blue, 48" x 23" Poncho-style with convenient front snaps and open sides. Generously sized for modesty.

(425-5731)50/case

TIDI Stretchable Cape, Nonwoven Dark blue, one size fits most. Latex-free elastic band at top provides adjustable sizing. Openings on both sides of cape allow for easy access.

(895-0055)25/case



**Ultimate
Gowns—Scrim
Reinforced**

Blue
Provide patients with comfort, hygiene, and modesty.

Sewn with Ties, 30" x 42"50/case
(859-5168)

Sewn with Ties, XX-large, 36" x 44"25/case
(859-7465)

Reinforced Material and Glued Seams, 30" x 42"50/case
(895-1526)

3X-large, 40" x 44"25/case
(422-3893)

Deluxe 3-Ply Glued, 30" x 42"50/case
(895-0580)



Bikini Briefs

Nonwoven, Blue
• One size fits waists 20"–52"
• No hip restrictions
• Fabric is durable and allows minimum see-through
• Benefits:

- Physician needs to carry only one size for most patients
 - Product stretchability provides comfort to patients
 - Patient privacy is maintained
 - Each bikini is individually packaged
- (895-0014)200/case



**Orthopedic Exam
Shorts**

Dark Blue
Durable, submersible, moisture-repellent, and modestly opaque in nonwoven material. Sewn elastic, latex-free waistband. Ideal for physical therapy, hydrotherapy, sports medicine, and lower-body X rays and exams. Smaller sizes are perfect for children and female patients.

Medium (28"–34")100/case
(967-8109)

Contains: 25 per bag, 4 bags per case.
Large (32"–40")100/case
(895-7343)

Contains: 25 per bag, 4 bags per case.
X-large (38"–48")100/case
(895-3575)

Contains: 25 per bag, 4 bags per case.
Plus Size (Up to 60")50/case
(895-9375)

Contains: 25 per bag, 2 bags per case.



Peace-of-Mind for You and Your Patients

**Dependable, convenient
and comfortable protection from TIDI**

Patient Gowns & Capes

ITEM NO.	DESCRIPTION	COLOR	SIZE	QTY.
(101-5598)	Duo-Ware Gown, Tissue/Poly	Sky Blue	30" x 42"	50/case
(116-0692)	SMS Non-Woven Gown	Blue	30" x 42"	50/case
(859-5168)	Scrim Gown, 2-Ply, XL	Blue	30" x 42"	50/case
(859-7465)	Scrim Gown, 2-Ply, XXL	Blue	36" x 44"	25/case
(422-3893)	Scrim Gown, 2-Ply, XXXL	Blue	40" x 44"	25/case
(425-5731)	Ultra Mammo Cape, Poncho-Style, Nonwoven	Dark Blue	48" x 23"	50/case
(108-1753)	Tissue/Poly Cape	Sky Blue	30" x 21"	100/case




X-Ray Filing Envelopes

ITEM NO.	DESCRIPTION	COLOR	SIZE	QTY.
(101-4868)	Storage Envelopes, Heavyweight	Golden Kraft	14.5" x 17.5"	500/case
(101-9137)	Storage Envelopes, Lightweight	Golden Kraft	14.5" x 17.5"	500/case
(115-3538)	Storage Envelopes, Lightweight	Golden Kraft	10.5" x 12.5"	100/box
(115-4284)	Storage Envelopes, Lightweight	Golden Kraft	8.5" x 10.5"	100/box



Available through
HENRY SCHEIN®
MEDICAL
To Order: 1-800-772-4346
To Fax: 1-800-329-9109
Online: www.henryschein.com/medical

www.tidiproducts.com



HENRY SCHEIN
FINANCIAL SERVICES
Business Solutions for Healthcare Practitioners

Join Our Community
Join Us on Facebook
[facebook](#)

Henry Schein Financial Services is not a bank, does not represent itself as such, and does not conduct banking activities. © 2012 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.

Discover Your Practice's Potential With Henry Schein Financial Services!

Find out how our low-cost, full-service equipment and practice-financing capabilities can help you work toward a secure future. Our specialists can facilitate by helping your practice achieve greater levels of productivity and profitability for the competitive edge you need in today's marketplace.

- Low rate, competitive financing
- Simple application process
- Same-day approvals
- Flexible programs
- Working capital

1-800-443-2756
or e-mail: hsfs@henryschein.com
www.henryschein.com/hsfs



Which Rewards Program Will You Choose?

Cashback Rebate Rewards

- Earn 2% **CASH BACK** on qualifying Henry Schein transactions up to \$100,000 a year†
- Earn UP TO 2% **CASH BACK** on other qualifying transactions†
- Get **AUTOMATIC REDEMPTIONS** in \$50 increments, posted to your account†

OR

Maximum Rewards®

- Earn 1.5 **POINTS** per \$1 in qualifying Henry Schein transactions*
- Earn 1 **POINT** per \$1 in other qualifying transactions*
- Redeem for **YOUR CHOICE** of exciting travel, merchandise, or gift cards

HENRY SCHEIN
FINANCIAL SERVICES
Business Solutions for Healthcare Practitioners

1-866-398-9296 or
www.henryschein.com/creditcard

The issuer of the card reserves the right to change the credit card account terms should your credit profile change or for other reasons described in the Cardmember Agreement. Cards are issued by First Bankcard, a division of First National Bank of Omaha. †You will receive 2% cash back on up to \$100,000 annually on qualifying Henry Schein transactions and 0.50% on qualifying Henry Schein transactions above \$100,000. *Complete details, including restrictions, limitations and exclusions will be available when you become a Cardmember. Henry Schein Financial Services is not a bank, does not represent itself as such, and does not conduct banking activities. © 2012 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.

HOLOGIC®

Henry Schein offers the complete line of Bone Densitometers from Hologic®.

Studies show 55% of women over 50 have osteoporosis or low bone mass. If you had the capability to test patients in a matter of minutes, how many patients could be tested per day?

Bone Density testing is safe, noninvasive and easy to do in your office. Great reimbursement rates make osteoporosis assessment very attractive and important for your patient care.



Your Henry Schein Sales Consultant can schedule our Imaging Equipment Specialist to help you evaluate bone density testing.

New from

 HENRY SCHEIN®
MEDICAL

SCANORA® 3D

The in-office Cone Beam CT system for ENT physicians and surgeons.

Short scan times from 10 to 26 seconds and effective exposure times from 2.25 to 6 seconds have the benefits of saving patient dose. Together with the integrated motorized seat and the rigid patient support, they increase the patient comfort and the diagnostic performance of the SCANORA® 3D system.

System highlights

- Stable seated patient
- Fully motorized chair
- 90 kV, 4–12.5 mA, pulsed X-ray generation
- Scan time 10 to 26 seconds
- Reconstruction time shorter than 1 minute
- Isotropic voxel size ranging from 0.133 to 0.35 mm
- 12" HD ClearTouch Panel
- Compatible with leading navigation systems
- DICOM/PACS compatibility
- Small footprint (2.24 m²)

Contact your Henry Schein Representative



Approved
Partner

SIEMENS

Ultrasound

Henry Schein brings the ACUSON family of Siemens products to your clinical practice



ACUSON X300™ ultrasound system, premium edition (PE)

Take full advantage of the broad range of capabilities this product offers, including Abdominal, Cardiac, Stress Echo, Vascular, OB/GYN, Urology, and Superficial/Small Parts imaging.

- The all-around performer, designed to address the diverse needs of any practice
- Ergonomic features facilitate ease of use in small spaces
- Highly configurable reporting package created specifically for the office environment



ACUSON X150™ ultrasound system

Ideal for general scanning including Abdominal, OB/GYN, Urology, and Superficial/Small Parts imaging. A scalable product designed to grow with your practice.

- Add powerful, advanced features at any time, as your clinical needs evolve
- Ultra-compact with an intuitive, easy-to-learn interface
- Highly configurable reporting package created specifically for the office environment



ACUSON Cypress™ cardiovascular system PLUS

A comprehensive cardiac and vascular ultrasound system in a small, hand-carried package ideal for the office environment.

- Supports adult and pediatric echocardiography, Stress Echo, and Peripheral and Abdominal Vascular imaging
- Easily carried from room to room or office to office
- Network connectivity with CypressViewer™ ultrasound software provides a cost-effective solution for offline review on a PC

SIEMENS



HENRY SCHEIN®

BRAND PRODUCTS



HENRY SCHEIN®
ULTRASOUND GELS
Blue

Multipurpose, hypoallergenic, medium-viscosity ultrasound gel that is bacteriostatic for patient satisfaction. Excellent conductive agent with superior lubricating and clinical features for the most common professional diagnostic and imaging procedures. Simple operation and disposal. Long shelf life.

- 8.5-fl-oz Bottle (900-4094)ea
- 5-liter Bottle (900-4095).....ea



HENRY SCHEIN®
CLEAR ULTRASOUND GELS

Multipurpose, hypoallergenic, medium-viscosity gel is clear, dye-free, and bacteriostatic for patient satisfaction. It is also an excellent conductive agent with superior lubricating and clinical features, making it suitable for the most common professional diagnostic and imaging procedures. For use in therapeutic ultrasound, diagnostic imaging, and other external procedures that require a thick medium with excellent lubricating and glide purposes. Will spread evenly without running or drying under standard clinical applications.

- 8.5-oz Bottle (900-4352)ea
- 5-liter Bottle (900-4353)ea



HENRY SCHEIN®
ULTRASOUND PAPER #110HD

For Use with 110mm Sony and Mitsubishi Printers
High-density video imaging media.

- 110mm x 20 meters (900-4305)5 rolls/box

HENRY SCHEIN®
BLACK AND WHITE THERMAL VIDEO PRINTER PAPER

Ultrasound Film Media

- 110mm x 20m Roll HD (High Density) (112-5630)5/box
- 110mm x 20m Roll HG (High Gloss) (112-5631)5/box
- 110mm x 20m Roll STD (Standard) (112-5632)5/box



Aquasonic® 100 Ultrasound Transmission Gels

Most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required.

- Acoustically correct for the broad range of frequencies used
- Hypoallergenic, bacteriostatic, and nonirritating
- No formaldehyde
- Not a spermicide
- 2-oz Tube, Clear (101-0693)ea
- 60-g Tube, Blue (100-7634)ea

- 0.25-liter Bottle, Blue (101-1042)ea
- 1-liter Bottle with Dispenser Cap, Clear (101-2417)ea
- 1-liter Bottle with Dispenser, Blue (100-5242)ea
- 5-liter Bottle with Dispenser, Clear (101-9866)ea
- 5-liter SONICPAC® with Refillable Dispenser, Blue (100-1833)ea
- 8.5-oz Dispensers (727-9534)12/box
- 5-liter Econopaks, Clear (426-2081)4/case



PARKER LABORATORIES, INC.

Transeptiq®
Cleansing Solution

Multipurpose spray for cleaning ultrasound transducer/probe surfaces. Does not contain chloride, phenols, or glutaraldehyde, and will not cause transducer membrane swelling.

- 250 mL (727-4217)ea



PARKER LABORATORIES, INC.

Aquagel®
Lubricating Gels

- Convenient flip cap
- For hospital, medical, and personal use
- Economical ½ gallon
- Non-runny, viscous gel
- Not a contraceptive
- Hypoallergenic, bacteriostatic, and nonirritating

- 142-g Tube (193-1722)20/box
- ½-gallon Bottle with Pump (101-4169)ea





X-ray Filing Storage Envelopes, Kraft Paper
 8½" x 10½"100/box
 (115-4284)
 10½" x 12", Lightweight100/box
 (115-3538)
 14½" x 17½", Lightweight500/case
 (101-9137)
 14½" x 17½", Heavyweight500/case
 (101-4868)



X-ray File Pocket Kraft Envelopes, Preprinted
 These envelopes safely store X-rays. Handy thumb notch for easy X-ray removal.
 • Preprinted for patient and doctor identification
 • 32-lb kraft paper
 8½" x 10½"250/box
 (366-0371)
 11½" x 14½"250/box
 (366-6527)
 14" x 17"100/box
 (366-8310)



X-RAY MAILER
 11" x 13"50/case
 (927-9654)
 15" x 18"100/pkg
 (101-0198)



X-Ray Film Filing Envelopes
 8½" x 10½"ea
 (557-5703)
 10½" x 12"ea
 (557-4429)
 14½" x 17"500/box
 (303-4837)



Film Filing Envelopes
 Manufactured from extra-heavy-duty envelope kraft, this wide range of envelope sizes features printed identification data spaces at both the upper left and right corners. They're tough enough to withstand use and prevent X-ray curling.
 8" x 10"500/box
 (108-3719)
 10" x 12"500/box
 (808-4214)
 14" x 17"500/box
 (226-3363)



Mailing Envelopes
 This line features individual cardboard stiffeners for film protection and a flap with a metal clasp. Manufactured from heavy-duty envelope kraft.
 10" x 12"100/box
 (782-6032)
 14" x 17"100/box
 (594-6094)



X-ray Utility File
 Easily wall mounted via 2 keyholes. The file's generous size will hold your largest cassettes, films, jackets, files, or whatever else you need to keep conveniently at hand. The file is manufactured of steel with a baked-on white finish. The individual bins have Plexiglas® front panels that allow you to see at a glance what is inside. The file has a unique rubber mat on the bottom of the bin that cushions anything placed on it.
 Singleea
 (989-8030)
 Doubleea
 (375-2065)



Film Filing Cart
 This versatile unit has hundreds of uses:
 • Shutting films from reading room to file area
 • Transporting film from processing to reading room
 • Transporting films from one unit to another
 • Handy for hospital rounds
 It's constructed of heavy satin anodized aluminum. Four hooded wheel casters provide total mobility, prevent tipping. Heavy duty wire dividers facilitate filing.
 Filing Carts
 Stand-Up (10 Dividers)
 40½"H x 21½"W x 32½"L
 (292-3250)ea
 Desk Height (6 Dividers)
 28¾"H x 17¼"W x 24¾"L
 (396-5639)ea



File Stack
 This ingenious filing system holds more, takes less floor space, and grows as your needs grow, providing total economy. File Stack takes only 19½" x 32" of valuable floor space and stacks as high as you want. Made of sturdy 20-gauge steel, individual units are secured to each other by bolts for total rigidity. Individual File Stack units are supplied with 3 adjustable dividers which can be placed where you want them for complete versatility of interior arrangement. Paneled doors slide easily at the touch of a finger for ready access. Plastic identification strips clip on easily to hold label inserts. Overall finish is in attractive, easy-to-clean white, baked-on polyester powder coating.
 File Stack Cabinet (50 lb)ea
 (594-5805)
 Base (12½ lb)ea
 (594-4623)
 Extra Dividersea
 (594-1652)



Color-Coded Chart Dividers

- Organize your patient information for quick and easy reference
- 2-hole punched
- Sets are collated, tabs 1-8
- Size: Letter, 8 1/2"W x 11 1/2"H with 3/8" Tab, 7 Pt. manila Dividers25/Pk

Specify:
(367-0084) History & Physical
(367-0085) Progress Notes
(367-0086) Correspondence
(367-0087) Lab/X-ray
(367-0088) Insurance
(367-0090) Miscellaneous
(367-0091) Blank



Color-Coded 11 Pt. End-Tab Folders

- Color code large groups of folders for easy return and access
- 10% postconsumer recycled content
- Generous 2" capacity prong "B" style fasteners
- 2-ply end-tab is embossed for easy label placement
- Size: Letter, 12 1/4"W x 9 1/2"H, scored for 3/4" expansion
- 11 Pt. Color-Coded Folders Fasteners Positions 1 & 3.....50/Box




Specify:
 Pink **(367-2135)**
 Orange **(367-1388)**
 Yellow **(367-0165)**
 Gray **(367-9377)**
 Blue **(368-0148)**
 Red **(367-8590)**
 Purple **(367-4353)**
 Green **(368-0149)**



HAZARD CLASS		
PROTECTIVE EQUIPMENT	COLOR/SHAPE CODES	DEGREE OF HAZARD
1	COMBUSTIBLE	0 - MINIMAL
2	REACTIVE	1 - SLIGHT
3	HEALTH	2 - MODERATE
4	PROTECTIVE	3 - SERIOUS
5	COMBUSTIBLE	4 - EXTREME
6	HEALTH	
7	HEALTH	
8	HEALTH	
9	HEALTH	
10	HEALTH	

683-5901
 Size: 11"W x 17"H



Allergy Label
 1 1/2"W x 7/8"H, Level A
(367-1244).....250/roll



Chart Label-X-Rays
 Red
(367-2598).....250/roll



Medical Alerts Label
 3 1/4"W x 1 1/4"H, Level B
(367-1255).....250/roll



Allergic to Label
 3 1/4"W x 1 1/4"H, Level B
(367-2612).....250/roll




Color-Coded X-ray Jackets-CC6 Style Standard Box Pack

- 11-point manila stock
- Scored for 1/2" expansion
- Full height side gussets
- Color bars printed on both ends
- Black loan bar printed on one end
- Offset thumbcut at top of both panels
- 10 distinctive colors match Smead numeric label colors
- 2"H x 3/4"W bars with reverse numerals printed on both sides of tab
- Two color blocks per jacket numbered 00 thru 99 (100 combinations)
- One color block per jacket in bottom position for use as aging year scale
- Front and back panels printed with medical information forms
- 10% recycled content, 10% post-consumer material
(975-4588).....100/box




End-Tab Manila Folders

- 11 Pt., Positions 3 & 5 **(367-5334)**50 per Box
- 11 Pt., Positions 1 & 3 **(367-8314)**50 per Box
- 11 Pt., Position 1 **(367-5129)**50 per Box
- 11 Pt., No Positions **(367-1105)**50 per Box



Fastener Positions



ALPHA-Z™ ACC SERIES LABELS

Size: 1"W x 1 1/4"H
 These durable labels can be used with either end-tab or top-tab products.

Starter Kit-File Box System

Includes file box, indexes, and 2,200 assorted labels.
(367-0876)Ea

Starter Kit-Roll System

Includes display unit and 6,500 assorted labels.
(366-1261)Ea

Description	Qty.
Pack	100/PK
Roll	250/RL

Pack Roll

A	367-0847	366-1034
B	367-0848	366-7536
C	367-0849	366-9064
D	367-0850	366-8712
E	367-0851	366-6876
F	367-0852	366-8970
G	367-0853	366-5868
H	367-0854	366-1979
I	367-0855	366-7306
J	367-0856	366-6519
K	367-0859	366-8001
L	367-0860	366-4463
M	367-0861	366-1798
N	367-0862	366-4198
O	367-0863	366-8869
P	367-0864	366-9522
Q	367-0865	366-3440
R	367-0867	366-9901
S	367-0868	366-3612
T	367-0869	366-1951
U	367-0870	366-7541
V	367-0871	366-9503
W	367-0872	366-1599
X	367-0873	366-5002
Y	367-0874	366-8713
Z	367-0875	366-7574

Call:
1-800-P-SCHEIN
(1-800-772-4346)
8am-9pm (et)

Fax:
1-800-329-9109
24 Hours


www.henryschein.com/medical

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products, excluding flu vaccine and Synagis®. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2012 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. †For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

 **HENRY SCHEIN®**
Henry Schein, Inc. • 135 Duryea Road • Melville, NY 11747

PRSR STD
U.S. POSTAGE
PAID
HENRY SCHEIN, INC.

12MS7244



Henry Schein DxRx Solutions

*We're there when you need us—
with our DxRx Solutions Hotline*

Do you have questions regarding the diagnostic, pharmaceutical and vaccine items that you purchase from Henry Schein? Need help finding a CPT code, package insert, or looking for brand to generic comparisons?

The Henry Schein DxRx Solutions Hotline provides information on the following topics for diagnostics, pharmaceuticals and vaccines:

- Indications and usage as provided on a package insert
- Package inserts
- Brand to generic comparisons
- CPT codes and J codes†

†For informational purposes only. Source: Reimbursement Codes.com and/or www.cms.gov <<http://www.cms.gov>>.
Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

*For informational purposes only; does not constitute legal advice. Source: www.immunize.org



- CLIA Product Status
- Storage Requirements
- ACIP recommendations for vaccinations
- State Law requirements for school and day care entry*
- Latest updates on supply issues, and much more

 **HENRY SCHEIN®**
MEDICAL

1.800.PSCHEIN
www.henryschein.com/dxrxsolutions

Rely on Us

DxRxSolutions™
DIAGNOSE • TREAT • PREVENT

1.877.523.SHOT 8:30am–5:00pm, est
1.877.523.7468
dxrxsolutions@henryschein.com