



# Product Catalog

[www.morrisonmed.com](http://www.morrisonmed.com) 800 438 6677 toll free 877 469 9696 toll free fax

2008 • 2009



Morrison Medical  
3735 Paragon Dr.  
Columbus, Ohio  
43228

It is the mission of Morrison Medical to:

- Provide our Customers with a wide range of quality emergency medical products, superior service and to conduct business with integrity as we have since 1970.
- Provide our Employees with a working atmosphere that offers stability and encouragement to realize the ambitions they have for themselves and their families.
- Remain independent and provide a competitive return on investment.

Hello,

We at Morrison Medical wish to extend warm wishes to each of you and thank you for your continued support and guidance.

It is hard to believe it is time for a new catalog.

You will notice several additions to our product line. I would like to introduce the new Morrison Medical Instant Cold Packs. They are made by using an environmentally safe chemical and water. Because of this, there is no danger of skin irritation and there are no special shipping requirements. Please contact us for samples and pricing.

We are dedicated to helping our distributors any way we can. Please feel free to call us for help on a special product, customized marketing literature, customer service assistance, contract manufacturing or a new product idea.

Thank you again for your support, guidance and friendship.

Sincerely,

Don Evans  
Vice President / General Manager



Don Evans  
Vice President/General Manager



Della Thompson  
Controller



Lena Smith  
Director of Sales/Marketing



Grace Taylor  
Plant Manger



Andrea Andrako  
Customer & Marketing Specialist



Danya Bindner  
Staff Accountant



Sherwin Meyers  
Sales/Marketing



Christina Balzer  
Assistant Plant Manger

<b>Cold Packs / Ice</b>		<b>Big Blue/Super Blue Head Immobilizers</b> .....	20
Cold Packs / Ice .....	4	Pediatric Head Immobilizer Blocks .....	20
<b>Silk Screening</b>		Head Vise III (Disposable) .....	21
Silk Screening .....	5	Improved Head Vise II (Disposable) .....	21
<b>New Product Preview</b>		Original Head Vise II (Disposable) .....	21
New Products .....	6	Improved Head Vise I .....	21
<b>Straps</b>		Original Head Vise I .....	22
Patho-Shield Straps .....	7	Head Blocks .....	22
Nylon Straps .....	8	Sticky Blocks .....	22
Polypropylene Straps .....	9	Foam EEZ Head Stabilization System .....	22
Impervious Straps .....	10	KED Head/Chin Straps .....	23
Disposable Straps .....	11	Head or Chin Straps (Disposable) .....	23
Disposable Stick & Cuts .....	11		
<b>Strap Systems</b>		<b>Splints/Carrying Cases/ IV Armboards</b>	
Foot Strap .....	11	Padded Wood Board Splints .....	24
Best Strap Systems .....	12	Splint Carrying Case with Webbing Handle .....	24
Gemini Strap System .....	12	Splint Carrying Case with Over-the-shoulder Strap .....	24
Loop End Gemini Strap System .....	13	Zipper Splint Carrying Case .....	24
Orion Strap System .....	13	Cold Splint .....	24
Shoulder Harness Strap System .....	14	VEL Splints .....	25
Bow Tie Strap .....	14	Padded Cardboard Splints .....	25
Hip Hugger System .....	14	Corrugated Splints .....	25
Hip Huggie System .....	14	Rugged Cardboard padded Board Splints .....	26
Spectrum System .....	15	Locking Leg Splints .....	26
		IV Armboards .....	27
<b>Limb Restraints</b>		<b>Accessories/Miscellaneous</b>	
Limb Lok .....	16	Toe Glides .....	28
Ankle or Wrist Straps with Velcro .....	16	Sand Bags .....	28
Economy "Deluxe" Limb Restraints .....	16	Umbilical Clamps .....	28
Quick Check Limb Holder .....	17	Umbi Cutter .....	28
Heavy Duty Wrist Restraints .....	17	Baby Bunting .....	29
Heavy Duty Ankle Restraints .....	17	Survival Bag .....	29
Traction Ankle Hitch .....	17	Disposable Pillows .....	29
		Staph-Check Pillows .....	29
<b>Gait Belts/holders</b>		Backboard Pad .....	29
Exam Glove Holder .....	18	<b>Kits</b>	
Extrication & Fire Glove Holder .....	18	First Aid Kits .....	30
Equipment Holding Strap .....	18	Clean Up Kits .....	32
Gait Belts .....	18	Isolation Kits .....	32
		Influenza Protection Kits .....	32
<b>Stretchers/Transfer Systems</b>		Emergency OB Kits .....	33
Perry Pouch .....	19	Pick Up Powder .....	33
Perry Pouch II .....	19		
Sit Stretcher .....	19	Terms & Conditions .....	34
Lift Assist .....	19	Association Application .....	35
Soft Stretchers .....	19		
<b>Head/Neck Immobilizers</b>			
Soft Cervical Collars .....	20		

# MORRISON MEDICAL

## Insta-Cold Compress



Shake contents to bottom of bag.



Squeeze firmly to activate.



Gently shake bag.



Wrap in towel or cloth and apply .

**Insta-Cold packs are UREA BASED!!**  
**ENVIRONMENTALLY SAFE!**

Item #6601 5" x 7"  
 24 Packs /Case

Item #6600 6" x 9"  
 24 Packs /Case

AMERICAN MEDICAL SUPPLY & EQUIPMENT CO. INC. PHOENIX, ARIZONA

## Cold Packs

## Instant Disposable Cold Packs

Environmentally friendly urea based instant ice packs

- Easily activated with one hand
- Urea formula is safe and environmentally friendly
- Does not require special packaging requirements when shipping
- Available in 2 sizes

6601 5" x 7", 24 packs / case

6600 6" x 9", 24 packs / case



## Silk Screening

## Our Screen-printing capabilities

Identify Your Emergency Equipment with Your Name or Logo

- Helps prevent accidental loss of equipment
- Improves recovery of equipment if left behind
- Customized screen printing available on:

- Straps
- Bags
- Soft Stretchers Bags
- Sit Stretchers Bags
- Perry Pouch Bags

Contact Morrison Medical at 800.438.6677 or your distributor for details and artwork requirements.





Instant Disposable Cold Packs



Foam EEZ Head Stabilization System



Patho-Shield Straps  
New Colors!!



Limb-Lok Strap System



Bow Tie Strap System



Lift Assist



Hip Huggie Strap System

## Patho-Shield Straps

Patho-Shield straps have a 1,600 pound break strength. 100% impervious to blood and other body fluids. Very flexible and easily adjusts. Orange **OR**, Black **BK**, Yellow **YL**

### Patho-Shield Strap

1310P

2 Piece, metal push button buckle with metal swivel speed clip ends, 5'



### Patho-Shield Strap

1308P

1 Piece, metal push button buckle, 7'



### Patho-Shield Strap

1360P

2 Piece, plastic side release buckle with metal swivel speed clip ends, 5'



### Patho-Shield Strap



1300P

2 Piece, metal push button buckle with loop buckle ends, 5'



## Nylon Webbing

PLEASE SPECIFY COLOR

Orange **OR**, Black **BK**, Yellow **YL**, Red **RD**, Blue **BL**, Maroon **MA**, Green **GR**

Metal 2 Piece Push Button Buckle (fig. A) Nylon Webbing

- 1200 Loop-Lok (fig. H), 2 piece, 5', each
- 1200R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1201 Loop-Lok (fig. H), 2 piece, 7', each
- 1201R Loop-Lok Roller Buckle (fig. K), 2 piece, 7', each
- 1208 1 piece, 7', each
- 1209 1 piece, 9', each
- 1210 Swivel Speed Clip (fig. I), 2 piece, 5', each
- 1211 Swivel Speed Clip (fig. I), 2 piece, 7', each
- 1214 Big Mouth Swivel Speed Clip (fig. L), 2 piece, 5', each
- 1215 Non-Swivel Speed Clip (fig. J), 2 piece, 5', each
- 1216 Non-Swivel Speed Clip (fig. J), 2 piece, 7', each
- 1217 Extension strap, 3', each
- 1218 Extension strap, 4', each

Metal Cam (Friction) Buckle (fig. B) Nylon Webbing

- 1228 1 piece, 7', each
- 1229 1 piece, 9', each

Metal Roller (Friction) Buckle (fig. C) Nylon Webbing

- 1234 1 piece, 7', each
- 1235 1 piece, 9', each

Metal Drop Jaw Buckle (fig. D) Nylon Webbing

- 1238 Loop-Lok (fig. H), 2 piece, 5', each
- 1238R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1241 1 piece, 9', each

Metal 2" Double "D" Rings (fig. E) Nylon Webbing

- 1245 1 piece, 9', each

Plastic 2 Piece Side Release Buckle (fig. F) Nylon Webbing

- 1250 Loop-Lok (fig. H), 2 piece, 5', each
- 1250R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1251 Loop-Lok (fig. H), 2 piece, 7', each
- 1251R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1258 1 piece, 7', each
- 1259 1 piece, 9', each
- 1260 Swivel Speed Clip (fig. I), 2 piece, 5', each
- 1261 Swivel Speed Clip (fig. I), 2 piece, 7', each
- 1264 Big Mouth Swivel Speed Clip (fig. L), 2 piece, 5', each
- 1265 Non-Swivel Speed Clip (fig. J), 2 piece, 5', each
- 1266 Non-Swivel Speed Clip (fig. J), 2 piece, 7', each

Plastic Cam (Friction) Buckle (fig. G) Nylon Webbing

- 1278 1 piece, 7', each
- 1279 1 piece, 9', each





Polypropylene Webbing

PLEASE SPECIFY COLOR

Orange **OR**, Black **BK**, Yellow **YL**, Red **RD**, Blue **BL**

Metal 2 Piece Push Button Buckle (fig. A) Polypropylene Webbing

- 1300 Loop-Lok (fig. H), 2 piece, 5', each
- 1300R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1301 Loop-Lok (fig. H), 2 piece, 7', each
- 1301R Loop-Lok Roller Buckle (fig. K), 2 piece, 7', each
- 1308 1 piece, 7', each
- 1309 1 piece, 9', each
- 1310 Swivel Speed Clip (fig. I), 2 piece, 5', each
- 1311 Swivel Speed Clip (fig. I), 2 piece, 7', each
- 1314 Big Mouth Swivel Speed Clip (fig. L) 2 pieces, 5', each
- 1315 Non-Swivel Speed Clip (fig. J), 2 piece, 5', each
- 1316 Non-Swivel Speed Clip (fig. J), 2 piece, 7', each
- 1317 Extension strap, 3', each
- 1318 Extension strap, 4', each



Metal Cam (Friction) Buckle (fig. B), Polypropylene Webbing

- 1328 1 piece, 7', each
- 1329 1 piece, 9', each

Metal Roller (Friction) Buckle (fig. C), Polypropylene Webbing

- 1334 1 piece, 7', each
- 1335 1 piece, 9', each

Metal Drop Jaw Buckle (fig. D), Polypropylene Webbing

- 1338 Loop-Lok (fig. H), 2 piece, 5', each
- 1338R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1341 1 piece, 9', each

Plastic 2 Piece, Side Release Buckle, (fig. F), Polypropylene Webbing

- 1350 Loop-Lok (fig. H), 2 piece, 5', each
- 1350R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1351 Loop-Lok (fig. H), 2 piece, 7', each
- 1351R Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1358 1 piece, 7', each
- 1359 1 piece, 9', each
- 1360 Swivel Speed Clip (fig. I), 2 piece, 5', each
- 1361 Swivel Speed Clip (fig. I), 2 piece, 7', each
- 1364 Big Mouth Swivel Speed Clip (fig. L), 2 piece, 5', each
- 1365 Non-Swivel Speed Clip (fig. J), 2 piece, 5', each
- 1366 Non-Swivel Speed Clip (fig. J), 2 piece, 7', each


Metal 2" Double "D" Rings (fig. E), Polypropylene Webbing

- 1345 1 Piece, 9', each

Plastic Cam (Friction) Buckle (fig. G), Polypropylene Webbing

- 1378 1 piece, 7', each
- 1379 1 piece, 9', each

## Impervious Webbing



PLEASE SPECIFY COLOR

Orange **OR**, Black **BK**, Yellow **YL**, Blue **BL**

Metal 2 Piece Push Button Buckle (fig. A), Impervious

- 1300X Loop-Lok (fig. H), 2 piece, 5', each
- 1300XR Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1301X Loop-Lok (fig. H), 2 piece, 7', each
- 1301XR Loop-Lok Roller Buckle (fig. K), 2 piece, 7', each
- 1308X 1 piece, 7', each
- 1309X 1 piece, 9', each
- 1310X Swivel Speed Clip (fig. I), 2 piece, 5', each
- 1311X Swivel Speed Clip (fig. I), 2 piece, 7', each
- 1314X Big Mouth Swivel Speed Clip (fig. L), 2 piece, 5', each
- 1315X Non-Swivel Speed Clip (fig. J), 2 piece, 5', each
- 1316X Non-Swivel Speed Clip (fig. J), 2 piece, 7', each
- 1317X Extension strap, 3', each
- 1318X Extension strap, 4', each

Metal Cam (Friction) Buckle (fig. B), Impervious

- 1328X 1 piece, 7', each
- 1329X 1 piece, 9', each

Metal Roller (Friction) Buckle (fig. C), Impervious

- 1334X 1 piece, 7', each
- 1335X 1 piece, 9', each

Metal Drop Jaw Buckle (fig. D), Impervious

- 1338X Loop-Lok (fig. H), 2 piece, 5', each
- 1338XR Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1341X 1 piece, 9', each

Metal 2" Double "D" Rings (fig. E), Impervious

- 1345X 1 piece, 9', each

Plastic 2 Piece Side Release Buckle (fig. F), Impervious

- 1350X Loop-Lok (fig. H), 2 piece, 5'. Each
- 1350XR Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1351X Loop-Lok (fig. H), 2 piece, 7', each
- 1351XR Loop-Lok Roller Buckle (fig. K), 2 piece, 5', each
- 1358X 1 Piece, 7', each
- 1359X 1 Piece, 9', each
- 1360X Swivel Speed Clip (fig. I), 2 piece, 5', each
- 1361X Swivel Speed Clip (fig. I), 2 piece, 7', each
- 1364X Big Mouth Swivel Speed Clip (fig. L), 2 piece, 5', each
- 1365X Non-Swivel Speed Clip (fig. J), 2 piece, 5', each
- 1366X Non-Swivel Speed Clip (fig. J), 2 piece, 7', each

Plastic Cam (Friction) Buckle (fig. G), Impervious

- 1378X 1 piece, 7', each
- 1379X 1 piece, 9', each

## Disposable Straps

### The Foot Strap

Disposable Poly, Foot Strap, 1 Piece, Plastic Buckle

- Use black strip as a guide
- Place black strip under feet. Criss Cross over the ankles. Loop ends thru hand holes of backboard and fasten buckles. Pull to tighten.
- Single use for backboard
- 2" polypropylene webbing

1000 1 Piece Disposable Foot Strap, each



### Disposable Poly, Strap, 2 Piece, Plastic Buckle

- Single use for backboard, not for transport
- 2" polypropylene, 5' or 9' with different style ends
- Colors for disposable: Orange **OR** or White **WH**

1390 Loop-Lok , 2 Piece, 5', each  
 1390-36 Loop-Lok , 2 Piece, 5', each, case of 36  
 1392 1 Piece, 9', each  
 1392-36 1 Piece, 9', case of 36  
 1394 Plastic Swivel Speed Clip , 2 Piece, 5', each  
 1394-36 Plastic Swivel Speed Clip , 2 Piece, case of 36



### Disposable Poly, Strap, Double Adjusting Buckle

- Replace strap when needed
- Single use for backboard, not for transport
- 2" polypropylene, 6', one piece strap or Loop-Lok™, 2 piece, 7', strap
- 2" plastic side release double adjusting buckle
- Colors Orange **OR** or White **WH** (specify color)

1396 Double adjust buckle, 1 piece, 6', each  
 1397 Replacement strap (for 1396) 6', each  
 1396L Double adjust buckle, Loop-Lok™, 7', each  
 1397L Replacement strap (for 1396L) 7', each



### Disposable Stick & Cut Strap

- Set of 3 Vinyl straps, 2 x 60", with 6" of adhesive on each end
- Ends feed through hand holds and stick back onto strap
- One time use, cut strap off to remove patient
- Adjusts to about any size patient
- X-ray translucent
- Available in Orange **OR**, Black **BK**, Yellow **YL**, Blue **BL**

1398 Disposable Stick & Cut Straps, set of 3



### Velcro® Keeper

1275 Nylon Webbing Loop-Lok , 2 piece, 5', each  
 1375 Polypropylene Webbing Loop-Lok , 2 piece, 5', each  
 1375X Impervious Loop-Lok 2 piece, 5', each



## Best Strap Systems

## Best Strap System, Super Deluxe

- Total strap system for fast, firm total body immobilization
- Two “suspender” type over-the-shoulder straps
- 5 adjustable torso/limb straps, with adjustable center strap
- Nylon webbing with Velcro® closures and EZ Peel ends
- All straps color-coded for easy application, 12 points of attachment
- No choking “v” by the patient’s neck

1280 Best Strap System

1691 Best Strap System carrying case (A)



## Economy Best Strap System With Velcro®

- Total strap system for fast, firm total body immobilization
- Two “suspender” type over-the-shoulder straps
- No choking “v” by patient’s neck
- 4 adjustable torso/limb straps, with adjustable center strap
- Black polypropylene webbing with Velcro® closures

1381 Economy Best Strap w/Velcro®



## Economy Best Strap System With Buckles

- Total strap system for fast, firm total body immobilization
- Two “suspender” type over-the-shoulder straps
- No choking “v” by patient’s neck
- 4 adjustable torso/limb straps, with adjustable center strap
- Black polypropylene webbing with Plastic Side Release Buckles

1382 Economy Best Strap w/buckles



## Gemini Strap Systems

## Gemini Strap System (patent #6, 092, 525)

- 4 alternating color-coded body straps
- Suspender type shoulder straps
- Snap hooks with metal spring closures
- Adjustable with the super lock buckles

1383X Gemini Strap System



## Loop End Gemini Strap Systems

### Gemini Strap System (A)

- Suspender type shoulder straps
- 4 body straps
- Made with blue polypropylene webbing
- Two end styles to choose from: Loop-Lok and Velcro® Loop

1383A-L Gemini Strap System A Loop-Lok  
1383A-V Gemini Strap System A Velcro® Loop



### Gemini Strap System (B)

- Suspender type shoulder straps
- 3 body straps
- Made with blue polypropylene webbing
- Two end styles to choose from: Loop-Lok and Velcro® Loop

1383B-L Gemini Strap System B Loop-Lok  
1383B-V Gemini Strap System B Velcro® Loop



### Gemini Strap System C


- 3 body straps
- Made with blue polypropylene webbing
- Two end styles to choose from: Loop-Lok and Velcro® Loop

1383C-L Gemini Strap System C Loop-Lok  
1383C-V Gemini Strap System C Velcro® Loop



### Replacement loop ends

1383L-F Loop-Lok End  
1383V-F Velcro® End



## Orion System

### Orion Strap System

- Three vinyl covered foam straps with Velcro® closures
- Loops through hand holes of backboard
- 60" chest strap, 64" waist strap and 34" ankle strap
- Great for water rescue

1384 Orion Strap System



## Shoulder Harness Strap System



- Two shoulder straps and two piece torso strap
- Adjustable points each strap, w/metal push button buckle
- Attach to elevating backrest on gurneys
- Chrome plated buckles and hardware
- 2" nylon, available in seven (7) colors: Orange **OR**, Black **BK**, Yellow **YL**, Red **RD**, Blue **BL**, Maroon **MA**, Green **GR**

1285 Shoulder Harness System, Complete

1286 Shoulder Harness System, lap strap only

1287 Shoulder Harness System, shoulder straps only

Also available in Polypropylene, Impervious or Plastic  
CALL FOR PRICING

## Bow Tie Strap



- 1" Black Polypropylene
- Looped on both ends to attach to shoulder straps of harness
- Prevents shoulder straps of harness from slipping off of patient

1355A Bow Tie Strap, each

## Hip Hugger / Huggie Systems



## The Hip Hugger System

- Intended for suspected pelvic/hip injury stabilization
- Eases discomfort in transport
- Made with easily cleaned orange vinyl
- Color coded straps for easy identification in hurried situations

3200 Hip Hugger System



## The Hip Huggie System

- Intended for suspected pelvic/hip injury stabilization
- Eases discomfort in transport
- Made with easily cleaned orange vinyl
- 2 Color coded & numbered straps for easy application
- Smaller more economical version of the original Hip Hugger
- Comes with a 15" x 17 1/2" carrying case

3201 Hip Huggie System

## Spectrum Strap System

## The Spectrum Strap with Swivel Speed Clip End

- 4" wide strap for improved patient stabilization and patient comfort
- Made with highly visible, durable orange vinyl, black 2" webbing and 4" black Velcro®
- X-ray translucent
- Metal Swivel Speed Clip Ends

1376XOR Spectrum Strap



## The Spectrum Strap with Loop-Lok™ End

- 4" wide strap for improved patient stabilization and patient comfort
- Made with highly visible, durable orange vinyl, black 2" webbing and 4" black Velcro®
- X-ray translucent
- Loop-Lok Ends

1377XOR Spectrum Strap



## The Spectrum Pelvic Immobilizer

- Extra wide, 12" pelvic strap
- Improves patient stabilization and comfort
- Made with highly visible, orange vinyl
- X-ray translucent
- 4 Velcro® Ends

1374XOR Spectrum Pelvic Immobilizer




800 438 6677 toll free

614 469 9696 fax

[www.morrisonmed.com](http://www.morrisonmed.com)

## Limb Lok



## Limb Lok

- Secures wrist of unconscious patient to backboard strap.
- Strong polypropylene 1" black strap
- Protects limbs during transport.
- Adjusts to fit various sizes
- 1" plastic side release buckle



1115 Limb Lok, each

## Ankle Or Wrist Straps With Velcro®

## 2" Nylon webbing Velcro® closures

- Wrist straps secures to top of cot rail (set of two individual straps)
- Ankle straps secures to main cot frame (pair attached to anchor)
- Seven (7) colors: Orange **OR**, Black **BK**, Yellow **YL**, Red **RD**, Blue **BL**, Maroon **MA**, Green **GR**

1290 Wrist Strap, pair

1291 Ankle Strap, each



## Economy "Deluxe" Limb Restraints

- Soft breathable foam
- Strong polypropylene 1" strap
- Velcro® closure holds around patient limb, allowing "hands-free" tying of straps
- Several styles: "tie down straps" or double "D" rings
- Available in cases

1103 Adult deluxe, tie down style, pair (fig. A)

1103-24 Adult deluxe, tie down style, cs/24 (fig. A)

1103-48 Adult deluxe, tie down style, cs/48 (fig. A)

1114 Adult double "D" ring style, pair (fig. B)

1114-48 Adult double "D" ring style, cs/48 (fig. B)

## Infant Sizes

Fits - 0-2 years old (6"x1 1/4" Cuff, 30" Strap)

1120 Infant Deluxe, tie down style, pair


1120-48 Infant Deluxe, tie down style, cs/48

## Pediatric Sizes

Fits - 2-5 years old (9"x2" Cuff, 40" Strap)

1124 Pediatric Deluxe, tie down style, pair

1124-48 Pediatric Deluxe, tie down style, cs/48






## Ankle Or Wrist Straps With Velcro®

## Quick Check Limb Holder

- Large 3 x 12½" with 60" strap
- Soft breathable foam
- Unique wrist adjustment with buckle
- Velcro® holds while you fasten buckle
- Doesn't "cinch" like other restraints, makes it quick to check patient
- Reusable with proper care
- Pre-adjustable
- Doesn't require resetting or retying
- Available in children sizes

1130 Quick Check limb restraint with buckle, 1 pair



## Heavy-Duty Wrist Restraints – Reusable

- Comfortable blue neoprene over nylon webbing
- Velcro® to assist restraint placement
- Quick release plastic buckles for fast on/off application
- Each side has two straps with double "D" ring closures for two-point anchoring
- 4 welded "D" rings for quick, strong application

1172 Heavy-Duty Wrist Restraint, pair



## Heavy-Duty Ankle Restraints – Reusable

- Comfortable red neoprene over nylon webbing
- Velcro® to assist restraint placement
- Quick release plastic buckles for fast on/off application
- Each side has two straps with double "D" ring closures for two-point anchoring
- 4 welded "D" rings for quick, strong application

1173 Heavy-Duty Ankle Restraint, pair



## Traction Ankle Hitch, Adjustable

- Totally wraps around ankle with Velcro® closure
- Made from padded orange vinyl with 1" polypropylene straps

4010C Ankle Hitch, each



## Holders



### Exam Glove Holder

- Made with orange vinyl with black binding around bag
- Polypropylene belt loop
- Holds up to 3 pair of exam gloves (Included)
- To reorder gloves, please contact a Morrison Medical customer services rep.
- 3½" X 3½"

7024 Exam Glove Holder, each



### Extrication & Fire Glove Holder

- Holder attaches to belt with Heavy Duty 1" metal swivel speed clip
- Gloves are easily accessible with Velcro® closure
- Available in 1" Black Polypropylene **BK** or 1" Gray Nylon **GY**

7020 Glove Holder, each



### Equipment Holding Straps (Black Only)

- Loop-Lok™ black polypropylene strap (1" x 36")
- 2 piece 1" plastic side release buckle
- 2 metal screw on brackets (screws not included)

1355 Equipment Holding Strap, 1 complete set

## Gait Belts



### Gait Belts

- Assist in safe transfer of patients
- 2" wide one piece belt available standard length 60" or extra large length 84"
- Buckle styles: Plastic Side Release and Metal Drop Jaw
- Durable washable belts in cotton: 2" cotton pastel stripe or 2" cotton plain natural

7038 Metal Drop Jaw 60", each

7038-7 Metal Drop Jaw 84", each

7050 Plastic Side Release 60", each

7050-7 Plastic Side Release 84", each



Perry Pouch

The Perry Pouch

- Stretcher cinches around patient when lifted
- Fast, safe patient transportation
- Compact and lightweight, 1 lb., 4 oz.
- Made with orange vinyl and 1" black polypropylene webbing
- Black vinyl handles for "easy on your hands" lifting
- Carrying Case Included
- Supports up to 500 pounds

3000 Perry Pouch, 28" x 31", each



The Perry Pouch II

- Stretcher cinches around patient when lifted
- Fast, safe patient transportation
- Compact and lightweight, 1 lb., 4 oz.
- Made with orange vinyl and 1" black polypropylene webbing
- Double black vinyl handles for easy "hands off" in confined areas
- Carrying Case Included
- Supports up to 500 pounds

3001 Perry Pouch II, 28" x 31", each




Stretchers

Sit Stretcher

- Compact Stretcher, can be used in confined areas
- Features 6 handholds
- Continuous straps
- Made from orange vinyl and 2" yellow nylon straps
- Carrying Case Included

3008 Sit Stretcher, 20" x 34", each



Lift Assist

- Used to transition patient from a sitting to a standing position.
- Made with durable orange vinyl with reinforced nylon webbing handles.
- Carrying Case Included
- Support up to 500 pounds.

3002 Lift Assist, 24 1/2" x 11 1/2", each



The Soft Stretchers, Regular to XXXL

- Specially designed for moving extra large patients
- Great for transporting patients out of confined areas or stairways.
- XXXL is 60 x 80" (size of a queen size bed)
- Each size supports over 2000 pounds
- Allows you to move several children or more than one average adult
- Continuous straps with vinyl handle for easier lifting
- Carrying Case Included
- Available in 3 sizes
- Orange **OR**, or Camouflage **CA** available


3010 Regular Soft Stretcher, 30" x 70" (A), each

3020 Large Soft Stretcher, 40" x 71" (B), each

3030 XXXL Soft Stretcher, 60" x 80" (C), each



## Soft Collars



## Cervical Collar – Universal Size

- Extended Velcro® closure for greater neck size range
- One size fits all – 4" x 21¾"
- Firm foam covered with washable stockinette

0420 4" x 21¾" Universal size

## Head Immobilizers



## Big Blue Head Immobilizer

- Vinyl coated foam for easy disinfecting
- Extra large pillows provide maximum support
- Ear holes for patient assessment
- Pillows firmly attach to base with wide Velcro®
- Position on backboard and secure with mounted straps with sewn-on "D" rings
- Big Blue includes base with attached straps, two large pillows and two head/chin straps
- Blue **BL**, Orange **OR**, Black **BK** (great for tactical use)

1420 Big Blue Head Immobilizer, each

1420-6 Case of 6




## Super Blue Head Immobilizer

- Vinyl coated foam for easy disinfecting
- Ear holes for patient assessment
- Head/chin straps adjust for firm comfortable hold
- Velcro® hook on base for pre-positioning on board
- Matching Velcro® loop for your board available #8912
- Attach with two 1" x 42" straps, included (removable)
- Complete Super Blue includes base, 2 foam pillows, 2 head/chin straps and two 1" x 42" securing straps

Blue **BL** only

1449 Super Blue Head Immobilizer, each

1449-6 Case of 6



## Pediatric Head Immobilizer Blocks

- Will fit Big Blue and most other popular head immobilizers
- Vinyl coated foam for easy disinfecting
- 2½" x 5½" x 3¾" pillow for smaller heads

Blue **BL** only


1432 Pediatric pillows, one pair

Head Immobilizers

Head Vise III (Disposable)


- All weather disposable head immobilizer
- Easy to apply
- Universal – will fit on most backboards and patients
- Easy access to ears
- Base measures 4 x 15"
- Extra wide wedges, 3", for greater stability
- Includes two blue all weather straps, 2 x 36" each

1470 Head Vise III  
 1470-24 Case of 24



Improved Head Vise II (Disposable)

- Made of sturdy, moisture resistant cardboard
- Includes comfortable foam head pad
- Quickly attaches to board with improved EZ peel tabs, easier to apply
- Improved access to ears
- X-ray translucent
- New compact size 5" x 15", stores flat
- Single use
- Adjustable for children and adults
- Can be written on



1465 Head Vise II, each  
 1465S Head Vise II w/#1475 strap, each

Original Head Vise II Head Immobilizer (Disposable)

- Made of sturdy, moisture resistant cardboard
- Includes comfortable foam head pad
- Quickly attaches to board with EZ peel tabs, use tape to secure patient, or purchase #1475 straps
- X-ray translucent, stores flat, for one time use
- Adjustable for children and adults
- Can be written on



1460 Original Head Vise II  
 1460S Original Head Vise II w/#1475 strap, each  
 1460-10 Case of 10  
 1460-50 Case of 50


Improved Head Vise I (Reusable)

- Made of sturdy, impervious corrugated plastic
- Includes comfortable foam head pad
- Quickly attaches to board with improved EZ peel tabs, easier to apply
- Improved access to ears
- X-ray translucent, easy to disinfect
- Compact size 5" x 15", stores flat
- Adjustable for children and adults, compatible with helmets
- Can be written on
- Sturdy-T construction, the larger surface gives more



1455 Head Vise I, each  
 1455S Head Vise I w/#1475 strap, each

## Head Immobilizers



## Original Head Vise I Head Immobilizer (Reusable)

- White, heavy duty, impervious corrugated plastic
- Completely waterproof, can be written on
- Adjustable for children and adults, compatible with helmets
- Attaches to board quickly with EZ pull tabs, use tape to secure patient or purchase #1475 straps
- Includes comfortable foam pad
- X-ray translucent, easy to disinfect, stores flat

1450	Original Head Vise I
1450S	Original Head Vise I w/#1475 strap, each
1450-18	Case of 18
1456-12	Head Vise I – Foam replacement, pack of 12 each



## Head Blocks (disposable) Foam Orange Blocks

- Rectangular (4" x 5" x 10") foam blocks sealed in highly visible orange vinyl, complete with handle
- Use anywhere for firm but lightweight support
- Blocks float, ideal for water rescue
- Can be written on
- Set includes 1 pair of head blocks plus 1 pair of head/chin straps



## Head Blocks Only

0320	Head Blocks, 1 set
0320-9	Case of 9 sets
0320-30	Case of 30 sets

## Head Blocks with Straps

0322	Set, blocks (0320) & straps (1475)
0322-9	Case of 9 sets
0322-30	Case of 30 sets



## Sticky Blocks Head Immobilizers

- Contoured for maximum head stabilization
- Universal size (3" x 5½" x 10<sup>3</sup>/<sub>8</sub>" )
- Comfortable foam, light blue
- Provides firm support
- Each block comes with Sticky tape affixed to the bottom; pull protective tape, position each Sticky Block and press to attach to any backboard
- Use tape to secure patient, or purchase Morrison Medical #1475 Head & Chin Straps
- Operational in all weather conditions

1400	Sticky Blocks, pair
1400-12	Case of 12 pair
1400-50	Case of 50 pair

Head Immobilizers

Foam-EEZ Head Stabilization System (Disposable)

- All weather disposable head immobilizer
- Easy to apply
- Universal – will fit on most backboards and patients
- Specially designed foam blocks to accommodate easy access to ears
- 4 each of 4 ½" x 4 ½" cervical alignment foam pads
- Includes 2 black velcro head & chin straps.
- Comfortable yet sturdy foam

1401 Foam-EEZ Head Stabilization System , each




Head and Chin Straps

KED Head/Chin Strap Disposable

- Tan
- Made of loop (soft) Velcro®-like material
- Apply over head or chin and press onto hook material
- 28" long overall
- Disposable

1478TAN KED head/chin strap replacement, pair



Head or Chin Straps, Disposable

- Orange self adhesive head/chin strap, 2" x 36"
- Peel backing off, and it sticks to itself
- Cut off and throw away
- EZ Peel back

1475 Self adhesive head/chin strap, disposable, each



Disposable Head & Chin Straps for Big Blue and Super Blue

- 2 pound White Crosslink Polyethylene Foam
- 2" x 7" x ½"
- 1" Black Velcro® Hook and Loop

1479 Disposable Head & Chin Straps for Big Blue and Super Blue, pair



Splints

Economy Padded Wood Board Splints

- Padded board splints are made of 1/2" x 3" solid wood core with patient safe rounded corners
- Orange reusable vinyl covered
- 1/2" foam on one side
- Large size range: 12", 15", 16", 18", 24", 32", 36", 42", 48", 54" & custom sizes

1800	Set of 2 each 15", 36", 54"
1812	12" Padded Board Splint
1815	15" Padded Board Splint
1815-150	15" Padded Board Splint, Case of 150
1816	16" Padded Board Splint
1818	18" Padded Board Splint
1824	24" Padded Board Splint
1832	32" Padded Board Splint
1836	36" Padded Board Splint
1836-75	36" Padded Board Splint, Case of 75
1842	42" Padded Board Splint
1848	48" Padded Board Splint
1854	54" Padded Board Splint
1854-50	54" Padded Board Splint, Case of 50
1875	1800 set with 1680 carrying case



Carrying Cases

Splint Carrying Case with Webbing Handle

- Has Velcro® closure and webbing handle
- Holds 6 or more splints plus accessories
- Orange nylon reinforced vinyl

1680 Splint Carrying Case



Splint Carrying Case with Over-the-shoulder Strap

- For hands-free carrying
- 9 1/4" x 56", Velcro® closure
- A webbing carrying handle plus over-the-shoulder (rifle-sling) style handle
- Holds 6 or more splints plus accessories
- Your choice of deluxe blue cordura **BL** or orange vinyl **OR**

1685 Over-the-shoulder strap carrying case



Zipper Splint Carrying Case

- 9 1/4" x 56", with full zipper closure, webbing handle
- Holds 6 or more splints & accessories

1681 Case with zipper closure and webbing handle



Cold Splint

Cold Splint

- Disposable ice pack adhered to center of splint is easily activated with mild force.
- Features easy, removable Velcro® straps
- Micro perf edges protect from paper cuts on cardboard

1700	Cold Splint, Cardboard, 12" each
1701	Cold Splint, Cardboard, 18" each
1702	Cold Splint, Cardboard, 24" each





VEL Splints / Cardboard & Plastic Splints

VEL Splints

Velcro® Enclosed Limb Cardboard or Corrugated Plastic Splints


- Removable Velcro® Straps
- Splints feature 5/8" punch-out holes for tie downs
- Micro Perf edges to protect from paper cuts on cardboard
- 4 sizes: 12" x 9", 18" x 9", 24" x 14" & 34" x 19"
- 3/8" Full Foam (FF)
- 24" and 34" lengths offer an extra center crease for firm support of thin limbs
- Available in cases – call for pricing



Cardboard		Corrugated Plastic	
1707	FF, 12" (A)	2707	FF, 12" (B)
1750	FF, 18" (A)	2750	FF, 18" (B)
1760	FF, 24" (A)	2760	FF, 24" (B)
1770	FF, 34" (A)	2770	FF, 34" (B)

Cardboard or Corrugated Plastic Folding Splints

- Splints feature 5/8" punch-out holes for tie downs
- Micro Perf edges to protect from paper cuts on cardboard
- 4 sizes: 12" x 9", 18" x 9", 24" x 14" & 34" x 19"
- 4 styles:
  - A plain (no foam)
  - B center strip 3/8" foam (CF)
  - C full 3/8" foam (FF)
  - D 3 strips of 1/2" foam (each side & center) (3PC)
- Cardboard is white on one side.
- 24" and 34" lengths offer an extra center crease for firm support of thin limbs
- Available in cases – call for pricing



Cardboard Folding Splints

1500	Plain, 12" (Fig. A)	1504	CF, 12" (Fig. B)
1510	Plain, 18" (Fig. A)	1530	CF, 18" (Fig. B)
1520	Plain, 24" (Fig. A)	1540	CF, 24" (Fig. B)
1570	Plain, 34" (Fig. A)	1573	CF, 34" (Fig. B)
1507	FF, 12" (Fig. C)	1548	3PC, 12" (Fig. D)
1550	FF, 18" (Fig. C)	1558	3PC, 18" (Fig. D)
1560	FF, 24" (Fig. C)	1568	3PC, 24" (Fig. D)
1576	FF, 34" (Fig. C)	1588	3PC, 34" (Fig. D)

Corrugated Plastic Folding Splints (Orange)

2500	Plain, 12" (Fig. A)	2504	CF, 12" (Fig. B)
2510	Plain, 18" (Fig. A)	2530	CF, 18" (Fig. B)
2520	Plain, 24" (Fig. A)	2540	CF, 24" (Fig. B)
2570	Plain, 34" (Fig. A)	2573	CF, 34" (Fig. B)
2507	FF, 12" (Fig. C)	2548	3PC, 12" (Fig. D)
2550	FF, 18" (Fig. C)	2558	3PC, 18" (Fig. D)
2560	FF, 24" (Fig. C)	2568	3PC, 24" (Fig. D)
2576	FF, 34" (Fig. C)	2588	3PC, 34" (Fig. D)


## Splints

## Rugged Cardboard Padded Board Splints



- 3 ply, cardboard combines the strength of wood and the light weight of cardboard
- Sealed in green vinyl and splints are padded on one side
- 3½" wide , 1" thick foam

2212	12" Splint
2215	15" Splint
2216	16" Splint
2218	18" Splint
2224	24" Splint
2232	32" Splint
2236	36" Splint



## Locking Leg Splint/Heavy Duty

- Double wall cardboard tapered to provide firm support for entire leg; multiple tie down holes
- Stores flat; assembles to 4½" x 14" x 34" with built-in double locking tabs to secure

1590	Locking Leg Splint
1590-49	Case of 49

## Custom Bags

Morrison Medical now offers custom size/style bags. Some styles can also be screen printed. Please call your distributor or Morrison Medical for a quote.



## IV Armboards

## Disposable IV / Armboards

- Strong double wall corrugated core
- Comfortable foam padding
- Vinyl meets FDA requirements
- Dielectrically sealed
- Washable satin finish

1002-50	1 x 4" case of 50
1003-50	1½ x 3" case of 50
1004-50	2 x 4" case of 50
1005-50	2 x 6" case of 50
1005-100	2 x 6" case of 100
1010-50	2 x 9" case of 50
1010-100	2 x 9" case of 100
1008-50	3 x 6" case of 50
1009-50	3 x 9" case of 50
1009-100	3 x 9" case of 100
1012-50	3 x 12" case of 50
1012-100	3 x 12" case of 100
1016-50	3 x 15" case of 50
1016-100	3 x 15" case of 100
1019-50	3 x 18" case of 50
1019-100	3 x 18" case of 100



## Reusable IV/ Armboards

- Strong plywood core
- Comfortable padding
- Vinyl meets FDA requirements
- Dielectrically sealed
- Washable satin finish
- Premium quality, lightweight and affordable
- "L" shaped (4" x 5") footboards available to prevent infant foot drop, padded on both sides

1032-50	1 x 4" case of 50
1020-50	1½ x 3" case of 50
1031-50	2 x 4" case of 50
1023-50	2 x 6" case of 50
1029-50	2 x 9" case of 50
1039-50	3 x 6" case of 50
1024-50	3 x 9" case of 50
1024-100	3 x 9" case of 100
1026-50	3 x 12" case of 50
1041-50	3 x 15" case of 50
1042-50	3 x 18" case of 50
1042-100	3 x 18" case of 100



#### Toe Glides

- Protection from blisters, calluses or turf toe
- Nearly frictionless space age material
- Worn next to skin under socks
- Can be folded over to fit
- Non-toxic, non-allergenic and flame resistant
- Good for walking, sports, breaking in new shoes
- Used by many professional and college sports teams
- Available in 3 sizes, size by width only  
Regular (narrow), Large (average), Extra Large (wide)

0680-48 Regular (narrow), case of 48

0681-48 Large (average), case of 48

0682-48 Extra Large (wide), case of 48



#### Handy Sandy Sand Bags

- Outer bag is nylon reinforced vinyl
- Dielectrically sealed to prevent leakage or dust
- Covers are color coded nylon reinforced vinyl
- Sewn in weight labels
- Complete with convenient handle

0360 Handy Sandy Sand Bag, 3 lbs., Yellow **YL**, 4½" x 9½"

0360-4 Case of 4 bags, 3 lbs. each

0370 Handy Sandy Sand Bag, 5 lbs., Orange **OR**, 5½" x 12"

0370-4 Case of 4 bags, 5 lbs. each

0380 Handy Sandy Sand Bag, 7 lbs., Red **RD**, 6½" x 15½"

0380-4 Case of 4 bags, 7 lbs. each

0390 Handy Sandy Sand Bag, 10 lbs., Blue **BL**, 6½" x 17"

0390-4 Case of 4 bags, 10 lbs.



#### Umbilical Clamps (A)

- Sterile umbilical clamps, pair

0982 Sterile umbilical clamps, pair



#### Umbi-Cutter (B)

- Clamps and seals the umbilical cord by means of two "V" shaped clamps
- Cuts the umbilical cord between the two clamps
- Sterile and ready to use, disposable

0984 Umbi-Cutter, each

### Survival Bag

- Pocket-size emergency blanket
- Has all the advantages of the space blanket
- It provides total body protection by reflecting heat back while keeping out rain, wind, or snow
- Used by all major life-raft manufacturers as essential emergency personal protection
- Has Icelandic Marine Inspectorate approval
- Useful item in any pack, car, boat, plane or great for recreational use

0980S Side Entry Survival Bag, each

0980T Top Entry Survival Bag, each



### Foil Baby Bunting

- Conserves heat in newborn infants
- One-piece bunting with hood
- Measures 19" x 26"

0981 Foil Baby Bunting, each



## Pillows

### Disposable Foam Pillows

- Foam pillow with comfortable sturdy blue cover
- Small – 12" x 16", Large – 16" x 21"
- Both sizes allow a standard pillowcase to be tucked under mattress to minimize movement

0210 Small Disposable, 12" x 16"

0216 Large Disposable, 16" x 21"



### Staph-Check Foam Reusable Pillows

- Blue Staph-Check Pillows are non-absorbent, anti-bacterial, anti-static & flame retardant
- Easily disinfected
- Durable and long lasting
- Use with or without pillowcases

0218 Small Staph-Check, 4" x 12" x 16"

0220 Large Staph-Check, 4" x 16" x 21"



### Backboard pad

- White Polyurethane Foam
- 52" x 9" x 1/2"
- 6" x 1 3/4" double stick tape at top of foam to secure to backboard
- Provides comfort to the back and spine of patient
- Individually packaged in a 4 1/2" x 5 1/2" bag for easy storage

3040 Backboard Pad



## First Aid Kits



## Personal First Aid Kit

- 12 each Adhesive Bandage Strip, ¾" x 3"
- 4 each Adhesive Bandage Strip, 1" x 3"
- 2 each Adhesive Bandage, Finger Tip, Large Flexible
- 2 each Adhesive Bandage, Knuckle, Flexible
- 2 each Vinyl Exam Gloves, Non-Sterile, Large
- 5 each Soothe-A-Sting, Swab
- 6 each Alcohol Prep Pad, Non-Sterile, 2-Ply
- 3 each Sani-Hands Towelette Wipes
- 1 each Plastic Case
- 3 each Benzalkonium Antiseptic Towelette Wipes
- 1 each First Aid Chart & Illustration
- 2 each First Aid Burn Cream
- 1 each Scissors, For General Use
- 1 each Green Plastic Forceps
- 3 each #2 Safety Pins
- 3 each 3" x 3" Gauze Pad, Sterile, 12-Ply
- 1 each Gauze Roll, Stretch Bandage, Non-Sterile, 2" x 5"
- 1 each ½" x 2" Tape

6500 Personal First Aid Kit, each



## 10 Person First Aid Kit

- 20 each Adhesive Bandage Strip, ¾" x 3"
- 20 each Adhesive Bandage Strip, 1" x 3"
- 10 each Adhesive Bandage, Finger Tip, Small Flexible
- 10 each Adhesive Bandage, Knuckle, Flexible
- 2 each Adhesive Bandage Strip, 2" x 4½" X-Large, Plastic
- 3 each 3" x 3" Gauze Pad, Sterile, 12-Ply
- 1 each Eye Pads, Sterile, 2 1/8" x 2 5/8"
- 3 each First Aid Burn Cream
- 1 each Ammonia Inhalant Pad
- 2 each Vinyl Exam Gloves, Non-Sterile, Large
- 6 each Benzalkonium Antiseptic Towelette Wipes
- 1 each Instant Cold Pack, Disposable
- 1 each Metal Scissors, Blunt, 4"
- 1 each Green Plastic Forceps
- 1 each Acetaminophen (Non-Aspirin), 325mg
- 1 each First Aid Chart & Illustration
- 1 each Plastic With Rubber Gasket #10 Case
- 1 each Abdominal Wound Pad, Sterile, 5" x 9"

6510 10 Person First Aid Kit, each

## First Aid Kits

## 25 Person First Aid Kit

- 3 each Triple Antibiotic
- 1 each First Aid Chart & Illustration
- 1 each Empty #2 or Plastic Container
- 50 each Bandage, Adhesive Strips,  $\frac{3}{4}$ " x 3", Sheer
- 1 each Gauze Roll, Stretch Bandage Non-Sterile 1 x 5"
- 1 each Green Plastic Forceps
- 10 each Povidone-Iodine Pad
- 3 each Bandage, Adhesive Strips, Finger Tips, Small Flexible
- 2 each Gloves, Exam, Latex, Non-Sterile Large
- 1 each Eye Wash, Irrigate, Generic,  $\frac{1}{2}$  oz.
- 3 each Alka-Seltzer® Tabs
- 6 each Aspirin, Generic, Regular Strength, 325mg
- 1 each Bandage, Elastic, (ACE®), Non-Sterile, 2"
- 1 each Gauze Roll, Stretch Bandage Non-Sterile 3 x 5"
- 3 each Gauze Pad, Sterile, 2" x 2", 12 Ply
- 6 each Gauze Pad, Sterile, 3" x 3", 12 Ply
- 1 each Bleached Tape,  $\frac{1}{2}$ " x 5"
- 10 each Bandage, Adhesive Strips,  $\frac{3}{8}$ " x  $\frac{1}{2}$ " Sheer JR
- 3 each Bandage, Adhesive Strips, 2" x 4 $\frac{1}{2}$ " X-L Plastic
- 3 each Bandage, Adhesive Strips, Knuckles, Flexible
- 1 each Cold Pack, Instant, Disp., Small Kit Size 1's in box
- 3 each Ammonia Inhalant Pad
- 2 each Eye Care, Eye Pads Sterile 1 –  $\frac{5}{8}$ " x 2 $\frac{5}{8}$ "
- 1 each Triangle Bandage, Non-Sterile, 40" x 40" x 51"
- 10 each Benzalkonium Antiseptic Wipes Towelette
- 1 each First Aid Burn Cream,  $\frac{1}{2}$  oz.
- 6 each Aspirin, Generic, Regular Strength, 325mg

6520 25 person first aid kit each

6550 50 person first aid kit also available



First Aid Kits available in plastic, metal or fanny pack



First Aid Kits are now available in fanny packs and can also be sold separately.

## Kits



## Universal Clean-Up System (Customized Kits Available)

Complete kit contains:

- 1 each Long sleeve gown
- 1 pair LATEX FREE gloves
- 1 each Mask with shield
- 1 each 4" scoop/scraper
- 1 each Biohazard Bag
- 1 each Disinfectant Cloth
- 1 each BZK Towelette
- 1 each 21 gram package Pick-Up Powder

6005 Universal Clean-Up System, each



## Clean-Up System II (Customized Kits Available)

Complete kit contains:

- 1 pair LATEX FREE Exam gloves
- 1 each Pick-Up Powder
- 1 each Super Scooper spatula
- 1 each EZ Zip bag
- 1 each BZK towelette
- 1 each Sani-Cloth germicidal cleaner
- 3 each Paper towels
- 1 each Trash bag with tie
- 1 each Red disposable Bag with biohazard warning and ties
- 1 each Easy-to-follow instructions
- 1 each 21 gram package Pick-Up Powder

6006 Universal Clean-Up System II, each



## Isolation Kit, Disposable (Customized Kits Available)

Complete kit contains:

- 1 each Long sleeve gown
- 1 pair LATEX FREE Gloves
- 1 each BZK Towelette
- 1 each Biohazard Bag
- 1 each Mask with shield
- 1 pair Shoe covers
- 1 each Bouffant cap

6300 Isolation Kit, each



## Influenza Protection Kit (Customized Kits Available)

Complete kit contains:

- Non-Latex Exam Gloves
- Biohazard Bag (s)
- PDI Super Sani Cloths
- BZK Towelettes
- Shields with Mask
- Instructions for Use

6400 Personal Influenza Protection Kit, each

6401 Basic Influenza Protection Kit, each

6402 Deluxe Influenza Protection Kit, each

6403 Family Influenza Protection Kit, each

6404 Office Influenza Protection Kit, each



## Kits

## Emergency OB Kits (Customized Kits Available)

- Everything necessary for an emergency birth
- All items are disposable
- Available in clear pouch or sturdy cardboard box

## Both kits contain:

- 1 Pair LATEX FREE sterile exam gloves
- 1 Disposable sterile scalpel
- 1 Plastic lined under pad
- 4 Sterile gauze sponges
- 2 Sterile Umbilical clamps
- 1 Plastic placenta bag with ties
- 1 Apgar scoring chart
- 1 OB pad, sterile
- 1 Receiving blanket
- 3 Non woven towels
- 1 Sterile Bulb syringe
- 1 Plastic apron
- 1 OB towelette
- 1 BZK towelette

0965NL OB Kit in clear plastic bag

0967NL OB Kit in cardboard box



## Pick-Up Powder – Featured in Clean-Up System II

- Absorbs 40 times its own weight
- Instantly “gels” body fluid and other liquids
- “Gelled” spills do not run, are easy to remove
- Sprinkle Pick-Up Powder in biohazard disposal bags to retard or prevent leakage
- Sanitizes and deodorizes as it works
- Less than 1% chlorine also assists in odor control

6101-12 21 gram package, case of 12

6101-100 21 gram package, case of 100



If you would like a quote for a customized kit,  
please call your distributor.



800 438 6677 toll free

614 469 9696 fax

www.morrisonmed.com



### TERMS, PRICING & PAYMENTS

NET THIRTY DAYS FOR ALL ESTABLISHED ACCOUNTS. Due to fluctuating costs and continual product improvements, prices are subject to change without notice. Past due accounts are subject to a charge of 1 ½% per month; plus a \$15.00 re-billing charge to cover postage and handling. When you place your order you accept all responsibility for proper and timely payment. Any Past-Due account will be turned over for collection, in which case the purchaser agrees to be responsible for all cost of collection and attorney's fees plus interest. Returned checks are subject to a \$30.00 processing charge plus \$15.00 re-billing charge. Wire transfer charges will be charged to buyer.

### RETURNS – NO CASH REFUNDS

No returns on privately labeled and screen printed products except where error has occurred by fault of Morrison Medical. No merchandise may be returned without prior written approval. Returns must include the Morrison return authorization number, date of original invoice and invoice number. Merchandise held longer than 30 days or customized items are non-refundable. No cash refunds are issued; returns will only be credited to your account. Authorized returns must be in new condition and are subject to inspection before credit issuance and are subject to a 20% restocking fee with minimum of \$20.00

### DROP SHIPMENTS

There is no charge for drop shipments valued over \$100.00. Drop shipments under \$100.00 will be subject to a \$15.00 service/handling charge.

### SPLIT CASES

Split cases are subject to a 10% charge.

### MINIMUM ORDER

Any order less than \$75.00 may be subject to a \$15.00 service/handling charge.

### INTERNATIONAL ORDER INFORMATION

All international customers are required to prepay their orders before they can be processed. Custom fees, local taxes, if applicable, freight, handling and all wire transfer fees are the responsibility of the customer. International orders are shipped complete unless otherwise requested. There is a standard fee of \$45.00 for preparing export documentation (except Canada). Payments for all orders must be made in U.S. Currency.

### FREIGHT POLICY

ALL SHIPMENTS ARE BEST WAY, FOB – COLUMBUS, OHIO unless otherwise specified at time of order. Normally shipments are sent via UPS with all charges added to the invoice. C.O.D. orders are \$10.00 per package. Despite our best efforts, items may occasionally be damaged in transit. Check all shipments for exterior damage or opened cartons BEFORE you sign for the goods, note any damage on the bill of lading or with driver. If the damage is obvious upon receipt, you may refuse the shipment. All claims for concealed loss or damage should be reported to our Customer Service Department immediately. YOU must file claims for shortages or damages directly with the common carrier. When reporting damage be sure to hold all containers and packing materials for inspection. All claims must be filed within 15 days of shipment.

### PRIVATE LABELING/SCREEN PRINTING

We invite inquiries into private labeling and screen printing. Please send specifications, quantities and other requirements to our main office. We will furnish you with a prompt and accurate quote. Minimum orders and lead times vary with each product.

### DISCLAIMERS

All products manufactured and/or sold by Morrison Medical are designed exclusively for use by trained professionals. Morrison Medical will not be liable for damages or injury to persons or property, resulting from the misuse of these products.

This catalog supercedes and nullifies all previous catalogs. We are not responsible for typographical errors in either our catalog or published price list.

Copyright© 2008 Morrison Medical. All rights reserved. Reproduction in whole or part by any means without prior written permission is prohibited.

DATE: \_\_\_\_\_

APPLICANT'S NAME \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

STATE \_\_\_\_\_

ZIP CODE \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

FAX NUMBER \_\_\_\_\_

EMAIL ADDRESS \_\_\_\_\_

 Corporation  Partnership  Proprietorship  Other \_\_\_\_\_**BILLING INFORMATION**

COMPANY NAME \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

FAX NUMBER \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

STATE \_\_\_\_\_

ZIP CODE \_\_\_\_\_

CONTACT PERSON OR DEPARTMENT \_\_\_\_\_

EMAIL ADDRESS \_\_\_\_\_

**SHIPPING INFORMATION**

COMPANY NAME \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

FAX NUMBER \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

STATE \_\_\_\_\_

ZIP CODE \_\_\_\_\_

**COMPANY OFFICERS**

NAME \_\_\_\_\_

TITLE \_\_\_\_\_

NAME \_\_\_\_\_

TITLE \_\_\_\_\_

**TRADE REFERENCES**

NAME \_\_\_\_\_

ADDRESS \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

FAX NUMBER \_\_\_\_\_

NAME \_\_\_\_\_

ADDRESS \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

FAX NUMBER \_\_\_\_\_

NAME \_\_\_\_\_

ADDRESS \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

FAX NUMBER \_\_\_\_\_

**BANK REFERENCE**

NAME \_\_\_\_\_

PHONE NUMBER \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

STATE \_\_\_\_\_

ZIP CODE \_\_\_\_\_

ACCOUNT NUMBER \_\_\_\_\_

ACCOUNT REPRESENTATIVE \_\_\_\_\_



3735 Paragon Drive  
 Columbus, Ohio 43228  
 (614) 461-4400 • Fax (614) 469-9696



Morrison Medical  
3735 Paragon Drive  
Columbus, Ohio 43228

614 461 4400 phone  
614 469 9696 fax

[www.morrisonmed.com](http://www.morrisonmed.com)

800 438 6677 toll free  
877 469 9696 toll free fax

## Product Catalog 2008 • 2009

