

Endo Catch™ Gold Disposable Specimen Pouch



Ethicon
EndoPouch Retriever™*

The image shows the Ethicon EndoPouch Retriever, a metal instrument with a hook at the top. A white specimen pouch is attached to the bottom of the instrument. The pouch is held in place by a thin, blue wire that forms an incomplete ring around the top edge of the pouch. The pouch is partially filled with a clear liquid.

INCOMPLETE RING



Covidien
Endo Catch™ Gold Device

The image shows the Covidien Endo Catch Gold Device, a metal instrument with a hook at the top. A white specimen pouch is attached to the bottom of the instrument. The pouch is held in place by a thick, continuous metal ring that completely encircles the top edge of the pouch. The pouch is partially filled with a clear liquid.

CONTINUOUS RING

IT'S IN THE BAG!

Endo Catch™ Gold Device

Disposable Specimen Pouch

Proprietary Continuous Ring


The bag will remain open and not slide off the ring during specimen insertion

Bag Strength

Polyurethane pouch allows for convenient removal of tissue specimens such as Appendix, Gallbladder, Ovaries. See IFU for details.


34.5 cm Shaft Length

For optimal reach




Bag Material Strength Comparison

Covidien Endo Catch™ Gold Device vs. Ethicon EndoPouch Retriever™



Weld Strength Comparison

Covidien Endo Catch™ Gold Device vs. Ethicon EndoPouch Retriever™



Data on file, August 2002.

“I use the Endo Catch™ Gold device to remove the appendix and gallbladder from the peritoneal cavity to help prevent wound contamination (appendix) and stone spillage (gallbladder).”

K.D.K., MD, General Surgeon

“I like that the continuous ring remains open when fully expanded to make it easier to insert the specimen in the bag. I use it for all my lap chole, appy, splenectomy cases and occasional lap gastric bypass when I have larger bowel specimens to remove.”

B.K.O., MD, General Surgeon

Order Codes	Description
173050G	Endo Catch™ Gold 10 mm Specimen Pouch
173049	Endo Catch™ II 15 mm Specimen Pouch

COVIDIEN, COVIDIEN with logo, COVIDIEN logo and positive results for life are U.S. and internationally registered trademarks of Covidien AG.

™ Trademark of its respective owner. Other brands are trademarks of a Covidien company. © 2012 Covidien 6.12 P120039