



Specimen Collection
1-800-962-9888



COVIDIEN

positive results for life™

Table of Contents

SPECIMEN COLLECTION

PRECISION™ SPECIMEN CONTAINERS

Premium,™ Sterile	4
Precision,™ Sterile & Non-Sterile	5

PRECISION™ STOOL COLLECTION SYSTEM

Stool Collector	6
-----------------------	---

PRECISION™ SPUTUM COLLECTION KITS

Sputum Kit with Tube	6
----------------------------	---

PRECISION™ DISPOSABLE TISSUE GRINDER SYSTEMS

Large and Small Tissue Grinder	6
--------------------------------------	---

PRECISION™ 24 HOUR COLLECTION SYSTEMS

3000mL Capacity	7
3000mL & 3700mL Capacity with Pour Spout	7

PRECISION™ STERILE MIDSTREAM COLLECTION SYSTEMS

Urine Specimen	8
Catheter Urine Specimen	8
Preservative Kits	8

PRECISION™ SOLUTION BOWLS

16oz. and 32oz. Capacity	8
--------------------------------	---

Precision™ Specimen Containers



2200SA



2210SA

Precision™ Premium 1.5 oz. & 4 oz. Specimen Containers - Sterile

- Clear plastic construction with graduations offers easy observation and measurement of contents
- Highly visible “positive seal” indicator clearly shows container is sealed when red bars are lined up
- Exclusive metal lid with plastic lining provides a leak-resistant seal to minimize leakage of specimen during transportation

Code	Size	Description	Ship Case
2200SA	4oz.	Sterile Specimen Container; 4 Boxes of 50 ea. Dimensions: 3 1/2” H x 2” W	200
2205SA	4oz.	Sterile Specimen Container w/out “Positive Seal” Indicator; Bulk Packed Dimensions: 3 1/2” H x 2” W	400
2210SA	1.5oz.	Sterile Specimen Container, Dimensions: 1 3/4” H x 1 3/4” W	200



Precision™ Premium 4 oz. Specimen Wide Mouth Container - Sterile

- Clear plastic construction with graduations offers easy observation and measurement of contents
- Highly visible “positive seal” indicator clearly shows container is sealed when red bars are lined up
- Exclusive metal lid with plastic lining provides a leak-resistant seal to minimize leakage of specimen during transportation

Code	Size	Description	Ship Case
2300SA	4oz.	Sterile Specimen Container with Wide Mouth, Dimensions: 2 3/4” H x 2 1/2” W	200



Precision™ Premium 4 oz. Operating Room Specimen Container - Sterile

- In blister packaging

Code	Size	Description	Ship Case
2600SA	4oz.	O.R. Sterile Specimen Container, in blister packaging, Dimensions: 2 3/4” H x 2 1/2” W	100

Precision™ Specimen Kits & Containers - LATEX FREE



19000 3.5 oz. Graduated Containers

Precision™ Specimen Container

- Complete line of specimen containers
- Allows for sterile specimen collection
- Protects against leakage
- Clearly marked gradients

Code	Size	Description	Ship Case
------	------	-------------	-----------

19000 Series Graduated Containers, Tamper-Evident Seal, Ideal for pneumatic tube transport systems. Dimensions: 3" H x 2" W

19000	3.5oz.	Sterile, Tamper-Evident Seal, 25 per bag	200
19500	3.5oz.	Non-sterile specimen containers with lined metal caps - no label	400



17000 Series 4oz. Graduated Containers

17000 Series 4 oz. Graduated Containers, Polyethylene Screw Caps.

Dimensions: 3 H x 2.5" W

17000	4oz.	Sterile, individually wrapped with screw-on caps - no label	200
17077	4oz.	Sterile, individually wrapped with screw-on caps, attached ID label	200
17099	4oz.	Sterile, O.R. packaged in blister pack	100
SU17500	4oz.	Non-sterile containers with full screw-on polyethylene caps	400
17400	4oz.	Non-sterile full turn clean containers without caps	500
17800	4oz.	Non-sterile full turn screw-on blue caps for 17400 containers	500

14000 Series 6 oz. Graduated Containers, Metal Screw-on Caps.

Dimensions: 3" H x 2 1/2" W

14000	6oz.	Sterile, individually wrapped with metal screw-on leakproof cap and attached ID label	100
14034	6oz.	Sterile, O.R. packaged with screw cap	100
14500	6oz.	Non-sterile containers with metal screw caps	250



14000 Series 6 oz. Graduated Containers

13000 Series 5 oz. Graduated Containers, Snap Caps.

Dimensions: 2 3/4" H x 2.5" W

13000	5oz.	Sterile, individually wrapped with snap caps	200
13400	5oz.	Sterile, O.R. packaged with snap caps	100
12500	5oz.	Non-sterile containers with blue snap caps	400
12400	5oz.	Non-sterile containers without caps	500
12800	5oz.	Non-sterile blue snap caps for 12400 container	500



13000 Series 5 oz. Graduated Containers

Precision™ Stool Collection System



2450SA

Precision™ Stool Collector

- Quick, easy way to collect and transport stool specimen
- Eliminates the need for gathering additional samples
- Large, rigid spoon for a more consistent sample size
- Sampling cap allows for safe, gentle release of gas pressure
- Small sample opening allows for minimal odor and easy access without exposure to full specimen
- Specially threaded twist-on cap prevents loss of specimen minimizing chances of leakage

Code	Description	Ship Case
2450SA	Stool Collector Dimensions: 3 1/2" H x 2" W	100

Precision™ Sputum Collection System



2650SA

Precision™ Sputum Collector Kit

- Easy to use, safe for transport through pneumatic tube systems
- Built-in funnel simplifies the collection of sputum
- Instructions are clearly printed on container for ease of use
- Built-in 50mL graduated centrifuge tube eliminates the need to transfer specimen from container to tube for processing
- Closed centrifuge tube fits in container after use for transport to the lab

Code	Description	Ship Case
2650SA	Sterile Sputum Kit with Tube Dimensions: 5 3/4" H x 2 1/2" W	50

Precision™ Tissue Grinder



3500SA

Precision™ Large & Small Disposable Tissue Grinder Systems

- Ready to use in a sterile peel pouch
- Completely disposable after grinding, no clean-up or sterilization required
- Specimen can be cooled during grinding
- Grinder tube doubles as a specimen container for transport and storage
- Protective sleeve helps to contain any aerosolization created during grinding
- Tissue is ground quickly due to fresh grinding surface

Code	Description	Ship Case
3500SA	Large Tissue Grinder with 50mL Tube	10
3505SA	Small Tissue Grinder with 15mL Tube	10



3505SA

Precision™ 24 Hour Specimen



Precision™ 24 Hour Specimen Collection Container

- 3000mL capacity
- Lined plastic screw cap minimizes leakage
- Easy-to-read graduations

Code	Size	Description	Ship Case
5000SA	3000mL	24 Hour Specimen Collection Container	40




Precision™ 24 Hour Specimen Collection Container

- Convenient pouring spout allows for quick, safe transfer of urine sample
- Self-venting when pouring
- Lined plastic screw cap minimizes leakage
- Easy-to-read graduations

Code	Size	Description	Ship Case
5010SA	3000mL	24 Hour Specimen Collection Container with Pour Spout	40
5011SA	3700mL	24 Hour Specimen Collection Container with Pour Spout	24

Precision™ Sterile Midstream Collection Systems with Preservative - LATEX FREE



7000SA

Precision™ Urine Specimen & Catheter Urine Collection Systems - Sterile

- Leak resistant.
- Effective for transport in pneumatic tube systems
- Maintain specimen integrity for up to 24 hours

Code	Description	Ship Case
7000SA	Urine Specimen System with Preservative 50 in dispenser	200
9500SA	Catheter Urine Collection Kit	100



9500SA



8000SA

Precision™ Midstream Preservative Kit with and without Funnel - Sterile

- Leak resistant
- Effective for transport in pneumatic tube systems
- Maintain specimen integrity for up to 24 hours
- Transparent blister package allows patient to see components and view instructions
- Metal twist cap with unique plastisol lining minimizes leakage of specimen during transportation

Code	Description	Ship Case
8000SA	Midstream Preservative Kit 40 in dispenser	200
9100SA	Midstream Preservative Kit with Funnel	100



9100SA

Kendall™ Solution Bowls - LATEX FREE



Kendall™ Solution Bowls

- A safe, sensible alternative to stainless steel
- Made of medical-grade polypropylene
- Available in two sizes: 16 oz. and 35 oz. capacity
- Sterile and non-sterile
- Ideal for operating/emergency rooms and general use
- Disposable

Code	Size	Description	Ship Case
61000	16oz.	Plastic Solution Bowl, Individually Sterile Packed	75
61200	35oz.	Plastic Solution Bowl, Individually Sterile Packed	50
65000	16oz.	Plastic Solution Bowl, Bulk Packed, Non-Sterile	500
65600	35oz.	Plastic Solution Bowl, Bulk Packed, Non-Sterile	250
61299	35oz	Plastic Solution Bowl, Double Wrapped Package, Sterile	36

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
12400	13000 Series 5 oz., Non-sterile	5	2450SA	Precision™ Stool Collector	6
12500	13000 Series 5 oz., Non-sterile	5	2600SA	Precision™ Premium 4 oz. Operating Room Specimen Container	4
12800	13000 Series 5 oz., Non-sterile	5	2650SA	Precision™ Sterile Sputum Kit with Tube	6
13000	13000 Series 5 oz., Sterile	5	3500SA	Precision™ Large Tissue Grinder with 50mL Tube	6
13400	13000 Series 5 oz., Sterile	5	3505SA	Precision™ Small Tissue Grinder with 15mL Tube	6
14000	14000 Series 6 oz. Graduated Container, Sterile	5	5000SA	Precision™ 3000mL 24 Hour Specimen Collection Container	7
14034	14000 Series 6 oz. Graduated Container, Sterile	5	5010SA	Precision™ 3000mL 24 Hour Specimen Collection Container	7
14500	14000 Series 6 oz. Graduated Container, Non-sterile	5	5011SA	Precision™ 3700mL 24 Hour Specimen Collection Container	7
17000	17000 Series 4 oz. Graduated Container, Sterile	5	61000	Kendall™ 16 oz. Plastic Solution Bowl	9
17077	17000 Series 4 oz. Graduated Container, Sterile	5	61200	Kendall™ 35 oz. Plastic Solution Bowl	9
17099	17000 Series 4 oz. Graduated Container, Sterile, O.R. packaged	5	61299	Kendall™ 35 oz. Plastic Solution Bowl	9
17400	17000 Series 4 oz. Graduated Container, Non-sterile	5	65000	Kendall™ 16 oz. Plastic Solution Bowl	9
17800	17000 Series 4 oz. Graduated Container, Non-sterile	5	65600	Kendall™ 35 oz. Plastic Solution Bowl	9
19000	19000 Series 3.5 oz. Graduated Container, Tamper-Evident Seal, Sterile	5	7000SA	Urine Specimen System with Preservative	8
19500	19000 Series 3.5 oz Graduated Container, Non-sterile	5	8000SA	Precision™ Midstream Preservative Kit	8
2200SA	4 oz. Precision™ Premium Specimen Container	4	9100SA	Precision™ Midstream Preservative Kit with Funnel	8
2205SA	4 oz. Precision™ Premium Specimen Container	4	9500SA	Catheter Urine Collection Kit	8
2210SA	1.5 oz. Precision™ Premium Specimen Container	4	SU17500	17000 Series 4 oz. Graduated Container, Non-sterile	5
2300SA	Precision™ Premium 4 oz. Specimen Wide Mouth Container	4			

NOTES



COVIDIEN, COVIDIEN with logo and TM marked brands are trademarks of Covidien AG or an affiliate.
© 2008 Covidien AG or affiliates. All rights reserved.

H-6200SCP XM 11/08

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

800-962-9888
508-261-8000

WWW.COVIDIEN.COM