

Medidyne MAKING A DIFFERENCE



FAST PACKS FOR MEDICS



# COMPANY

## FAST PACKS FOR MEDICS

And we're getting faster. Innovating new specialized packs and modules. Expanding mission-specific product lines, for professionals from urban crews to tactical medics. Refining our right-size philosophy so you get a pack that fits. And earning recognition from the industry press.

StatPacks builds a range of high-performance packs and modules. Every one designed to work with your body. Advanced ergonomics. Smart weight distribution. Comfortable details to limit wear and tear. And hands-free designs to keep your body agile.

Our modular system lets you build a set-up to suit your mission, your body, your needs. When the call requires, add specialized extras or switch out internal modules. Mesh and urethane pockets let you identify contents at a glance. Modular organization — click, zip, roll — makes good medics faster.





*"I am using the "Back Up" for My Local Volunteer Fire Department as well as Paid EMS. We are a rural service, therefore I need a product that I can carry a wide range of equipment. I have used several "No-name brand" bags and none of them have proven useful for what I need. I have found this bag to be an excellent product. It is essential in the care of our patients prior to arrival at our unit."*

*~ Tyler K., Basic EMT from Palmer, TN*

*"I can truly appreciate my pack when the days are long and calls are far. Being able to have good back support and easy access make my job so much easier, plus the quality of material is leaps and bounds from other bags I've used."*

*- John M., from Winter Park, FL*

*"Firstly let me congratulate you on a fantastic product. I have been searching high and low for a pack which not only provides a robust and efficient level of service for me to conduct my EMS activity, but also looks the part. I have finally found it and now widely promote it and use it within our organization. Thanks guys. We use the gear all the time. Great stuff."*

*~ Bill B., Paramedic from Ivanhoe, Australia*

*"I use my code 7 for everyday use on the ambulance. I love it and all the room. The comfort is great. I can have it loaded with gear and a laptop and it's not digging into my shoulders or back. I wouldn't change a thing. Keep up the great work!"*

*~ Alex D., Basic EMT from Merced, CA*

*"I have used StatPacks in the air medical environment for several years. We recently chose the Med Slinger for use as our anesthesia code bag and we are quite*

*impressed. We've found the pack to be durable and well organized. It improves our ability to respond to crisis situations within our facility. Keep up the great work."*

*~ B.J., Paramedic from Wyoming, PA*

*"The Perfusion is used on a daily basis as a primary response pack on a moderately busy paramedic unit. The pack fits snugly to anyone's back. It also has plenty of room to store everything you need in a primary bag. We have used other backpack style bags, but the Perfusion is top of the line!"*

*~ Trevor K., Paramedic from Lenoir, NC*

*"I am using the "Perfusion" primarily for volunteering in an EMS agency, and working with a 911 agency. This bag is one of the easiest I have ever found to carry on my back, and even on my shoulders. If I have it on my back, my partner is able to work off of me as I am tending to the patient (or even at times, the other patient). I find that the quick zip access to the main compartment is my favorite feature, because everything that I need (airway control, IV access, trauma dressings) is in one compartment. This is one of the best packs I have ever owned."*

*~ Danielle S., Paramedic from Glen Head, NY*

*"We are taking in only what we need – not the entire ambulance with your Tidal Volume. We use it for the BLS first in kit. The pack is stocked with what you might need in the first 5 min of a call. We will continue to purchase your products."*

*~ Mike L., Paramedic from Elma, WA*

*"I made the Whiner into a small trauma bag for when I have to get to a site in a hurry. There are no belt packs to compare. I love this pack!!"*

*~ Cynthia W., First Responder from Helena, MT*

*"The Traverse is an impressive invention. It is very sturdy, very practical. I've snagged it on tree limbs and I've even managed to close it in a vehicle door, yet there is only some cosmetic damage, all the wear has had no effect on the actual structure of the bag. What a wonderful design."*

*~ Mathew C., Basic EMT from Medford, OR*

*"Hi, my name is Dr. G a missionary doctor working in Costa Rica but also an Emergency Doctor for the National Emergencies Commission, CNE, and a Fire Department Doctor. I have many years using the StatPacks Kits, In the city, in the mountains, in a helicopter emergency flight or in a disaster area, the design of the kit lets me keep everything handy, without the worry of having a different kit for a different situation. I have Pigmentary Retinitis disease – a chronic and degenerative disease of the retina, so I am losing my vision a little every day, StatPacks help me because the design lets me keep each tool or instrument in a specific place so when I need something I already know where to find it and most of the time I don't have to use my eyes to find what I am looking for. Thanks StatPacks for offering the world such smart designs."*

*~ Dr. G. MD from Costa Rica*

## TABLE OF CONTENTS

**4 G2 EMS**

**10 G1 EMS**

**18 SPORT**

**20 OFF-DUTY**

**22 FIRE/SAR**

**24 TACTICAL**


# G2 OVERVIEW

## APPLIED LESSONS IN ADVANCE


The second generation of StatPacks is born from the proven products we have made for years. The StatPacks G2 line embodies a refined approach to our core values that all focus on making you faster and better at saving lives.

### CELL SPACING : THE BUILDING BLOCKS OF ORGANIZATION


Each product is built around a 6" x 12" x 6" Cell. Each Cell is specialized to hold supplies for a specific purpose. The G2 packs are designed to house a given number of cells. The I.V. cell is half the size of a full Cell, The O2 Cell is twice the size of a full cell. All of them work together to provide a customizable system that is orderly and effective at keeping supplies and tools quickly at hand.



TECHNICIAN 6 CELL



RESPONDER 4 CELL



CLINICIAN 3 CELL



OXYGEN CELL



MED CELL, AIR CELL  
& UNIVERSAL CELLS



I.V. CELL

### INTERNAL ANATOMY: THE G2 APPROACH TO SPEED

Each Cell sets the stage for clear thinking and quick reaction. They are designed to reveal your supplies in a way unique to StatPacks – intuitively grouped and securely protected. G2 Cells unfold to reveal all your vital supplies at once, and puts away just as quickly.





# ADVANCED FASTOLOGY

## DEFENSE MECHANISM: EVOLVED DURABILITY & SAFETY

A no-holds-barred approach to durability features reinforced wear areas, over-sized zippers, and an overall heightened standard of quality right down to the last stitch. Foam molded back-panels provide comfortable, balanced transport while resisting wear. Additionally, each G2 product is highly reflective to keep you and your gear safe at night.



# G2 EMS CELLS



## G2 AIRWAY CELL # G24006GN

Foam molded quad-fold design helps keep basic and advanced airway adjuncts organized and at your fingertips. Laryngoscope blades, OPA's, and NPA's are neatly arranged for fast access. Endotracheal tubes, suction catheters, stylets, and ETCO2 meters can be organized in the middle section. Considered to be one full cell in the Cellular System.

- Unique foam molded design keeps medications protected
- Innovative fold-out design is compact, yet quickly reveals well-organized airway supplies and tools
- Specifically designed to stack neatly in StatPacks G2 EMS packs and also fits in the G1 Perfusion and Backup
- Color-coded for easy identification
- Two handles with "AIRWAY" written to identify
- ID window to label any way you want
- Available color: Green

| CAPACITY: 470 in<sup>3</sup> | WEIGHT: 1.5 lbs | DIMENSIONS: H: 12" x W: 6" x D: 6.5" |



## G2 MEDICINE CELL # G24005RE

Foam molded quad-fold design helps keep medications organized and at your fingertips. IV catheters, needles, syringes are neatly arranged for simple identification. The "Vial Strand" holds multi dose meds tightly in place with pull tabs for easy removal. Prefilled doses can be organized in the middle section. Considered to be one full cell in the Cellular System.

- Unique foam molded design keeps medications protected
- Innovative fold-out design is compact, yet quickly reveals well-organized medication supplies and tools
- Specifically designed to stack neatly in StatPacks G2 EMS packs and also fits in the G1 Perfusion and Backup
- Color-coded for easy identification
- Two handles with "MEDICINE" written to identify
- ID window to label any way you want
- 4 Vial Strands that holds multi dose vials securely in place
- Accommodates a Sharps Shuttle™
- Available color: Red

| CAPACITY: 470 in<sup>3</sup> | WEIGHT: 1.5 lbs | DIMENSIONS: H: 12" x W: 6" x D: 6.5" |



## VIAL STRAND 4 PACK # G25000BK

An ingenious way to securely hold drug vials while keeping them quickly identifiable and easily retrievable. Military spec Velcro hook on back attaches to any compartment with Velcro loop. Each strand holds 5 vials of any size.

| DIMENSIONS: H: 1.5" x W: 5" |



## G2 I.V. CELL # G24004BU

Foam molded book style design helps keep IV therapies organized and at your fingertips. Drip sets, IV bag(s), catheters, needles and Start Kits are neatly arranged for fast access. It doubles as a pressure bag for rapid fluid infusion. Considered to be a half cell in the Cellular System.

- Unique foam molded design keeps contents protected
- Specifically designed to stack neatly in StatPacks G2 EMS packs and also fits in the G1 EMS packs
- Two handles with "Intravenous" written to identify
- Clear plastic window to identify IV solution
- Area for pressure gauge for rapid fluid infusion
- Hot/Cold packs pocket to warm or cool fluid(s)
- Accommodates a Sharps Shuttle™
- Color-coded for easy identification
- Available color: Blue

| CAPACITY: 210 in<sup>3</sup> | WEIGHT: 0.6 lbs | DIMENSIONS: H: 12" x W: 6" x D: 3" |



## G2 UNIVERSAL CELL # G24008BK

Keep your random supplies together and quickly accessible. Foam molded design helps keep inside contents organized and at your fingertips.

- Unique foam molded design keeps contents protected
- Clear plastic window to identify inside contents
- Specifically designed to stack neatly in StatPacks G2 EMS packs and also fits in the G1 Perfusion and Backup
- Perfect for keeping random supplies together and organized holds a BVM nicely
- Color-coded for easy identification
- Two handles for easy grab
- Available color: Black

| CAPACITY: 470 in<sup>3</sup> | WEIGHT: 1.5 lbs | DIMENSIONS: H: 12" x W: 6" x D: 6.5" |



## I.D. CARD # G25001YE

This high-vis card folds to fit in the I.D. window of any StatPack. Use the pre-printed Pediatric and Adult tags or make your own using the reverse side. Made from recyclable plastic and will hold up to water and endless use.



# G2 EMS PACKS



## G2 TECHNICIAN 6 CELL

# G2400RE

*This simple, box-style carry-all pack is stackable, with shoulder straps for safe, hands-free hauling. This unique design keeps straps out of the dirt and muck when you're accessing contents. Internal Cells stack neatly and securely. Accommodates the QuickRoll Intubation Module, IV Module, and small Drug Module, as well as the stackable Cellular System. The removable Roller Cart allows rolling assisted transport.*



- Unique design keeps shoulder straps off the ground and out of the dirt when accessing gear
- Specifically designed for the StatPacks Cell System and accommodates all StatPacks modules in many different configurations
- Skids protect high wear areas
- Tarpaulin bottom panel for protection from abrasion
- Reflective taping for extra nighttime safety
- Removable Roller Cart
- Internal Velcro lining strips to secure modules
- LED light Attachment point (lights sold separately)
- Available Color: Red

| CAPACITY: 6500 in<sup>3</sup> | WEIGHT: 14.0 lbs | DIMENSIONS: H: 30" x W: 24" x D: 9" |



## G2 RESPONDER 4 CELL

# G2401RE

*A medium-sized street pack for medics who carry less. Boxy and compact, with ergonomic hauling. Unique design keeps your shoulder straps off the ground when you're accessing contents. Accommodates the QuickRoll Intubation Module, IV Module, and small Drug Module, as well as the stackable Cellular System.*



- Unique design keeps shoulder straps off the ground and out of the dirt when accessing gear
- Specifically designed for the StatPacks Cell System and accommodates all StatPacks modules in many different configurations
- Skids protect high wear areas
- Tarpaulin bottom panel for protection from abrasion
- Reflective taping for extra nighttime safety
- Internal Velcro lining strips to secure modules
- LED light Attachment point (lights sold separately)
- Available Color: Red

| CAPACITY: 3600 in<sup>3</sup> | WEIGHT: 6.0 lbs | DIMENSIONS: H: 29" x W: 18" x D: 7" |



## G2 CLINICIAN 3 CELL

# G24002RE



Designed as a quick first in pack for Medics who only need the basic ABC supplies. Accommodates the QuickRoll Intubation Module, IV Module, and small Drug Module, as well as the stackable Cellular System.

- Specifically designed for the StatPacks Cell System and accommodates all StatPacks modules in many different configurations
- Tarpaulin bottom panel for protection from abrasion
- Reflective taping for extra nighttime safety
- Internal removable divider shelf to separate Cells
- LED light Attachment point (lights sold separately)
- Available Color: Red



| CAPACITY: 3000 in<sup>3</sup> | WEIGHT: 3.0 lbs | DIMENSIONS: H: 22" x W: 19" x D: 7" |







## G2 TIDAL VOLUME # G24003RE

The best way to carry oxygen. Securely hold a regular or Jumbo D size oxygen cylinder and breathing equipment. Carry your oxygen on your back, or hook it to a cot. Opens quickly to reveal breathing equipment.

- Foam-lined construction maintains pack shape and protects contents
- Outside accessible urethane window to see regulator and access tank
- Main compartment securely holds any regular D or Jumbo D size cylinder
- Secure Velcro/Buckle attachments hook to a cot
- G1 Oxygen Module included and Velcro's securely in place
- Padded shoulder straps double as shoulder sling, offering hands-free comfort during transport
- Load handles for quick/easy carry
- O<sub>2</sub> wrench lanyard
- Reflective taping and tie-down loops for nighttime safety
- LED light Attachment point (lights sold separately)
- Available Color: Red

| CAPACITY: 1500 in<sup>3</sup> | WEIGHT: 3.5 lbs | DIMENSIONS: H: 22" x W: 8" x D: 9" |



## G2 OXYGEN CELL # G24007BK

Built to accommodate most D and Jumbo D tanks. Considered to be a Double Cell In the Cellular System.

- Foam-lined construction maintains rigidity
- Two adjustable locking cam buckles with strong hypalon material hold the tank firmly to the frame
- Hook Velcro on back panel attaches securely in packs
- Top pocket holds tubing, Nasal cannula's, and masks
- Secure Velcro/Buckle attachments hook to a cot
- Load handles for quick/easy carry
- O<sub>2</sub> wrench lanyard
- Available Color: Normal Black

| CAPACITY: 250 in<sup>3</sup> | WEIGHT: 1.0 lbs | DIMENSIONS\*: H: 20" x W: 6" |  
\* FITS "D" & "JUMBO D" CYLINDERS



## G2 MCI VEST/PACK # G24009YE

Keeps personal and medical equipment close at hand. Designed for MCI's and medics who require both a pack and a highway safety vest. One adjustable size fits everyone on staff comfortably. The Vest/Pack is also an ergonomic alternative for medics who cannot (or choose not to) wear a backpack.

- One-size-fits-all with adjustments on the shoulders and sides
- Two main organizer pockets for easy access to equipment
- Radio/cell phone webbing attachment points
- Pen and scissor holders
- Removable name plates on front and back for identification of personnel
- Front/Back reflective taping for extra night time safety
- LED light Attachment point (lights sold separately)
- ANSI material for daytime safety (meets ANSI 207 / ISEA 2006 requirements)
- Available Color: ANSI Yellow



| CAPACITY: 150 in<sup>3</sup> (Front Pockets) | WEIGHT: 1.0 lbs | DIMENSIONS: H: 21" x W: 17" x D: 8" | (Adjustable SM - XXL)

## REPLACEMENT NAME PLATE

| EMS # G25002E | FIRE # G25003F | POLICE # G25004P |

Quickly and easily change the identification panel with EMS, FIRE, and POLICE nameplates held in place with sturdy hook and loop strips. Call for custom options.





# G1 OVERVIEW

## APPLIED LESSONS IN AD

*The first generation of StatPacks gets better. The same great products that have earned your rave reviews over the past 6 years also responds to feedback from the most demanding and abusive crews in the field. Fast is what doesn't slow you down. When it comes to durability, there is nothing normal about the wear and tear a StatPack goes through – so we've made the G1 products abnormally strong. Our message is the same: Two straps are better than one to keep you balanced, hands free, and safe. Speed is king and effective management of your supplies keeps you thinking clearly and responding quickly even under the most stressful circumstances.*



### WEAR POINTS:

- **ZIPPERS**

Beefy, #10 YKK zippers are used throughout; Custom zipper-pulls keep stress off the chain; Hypalon "garages" keep the wear areas of zippers protected; Reinforced stitching resists blow-out.

- **SUPER STRONG FABRICS**

Nearly indestructible coated EPO fabric used in areas of most wear.

- **IMPROVED SIDE POCKETS**

Newly engineered side pockets hold more, open easier, and look better than ever before.

- **I.D. WINDOWS**

Make it easy to keep your gear with the right rig, and identify different set-ups you create.

- **REFLECTIVE**

A standard, 2" strip of reflective on each pack keep you and your pack more visible at night



# ADVANCED FASTOLOGY



# G1 MODULES

## G1 QUICKROLL INTUBATION MODULE

# G11028RE / G11028TK

Designed to organize basic and advanced airway equipment. The single-buckle closure allows fast, one handed access. The perfect accessory for any StatPacks® pack.

- Clear internal durable urethane pockets allow quick identification of contents
- Organizes laryngoscope, blades, as well as intubation tubes, syringes and forceps
- Holds OPA's / NPA's and gloves
- Single-buckle closure opens quickly, revealing all airway control tools
- Fits nicely in all StatPacks® packs
- Available Colors: Red and Tactical Black (See page 28)

| CAPACITY: 400 in<sup>3</sup> | WEIGHT: 11 oz. | DIMENSIONS: H: 8" x W: 10" x D: 5" (Rolled) |



## G1 DRUG MODULE

# G11027RE / G11027TK

Compact and organized, this module brings all your meds together while offering a quick way for central supply or pharmacies to re-stock. Lockable kiss-zippers help you keep track of used modules. Save time by checking in used modules for restocked sealed ones.

- Transparent durable urethane outer window for quick medication identification
- External pocket window for medical drug listing or unit identification
- Lockable kiss-zippers for pharmacy seal
- Clear durable urethane pockets accommodate extra needles and syringes
- Attaches to the inside of larger StatPacks products with hook-&-loop panel
- Available Colors: Red and Tactical Black (See page 28)

| CAPACITY: 240 in<sup>3</sup> | WEIGHT: 9 oz. | DIMENSIONS: H: 6" x W: 8" x D: 5" |



## G1 IV MODULE

# G11026RE / G11026TK

Holds everything you need for a quick IV setup in one convenient module. The IV Module makes it easy to share on multi victim accident scenes.

- Accommodates all IV bag sizes 1 L, 500 mL, 250 mL
- Keeps IV Start Pack™ ready for use
- Front and inside pockets organize IV tubing and catheters together with the IV bag
- Transparent durable urethane window for quick solution identification
- Single-buckle closure securely holds contents in place and allows quick, one-handed access
- Top loop hangs the entire module with the IV solution
- Adjustable cam-lock buckle allows you to pressurize contents for rapid fluid infusion
- Available Colors: Red and Tactical Black (See page 28)

| CAPACITY: 210 in<sup>3</sup> | WEIGHT: 6 oz. | DIMENSIONS: H: 10" x W: 6" x D: 3.5" (Rolled) |





## G1 SMALL UNIVERSAL MODULE

# G11080BK

- Transparent durable urethane outer window for quick equipment identification
- Foam lined for equipment protection
- Stackable in Technician, Responder, Clinician, BackUp and Stealth packs
- Can be used as a standalone module
- Available Color: Black

| CAPACITY: 300 in<sup>3</sup> | WEIGHT: 7 oz | DIMENSIONS: H: 6" x W: 9.5" x D: 5.5"



## G1 FOLDAWAY PHARMACY MODULE

# G11082BK

- Mesh pockets that fit pre-filled syringes and vials
- Foam lined for equipment protection
- "MED" printed on top flap for identification
- Stackable in Technician, Responder, Clinician, BackUp and Stealth packs
- Can be used as a standalone module
- Available Color: Black with Red Handle

| CAPACITY: 500 in<sup>3</sup> | WEIGHT: 1 lb | DIMENSIONS: H: 6" x W: 8" x D: 10.5"



## G1 FOLDAWAY AIRWAY MODULE

# G11083BK

- Mesh pockets that fit laryngoscopes, NPA's, OPA's, and ET Tubes
- Foam lined for equipment protection
- "AIR" printed on top flap for identification
- Stackable in Technician, Responder, Clinician, BackUp and Stealth packs
- Can be used as a standalone module
- Available Color: Black with Green Handle

| CAPACITY: 500 in<sup>3</sup> | WEIGHT: 1 lb | DIMENSIONS: H: 6" x W: 8" x D: 10.5"



## G1 OXYGEN MODULE

# G11084GR

- Locking buckles for secure fit with tank
- Fits D or Jumbo-D size O2 Cylinder
- Aggressive Velcro on back panel for secure attachment in main packs
- Secure Velcro/Buckle attachments hook to a cot
- Load handles for quick/easy carry
- Fits in Technician, Responder, Clinician, Tidal Volume, BackUp and Stealth packs
- Available Color: Gray

| CAPACITY: D or Jumbo-D cylinder | WEIGHT: 12 oz | DIMENSIONS: H: 16.5" x W: 5" x D: 1" |



# G1 PACKS

## G1 GOLDEN HOUR # G11021RE / G11021TK

The perfect pack for use in close-quarter areas, such as helicopters. Compact and comfortable, built to organize and carry essential equipment. The unique wide-mouth main compartment opening allows easy access to contents and provides convenient, upright availability from the captain's chair. Ergonomic shoulder straps and a breathable padded rear panel.

- Foam-lined construction maintains pack shape and stays upright for efficient access to contents
- Clear, durable urethane window on top inside pouch allows quick view of contents
- Padded, ergonomic shoulder straps for maximum comfort
- Removable adjustable inside shelf for customized organization
- Designed to carry StatPacks® drug module, intubation module, IV module and accessories
- Large side pockets allow quick access to vital sign tools independent of the main compartment
- Available colors: Red and Tactical Black (see page 26)

| CAPACITY: 2400 in<sup>3</sup> | WEIGHT: 3.1 lbs | DIMENSIONS: H: 18" x W: 17" x D: 8" |

## G1 LOAD N' GO # G11022RE / G11022TK

The ideal pack for fixed wing and critical care transport, ambulance and fire engines. No pack keeps equipment more organized and quickly accessible thanks to thoughtfully designed transparent sleeves, mesh pockets, and elastic holsters that are arranged in three layers.

- Foam-lined construction maintains pack shape and protects contents
- Clear, durable urethane windows, mesh pockets, and elastic holsters keep contents in place and well-organized
- Ergonomic, padded shoulder straps and back panel for comfort during transport
- QuickZip access to front compartment
- Large side pockets allow quick access to vital sign tools
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, and IV module
- LED light Attachment point (lights sold separately)
- Available in Blood Borne Pathogen Resistant material (red)
- Available Colors: Red and Tactical Black (see page 26)

| CAPACITY: 2600 in<sup>3</sup> | WEIGHT: 4.5 lbs | DIMENSIONS: H: 20" x W: 19" x D: 7" |



OPTIONAL  
BBP RESISTANT  
MATERIAL



# G1 PERFUSION

# G11034RE / G11034BK

Great for all EMS Clinicians. Two main compartments and three external pockets allow superior organization of large amounts of BLS and ALS equipment. The hybrid design can be worn as a back pack or side sling, or carried by hand.

- Foam-lined construction maintains pack shape and protects contents
- Fully adjustable hook-&-loop shelving can be custom-configured to individual preference and equipment needs
- Stow-away padded shoulder straps offer back-pack and side-carry options
- QuickZip access to main compartment
- Large side pockets allow quick access to vital sign tools, IV and airway equipment
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, and IV module
- LED light Attachment point (lights sold separately)
- Available in Blood Borne Pathogen Resistant material (red)
- Available Colors: Red and Black

| CAPACITY: 3000 in<sup>3</sup> | WEIGHT: 5.7 lbs | DIMENSIONS: H: 22" x W: 17" x D: 8" |



OPTIONAL  
BBP RESISTANT  
MATERIAL

# G1 BACKUP

# G11023RE / G11023TK

Exceptional in urban or rural areas where "stay and play" medical protocol is used. Customizable shelving allows superior organization of virtually any combination of equipment and supplies. The advanced harness system helps you tackle daunting flights of stairs and hard-to-reach accident scenes safely and with free hands.

- Foam-lined construction maintains pack shape and protects inside contents
- Fully adjustable Velcro shelving can be custom configured to individual preference
- Ergonomic, padded shoulder straps and back panel offer comfort during transport
- Removable padded waistbelt
- QuickZip access to main compartment
- Large side pockets allow quick access to vital sign tools
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, and IV module
- LED light Attachment point (lights sold separately)
- Available colors: Red and Tactical Black (See page 27)

| CAPACITY: 3800 in<sup>3</sup> | WEIGHT: 6.0 lbs | DIMENSIONS: H: 25" x W: 18" x D: 8.5" |



# SPECIALIZED PACKS

## JOULE AED # G11040RE

Any Size of Automated External Defibrillators will fit. Open the main compartment and the pack will lie open for access to the AED and supplies at the same time. Great for any location that requires an on-hand AED and supplies.

- Foam-lined construction maintains pack shape and protects contents
- Clear durable plastic window allows quick view of AED status
- Main compartment features many mesh pockets that give visibility and organization to contents
- Ergonomic, padded shoulder straps and back panel for comfort during transport
- QuickZip access to main compartment for fast access to contents and AED
- Designed to carry StatPacks® drug module, intubation module, and IV module
- Available Color: Red

| CAPACITY: 1700 in<sup>3</sup> | WEIGHT: 2.5 lbs | DIMENSIONS: H: 18" x W: 13.5" x D: 7" |

## QUICKLOOK AED # G11039RE

Built to accommodate most Automated External Defibrillators for extremely prompt access. Perfect for use in public areas such as airports, coliseums, shopping malls, or homes as an all-in-one first responder bag. Simple design for the layperson as well as experienced clinicians.

- Foam-lined construction maintains pack shape and protects contents
- Single-buckle AED pocket offers quick access to the Defibrillator
- Clear durable urethane window allows quick view of AED status
- Main compartment features three removable stick-in-place pouches that give visibility and organization to contents
- Ergonomic, padded shoulder straps and back panel for comfort during transport
- QuickZip access to main compartment through the back panel, maintaining access to contents and AED
- Large side pockets allow quick access to vital sign tools
- Designed to carry StatPacks® drug module, intubation module, and IV module
- Available Color: Red

| CAPACITY: 1800 in<sup>3</sup> | WEIGHT: 3.2 lbs | DIMENSIONS: H: 18" x W: 14" x D: 7.5" |





## G1 BREATHER # G12024RE

An all inclusive pack for airway control and oxygen therapy. The Breather is designed to hold the O2 module with a "D" or Jumbo "D" size tank. It fits many basic and advanced airway adjuncts and oxygen delivery devices. The Breather is a versatile pack that can be used in many different departments and setting, meeting local protocol for oxygen and airway requirements.

- Foam-lined construction maintains pack shape and protects contents
- Clear, durable urethane windows, mesh pockets, and elastic holsters keep contents in place and well-organized
- Ergonomic, padded shoulder straps and back panel for comfort during transport
- QuickZip access to front compartment
- Large side pockets allow quick access to vital sign tools
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, IV module, and O2 Module
- LED light Attachment point (lights sold separately)
- Available Color: Red

| CAPACITY: 2600 in<sup>3</sup> | WEIGHT: 4.5 lbs | DIMENSIONS: H: 20" x W: 19" x D: 7" |



## G1 BOLUS # G12025RE

A Medication pack that works in cooperation with the Breather. It holds a large supply of oral and intravenous medication for basic and advanced life support. Medication vials and pre-filled syringes are easily organized as well as IV set-ups and fluids. The Bolus is a versatile pack that can be used in many different departments and settings, meeting local protocol for medication administration.

- Foam-lined construction maintains pack shape and protects contents
- Clear, durable urethane windows, mesh pockets, and elastic holsters keep contents in place and well-organized
- Ergonomic, padded shoulder straps and back panel for comfort during transport
- QuickZip access to front compartment
- Removable center divider that holds vials, syringes, and IV therapy equipment
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, and IV module
- LED light Attachment point (lights sold separately)
- Available Color: Red

| CAPACITY: 1800 in<sup>3</sup> | WEIGHT: 5.0 lbs | DIMENSIONS: H: 16" x W: 15.5" x D: 7.5" |





# SPORT



## MED SLINGER

# G11036RE / G11036BK

A soft sling pack that can work as a personal gear bag or a medical equipment pack. Designed for medics who do not like to wear double shoulder-strap packs. Still ergonomic, it slings quickly off and on for easy wearing. Plenty of transparent and mesh pockets keep you organized and efficient.

- Single, ergonomic shoulder strap and padded back make the pack comfortable to wear
- Sling design allows wearer to swing the pack to the front for quick access to contents without removing pack
- Quick-release adjustable buckle on the shoulder strap provides a custom fit and quick in and out.
- Two main compartments are efficiently organized with transparent and mesh pockets for quick identification and access to contents
- Available Colors: Red and Black

| CAPACITY: 1400 in<sup>3</sup> | WEIGHT: 2.0 lbs | DIMENSIONS: H: 19" x W: 13" x D: 6" 14 |



## COMPETITOR

# G11031RE / G11031BK

Our largest waist pack. The low profile design helps the pack to ride tight in your lower back. This makes for a comfortable manner to carry your supplies. Athletic trainers to Paramedics will be able to carry the basics with no problem.

- Two removable accessory belt pockets provide quick access to supplies
- Main body of pack features two zippered compartments and two side pockets to organize gear and supplies
- Padded waist belt for hands-free convenience and comfort
- Available Colors: Red and Normal Black

| CAPACITY: 900 in<sup>3</sup> | WEIGHT: 1.5 lbs | DIMENSIONS: H: 7" x W: 18" x D: 7" |





# ELEVATE

# G11030RE / G11030BK

A mid-sized waist pack designed to carry first responder gear comfortably and efficiently. Great for Sports Trainers, Ski Patrol, or First Responders who need to carry just enough equipment for fast patient stabilization. Also handy as a personal pack.

- Two removable accessory belt pockets provide quick access to supplies
- Large main compartment accommodates plenty of first responder supplies
- Easy to wear in back for out-of-the-way comfort, or in front for quick access
- Buckled gear loop gives fast access to tape
- Padded waist belt for hands-free convenience and comfort
- Available Colors: Red and Black

| CAPACITY: 400 in<sup>3</sup> | WEIGHT: 1.2 lbs | DIMENSIONS: H: 8" x W: 14" x D: 3" |



# WHINER

# G11029RE / G11029BK

A lightweight waist pack specifically designed to organize provisions needed for basic complaints. Makes it easy to assist with minor cuts, bruises, aches and pains. Great for Lifeguards, Beach, Special Event, and Bike Medics. Also handy for carrying personal items, i.e. pen light, gloves, wallet, or keys.

- Two removable accessory belt pockets provide quick access to supplies
- Main body of pack features two zippered compartments to organize gear and supplies
- Padded waist belt for hands-free convenience and comfort
- Available Colors: Red and Normal Black

| CAPACITY: 200 in<sup>3</sup> | WEIGHT: 14 oz | DIMENSIONS: H: 6" x W: 13" x D: 3" |



# OFF-DUTY PACKS



## CODE 7

# G11038RE / G11038BK

*Style meets function in this personal gear pack that keeps your personal supplies organized and close at hand for overnight shifts or long on-call days. Fits nicely between the seats in the cab of an ambulance.*

- Ergonomic shoulder straps and padded back make the pack very comfortable to wear
- Main compartments are efficiently organized to hold books, clothes, pens, CD's, and keys
- Padded inside sleeve accommodates most laptop computers
- Side pockets keep items organized and quickly accessible
- Front pocket can be easily personalized with company logo or department patch with standard embroidery machine
- Available Colors: Red and Normal Black

| CAPACITY: 1900 in<sup>3</sup> | WEIGHT: 1.0 lbs | DIMENSIONS: H: 18" x W: 13.5" x D: 9" |



## TRAVELER

# G11041RE / G11041BK

*This personal gear pack doubles as a back pack and a side messenger bag that keeps your personal gear structured and in one place. Great for holding clothes, sundries, and items for 24 hour shifts. Features lots of organization for computer gear and personal supplies.*

- Tuck away shoulder straps and padded back make the pack very comfortable to wear
- Removable side carry strap
- Main compartments are efficiently organized to hold books, clothes, pens, CD's, and keys
- Padded inside sleeve protects even the largest laptop computers
- Side pockets keep items organized and quickly accessible
- Available Colors: Red and Normal Black

| CAPACITY: 1300 in<sup>3</sup> | WEIGHT: 3.0 lbs | DIMENSIONS: H: 16" x W: 14" x D: 6" |



# ACCESSORIES

## LED LIGHTS

| RED # G11097RE | WHITE # G11097WH |

- Attaches to any Statpack belt or fixed webbing loop
- Steady on, and blinking modes
- Easily replaceable standard battery
- Red or white color



## VIAL STRAND 4 PACK

# G25000BK

An ingenious way to securely hold drug vials while keeping them quickly identifiable and easily retrievable. Military spec Velcro hook on back attaches to any compartment with Velcro loop. Each strand holds 5 vials of any size.

| DIMENSIONS: H: 2" x W: 5" |



## I.D. CARD

# G25001YE

This high-vis card folds to fit in the I.D. window of any StatPack. Use the pre-printed Pediatric and Adult tags or make your own using the reverse side. Made from recyclable plastic and will hold up to water and endless use.



## TEE SHIRT

| SMALL # G11059S | MEDIUM # G11059M | LARGE # G11059L |  
| XLARGE # G11059XL | XXLARGE # G11059XXL |

- SM - XXL available
- Black shirt with color logo
- 50/50 Heavy Cotton/Poly fabric



## TRAVEL MUG

# G11060GR

- 16 oz. double insulated mug
- Stainless steel outer, plastic inner.
- Secure sipper lid



Visit [StatPacksStore.com](http://StatPacksStore.com)  
for additional accessories and specials

# WILDLAND FIRE/SAR

## FIRELINE # G24020GR

New modular design that starts with a web gear belt and can interchange with holsters to create different configurations. The main pack can detach along with the side holsters that hold fuzee's and water. Reattach anywhere there is webbing. The fire shelter holster is removable and can be worn separately and in different set ups.



- Wide-mouth top pocket flap opens to a top-load main compartment
- Accommodates standard fire helmet under top pocket
- Internal mesh divider sleeves keep inside contents separate
- Fire-shelter stash bag gives fast, one-handed access to new-generation as well as old style fire-shelters
- Fire-shelter stash bag is removable and adjustable to prevent "pack rash" (can be worn separately)
- Bottle holsters hold up to 4 standard-issue water bottles, and feature foam construction to maintain shape, cinch cords and grommets drain holes to keep dirt and water out
- Ergonomic, adjustable, padded shoulder straps provide custom fit and superior freedom of movement
- Padded, back panel equipped with hydration bladder stash — holds up to 100-fluid oz. bladder
- Padded waist belt with bidirectional compression for optimum weight distribution and stability

| CAPACITY: 4000 in<sup>3</sup> | WEIGHT: 5.0 lbs | DIMENSIONS: H: 26" x W: 22" x D: 8" |  
(With all Holsters attached)





# QUICKDRAW

# G11024RE / G11024TK

Combines the technical features of the most advanced sport climbing packs with the practical and organizational needs of search and rescue. Designed for a 3:1 haul system, the QuickDraw allows for efficient access to gear while providing comfort and durability.

- Wide-mouth top for easy access to contents
- Web- helmet holder for quick, secure attachment
- Internal divider shelf zips in to divide off rope section or zips out for one large compartment area
- Removable rope bag makes storing up to 165 feet of climbing rope fast and easy; rolls out flat to protect rope from harmful dirt during use
- Mesh side pockets provide viewable, breathable storage for extra gear
- Ergonomic, adjustable padded shoulder straps for custom fit and superior freedom of movement
- Padded, back panel equipped with hydration bladder stash — holds up to 100-fluid oz. bladder
- Padded waist belt with bidirectional compression for optimum weight distribution and stability
- Large gear loops for climbing protection gear
- Available Colors: Red and Tactical Black

| CAPACITY: 1900 in<sup>3</sup> | WEIGHT: 3.6 lbs | DIMENSIONS: H: 21" x W: 12" x D: 7.5" |



# SPOTFIRE

# G11035BK

Wildland firefighters will appreciate all the industrial features of a sophisticated outdoor trail hiking pack combined with simplicity and sensible organization. Helps reduce fatigue by positioning the weight low on the back while distributing the load on hips and shoulders.

- Wide-mouth top pocket flap opens to a top-load main compartment
- Accommodates standard fire helmet under top pocket
- Internal mesh divider sleeves keep inside contents separate
- Fire-shelter stash bag gives fast, one-handed access to new- generation as well as old style-shelters
- Fire-shelter stash bag is removable and adjustable to prevent "pack rash" (can be worn separately)
- Bottle holsters hold up to 4 standard-issue water bottles, and feature foam construction to maintain shape; cinch cords and drain holes to keep dirt and water out
- Ergonomic, adjustable, padded shoulder straps provide custom fit and superior freedom of movement
- Padded, back panel equipped with hydration bladder stash — holds up to 100-fluid oz. bladder
- Padded waist belt with bidirectional compression for optimum weight distribution and stability
- Available Color: Normal Black

| CAPACITY: 1200 in<sup>3</sup> | WEIGHT: 3.5 lbs | DIMENSIONS: H: 14" x W: 10.5" x D: 8" |





# TACTICAL OVERVIEW

## COMBAT FASTOLOGY

*No need to preach the benefits of fast here. Combat medics have always understood the idea of carrying only what you need while still being prepared for the worst. StatPacks Tactical products feature the same upgrades as the EMS line, including reinforced wear areas, oversized zippers, and other performance-enhancing tweaks thanks to generous feedback from medics across the world. Tactical packs come standard in black but special orders are available in other military color patterns – call for more information regarding special orders.*

### TACTICAL ENHANCEMENTS:

- **ZIPPERS**

Beefy, #10 YKK zippers are used throughout; Custom zipper-pulls keep stress off the chain; Hypalon “garages” keep the wear areas of zippers protected; Reinforced stitching resists blow-out.

- **SUPER STRONG FABRICS**

Nearly indestructible coated EPO fabric used in areas of most wear.

- **IMPROVED SIDE POCKETS**

Newly engineered side pockets hold more, open easier, and look better than ever before.

- **SPECIAL MAKE AVAILABLE**

Call for special orders in military spec camouflage.







AUTO

SEMI





# TACTICAL PACKS



## GOLDEN HOUR # G11021TK

*The perfect pack for use in close-quarter areas, such as helicopters. Compact and comfortable, built to organize and carry essential equipment. The unique wide-mouth main compartment opening allows easy access to contents and provides convenient, upright availability from the captain's chair. Ergonomic shoulder straps and a breathable padded rear panel.*

- Foam-lined construction maintains pack shape and stays upright for efficient access to contents
- Clear, durable urethane window on top inside pouch allows quick view of contents
- Padded, ergonomic shoulder straps for maximum comfort
- Removable adjustable inside shelf for customized organization
- Designed to carry StatPacks® drug module, intubation module, IV module and accessories
- Large side pockets allow quick access to vital sign tools independent of the main compartment
- Special make Camouflage colors are available – please contact us for details
- Available Colors: Tactical Black and Red (see page 14)

| CAPACITY: 2400 in<sup>3</sup> | WEIGHT: 3.1 lbs | DIMENSIONS: H: 18" x W: 17" x D: 8" |



## LOAD N' GO # G11022TK

*The ideal pack for fixed wing and critical care transport, ambulance and fire engines. No pack keeps equipment more organized and quickly accessible thanks to thoughtfully designed transparent sleeves, mesh pockets, and elastic holsters that are arranged in three layers.*

- Foam-lined construction maintains pack shape and protects contents
- Clear, durable urethane windows, mesh pockets, and elastic holsters keep contents in place and well-organized
- Ergonomic, padded shoulder straps and back panel for comfort during transport
- QuickZip access to front compartment
- Large side pockets allow quick access to vital sign tools
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, and IV module
- Special make Camouflage colors are available – please contact us for details
- Available Colors: Tactical Black and Red (see page 14)

| CAPACITY: 2600 in<sup>3</sup> | WEIGHT: 4.5 lbs | DIMENSIONS: H: 20" x W: 19" x D: 7" |



## STEALTH # G11033TK

Incorporates advanced features of modern outdoor long-haul hiking packs with the realistic needs and governmental requirements of Military, SWAT, and Special Ops medics. Created to help Tactical Medics carry necessary medical equipment in special operation situations.

- Foam-lined construction maintains pack shape and protects contents
- Main compartment houses three removable sub-packs that are held in place with hook-&-loop backing
- Allows a customizable set-up with any combination of StatPacks modules including the QuickRoll Intubation, Drug and IV modules
- Quick-zip access to main compartment to offer fast, smooth action
- Single-buckle top closure allows quick in and out
- Ergonomic, padded shoulder straps, back panel and waist belt offer comfort during long transport
- Padded, back panel holds hydration bladder
- Large side pockets allow quick access to equipment
- Webbing loops on outside of pack for various attachments
- Load handle for quick/easy carry
- Silent zipper-pulls for quiet situations
- Special make Camouflage colors are available – please contact us for details
- Available Color: Tactical Black



| CAPACITY: 3300 in<sup>3</sup> | WEIGHT: 6.0 lbs | DIMENSIONS: H: 23" x W: 18" x D: 8" |

## BACKUP # G11023TK

Exceptional in urban or rural areas where "stay and play" medical protocol is used. Customizable shelving allows superior organization of virtually any combination of equipment. The advanced harness system helps you tackle daunting flights of stairs and hard-to-reach accident scenes with ease and safety.

- Foam-lined construction maintains pack shape and protects inside contents
- Fully adjustable Velcro shelving can be custom configured to individual preference
- Ergonomic, padded shoulder straps and back panel offer comfort during transport
- Removable padded waistbelt
- QuickZip access to main compartment
- Large side pockets allow quick access to vital sign tools
- Designed to carry StatPacks® drug module, QuickRoll intubation kit, and IV module
- Special make Camouflage colors are available – please contact us for details
- Available Colors: Tactical Black and Red (see page 15)



| CAPACITY: 3800 in<sup>3</sup> | WEIGHT: 6.0 lbs | DIMENSIONS: H: 25" x W: 18" x D: 8.5" |

## QUICKROLL INTUBATION MODULE # G11028TK

Designed to organize basic and advanced airway equipment. The single-buckle closure allows fast, one handed access. The perfect accessory for any StatPacks® pack.

- Clear internal durable urethane pockets allow quick identification of contents
- Organizes laryngoscope, blades, and ts as well as intubation tubes, syringes and forceps
- Holds OPA's / NPA's and gloves
- Single-buckle closure opens quickly, revealing all airway control tools
- Fits nicely in all StatPacks® packs
- Available Colors: Red and Tactical Black

| CAPACITY: 400 in<sup>3</sup> | WEIGHT: 11 oz. | DIMENSIONS: H: 8" x W: 10" x D: 5" (Rolled)



## DRUG MODULE # G11027TK

Compact and organized, this module brings all your meds together while offering a quick way for central supply or pharmacies to re-stock. Lockable kiss-zippers help you keep track of used modules. Save time by checking in used modules for restocked sealed ones.

- Transparent durable urethane outer window for quick medication identification
- External pocket window for medical drug listing or unit identification
- Lockable kiss-zippers for pharmacy seal
- Clear durable urethane pockets accommodate extra needles and syringes
- Attaches to the inside of larger StatPacks products with hook-&-loop panel
- Available Colors: Red and Tactical Black

| CAPACITY: 240 in<sup>3</sup> | WEIGHT: 9 oz. | DIMENSIONS: H: 6" x W: 8" x D: 5" |



## IV MODULE # G11026TK

Holds everything you need for a quick IV setup in one convenient module. The IV Module makes it easy to share on multi victim accident scenes.

- Accommodates all IV bag sizes 1 L, 500 mL, 250 mL
- Keeps IV Start Pack™ ready for use
- Front and inside pockets organize IV tubing and catheters together with the IV bag
- Transparent durable urethane window for quick solution identification
- Single-buckle closure securely holds contents in place and allows quick, one-handed access
- Top loop hangs the entire module with the IV solution
- Adjustable cam-lock buckle allows you to pressurize contents for rapid fluid infusion
- Available Colors: Red and Tactical Black

| CAPACITY: 210 in<sup>3</sup> | WEIGHT: 6 oz. | DIMENSIONS: H: 10" x W: 6" x D: 3.5" (Rolled)





# G1 TRAVERSE #G11032TK

A straightforward leg-pack provides quick access to small amounts of BLS and ALS equipment. Putting supplies literally at your fingertips, access in tight areas is swift and easy. Designed specifically for Tactical Medics in Military, SWAT, and Special Ops. Smart design features a shelf effect to prevent content spillage.

- Easily attaches to any belt and secures comfortably to either leg
- Shelf-style design opens to 90° for easy access from a standing position
- Accommodates IV fluids, some airway equipment and trauma dressings; also great for non-medical gear for the entire team
- Silent Zipper pulls for quiet situations
- Special make camouflage colors are available – please contact us for details
- Available color: Tactical Black

| CAPACITY: 250 in<sup>3</sup> | WEIGHT: 7 oz | DIMENSIONS: H: 16.5" x W: 10.5" x D: 2" |



## ORDER INFO:

Call us at (435) 627.2265

or fax us your P.O. at (435) 214.0152

Place your order and learn more about StatPacks at [www.statpacks.com](http://www.statpacks.com). Be sure to check out our special promotions and free giveaways. We invite you to share your stories and give us your feedback on our products. See what the press is saying and get our quarterly news letter to keep you posted on how the industry is changing.

### Distributors:

For a current list of our domestic and International distributors please call or visit our website.

### Government or Military Sales:

Please contact us at [customerservice@statpacks.com](mailto:customerservice@statpacks.com) or call 435-627-2265 to get information on our GSA Schedule. We have all the information needed:

- Name of product or accessory.
- NSN (if available)
- Product #
- GSA Contract #
- Cage Code
- Tax ID #
- DUNS #



### Special Order Policy:

Don't see exactly what you're looking for? Give us a call. Our talented R & D would love to hear new ideas. If you have a special need we would be happy to give you a quote for a custom design. Minimum quantities will apply. Drop us an e-mail, or give us a call.

### Warranty Information:

StatPacks offers a lifetime warranty against material and workmanship defects. If a product fails due to a manufacturing flaw, we will repair or replace it free of charge. Damage due to normal wear, accident, improper care or negligence where we are not at fault will be repaired for a reasonable charge.

[customerservice@statpacks.com](mailto:customerservice@statpacks.com)



Special Thanks to Dixi Ambulance and Santa Clara Fire



690 N Industrial Drive  
St George, UT 84770  
Phone: 435.627.2265  
Fax: 435.214.0152  
[www.statpacks.com](http://www.statpacks.com)

**Danmark**

Medidyne ApS  
Skodsborgvej 311  
2850 Nærum  
  
Tel: +45 35 25 12 48  
info@medidyne.dk  
www.medidyne.dk

**Sverige**

Medidyne AB  
Rönnbärsvägen 8  
260 83 Vejbystrand  
  
Tel: + 46 431 448 050  
info@medidyne.se  
www.medidyne.se

**Norge**

Medidyne AS  
Hovfaret 13  
0275 Oslo  
  
Tel: +47 67 53 01 40  
info@medidyne.no  
www.medidyne.no

**Finland**

Medidyne Oy  
Perintötie 2b  
01510 Vantaa  
  
Tel: +358 20 734 4350  
info@medidyne.fi  
www.medidyne.fi