




RUSCH® NASOPHARYNGEAL AIRWAYS

Clinician-inspired materials and designs



AIRWAYS



LATEX-FREE ROBERTAZZI AIRWAY

- MEDIPRENE™ construction
- Flexible for maximum patient comfort
- Soft material to avoid damage to nasal passageway
- Rounded tip allows for gentle insertion
- Trumpet design for secure placement
- Emergency kits available – include foil packages of lubricating jelly
- Latex free, single use product



AMBER ROBERTAZZI AIRWAY

- Amber latex construction
- Soft and flexible for maximum patient comfort
- Rounded tip allows for gentle insertion
- Trumpet design for secure placement



PVC AIRWAY

- Soft PVC construction
- Inner striations for smooth suction catheter insertion
- Beveled tip for comfortable, gentle insertion
- Available in a full range of sizes, including pediatric
- Emergency kits available—include foil packages of lubricating jelly

LATEX-FREE ROBERTAZZI AIRWAY

PRODUCT #	SIZE	LENGTH	BOX QTY
123120	20 FR	100 MM	10
123122	22 FR	105 MM	10
123124	24 FR	110 MM	10
123126	26 FR	118 MM	10
123128	28 FR	125 MM	10
123130	30 FR	127 MM	10
123132	32 FR	120 MM	10
123134	34 FR	155 MM	10
123136	36 FR	170 MM	10

PRODUCT #	DESCRIPTION	SIZES	BOX QTY
1231	Latex-free Robertazzi Airway Emergency Pack	22 FR – 32 FR (6 AIRWAYS)	1
1232	Latex-free Robertazzi Airway Emergency Pack	20 FR – 36 FR (9 AIRWAYS)	1

AMBER ROBERTAZZI AIRWAY

PRODUCT #	SIZE	LENGTH	BOX QTY
123020	20 FR	110 MM	10
123022	22 FR	110 MM	10
123024	24 FR	118 MM	10
123026	26 FR	118 MM	10
123028	28 FR	128 MM	10
123030	30 FR	138 MM	10
123032	32 FR	140 MM	10
123034	34 FR	155 MM	10
123036	36 FR	170 MM	10

PVC AIRWAY

PRODUCT #	SIZE	LENGTH	BOX QTY
123312	12 FR	60 MM	10
123314	14 FR	72 MM	10
123316	16 FR	85 MM	10
123318	18 FR	105 MM	10
123320	20 FR	115 MM	10
123322	22 FR	125 MM	10
123324	24 FR	130 MM	10
123326	26 FR	140 MM	10
123328	28 FR	155 MM	10
123330	30 FR	165 MM	10
123332	32 FR	170 MM	10
123334	34 FR	175 MM	10
123336	36 FR	180 MM	10

PRODUCT #	DESCRIPTION	SIZES	BOX QTY
1233	PVC Airway Emergency Pack	12 FR – 18 FR (4 AIRWAYS)	1
1234	PVC Airway Emergency Pack	20 FR – 36 FR (9 AIRWAYS)	1



ADJUSTABLE FLANGE LATEX-FREE AIRWAY

- MEDIPRENE construction
- Adjustable flange for secure placement
- Soft, resilient material to avoid damage to nasal passageway
- Rounded tip allows for gentle insertion
- Available in a full range of sizes, including pediatric
- Latex free, single use product

ADJUSTABLE FLANGE LATEX-FREE AIRWAY

PRODUCT #	SIZE	LENGTH	BOX QTY
18542012	12 FR	95 MM	5
18542014	14 FR	95 MM	5
18542016	16 FR	95 MM	5
18542018	18 FR	95 MM	5
18542020	20 FR	125 MM	5
18542022	22 FR	170 MM	5
18542024	24 FR	170 MM	5
18542026	26 FR	170 MM	5
18542028	28 FR	170 MM	5
18542030	30 FR	170 MM	5
18542032	32 FR	170 MM	5
18542034	34 FR	170 MM	5
18542036	36 FR	170 MM	5



ADJUSTABLE FLANGE SILKOLATEX™ AIRWAY

- Silicone coated latex material
- Adjustable flange for secure placement
- Soft material to avoid damage to nasal passageway
- Rounded tip allows for gentle insertion
- Latex free, single use product

ADJUSTABLE FLANGE SILKOLATEX™ AIRWAY

PRODUCT #	SIZE	LENGTH	BOX QTY
185200200	20 FR	125 MM	10
185200220	22 FR	170 MM	10
185200220	24 FR	170 MM	10
185200260	26 FR	170 MM	10
185200280	28 FR	170 MM	10
185200300	30 FR	170 MM	10
185200320	32 FR	170 MM	10
185200340	34 FR	170 MM	10
185200360	36 FR	170 MM	10



ADJUSTABLE FLANGE RED RUBBER AIRWAY

- Soft red rubber construction
- Adjustable flange for secure placement
- Rounded tip allows for gentle insertion
- Available in a full range of sizes, including pediatric
- Reusable latex product

ADJUSTABLE FLANGE RED RUBBER AIRWAY

PRODUCT #	SIZE	LENGTH	BOX QTY
125200120	12 FR	95 MM	5
125200140	14 FR	95 MM	5
125200160	16 FR	95 MM	5
125200180	18 FR	95 MM	5
125200200	20 FR	115 MM	5
125200220	22 FR	125 MM	5
125200240	24 FR	170 MM	5
125200260	26 FR	170 MM	5
125200280	28 FR	170 MM	5
125200300	30 FR	170 MM	5
125200320	32 FR	170 MM	5
125200340	34 FR	170 MM	5
125200360	36 FR	170 MM	5

Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specializes in products and services for vascular access, respiratory, general and regional anesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

The Teleflex family of brands includes ARROW®, BEERE MEDICAL®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, PILLING®, PLEUR-EVAC®, RUSCH®, SHERIDAN®, SMD®, TAUT®, TFX OEM®, VASONOVA™ and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. Teleflex Incorporated (NYSE:TFX) has annual revenues of approximately \$1.9 billion and serves healthcare providers in 140 countries.

For detailed information, see teleflex.com

TELEFLEX MEDICAL

PO Box 12600

Research Triangle Park, NC 27709

Toll Free 866.246.6990

Phone 919.544.8000

Intl 919.361.3961

MEDIPRENE is a registered trademark of VTC Elatoteknik. Silkolatex and Teleflex are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.
© 2011 Teleflex Incorporated. All rights reserved. 2010-0509 v1