



Operating Room Positioning Products



COVIDIEN: Positioning & Protection by Design

WHY DEVON™ POSITIONING PRODUCTS?

Positioning patients properly in the Operating Room is critical to providing the best care possible and protecting against complications such as pressure ulcers and nerve damage. These preventable injuries can be devastating to the patient and have significant consequences. Selecting the appropriate positioning devices for each procedure is a collaborative process for the perioperative registered nurse, surgeon and anesthesia provider.

Designed to aid in proper positioning and pressure distribution, Devon Positioning Products provide a complete offering to help position and protect patients in the Operating Room. You can trust Covidien to provide quality products, personalized service, education and support.

Raspberry Swirl Foam is 1.35 pounds per cubic foot for optimal support. All items are latex free, single patient use and meet California Technical Bulletin 117 fire safety standards.

Covidien has a direct U.S. sales force of service-driven professionals focused on Operating Room supplies in both the Acute Care and Ambulatory Surgery Center markets. Devon Positioning Products are also available through all major medical/surgical distributors.

Partner with Covidien to develop the optimal Operating Room positioning portfolio to care for your patients.



Head Positioning



DEVON™ HEAD POSITIONERS

- Provides support, stability, elevation and safe alignment of the head and neck

Devon™ Prone and Lateral Head Positioning

- Supports and protects facial features in the prone position
- Designed to accommodate endotracheal tube
- Soft Touch positioner provides an additional soft foam layer for improved patient comfort



Slotted Head Positioner
31143160



Soft Touch Head Positioner
31163184




Frame Head Positioner
31143020

Reorder	Description	Case Quantity
31143160	Slotted Head Positioner	24
31163184	Soft Touch Head Positioner	28
31143020	Frame Head Positioner	28

COVIDIEN Operating Room Positioning Products Head Positioners

Devon™ Supine Head Positioning

- Provides balanced pressure distribution to back of head
- Available in various sizes for proper fit of all patients



Adult Head Positioner
31143129

Child Head Positioner
31143194



Head Rest Whole-in-One
31143152




Head Rest - 7" Bagel
31143145



Head Rest - 9" Bagel
31143137

Reorder	Description	Case Quantity
31143129	Adult Head Positioner	24
31143194	Child Head Positioner	24
31143152	Head Rest – Whole-in-One	36
31143145	Head Rest - 7" Bagel	72
31143137	Head Rest - 9" Bagel	36





Body Positioning

DEVON™ BODY POSITIONERS

- Provides support and stability required for safe positioning
- Protects bony prominences and other pressure points
- Designed to safely position patients for a wide range of procedures



Body Alignment Wedge
31143392

OR Table Wedge
31143103



Convoluted Table Pad
31163457



Foam Squares
572304



Convoluted Utility Pad
31143384



Chest Roll
31143418



Pole Protector
31143509

Reorder	Description	Case Quantity
31143103	45° OR Table Wedge, 18" x 8" x 8"	8
31143392	30° Body Alignment Wedge, 22" x 12" x 9"	8
31163457	Operating Room Table Pad, 72" x 20" x 2"	12
31143384	Utility Pad, 20" x 12" x 2"	12 pair
31143418	Chest Roll, 17" x 5"	12
31143434	Bolster, Small, 18" x 7"	12
31143442	Bolster, Large, 30" x 7"	10
31143509	Pole Protector	36
572304	Foam Square Positioner, 7" x 7" x 1"	288

COVIDIEN Operating Room Positioning Products **Body Positioners**

DEVON™ POSITIONING STRAPS

- Straps are specifically designed for OR tables and gurneys
- Soft, yet strong to ensure proper positioning for security and comfort
- Latex free



Disposable straps with adjustable hook & loop fasteners and durable rings for attachment



Knee & Body Strap
31143004



Knee & Body Strap
31142964



Arm Board Strap
31142980

Reorder	Description	Case Quantity
Reusable Straps		
31142964	Knee & Body Strap, 60" x 4", Non-conductive	6
31142931	Knee & Body Strap, 60" x 4", Conductive	6
31142998	Arm Board Strap, 26" x 1.5"	25
Disposable Straps		
31143004	Knee & Body Strap, 60" x 3"	12
31142980	Arm Board Strap, 32" x 1.5"	48
31143005	Knee & Body Strap, 60" x 4"	8
31143006	Knee & Body Strap, 80" x 4"	6



Extremities



DEVON™ EXTREMITY POSITIONERS

- Provides support and stability required for safe positioning
- Protects bony prominences and other pressure points
- Designed to safely position patients for a wide range of procedures

Devon™ Upper Extremity Positioners


- Variety of positioners provide support, stability and protection for shoulders, elbows, ulnar nerve, and bony prominences
- Designed to help maintain proper alignment
- Aids in pressure distribution of critical pressure points during the surgical procedure



Convoluted Arm Board Pad
31143467



Convoluted Ulnar Nerve Protector
31143095



Arm Cradle - Laminectomy
31143491

Reorder	Description	Case Quantity
31143095	Ulnar Nerve Protector	36 pair
31143467	Convoluted Arm Board Pad	12 pair
31143491	Arm Cradles – Laminectomy	12 pair

COVIDIEN Operating Room Positioning Products **Extremities**

Devon™ Lower Extremity Positioners

- Variety of positioners provide support, stability and protection for legs, knees, heels and ankles
- Designed to help maintain proper lower body alignment
- Aids in pressure distribution of lower body during the surgical procedure



Abduction Pillows (3 sizes)
31143087, 31143061 & 31143053



Cradle Positioner
31157988



Frog Leg Positioner
31143376



Convoluted Foot & Heel Protectors
31143178



Arthroscopic Well Leg Holder
31143277



Arthroscopic Knee Holder
31143046

Reorder	Description	Case Quantity
31143087	Abduction Pillow – Small, 18" x 12" x 6"	6
31143061	Abduction Pillow – Medium, 22" x 15" x 6"	6
31143053	Abduction Pillow – Large, 25" x 18" x 6"	6
31157988	Cradle Positioner	12
31143376	Frog Leg Positioner	5 pair
31143178	Foot & Heel Protectors	36 pair
31143277	Arthroscopic Well Leg Holder	18
31143046	Arthroscopic Knee Holder	10



Procedural Kits



COVIDIEN Operating Room Positioning Products Procedural Kits

DEVON™ PROCEDURAL KITS

- Designed for specific OR tables and procedures
- Provides all necessary components in one convenient and compact kit
- Reduces OR set-up time

Jackson Table Kit

Components:

1 Soft-Touch Headrest Pillow

1 Chest Pillow

1 Pair, Arm Cradles

4 Hip & Thigh Pad Covers




Reorder	Description	Case Quantity
31151010	Jackson Table Kit	6
31151011	Jackson Table Kit - without arm cradles	6
31151012	Jackson Table Kit - without headrest or arm cradles	6
31151020	Jackson Table Kit - without headrest	6

Supine Lithotomy Trendelenburg Kit

SLT Kit is rated for patients up to 450 pounds.

Components:

1 Graduated Foam Table Pad	1 Chest Foam
1 Drawsheet	1 Body Strap



Reorder	Description	Case Quantity
31151090	Supine Lithotomy Trendelenburg Kit	6

Andrews Table Kit

Components:

1 Soft-Touch Headrest Pillow	1 Pair, Arm Cradles
1 Chest Pillow	1 Seat Cover
1 Thigh Strap	1 Pair, Foam Heel Protectors



Reorder	Description	Case Quantity
31151030	Andrews Table Kit	6

Andrews Frame Kit

Components:

1 Soft-Touch Headrest Pillow	1 Seat Drape Cover
1 Thigh Strap	1 Pair, Heel Cup Protectors



Reorder	Description	Case Quantity
31151040	Andrews Frame Kit	3

Schlein Beach Chair Kit

Components:

- 1 Back Support Pad
- 1 Torso Support Strap
- 1 Chin Strap
- 1 Optional Padding for Patient's Head
- 1 Head Support
- 1 Arm Support Strap
- 1 Forehead Strap



Reorder	Description	Case Quantity
31151050	Schlein Beach Chair Kit	3

Wilson Frame Kit

Components:

- 1 Soft-Touch Headrest Pillow
- 2 Pad Sleeve Covers
- 1 Bar Protector
- 1 Pair, Arm Cradles
- 1 Drape Protector



Reorder	Description	Case Quantity
31151060	Wilson Frame Kit	6
31151061	Wilson Frame Kit - without arm cradles	6
31151062	Wilson Frame Kit - without headrest	6


Knee Positioning Kit

Components:

- 1 Arthroscopic Well Leg Holder
- Arthroscopic Knee Holder



Reorder	Description	Case Quantity
31143293	Knee Positioning Kit	10



Larger Patients

DEVON™ POSITIONERS FOR LARGER PATIENTS

- Safe, consistent positioning for larger patients
- Provides additional protection for patients vulnerable to skin integrity concerns
- Wider, thicker pieces designed to distribute weight and maintain circulation



Round Head Positioning
777001



Supine Positioning Kit
777004



Arm Positioner
777002



Leg Positioner
777003

Reorder	Description	Case Quantity
777001	Round Head Positioner	36
777002	Arm Positioners	6 pair
777003	Leg Positioners	6 pair
777004	Supine Positioning Kit	4

COVIDIEN Operating Room Gel Positioning Products
Pressure Relief & Patient Positioning





OASIS ELITE™* - PRONE HEAD REST

Developed in conjunction with the Operating Staff at a leading hospital, the Oasis Elite™* Prone Head Rest is anatomically designed so that the patient's head can comfortably rest face down with respiratory tubes in place. The unique design provides excellent support, stability, comfort and protection during surgery in the prone position, with openings on both sides. The product is lightweight, durable and easy to handle.

Reorder	Description	Dimensions
50003141	EL141 Prone Head Rest, Adult	11" x 9" x 6"

Surgical Positions: Prone



OASIS ELITE™* - SUPINE HEAD REST

This Supine Head Rest is ideal for positioning the patient during eye surgery and other supine procedures. Ergonomically shaped to position the patient's head, while effectively protecting against the onset of pressure sores.

Reorder	Description	Dimensions
50003151	EL151 Supine Head Rest, Adult	10.75" x 9.25" x 3.5"

Surgical Positions: Supine



OASIS ELITE™* - CLOSED HEAD RING

Unique gel & foam combination conforms comfortably to the patient's head, this head ring is suitable for use in many surgical procedures.

Reorder	Description	Dimensions
50003011*	EL011 Closed Head Ring, Adolescent	5.5" x 1.5"
50003012	EL012 Closed Head Ring, Adult	8" x 2"

Surgical Positions: Supine, Lithotomy, Lateral
**Specific for Infant and Adolescent use only*



OASIS ELITE™* - OPEN HEAD RING

Combination of gel and foam provides stability and protection to the patient's face during surgery. An opening allows for intubation.

Reorder	Description	Dimensions
50003021*	EL021 Open Head Ring, Infant/Adolescent	5.5" x 1.5"
50003022	EL022 Open Head Ring, Adult	8" x 2"

Surgical Positions: Prone, Lateral
**Specific for Infant and Adolescent use only*

COVIDIEN Operating Room Gel Positioning Products **Body**



AZURE™* - SACRAL PROTECTOR

The Sacral area is an area of serious risk if left unprotected. The Azure™ Sacral protector has been designed to provide support in the supine and lithotomy positions by relieving pressure and focusing on peripheral support.

Reorder	Description	Dimensions
50009300*	AZ300 Azure™ Sacral Protector, Small	15.75" x 15" x 1"
50009301	AZ301 Azure™ Sacral Protector, Large	17.75" x 17.5" x 1"

Surgical Positions: Supine, Lithotomy

**Specific for Adolescent use only*



OASIS ELITE™* - PRONE POSITIONER

The Prone Positioner gel and foam combination has been designed to provide effective positioning to the patient during surgical procedure in the prone position. For best protection, it is recommended to be used in conjunction with the Prone Head Rest.

Reorder	Description	Dimensions
50003903	EL903 Prone Positioner, Adult	24" x 21" x 5"

Surgical Positions: Prone



OASIS ELITE™* - DOME POSITIONER

This lightweight positioner can be used to position and protect the patient during various surgical procedures in all positions. It can be used in many different applications.

Reorder	Description	Dimensions
50003210*	EL210 Dome Positioner, Adolescent/Adult	14" x 6" x 3"

Surgical Positions: Prone, Supine

**Specific Adolescent use only for Prone Position*



COVIDIEN Operating Room Gel Positioning Products **Body**





OASIS™* - FLAT BOTTOM CHEST ROLLS

These 100% silicone gel chest rolls provide excellent positioning and stability and are available in a wide variety of sizes.

Reorder	Description	Dimensions
50008209*	OA209 Flat Bottom Chest Roll, Infant	8" x 3" x 2"
50008211	OA211 Flat Bottom Chest Roll, Adult	12" x 3.75" x 3"
50008213	OA213 Flat Bottom Chest Roll, Adult	20" x 3.5" x 3"

Surgical Positions: Supine, Prone, Lithotomy, Lateral

**Specific for Pediatrics*



OASIS™* - UNIVERSAL PATIENT POSITIONERS

These 100% silicone positioners are suitable for safe and effective positioning during all surgical procedures.

Reorder	Description	Dimensions
50008051	OA051 Universal Patient Positioners, Infant/Adolescent/Adult	13" x 2.75" x 2"
50008052	OA052 Universal Patient Positioners, Adolescent/Adult	11.5" x 4" x 2.5"
50008053	OA053 Universal Patient Positioners, Adolescent/Adult	13.5" x 4.75" x 2.75"

Surgical Positions: Supine, Prone, Lithotomy, Lateral, Pediatric



OASIS Elite™* - LATERAL POSITIONER

This positioning device permits a safe, quick and easy transition from supine to lateral. The patient's lower shoulder fits into the integral channel to provide exceptional patient comfort and nerve protection. The Lateral Positioner is designed for use in hip replacement, shoulder, kidney and other procedures. The built-in headrest and auxiliary support assists in proper spinal alignment and surgical site access.

Reorder	Description	Dimensions
50003900*	EL900 Lateral Positioner, Adolescent/Adult	28" x 20" x 6.5"

Surgical Positions: Lateral

**Adolescent use only for Lateral Position*

COVIDIEN Operating Room Gel Positioning Products **Table Pad**



OASIS™ - STANDARD OPERATING TABLE PADS

These Oasis™ table pads provide protection over the entire surface area of the operating table and are recommended for use with patients in all positions. The pads are flexible, soft, comfortable and versatile and distribute weight evenly over the entire surface area.

Reorder	Description	Dimensions
50008030	OA030 Standard Operating Universal Table Pad, Small	20.5" x 20.5" x 0.5"
50008032	OA032 Standard Operating Universal 3/4 Length Table Pad	71" x 20.5" x 0.5"

Surgical Positions: Supine, Prone, Lithotomy, Lateral, Pediatric




OASIS™ - THREE PIECE OPERATING TABLE PAD

The Oasis™ 3 Piece Operating Table Pad is designed for use on sectioned tables. The set consists of a head section, torso section and a foot section, with the latter two having a perineal cutout. The product provides a high degree of flexibility, versatility and protection during a variety of surgical procedures, particularly in the areas of Orthopaedics and Urology. Can also be used as a standard full length table pad.

Reorder	Description	Dimensions
50008039	OA039 3 Piece Operating Table Pad Set	68.5" x 20" x 0.5"
50008037	OA0337 Table Pad Torso Section	40" x 20" x 0.5"

Surgical Positions: Supine, Prone, Lithotomy, Lateral, Pediatric



OASIS™* - CONTOURED KNEE/ARM PAD

This contoured pad has been designed for use on a standard table extension. Its contoured shape provides excellent protection as well as keeping the arm or leg in place.

Reorder	Description	Dimensions
50008101	OA101 Contoured Knee/Arm Pad <i>Surgical Positions: Supine, Prone, Lithotomy, Lateral</i>	20" x 6" x 1.5"



Oasis™* - ARMBOARD PADS

Oasis™* Armboard Pads provide the necessary protection whenever the arm is extended on an armboard. Suitable for all types of armboards.


Reorder	Description	Dimensions
50008070	OA070 Armboard Pads, Infant/Adolescent/Adult	15.5" x 4.5" x 0.5"
50008071	OA071 Armboard Pads, Infant/Adolescent/Adult <i>Surgical Positions: Supine, Prone, Lithotomy, Lateral</i>	23.5" x 8" x 0.5"



OASIS™* ULNAR - ULNAR/BRACHIAL NERVE PROTECTORS

The Oasis™* Ulnar/Brachial Nerve Protectors are essential for providing protection to the sensitive areas of the elbows and arms.

Reorder	Description	Dimensions
50008060	OA060 Ulnar/Brachial Nerve Protector, Adolescent/Adult	15" x 4.5" x 0.25"
50008061	OA061 Ulnar/Brachial Nerve Protector, Adolescent/Adult <i>Surgical Positions: Supine, Lithotomy, Lateral</i>	18" x 8.25" x 0.25"



OASIS™* - ARM RETAINER PAD

The Oasis™* Arm Retainer Pad has been specially designed to fit perspex/acrylic arm supports. It can be easily fitted with Velcro™* attachments.

Reorder	Description	Dimensions
50008102	OA102 Arm Retainer Pad <i>Surgical Positions: Supine, Pediatric</i>	6" x 5" x 0.5"

COVIDIEN Operating Room Gel Positioning Products **Leg and Knee**





OASIS Elite™* - LATERAL LEG POSITIONER

This product is used to position the legs in a lateral position for kidney surgery and hip surgery. In addition to comfort, it facilitates optimal surgical positioning. This positioner also allows for knee and hip flexion. To use: Place patient on side, place upper limb on top (gel part of the product), and place lower limb underneath ledge.

Reorder	Description	Dimensions
50003800	EL800 Lateral Leg Positioner, Adult	23.5" x 8" x 0.5"

Surgical Positions: Lateral



OASIS™* - CRUTCH STIRRUP PADS

The Oasis™* Crutch Stirrup Pads are designed to fit most types of stirrups. Supplied in pairs, they come in a variety of sizes.

Reorder	Description	Dimensions
50008091	OA091 Crutch Stirrup Pads, Adolescent/Adult	15.5" x 9.5" x 0.5"
50008092	OA092 Crutch Stirrup Pads, Adolescent/Adult	20.5" x 12" x 0.5"

Surgical Positions: Lithotomy



OASIS Elite™* - HEEL PADS

The Oasis™* Elite Heel Pads provide support and maximum protection for the heels, even during the longest surgical procedures. They support the Achilles and effectively remove pressure from the heels. Supplied in pairs, the products are lightweight and easy to use.

Reorder	Description	Dimensions
50008130	EL130 Heel Pads, Adolescent <i>Surgical Positions: Supine</i>	7" x 4" x 2.75"







OASIS™ - HEEL AND ANKLE PROTECTOR

These unique pads are designed to fit around the heel and ankle like a boot with Velcro™ attachment. They provide protective padding for the heel when the patient is in the supine or lateral positions. Can also be used as a protector for the ankle and foot when using conventional ankle strap stirrups and lithotomy poles.

Reorder	Description	Dimensions
50008200	OA200 Heel & Ankle Protector, Adolescent/Adult	10" x 4.75" x 0.5"

Surgical Positions: Supine, Lithotomy, Lateral Adolescent Only: For Lateral Use



OASIS™ - CONTOURED KNEE/ARM PAD

This contoured pad has been designed for use on a standard table extension. Its contoured shape provides excellent protection as well as keeping the arm or leg in place.

Reorder	Description	Dimensions
50008101	OA101 Contoured Knee/Arm Pad	20" x 6" x 1.5"

Surgical Positions: Supine, Prone, Lithotomy, Lateral

Devon™ OR
Patient Positioning

- Trust the proven quality of Devon™ OR products
- Comprehensive portfolio of foam and gel devices, including table specific kits
- Personalized service, education and support
- Devon™ positioning foam is clinician preferred

*Positioning and
Protection
by Design*





COVIDIEN, COVIDIEN with logo, Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien Company
©2015 Covidien. All rights reserved.

H8933 xM 06/15

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

800-962-9888
508-261-8000

WWW.COVIDIEN.COM