

Lightweight Folding Wheelchairs

Traveler® L3

The Traveler L3 lightweight frame weighs less than 36 pounds making it the ideal K0003 Medicare coded wheelchair built for short and long term use. Designed with dual axles and multi-position caster forks, this wheelchair provides an adjustable seat height of 18 or 20 inches. The wheelchair comes standard with non-marking polyurethane tires and industry compatible front rigging.



Standard Features

- Black nylon upholstery with chart pocket on back
- Non-marking polyurethane rear tires
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003
- Latex Free
- Limited lifetime warranty

Traveler® L4

The Traveler L4 is a high-strength lightweight K0004 Medicare coded wheelchair. The Traveler L4 has dual axles and multi-position caster forks for an adjustable seat height of 18 or 20 inches. The wheelchair comes standard with non-marking polyurethane tires, anti-tippers and a chart pocket.



Standard Features

- Black nylon upholstery with chart pocket on back
- Hammer tone paint provides a long lasting and durable finish
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0004
- Latex Free
- One year limited warranty

E&J Navigator

The E&J Navigator easily converts from a classic transport chair to a lightweight self propelled wheelchair, using the included 24" quick-release axle rear wheels. The rear wheels are easily attached or removed with the push of a button. No additional tools are required.



Other features include convenient flip back desk arms and swing away footrests for safe patient transfer. The back of the Navigator folds down making transport and storage easy. The compact size, flexible configurations, lightweight and high strength frame make this the perfect companion chair.

Standard Features

- Transport chair and lightweight wheelchair all in one
- Available in 3 stylish colors – Black, Blue and Red
- 18" seat width
- Aluminum frame
- Black padded nylon upholstery
- Convenient flip-back desk arms
- Fold down back makes transport and storage easy
- Non marking polyurethane tires
- Weighs 33.5 lb with 24" rear wheels installed, 23 lb without 24" rear wheels (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code K0004
- Rear mount Anti-tippers
- Two year limited warranty

Transport Chairs

Aluminum Transport

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device and known for its durability. Its compact size and lighter weight allow a companion to lift and store the chair with ease.



Rear mounted attendant locks make it easy for a companion to engage. Available in 5 stylish colors options of Black, Red, Blue, Silver and Pink.

Standard Features

- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 17" and 19" seat widths
- Footrests are industry compatible and offer easy tool free adjustment
- Full length padded arms provide user comfort
- Weighs 20 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E10385
- One year limited warranty

Steel Transport

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device and known for its durability. Its compact size and lighter weight allow a companion to lift and store the chair with ease.



Rear mounted attendant locks make it easy for a companion to engage.

Standard Features

- Frame is finished in a durable hammertone coating
- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 17" and 19" seat widths
- Footrests are industry compatible and offer easy tool free adjustment
- Full length padded arms provide user comfort
- Weighs 24 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E10385
- One year limited warranty

Bariatric Transport

The Everest & Jennings Bariatric Transport Chair is designed as an affordable bariatric patient transfer device and known for its extra width. The generous 24" width and 400 lb weight capacity make it the right solution when fitting patients that require extra capacity in addition to a compact transport need.



Handle mounted attendant locks make it easy for a companion to engage.

Standard Features

- Lightweight aluminum frame with silver finish
- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 24" seat width
- Carry straps on seat allow for ease of lifting and folding
- Padded desk arms provide user comfort
- Weighs 31 lb (w/o front rigging)
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E1039
- Three year limited warranty

item No.	Description	item No.	Description
advantage (K0001)			
3h010220	Advantage 16X16 Desk Arms, Footrest	3h010230	Advantage 16X16 Desk Arms, Elevating Legrest
3h010120	Advantage 18X16 Desk Arms, Footrest	3h010130	Advantage 18X16 Desk Arms, Elevating Legrest
3h010320	Advantage 20X16 Desk Arms, Footrest	3h010330	Advantage 20X16 Desk Arms, Elevating Legrest
Traveler se (K0001/K0002)			
3e010220	Traveler SE 16X16 Desk Arms, Footrest	3e010230	Traveler SE 16X16 Desk Arms, Elevating Legrest
3e010120	Traveler SE 18X16 Desk Arms, Footrest	3e010130	Traveler SE 18X16 Desk Arms, Elevating Legrest
3e010320	Traveler SE 20X16 Desk Arms, Footrest	3e010330	Traveler SE 20X16 Desk Arms, Elevating Legrest
Traveler l X (K0001/K0002)			
3e020220	Traveler LX 16X16 Desk Arms, Footrest	3e020230	Traveler LX 16X16 Desk Arms, Elevating Legrest
3e020120	Traveler LX 18X16 Desk Arms, Footrest	3e020130	Traveler LX 18X16 Desk Arms, Elevating Legrest
3e020320	Traveler LX 20X16 Desk Arms, Footrest	3e020330	Traveler LX 20X16 Desk Arms, Elevating Legrest
Traveler l 3 (K0003)			
3F010220	Traveler L3 16X16 Desk Arms, Footrest	3F010230	Traveler L3 16X16 Desk Arms, Elevating Legrest
3F010120	Traveler L3 18X16 Desk Arms, Footrest	3F010130	Traveler L3 18X16 Desk Arms, Elevating Legrest
3F010320	Traveler L3 20X16 Desk Arms, Footrest	3F010330	Traveler L3 20X16 Desk Arms, Elevating Legrest
Traveler l 4 (K0004)			
3F020220	Traveler L4 16X16 Desk Arms, Footrest	3F020230	Traveler L4 16X16 Desk Arms, Elevating Legrest
3F020120	Traveler L4 18X16 Desk Arms, Footrest	3F020130	Traveler L4 18X16 Desk Arms, Elevating Legrest
3F020320	Traveler L4 20X16 Desk Arms, Footrest	3F020330	Traveler L4 20X16 Desk Arms, Elevating Legrest
Navigator Wheelchair & Transport chair			
eJNaV-18BK	E&J NAVIGATOR 18W BLACK	eJNaV-18r D	E&J NAVIGATOR 18W RED
eJNaV-18BI	E&J NAVIGATOR 18W BLUE		
Traveler hD heavy Duty (K0007)			
3G010320	Traveler HD 20X18 Desk Arms, Footrest	3G010330	Traveler HD 20X18 Desk Arms, Elevating Legrest
3G010420	Traveler HD 22X18 Desk Arms, Footrest	3G010430	Traveler HD 22X18 Desk Arms, Elevating Legrest
3G010520	Traveler HD 24X18 Desk Arms, Footrest	3G010530	Traveler HD 24X18 Desk Arms, Elevating Legrest
Traveler hospital Transport chair (hTc)			
3J010330	Traveler HTC 20X18 Desk Arms, Elevating Legrest	3J010530	Traveler HTC 24X18 Desk Arms, Elevating Legrest
3J010430	Traveler HTC 22X18 Desk Arms, Elevating Legrest		
Paramount XD heavy Duty (K0007)			
5PX10620	Paramount XD 26W Desk Arms, Footrest	5PX10630	Paramount XD 26W Desk Arms, Elevating Legrest
5PX10820	Paramount XD 30W Desk Arms, Footrest	5PX10830	Paramount XD 30W Desk Arms, Elevating Legrest
aluminum Transport chairs			
eJ789-1/eJ780-1	Transport Chair Green 17" / 19"	eJ785-1/eJ786-1	Transport Chair Black 17" / 19"
eJ781-1/eJ782-1	Transport Chair Red 17" / 19"	eJ787-1/ eJ788-1	Transport Chair Blue 17" / 19"
eJ783-1/eJ784-1	Transport Chair Silver 17" / 19"	eJ790-1 / eJ791-1	Transport Chair Pink 17" / 19"
steel Transport chair		Bariatric Transport	
eJ796-1	Transport Chair Silver Vein Steel 19"	eJ777-2	Transport Bariatric Chair 24W
advantage reclining Wheelchair			
3K010230	Advantage Recliner 16X17 Desk Arms, Elevating Legrest	3K010250	Advantage Recliner 16X17 Full Arms, Elevating Legrest
3K010130	Advantage Recliner 18X17 Desk Arms, Elevating Legrest	3K010150	Advantage Recliner 18X17 Full Arms, Elevating Legrest
3K010330	Advantage Recliner 20X17 Desk Arms, Elevating Legrest	3K010350	Advantage Recliner 20X17 Full Arms, Elevating Legrest

*For a complete list of sizes and options available, please visit www.grahamfield.com

Everest & Jennings®

Wheelchairs

Setting standards in quality
Designed with precision
Built to last





Everest & Jennings®

Value:

Graham-Field is known for a broad product offering at affordable pricing; without sacrifice to quality. Like you, we understand that the cost of a product is more than a price tag. The cost also includes the time and effort required to sell and service your customers. At Graham-Field we are proud to offer the value found in the Everest & Jennings line of wheelchairs.

Quality:

E&J has maintained a commitment to provide our customers quality manual folding wheelchairs since 1932. In a time when many manufacturers are producing products at the lowest possible cost, often cutting corners on quality, E&J is standing by its dedication to quality through product design, specification control, testing and inspection. It is easy to see why the E&J Traveler® name is synonymous with quality.

Service:

Our products are often only a day away with centralized distribution located in the metro areas of Atlanta, Ga.; St. Louis, Mo; and Los Angeles, Ca. E&J also carries a full line of parts and upholstery to meet your service needs. In addition, our customer service department and technical service teams are always ready to help with your product needs.



Manual Folding Wheelchairs

Advantage

The Advantage Wheelchair provides comfort and durability at an affordable price. Designed as a low maintenance wheelchair, it offers an attractive silver vein durable hammer tone steel finish and comes standard with padded nylon upholstery.

Standard Features

- Black, padded nylon upholstery with chart pocket on back
- Carbon Steel frame provides strength and durability
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Hammer tone paint provides a long lasting and durable finish
- Comfortable padded armrests
- Maximum Weight Capacity – 250 LB – EVENLY DISTRIBUTED
- Meets Medicare code K0001
- Latex Free
- One year limited warranty



Traveler® SE

The Traveler SE incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction and the latest industry standard front rigging. This wheelchair comes in a hammer tone finish and provides the most effective combination K0001/K0002 wheelchair available.

Standard Features

- Black vinyl upholstery with chart pocket on back
- Hammer tone paint provides a long lasting and durable finish
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002
- Latex Free
- One year limited warranty



Traveler® LX

The Traveler LX incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction and industry standard front rigging. The premium vinyl upholstery is reinforced with an extra heavy duty inner liner. This wheelchair comes in a chrome finish and provides the highest quality combination K0001/K0002 wheelchair available.

Standard Features

- Blue premium vinyl upholstery with chart pocket on back
- Chrome finish provides a long lasting and easy to clean finish
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Gray non-marking tires
- Two position axle and caster fork provide 2" seat height and adjustment
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002
- Two year limited warranty



Bariatric Folding Wheelchairs

Traveler® HD

The Traveler HD provides the ideal solution to the need for a heavy duty frame designed to support up to 500 lb users.

The Traveler HD comes with double crossbraces, heavy-duty 5/8" diameter axles, and a long lasting, embossed back and seat upholstery with a double inner liner for durability.

Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard). Meets Medicare K0007 specifications.

Standard Features

- Embossed black leatherette upholstery with double inner lining for durability
- Hammer tone paint provides a long lasting and durable finish
- Reinforced frame at all critical load-bearing points
- Chart pocket standard on upholstery
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007
- 5 year limited warranty



Paramount XD

The Paramount XD is ideal for the bariatric wheelchair user. It features 26" and 30" seat widths as well as a weight capacity of 650 lb. The Paramount XD is constructed of 1" tubing to meet high weight capacity demands.

Standard Features

- 19 1/2" seat depth
- Swingaway footrest and Elevating legrest options
- Maximum weight capacity of 650 lbs.
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Reinforced welds and gussets on frame
- 8" x 2" Solid casters
- 24" Solid rear wheels
- Meets Medicare code K0007
- Limited Lifetime Warranty on Frames and Crossbraces



Traveler® HTC

The Traveler HTC (Hospital Transport Chair) provides the ideal solution to the hospital/clinical transport need by incorporating special attendant and inventory control features with a heavy duty frame reinforced at all critical load-bearing points.

The Traveler HTC comes with an attendant push handle and rear mounted brakes that offer easier transport for the caregiver. In addition, the full width push handle acts as an anti-fold bar.

Standard Features

- Fixed attendant push handle also functions as an anti-fold bar
- Includes a fixed anti-theft post ideal in hospital/clinical settings
- Available in 20", 22" and 24" seat widths, 18" seat depth (1" longer than standard)
- Rear mounted attendant wheel locks
- Fixed desk arms and lock nuts on the front rigging ensure inventory controls
- Anti-tippers standard
- Industry compatible "Tool-free" front rigging with heel loops standard on swing away footrests
- Black vinyl upholstery with chart pocket on back
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007
- 5 year limited warranty



Specialty Folding Wheelchairs

Advantage Recliner

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. Anti-tippers and an anti-fold bar come standard. The durable hammer tone frame is designed to support up to 300 lb users.

Available in 16", 18" and 20" seat widths, with 17" seat depth.

Standard Features

- Manual full reclining back
- Removable full support headrest
- Detachable desk or full arms
- Hammer tone finish provides durability and style
- Low maintenance black vinyl upholstery
- Anti-tippers standard
- Anti-fold bar standard
- "Tool-free" adjustable elevating front rigging
- Chart pocket standard on back
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001/ E1226 (Recliner Back) / E0966 (Head Rest Extension)
- 5 year limited warranty



Accessories

high Pivot legrest		legrest		high Pivot legrest	
	90763030		90763031		90763435
Footrest with Composite Footplate		Elevating Legrest with Composite Footplate		High Pivot Elevating Legrest with Aluminum Footplate, Hammer tone Finish	
	90763430		90763431		90763135
Footrest with Aluminum Footplate		Elevating Legrest with Aluminum Footplate		High Pivot Elevating Legrest with Aluminum Footplate, Chrome Finish	
iV/O2 holders		Positioning Belts		anti-tippers	
	B10320F		92133040		90763565
O2/IV Holder Combination		Positioning Belt, Metal Buckle		Anti-tipper, Advantage	
	GF1844W		300-505		90763065
IV Holder		Positioning Belt, Auto Style		Anti-tipper, Traveler SE & HD	
	GFM905-1		eJ76-Bel T		90763065
Cylinder Carrier		Positioning Belt, Plastic Buckle		Anti-tipper, Traveler L3	