


2011 Equipment Catalog



- Monitoring Equipment
- Temperature Management
- Anesthesia Equipment
- MR Conditional Equipment
- Ventilators
- Vaporizers
- Gas Delivery Management
- Waste Gas Management
- Rescue & Resuscitation



Veterinary

Equipment Products

Table of Contents

Monitoring Equipment	p2
Temperature Management	p15
Anesthesia Equipment	p17
MR Conditional Equipment	p26
Ventilators	p28
Vaporizers	p30
Gas Delivery Management	p33
Waste Gas Management	p43
Rescue & Resuscitation	p45
Product Index	p47



SurgiVet® Advisor® Vital Signs Monitor



SurgiVet® Equator® Convective Warmer



SurgiVet® Universal CDS 9000
Small Animal Anesthesia Machine



SurgiVet® Classic T3™ Vaporizer

Monitoring Equipment

Common Features

PULSE OXIMETRY MONITORS

- Digital Pulse Oximetry with SAC technology-perfect for low perfusion and/or higher motion situations

ALL MONITORS

- Veterinary specific software



V1030

SurgiVet® V1030 Hand-Held Pulse Oximeter



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V1030	Hand-Held Pulse Oximeter	<ul style="list-style-type: none"> • 3311 5 ft. oximetry extension cable • V1700 Reflectance oximetry sensor • V1703 Universal "Y" oximetry sensor • Four (4) "AA" (LR6) alkaline batteries • Blue Protective Glove • A CD that contains the Operation Manual, Service Manual, and Quick Reference Guide. A printed Quick Reference Guide is also included. 	<ul style="list-style-type: none"> • Low perfusion performance • Motion tolerance • Pulse Amplitude Index (PI) bar-graph indicator to help identify and maintain a good test site • Memory storage for up to 99 patients • Quick selection of default alarm limit values by patient type: Cat, Dog, Horse, Factory and Custom • Audible and visual alarms • Multiple power choices (see p. 13) • 3 year warranty • Dimensions: 6.1" H x 3.3" W x 1.7" D (15.5 cm H x 8.4 cm W x 4.3 cm D) • Weight: 12 oz. (340 grams)



V3304

SurgiVet® V3304 Tabletop Pulse Oximeter



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V3304	Tabletop Pulse Oximeter	<ul style="list-style-type: none"> • 1613 Universal Power Supply with US plug AC/DC POWER SUPPLY 100-240VAC/12DC (US, JAPAN) • 1613EU Universal Power Supply with EU plug AC/DC POWER SUPPLY 100-240VAC/12DC (EURO) [OPTIONAL] • 1613UK Universal Power Supply with UK plug AC/DC POWER SUPPLY 100-240VAC/12DC (UK) [OPTIONAL] • 1613AU Universal Power Supply with AU plug AC/DC POWER SUPPLY 100-240VAC/12DC (Australia) [OPTIONAL] • V1700 Reflectance oximetry sensor • V1703 Universal "Y" oximetry sensor • 3111 5 ft Oximetry Extension Cable • V1857P Operations Manual 	<ul style="list-style-type: none"> • Digital SpO₂ and heart rate • Adjustable alarm limits • Visual and audible alarm signals • 99 patient memory • Ability to transfer real time data or trend data via serial cable connections • Optional printer is available (V3408) • Optional PC interface cable available (3362) • Dimensions: 3.24" H x 8.5" W x 5.5" D (21.6 cm H x 8.2 cm W x 14 cm D) • Weight: 30 oz. (850 grams)



V3395

SurgiVET® V3395 TPR MONITOR

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V3395	TPR Monitor	<ul style="list-style-type: none"> • V3396 Power Supply 120 VAC 9 VDC 1 A • V3397 Power Supply 220 VAC 9 VDC 1 A [Optional] • V1703 Universal "Y" Sensor • V1700 Reflectance Sensor • 3311 5 ft Oximetry Extension Cable • V3413 Temperature Cable • V3417 Temperature Probe • V3422 Respiration Sensor • WW1100 Airway Adapter, Straight • 20690B3 Protective Boot • V1907P Operations Manual 	<ul style="list-style-type: none"> • Digital SpO₂, heart rate, temperature, and respiration • Plethysmographic waveform display • Adjustable alarm limits • Visual and audible alarm signals • Optional infrared printer [8411] • Optional Reusable Esophageal/Rectal Temperature Probe [WWV3418] • Dimensions: 5.0" H x 4.38" W x 2.9" D [12.2 cm H x 11.1 cm W x 7.37 cm D] • Weight: 22 oz. (0.63 kg)



V60046

SurgiVet® V6004 Series Non-Invasive Blood Pressure Monitors



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH V6004 MONITOR	FEATURES
V60040	NIBP		
V60041	NIBP and Digital SpO ₂	<ul style="list-style-type: none"> • 1616 Universal Power Supply with US plug AC/DC POWER SUPPLY 100-240VAC/12DC (US, JAPAN) 	<ul style="list-style-type: none"> • Oscillometric non-invasive blood pressure technology • Systolic, Diastolic, and Mean NIBP values • Four available cuff sizes will comfortably fit a variety of patients
V60045	NIBP and Built-In Printer	<ul style="list-style-type: none"> • 1616EU Universal Power Supply with EU plug AC/DC POWER SUPPLY 100-240VAC/12DC (EURO) [OPTIONAL] 	<ul style="list-style-type: none"> • Two modes of use: SA mode for patients under 20 lbs and LA mode for patients over 20 lbs • Ability to transfer real time data or trend data via serial cable connections
V60046	NIBP, Digital SpO ₂ , and Built-In Printer	<ul style="list-style-type: none"> • 1616UK Universal Power Supply with UK plug AC/DC POWER SUPPLY 100-240VAC/12DC (UNITED KINGDOM) [OPTIONAL] • 1616AU Universal Power Supply with AU plug AC/DC POWER SUPPLY 100-240VAC/12DC (Australia) [OPTIONAL] • V1703 Universal "Y" Sensor (V60041 & V60046 Only) • V1700 Reflectance Sensor (V60041 & V60046 Only) • 3311 5 ft Oximetry Extension Cable (V60041 & V60046 Only) • V6148SA NIBP Cuff Kit [See pg. 9] • Built-in Printer with One Roll of Paper (V60045 & V60046 Only) • V1876P Operations Manual 	<ul style="list-style-type: none"> • Built-in printer or external printer (V3408) options available • Optional PC interface cable available [3362] • Base unit dimensions: 3.24" H x 8.5" W x 5.5" D [8.2 cm H x 21.6 cm W x 14.0 cm D] • Base unit weight: 3.5 lbs. (1.6 kg)

Monitoring Equipment



V8401

SurgiVet® V8401 Hand-Held Capnograph

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V8401	Hand-Held Capnograph	<ul style="list-style-type: none"> • 8403 Battery Eliminator 105-125 VAC, 60 Hz 	<ul style="list-style-type: none"> • Sidestream technology for application with intubated and non-intubated patients
V8400D	Hand-Held Capnograph w/ SpO ₂ Digital Pulse Oximetry	<ul style="list-style-type: none"> • 8406 Battery Eliminator 208-252 V (Optional) • 8404 Battery Charger 100-240 VAC 9 VDC 1.7 mA (Optional) • 8408 Battery (Optional) • 20690B2 Protective Boot • V1175 CO₂ Monitor Exhaust Kit • WW8214 Patient Attachment Kit Includes: <ul style="list-style-type: none"> • 8208 10-Gas Manifold Filters • 8211 10-4 ft Sample Lines • WW1100 10-Straight Airway Adapters • V1894P Operations Manual • V1700 Reflectance Sensor (V8400D Only) • V1703 Universal "Y" Sensor (V8400D Only) • V1894D Operations Manual (V8400D Only) 	<ul style="list-style-type: none"> • End tidal CO₂ numeric value, end tidal CO₂ waveform and respiration rate • Visual and audible alarm signals • Adjustable alarm limits • Optional infrared printer (8411) • Dimensions: 5.0" H x 4.38" W x 2.9" D (12.2 cm H x 11.1 cm W x 7.37 cm D) • Weight: 22 oz. (0.63 kg)
			RECOMMENDED ITEMS FOR MAINTENANCE <ul style="list-style-type: none"> • 8217 CO₂ Calibration Kit (5093, 8061, 8211, 8223) • 8208 10 Gas Manifold Filters



SurgiVet® V9004 Series Capnograph Monitor



V90041

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH 9004 MONITOR	FEATURES
V90040	Capnograph	<ul style="list-style-type: none"> • 1616 Universal Power Supply with US plug AC/DC POWER SUPPLY 100-240VAC/12DC (US, JAPAN) 	<ul style="list-style-type: none"> • Incorporates sidestream technology for application with intubated and non-intubated patients • End tidal CO₂ numeric value, end tidal CO₂ waveform, inspired CO₂, and respiration rate
V90041	Capnograph and Digital Pulse Oximetry	<ul style="list-style-type: none"> • 1616EU Universal Power Supply with EU plug AC/DC POWER SUPPLY 100-240VAC/12DC (EURO) [OPTIONAL] 	<ul style="list-style-type: none"> • Ability to transfer real time data or trend data via serial cable connections • Optional printer [V3408] • Optional PC interface cable available [3366]
V90043	Capnograph, Digital Pulse Oximetry, and FiO ₂ (Fractionally Inspired O ₂)	<ul style="list-style-type: none"> • 1616UK Universal Power Supply with UK plug AC/DC POWER SUPPLY 100-240VAC/12DC (UNITED KINGDOM) [OPTIONAL] • 1616AU Universal Power Supply with AU plug AC/DC POWER SUPPLY 100-240VAC/12DC (Australia) [OPTIONAL] • V1700 Reflectance Sensor [V90041 & V90043 Only] • V1703 Universal "Y" Sensor [V90041 & V90043 Only] • 3311 5 ft Oximetry Extension Cable [V90041 & V90043 Only] • V1175 CO₂ Monitor Exhaust Kit • 9193T FiO₂ Kit [V90043 Only] • Start-Up Kit Includes: <ul style="list-style-type: none"> • WW1100 6-Straight Airway Adapters • V1186 3-Nasal Sample Lines • 1178N 2-Aquaknot Moisture Traps • 8044 6-8 ft Sample Dead Space Adapters • V1182 2-3.5 mm Dead Space Adapters • 9048 1-13 mm CO₂ Disk Filter • V1872P Operations Manual 	<ul style="list-style-type: none"> • Dimensions: 3.5" H x 10" W x 5.5" D [8 cm H x 25.4 cm W x 13.9 cm D] • Weight: 4.8 pounds [2.2 kg]
V90041LF	Low Flow Capnograph and Digital Pulse Oximetry		<p>RECOMMENDED ITEMS FOR MAINTENANCE</p> <ul style="list-style-type: none"> • 8217 CO₂ Calibration Kit [5093, 8061, 8211, 8223] • 1178N 5-Aquaknot Moisture Traps • 9048 2-13 mm CO₂ Disk Filters

Monitoring Equipment



V9204

SurgiVet® Advisor® Vital Signs Monitor



- The Advisor® Vital Signs Monitor is the leading choice in multi-parameter monitoring. With its large, easy-to-read display and intuitive design, the Advisor® vital signs monitor is available in multiple configurations. The standard parameters include ECG, digital pulse oximetry, heart rate, non-invasive blood pressure and a built-in battery. Optional parameters include end tidal CO₂, respiratory rate, 2-channel temperature, 2-channel invasive blood pressure, and an integrated printer. The SurgiVet® Advisor® vital signs monitor features three monitoring modes that adjust ECG sensitivity, cuff inflation pressure and alarm limits that allow the monitor to be tailored to the physical requirement of the patient.

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH THESE ADVISOR® MODELS	FEATURES
V9200	3 Lead ECG, Digital SpO ₂ , Heart Rate, and NIBP	<ul style="list-style-type: none"> • WW3005 Power Cord 110 V 	<ul style="list-style-type: none"> • Four waveform color display screen
V9200IBP/Temp	3 Lead ECG, Digital SpO ₂ , Heart Rate, NIBP, IBP, and Temperature	<ul style="list-style-type: none"> • 3005A Power Cord 2 Prong-Europe (Optional) • V1700 Reflectance Sensor 	<ul style="list-style-type: none"> • Intuitive point and click user-interface • Ability to view all parameters trending information on screen
V9200P	3 Lead ECG, Digital SpO ₂ , Heart Rate, NIBP, and Built-In Printer	<ul style="list-style-type: none"> • V1703 Universal "Y" Sensor • 3311 Oximetry Extension Cable 	<ul style="list-style-type: none"> • Ability to transfer real time or trend data via serial cable connections
V9201	3 Lead ECG, Digital SpO ₂ , Heart Rate, NIBP, End Tidal CO ₂ , and Respiration	<ul style="list-style-type: none"> • V3110TL/V3110TLS 3 Lead ECG 3 Lead ECG • V6148SA NIBP Cuff Set (See pg. 9) 	<ul style="list-style-type: none"> • Optional PC interface cable available (3362) • Optional Veterinary Data Logger (WWW9220SYS) (See pg. 7)
V9201P	3 Lead ECG, Digital SpO ₂ , Heart Rate, NIBP, End Tidal CO ₂ , Respiration, and Built-In Printer	<ul style="list-style-type: none"> • 6111 Printer Paper (If a Printer is Ordered) • Operations Manual 	<ul style="list-style-type: none"> • Dimensions: 12" H x 13.25" W x 7" D (30.5 cm H x 34 cm W x 18 cm D) • Weight: 9-11 lbs. (4.1-5 kg) • Two year warranty
V9203	3 Lead ECG, Digital SpO ₂ , Heart Rate, NIBP, 2-Channel IBP, End Tidal CO ₂ , 2-Channel Temperature, and Respiration		
V9204	3 Lead ECG, Digital SpO ₂ , Heart Rate, NIBP, 2-Channel IBP, End Tidal CO ₂ , 2-Channel Temperature, Respiration, and Built-In Printer	<p>ADDITIONAL ITEMS INCLUDED WITH V9203 & V9204 MODELS</p> <ul style="list-style-type: none"> • V3417 Temperature Probe • V3413 Temperature Cable • WWV3418 Reusable Esophageal/Rectal Temperature Probe (Optional) • V6401 IBP Transducer Cable • V6402 Invasive Pressure Transducer • Start-Up Kit Includes: <ul style="list-style-type: none"> • 1178N 2-Aquaknot Moisture Traps (V9201 also) • 8044 6-8 ft Sample Lines (V9201 also) • V1186 3-Nasal Sample Lines (V9201 also) • WW1100 6-Straight Airway Adapters (V9201 also) • V1175 CO₂ Monitor Exhaust Kit (V9201 Also) • 9048 1-13 mm CO₂ Disk Filter • V1182 1-3.5 mm Dead Space Adapters 	<p>RECOMMENDED ITEMS FOR MAINTENANCE</p> <ul style="list-style-type: none"> • 8217 CO₂ Calibration Kit (5093, 8061, 8211, 8223) • 1179 2-CO₂ Absorber • 9048 2-13 mm CO₂ Disk Filter • 1178N 1-Moisture Trap <p>ADDITIONAL ITEMS INCLUDED WITH V9200 IBP/TEMP</p> <ul style="list-style-type: none"> • V3413 Temperature Cable • V3417 Temperature Probe • V6401 Transducer Cable • V6402 Invasive Pressure Transducer



WWV9230

SurgiVet® Advisor® 3 Parameter Vital Signs Monitor



- The Advisor® Multi-Parameter Vital Signs Monitor is also available in a basic three parameter configuration, providing veterinary specific vital signs measurements for non-invasive blood pressure, pulse oximetry, and ECG. With its large easy to read display and intuitive user interface, this configuration of the Advisor® vital signs monitor is ideal for the budget conscious buyer, offering the same high quality and reliability you've come to rely on from the SurgiVet® Advisor® vital signs monitor product family.

WWV9230	Advisor® 3 Parameter Monitor	<ul style="list-style-type: none"> • V1703 Universal "Y" Sensor • 3311 SpO₂ Extension Cable • V3110TL 3-Lead ECG Cable • V3110TLS ECG Lead Wires • 31543B2 Medium NIBP Cuff (5-15 cm inflation size) • 31594B1 6' NIBP Supply Hose • Power Cord • Operations Manual 	<ul style="list-style-type: none"> • Pulse oximetry with pulse rate • Non-invasive blood pressure • ECG with heart rate • 10.4" LCD with color monitor display • 2 year warranty on monitor • Dimensions: 12" H x 7" W x 13.25" D (30.5 cm H x 18 cm W x 34 cm D) • Weight: 9-11 lbs. (4.1-5 kg)
----------------	------------------------------	--	---

*Not Data Logger compatible



WWV9220

SurgiVet® Advisor® Veterinary Data Logger



PRODUCT CODE	DESCRIPTION	FEATURES
WWV9220SYS	Advisor® Veterinary Data Logger	<ul style="list-style-type: none"> • Record real time trend data at regular intervals • Capture a 30-second primary ECG waveform image • Capture a screen image • Upload data to a PC (not MACs); supports various formats (pdf, xls) • Supports Microsoft® Excel® 2002, 2003, 2007 • Choice of basic or detailed anesthesia record • SD (Secure Digital) flash card and card reader included • Dimensions: .89" H x 2.01" W x 2.30" D (22.5 mm H x 51.1 mm W x 58.5 mm D) • Weight: 38 grams (1.34 oz)

*Not compatible with WWV9230



V1703

Pulse Oximetry Sensors



- Pulse oximetry sensors work exceptionally well on pink perfused tissue such as a patient's tongue, toe webbing, vulva, prepuce, or ventral tail base.

PRODUCT CODE	DESCRIPTION	FEATURES
V1700	Reflectance Sensor	<ul style="list-style-type: none"> • The Reflectance Sensor can be applied to the ventral tail base of a patient and secured with a cohesive bandage. For cats and dogs, if the sensor is inserted rectally the red light should be shining on the sphincter muscle.
V1703	Universal "Y" Sensor	<ul style="list-style-type: none"> • The Universal "Y" Sensor is designed to be used on areas of thicker tissue. It is the standard clip used in most small animal practices.
V1707	"C" Sensor	<ul style="list-style-type: none"> • The "C" Sensor is perfect for cats, dogs, and horses. The sensor should be placed on the thigh, metatarsal, or metacarpal for smaller patients. The sensor should be placed over the Achilles tendon, tongue, prepuce, vulva, or through the toe webbing on large patients.
V3078	Mini Clip Sensor	<ul style="list-style-type: none"> • The Mini Clip Sensor can be used on a variety of animals ranging in size, but is ideal for small dogs, cats and some exotics.



V1700



V1707



V3078

Monitoring Equipment

Replacement Clips



PRODUCT CODE	DESCRIPTION	NOTE
V1709	10-Replacement "C" Clips	• For use with the V1707 sensor.
V1720	1-Replacement Universal "Y" clip	• For use with the V1703 sensor.



3311

5 Ft Oximetry Extension Cable



PRODUCT CODE	DESCRIPTION	NOTE
3311	5 ft Oximetry Extension Cable	• Provides the ability to use pulse oximeter sensors at farther distances from monitors.
3311L	15 ft Oximetry Extension Cable	• Provides the ability to use pulse oximeter sensors at farther distances from monitors.




V3110TLS

ECG Leads and Cables



- For use with V3404, V3404P, and Advisor® vital signs monitor.

PRODUCT CODE	DESCRIPTION
V3110TLS	3 Lead ECG Wires with Alligator Clips, 40"
V3110TLL	3 Lead ECG Wires with Alligator Clips, 60"
V3110TL	3 Lead ECG Cable (can be used with V3110TLS, V3110TLL)
V3406FL	5 Lead ECG Cable (can be used with V3110TLS, V3110TLL, V3406FLS, 3107)
V3406FLS	5 Lead ECG Wires with Alligator Clips, 40"
21-0460-51	ECG Cable 3-Lead Shielded SNAP
21-0461-51	ECG Wires 3-Lead 40" Shielded SNAP
21-0462-51	ECG Cable 5-Lead Shielded SNAP
21-0463-51	ECG Wires 5-Lead 40" Shielded SNAP



Cardio Companion

Cardio Companion

- The Cardio Companion is a device used to monitor core body temperature and ECG via an esophageal lead. Compatible with most ECG lead wires.

PRODUCT CODE	DESCRIPTION
V706000	Control Box
32491B2	Small Green Esophageal Lead
32491B3	Medium Blue Esophageal Lead
V706001	Cardio Companion with Small and Medium Lead
V706004	Cardio Companion with Two Small Leads
V706005	Cardio Companion with Two Medium Leads

Temperature and Respiration Cables



- For use with V3404P, V3395, and Advisor® Vital Signs Monitor.

PRODUCT CODE	DESCRIPTION
V3413	Temperature Cable
V3417	Temperature Probe (Disposable)
WWV3418	Reusable Esophageal/Rectal Temperature Probe
V3422	Respiratory Cable and Sensor (V3404P and V3395 Only)



V3413



V3417



WWV3418



V3422



V6148SA

Non-Invasive Blood Pressure Cuffs



- Fits patients comfortably because of a 360° bladder that eliminates pinching. The 360° bladder also eliminates the need to align an arrow/mark over an artery. The durable cuffs can withstand 10,000 inflations to 200 mmHg and are machine washable. A reinforced closure system ensures that the cuffs will stay on and can handle wear and tear.

PRODUCT CODE	DESCRIPTION	SIZE
31543B1	Small NIBP Cuff	3-9 cm
31543B2	Medium NIBP Cuff	5-15 cm
31543B3	Large NIBP Cuff	9-25 cm
31543B4	Extra Large NIBP Cuff	17-45 cm
31544B1	6 ft Hose	
31545B1	Cap	
V6148SA	NIBP Cuff Kit (Includes One of Each Listed Above)	

Optional Accessories For Non-Invasive Blood Pressure



PRODUCT CODE	DESCRIPTION
31544B2	3 ft Hose
31544B3	10 ft Hose

Invasive Pressure Accessories



- For use with V6400 and Advisor® Vital Signs Monitor.

PRODUCT CODE	DESCRIPTION
V6401	Transducer Cable
V6402	Invasive Pressure Transducer



V6402

Monitoring Equipment



8217

CO₂ Calibration



PRODUCT CODE	DESCRIPTION
8217	CO ₂ Calibration Kit: 19 Liter Tank of Calibration Gas, Regulator, and Sample Line Adapter (For use with V9004, V8401, and Advisor® Vital Signs Monitor)
5093	CO ₂ Calibration Gas Canister Only
8061	Calibration Gas Regulator Only
8223	Calibration Adapter (V9004 Series)

Agent Calibration



PRODUCT CODE	DESCRIPTION
8224	Agent Calibration Kit: 19 Liter Tank of Calibration Gas, Regulator, and Sample Line Adapter (V9400 Monitor)
1502	Agent Calibration Gas Canister (V9400 Monitor Only)
8061	Calibration Gas Regulator Only
8223	Calibration Adapter

CO₂ Accessories



PRODUCT CODE	DESCRIPTION
1178N	(5) Aquaknot Moisture Traps
1187	Aquaknot Adapter (V9004 and V9400 Monitors)
1179	(2) CO ₂ Absorber (Advisor® Vital Signs Monitor)
8208	10-Gas Manifold Filters (V8401 Monitor)
8075	Moisture Trap (Advisor® Vital Signs Monitor)
9048	(2) CO ₂ Disk Filter 13 mm (Rear Filter for V9004, V9400, and Advisor® Vital Signs Monitors)
9193T	FiO ₂ Sensor Kit: (1) 9190 FiO ₂ Sensor, 1-9181 Adapter "T", and (1) 9191 Cable (V90043 Monitor)



1178N



1187



1179



8208



8075



9048



9193T



V1175

CO₂ Monitor Exhaust Kit

PRODUCT CODE	DESCRIPTION
V1175	This kit will allow waste gas scavenging for CO ₂ monitors. Can be used with the Advisor® vital signs monitor, V8401, V9004, and V9400.

Sample Lines



PRODUCT CODE	DESCRIPTION
V1186	Nasal Sample Lines (qty 6)
8211	4 ft Sample Lines (qty 10)
8044	8 ft Sample Lines (qty 10)
WW1140	15 ft Extension Sample Lines (qty 10)



V1186



8211

Airway Adapters



- For use with CO₂ and Respiration Monitors.

PRODUCT CODE	DESCRIPTION
WW1100	Straight Airway Adapters 13 ID x 15 OD mm (qty 10)
1151	Pediatric Airway Adapters 13 ID x 15 OD mm (qty 10)
V1180	(2) 2.5 mm Dead Space Adapters
V1181	(2) 3.0 mm Dead Space Adapters
V1182	(2) 3.5 mm Dead Space Adapters



WW1100



1151



V1180

Monitoring Equipment



8411

Printers

PRODUCT CODE	DESCRIPTION
8411	Infrared Printer (Monitors: V3395, V3402, V8401, and V6400) *Power supply not included
WW1026SYS	Attachable Thermal Printer (attaches to the docking station) (Monitor: V1030)



WW1025SYS

Docking Stations

PRODUCT CODE	DESCRIPTION
WW1025SYS	Docking Station
WW10250US	Docking Station (OUS)
WW1028SYS	Docking Station and Printer System
WW10280US	Docking Station and Printer System (OUS)



Printer Paper

PRODUCT CODE	DESCRIPTION
3409	4 Pack (Printer: V3408) [2]
6012	4 Pack (Monitors: V60045, V60046) [1]
8416	4 Pack (Printer: 8411) (Not Pictured)
6111	4 Pack (Advisor® Vital Signs Monitor) [3]
WW1027	Printer Paper, 5 Pack (Not Pictured) (Printer: WW1026SYS)

Power Supplies

PRODUCT CODE	DESCRIPTION
1613	Universal Power Supply with US plug AC/DC POWER SUPPLY 100-240VAC/12DC (US, JAPAN) (V3304)
1613EU	Universal Power Supply with EU plug AC/DC POWER SUPPLY 100-240VAC/12DC (EURO) (V3304)
1613UK	Universal Power Supply with UK plug AC/DC POWER SUPPLY 100-240VAC/12DC (UNITED KINGDOM) (V3304)
1613AU	Universal Power Supply with AU plug AC/DC POWER SUPPLY 100-240VAC/12DC (Australia) (V3304)
1616	Universal Power Supply with US plug AC/DC POWER SUPPLY 100-240VAC/12DC (US, JAPAN) (V6004, V9004)
1616EU	Universal Power Supply with EU plug AC/DC POWER SUPPLY 100-240VAC/12DC (EURO) (V6004, V9004)
1616UK	Universal Power Supply with UK plug AC/DC POWER SUPPLY 100-240VAC/12DC (UNITED KINGDOM) (V6004, V9004)
1616AU	Universal Power Supply with AU plug AC/DC POWER SUPPLY 100-240VAC/12DC (Australia) (V6004, V9004)
WW1090	WW1090 Custom Lithium Ion Rechargeable Battery (30 hrs continuous use; 3hr recharge†) (V1030)
WW1095	Power Supply Universal with US/Japan Adapter 30 WATT, 100-240V, 50Hz (V1030)
WW1095OUS	Power Supply Universal with International Blades 30 WATT, 100-240V, 50Hz (V1030)
WW1096SYS	AC Power and Rechargeable Lithium Ion Battery System (V1030)
WW1096OUS	AC Power and Rechargeable Lithium Ion Battery System (OUS) (V1030)
WW1618	Power Supply 208-252 VAC, 13.2 VAC, 1 A (Monitors: V3404, V3404P, or V3408)
V3396	Power Supply 120 VAC, 9 VDC, 1 A (Monitor: V3395)
V3397	Power Supply 220 V, 9 VDC, 1 A (Monitor: V3395)
8404	Battery Charger 100-240 VAC, 9 VDC, 1.7 mA (Monitor: V8401) (For use with 8408 only)
8404A	Battery Charger 220 V, 9 VDC, 1.7 mA (Monitor: V8401) (For use with 8408 only)
8408	Rechargeable Battery (Monitor: V8401) (For use with 8404 or 8404A only)



8403

Battery Eliminator

PRODUCT CODE	DESCRIPTION
8403	Battery Eliminator 105-125 VAC, 60 Hz (Monitor: V8401)
8406	220 V Battery Eliminator 208-252 V (Monitor: V8401)



3005

Power Cords

PRODUCT CODE	DESCRIPTION
WW3005	Power Cord 110 V (Monitors: V6400 or Advisor® Vital Signs Monitor)
3005A	Power Cord 2 Prong-Europe (Monitors: V6400 or Advisor® Vital Signs Monitor)

Monitoring Equipment

Monitor Wall Mounting Options

- For use in areas where space is limited.

PRODUCT CODE	DESCRIPTION
3393	Advisor® Vital Signs Monitor Wall Mount
3394	Wall Mount Bracket (Monitor: V3304, V3404, V6004, V9004)
WW1098	Universal Mounting Bracket (Monitor: V1030)
WW1098C	Universal Mounting Bracket with IV Pole Clamp (Monitor: V1030)



3393



3394



8409

Monitor IV Pole Mounting Options

- Allows a monitor to be mounted to an IV pole.

PRODUCT CODE	DESCRIPTION
WW3356	IV Pole Bracket (Monitor: V3304, V3404, V6004, V9004)
8409	Pole Mount Bracket (Monitor: V8401, V3395)

Roll Stands

PRODUCT CODE	DESCRIPTION	
3387	SurgiVet® Advisor® Vital Signs Monitor Roll Stand	The Advisor® Vital Signs Monitor Roll Stand has two locking casters and a ten pound base for increased stability. A mounting plate and locking device allows for the easy attachment or removal of the monitor. A utility basket provides convenient storage. Dimensions: 45" H x 23" W.
3388	SurgiVet® Tabletop Monitor Roll Stand	This roll stand is similar to the SurgiVet® Advisor® Vital Signs Monitor Roll Stand, but it excludes the 10 pound base. Compatible with V3304, V3404, V3404P, V6004, and V9004 monitors.
7108	SurgiVet® Monitor Roll Stand	The 7108 Monitor Roll Stand provides a secure and inexpensive solution for transporting tabletop monitors. A plastic utility basket provides convenient storage and has a power cord access hole. Dimensions: 41" H x 23" W. Compatible with V3304, V3404, V3404P, V6004, and V9004 monitors.



3387



3388



7108



V707000

HOTLINE® Fluid Warmer



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V707000	HOTLINE® Fluid Warmer	<ul style="list-style-type: none"> • 32491B9 5 Warming Sets • 1 Operations Manual 	<ul style="list-style-type: none"> • Allows fluid or blood to be warmed to the appropriate range of 95-105° Fahrenheit • A patented patient line with triple lumen tubing maintains a layer of 107.6° Fahrenheit solution circulating around IV line all the way to the patient • Warm and ready to use in four minutes (from a cold start) • Accommodates a slow drip to a bolus • Can be mounted to an IV pole • Remanufactured • One-year warranty • Dimensions: 9.5" H x 8.3" W x 7.0" D • Weight: 7.95 lbs.



WWVEQ5000RS

SurgiVet® Equator® Convective Warming System



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
WWVEQ5000RS	Equator® Convective Warming System (Remanufactured)	<ul style="list-style-type: none"> • 32491B12 1 7 ft Large Air Delivery Hose • SW2001 1 Full Body Blanket (40" W x 80" L) • SW2002 1 Pediatric Full Body Blanket (40" W x 57.5" L) • SW2004 1 Lower Body Blanket (40" W x 64" L) • SW2009 1 Small Pediatric Underbody Blanket (26" W x 37" L) • SW2011 1 Large Pediatric Underbody Blanket (40" W x 54" L) • 1 Operations Manual 	<ul style="list-style-type: none"> • Lightweight • Compact size • Quiet air delivery • Can be mounted to an IV pole • Features a digital display and push button control for four temperature settings • Remanufactured unit: One year warranty • New unit: Three-year warranty • Dimensions: 11.75" H x 9.5" W x 7.5" D • Weight: 15 lbs.
WWVEQ5000	Equator® Convective Warming System (New)		

*Level 1® Snuggle Warm® Convective Warming Blankets-pg. 16

Accessories

PRODUCT CODE	DESCRIPTION
SW5-HOSE-N	5 ft Small Air Delivery Hose (Blue)
SW5-HOSE-7	7 ft Large Air Delivery Hose
F3-5000	Filter Kit

Temperature Management

Level 1® Snuggle Warm® Convective Warming Blankets



- The disposable blankets are constructed of two layers of non-woven polypropylene fabric coated with a layer of polyethylene. The blankets are designed to conform to the patient. Sold in ten packs only.

PRODUCT CODE	DESCRIPTION	SIZE	COMPATIBLE WITH	UNITS
SW2001*	Full Body Blanket	40" W x 80" L	Large Hose	10/pk
SW2002*	Pediatric Full Body Blanket	40" W x 57.5" L	Large Hose	10/pk
SW2003	Upper Body Blanket	80" W x 40" L	Large Hose	10/pk
SW2004*	Lower Body Blanket	40" W x 64" L	Large Hose	10/pk
SW2005	Neonate to Small Blanket	51" W x 28" L	Small Hose (Blue)	10/pk
SW2006	Preemie to Neonate Blanket	40" W x 22" L	Small Hose (Blue)	10/pk
SW2009*	Small Pediatric Underbody Blanket	26" W x 37" L	Large Hose	10/pk
SW2011*	Large Pediatric Underbody Blanket	40" W x 54" L	Large Hose	10/pk


*One sample blanket included with the SurgiVet® Equator® Convective Warming System.



SW2001*




SW2002*



SW2003



SW2004*



SW2005



SW2006



SW2009*



SW2011*

Common Features of All SurgiVet® CDS 9000 Small Animal Anesthesia Machines


- Pop-Off Valve designed to interface with an active waste gas evacuation system and maintain a passively full rebreathing bag
- Rebreathing and non-rebreathing capability
- Easy, quick change absorber canister
- Tapered 22 mm port on the inhalation/exhalation valve and bag port allow easy changing of the breathing circuit or bag
- Fresh gas port
- Adaptable to all styles of vaporizers



*Vaporizer not included WWV9000

SurgiVet® Universal CDS 9000 Small Animal Anesthesia Machine Pole Mount

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
WWV9000	Single Flowmeter (0.2-4 LPM)	<ul style="list-style-type: none"> • 32359B3 F-Circuit with 2 Liter Bag • 32359B6 Modified Jackson Rees with 1/2 Liter Bag • CO₂ Absorbent Refill Bag (2.2 lbs) • V7309 Pop-Off Valve • WWV1932 Operations Manual 	<ul style="list-style-type: none"> • Optional conversion kits available for flexibility • Pop-Off Occlusion Valve included • 3 piece assembly • 5 year warranty • Dimensions: 53.25" H x 26" W x 26" D • Weight: 23 lbs.
WWV9001	Dual Flowmeter (0-1,000 cc, 0.2-4 LPM)		
WWV9002	N ₂ O/O ₂ (N ₂ O: 0.2-4.5 LPM, O ₂ : 0.2-4 LPM)		
OPTIONAL ACCESSORIES			
<ul style="list-style-type: none"> • WWV90TT Tabletop Conversion Kit • WWV90WM Wallmount Conversion Kit • WWV90MB Monitor Bracket 			



WWV9010

SurgiVet® Universal CDS 9000 Small Animal Anesthesia Machine Tabletop

WWV9010	Single Flowmeter (0.2-4 LPM)	<ul style="list-style-type: none"> • 32359B3 F-Circuit with 2 Liter Bag • 32359B6 Modified Jackson Rees with 1/2 Liter Bag • CO₂ Absorbent Refill Bag (2.2 lbs) • V7309 Pop-Off Valve • WWV1932 Operations Manual 	<ul style="list-style-type: none"> • Optional conversion kits available for flexibility • Pop-Off Occlusion Valve included • Two built-in handles for ease of mobility • 5 year warranty • Dimensions: 14.5" H x 14.5" W x 10.5" D • Weight: 15 lbs.
WWV9011	Dual Flowmeter (0-1,000 cc, 0.2-4 LPM)		
WWV9012	N ₂ O/O ₂ (N ₂ O: 0.2-4.5 LPM, O ₂ : 0.2-4 LPM)		
OPTIONAL ACCESSORIES			
<ul style="list-style-type: none"> • WWV90PK Pole Mount Conversion Kit • WWV90WM Wall Mount Conversion Kit • WWV90MB Monitor Bracket 			

SurgiVet® Universal CDS 9000 Small Animal Anesthesia Machine Conversion Kit

PRODUCT CODE	DESCRIPTION	FOR USE WITH
WWV90TT	Tabletop Kit	Converts Pole Mount machines to Tabletop
WWV90WM	Wall Mount Kit	Converts Pole Mount machines to Tabletop
WWV90PK	Pole Mount Kit	Converts Tabletop machines to Pole Mount
WWV90MB	Monitor Kit	SurgiVet® Tabletop & Advisor® Vital Signs Monitors (V3304, V3404, V6004, V9230, V9200, V9400)

Anesthesia Equipment



V7025

Shown with optional vaporizer, V7275, V7276, V7277

SurgiVet® LDS 3000 Large Animal Anesthesia Machine

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V7025	LDS 3000	<ul style="list-style-type: none"> • 32387B3 2 64" Black Hose (2" diameter) • V7142 Stainless Steel "Y" Piece with Luer Fitting for End Tidal CO₂ Monitor Connection • 32346B7 30 Liter Bag • 5 CO₂ Absorbent Refill Bags (2.2 lbs) • Operations Manual 	<ul style="list-style-type: none"> • Easy to clean, chip proof powder coat finish • Dimensions: 53" H and 24" W diameter base with wheels • Four wheel base for maximum stability • Mobile
OPTIONAL ACCESSORIES			
		<ul style="list-style-type: none"> • V7140 Foal Absorber Canister • 32364B1 Large Animal Endotracheal Tube Funnel • V7114 Stand Alone Absorber Bracket Mount • V7275 DHV 1000 Ventilator • V7276 Cart (Pictured) • V7277 Shelf (Pictured) • V7138 Foal Circuit Set 	



V7276

Anesthesia and Ventilator Cart

PRODUCT CODE	DESCRIPTION	FEATURES
V7276*	Cart	The durable Anesthesia and Ventilator Cart is constructed of heavy gauge black powder coated steel. Five-inch caster wheels make the cart balanced and mobile. The cart is ideal for the LDS 3000 Large Animal Anesthesia Machine and DHV 1000 Large Animal Ventilator. Optional shelf is available.
V7277*	Shelf	

*V7025 Large Animal Anesthesia Machine and V7275 Large Animal Ventilator sold separately.

Common Features of All SurgiVet® Lab/Small Animal Anesthesia Machines

- Non-rebreathing system (ideal for patients 15 lbs. and under)

SurgiVet® Tabletop Research Anesthesia Machine



*Vaporizer not included V702002

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V702001	0.2-4 LPM O ₂ Flowmeter	<ul style="list-style-type: none"> • V1703 Rodent Circuit Set • 32359B2 Avian Research Circuit and 1/4 Liter Bag • Operations Manual 	<ul style="list-style-type: none"> • Great for feline only practices or research settings • Lightweight and portable
V702002	Zero-1,000 cc O ₂ Flowmeter		
V702003	Dual O ₂ Flowmeter (0-1,000 cc, 0.2-4 LPM)		
V702004	O ₂ , N ₂ O, and Fail-Safe (N ₂ O: 0-4.5 LPM, O ₂ : 0.2-4 LPM)		
V702005	Air and O ₂ Flowmeter (Air: 0-5 LPM, O ₂ : 0.2-4 LPM)		
V702006	Air Only Flowmeter (0-5 LPM)		

SurgiVet® Multi-Station Lab Animal Anesthesia Machine



V703001

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V704001	2 Station (Primary flowmeter 0-3 LPM, Station flowmeter 0-3 LPM)	<ul style="list-style-type: none"> • V7103 Rodent Circuit Set per Station • 1 36" Silicone Hose per Station • Operations Manual 	<ul style="list-style-type: none"> • Ideal for research • Multiple patients can be anesthetized at the same time • Each station can be controlled individually • Each station can deliver a different percentage of anesthetic without affecting any other station • Mobile
V703501	4 Station (Primary flowmeter 0-5 LPM, Station flowmeter 0-3 LPM)		
V703001	6 Station (Primary flowmeter 0-8 LPM, Station flowmeter 0-3 LPM)		

SurgiVet® Vaporstick Anesthesia Machine



V7015

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V7015	Single O ₂ Flowmeter (0.2-4 LPM)	<ul style="list-style-type: none"> • 32359B6 Modified Jackson Rees Non-Rebreathing System with 1/2 Liter Bag • 32359B2 Avian Research Circuit and 1/4 Liter Bag • V7121 "Z" Bracket for vaporizer • Operations Manual 	<ul style="list-style-type: none"> • Great for feline only practices or research settings • Can be used to provide supplemental oxygen for intensive care units in emergency or critical care situations • Compact design, mobile
V701501	Air and O ₂ Flowmeter (0-1,000 cc, 0.2-4 LPM)		
V701502	O ₂ , N ₂ O, and Fail-Safe (N ₂ O: 0-4.5 LPM, O ₂ : 0.2-4 LPM)		
V701504	Single O ₂ Flowmeter (0-1,000 cc)		

Anesthesia Equipment



V7321

V7321 Low Oxygen Pressure Alarm

PRODUCT CODE	DESCRIPTION	FEATURES
V7321	Low Oxygen Pressure Alarm	<ul style="list-style-type: none"> • Provides an audible alert when the oxygen supply is running low • Simple DISS fittings on the included oxygen hose make the alarm adaptable to any location such as a primary oxygen supply • Preset to alarm at 35 psi • Equipped with a test button • Utilizes a 9-volt battery that is included • Dimensions: 6.25" H x 3.75" W x 1.5" D



V7320

V7320 High Patient Pressure Alarm

PRODUCT CODE	DESCRIPTION	FEATURES
V7320	High Patient Pressure Alarm	<ul style="list-style-type: none"> • Provides an audible alert when a patient circuit reaches a dangerous pressure level • The alarm is adjustable, but is preset at 20 cm H₂O pressure • Can be used with rebreathing and non-rebreathing patient circuits • Equipped with a test button • Utilizes a 9-volt battery that is included • Dimensions: 6.1" H x 3.5" W x 1.2" D



Induction Chambers

Induction Chambers

PRODUCT CODE	DESCRIPTION	SIZE/MATERIAL	FEATURE
V711802	Small Induction Chamber	4" W x 4" H x 6" L (Acrylic)	• Calmly anesthetize patients in a stable container
V711801	Medium Induction Chamber	4" W x 4" H x 12" L (Acrylic)	



Breathing Bags

Breathing Bags

PRODUCT CODE	DESCRIPTION
32346B4	1/4 Liter Latex
32346B3	1/2 Liter Non-Latex
32346B1	1 Liter Non-Latex
32346B5	2 Liter Non-Latex
32346B6	3 Liter Non-Latex
32346B8	5 Liter Latex
32346B9	6 Liter Latex
32346B2	15 Liter Latex
32346B7	30 Liter Latex

Silicone Tubing

PRODUCT CODE	DESCRIPTION
32372B1	1/4" inner diameter medical grade silicone tubing is sold by the foot
32372B2	1/8" inner diameter medical grade silicone tubing is sold by the foot

Circuit Kit

PRODUCT CODE	DESCRIPTION	NOTE	ITEMS INCLUDED
V7130	Circuit Kit	The Circuit Kit is designed for new clinics	<ul style="list-style-type: none"> • 32359B3 1 Universal F-Circuit • 32359B3 1 Modified Jackson Rees Circuit • Breathing Bags (1 of each) 1/2 L, 1 L, 2 L, 3 L




32359B3

Universal F-Circuit

- The Universal F-Circuit is designed to maintain humidity and prevent heat loss. The F-Circuit fits almost any anesthesia machine designed to deliver rebreathing anesthesia.

PRODUCT CODE	LENGTH	INNER DIAMETER AT PATIENT END	OUTER DIAMETER AT PATIENT END	INNER DIAMETER AT VALVE END	OUTER DIAMETER AT VALVE END	ITEMS INCLUDED
32359B3	40"	15 mm	22 mm	22 mm	28 mm	2 Liter Bag
32359B4	60"	15 mm	22 mm	22 mm	28 mm	
32359B10	120"	15 mm	22 mm	22 mm	28 mm	

*[1] Flex Tubing-V7133, [2] F-Circuit Elbow-32343B10, [3] Adapter 22 mm Male x Female (Cuff)-32343B6, [4] 2 Liter Bag-32346B5




32235B1

Standard "Y" Circuit

- The Standard "Y" Circuit fits almost any anesthesia machine designed to deliver rebreathing anesthesia.

PRODUCT CODE	DESCRIPTION	LENGTH	INNER DIAMETER AT PATIENT END	OUTER DIAMETER AT PATIENT END	INNER DIAMETER AT VALVE END	OUTER DIAMETER AT VALVE END
32235B1	Pediatric	36"	15 mm	22 mm	22 mm	27 mm



32359B6

Modified Jackson Rees

PRODUCT CODE	DESCRIPTION	FEATURES
32359B6	Modified Jackson Rees	The Modified Jackson Rees is designed for a non-rebreathing system. The control valve allows you to close the system for manual ventilation and then reopen for spontaneous breathing. The waste gas is evacuated from the circuit by an evacuation hose with a 19 mm fitting.

*[1] Adapter Supply Line-32343B15, [2] 1/2 Liter Bag-32346B3, [3] 13 mm Male x 19 mm Female Adapter-32343B27



V710401

C-Pram Breathing Circuit

- The C-Pram is a non-rebreathing circuit designed to maintain humidity and reduce patient heat loss.

PRODUCT CODE	LENGTH	INNER DIAMETER AT PATIENT END	OUTER DIAMETER AT PATIENT END	INNER DIAMETER AT WASTE GAS END	OUTER DIAMETER AT WASTE GAS END	ITEMS INCLUDED
32135B1	60"	16 mm	22 mm	15 mm	19 mm	NA
V710401	120"	16 mm	22 mm	15 mm	19 mm	Supply adapter with hose to connect to fresh gas outlet on anesthesia machine

*[1] Adapter Supply Line-32343B15, [2] Silicone Tubing-32372B1, [3] 22 mm Male x 22 mm Male Adapter-32001B1, [4] 22 mm Female x 19 mm Adapter-32343B23

Anesthesia Equipment



32359B2

Avian Circuit

PRODUCT CODE	DESCRIPTION	FEATURES
32359B2	Avian Circuit	The Avian Circuit is a non-rebreathing circuit designed for small patients where dead air space is an issue. The control valve allows you to close the manual ventilation and then reopen for spontaneous breathing. Waste gas is evacuated from the circuit by an evacuation hose with a 19 mm fitting.

*[1] Adapter Supply Line-32343B15, [2] 1/4 Liter Bag-32346B4, [3] 13 mm Male x 19 mm Female Adapter-32343B27



V7103

Rodent Circuit Set

PRODUCT CODE	DESCRIPTION	FEATURES
V7103	Rodent Circuit Set	The Rodent Circuit Set includes a mask, diaphragm without a hole, breathing circuit, supply hose, adapters, supply line adapter, 3 ft silicone tubing, 22 mm male x 22 mm male adapter (32001B1), and a 22 mm female x 22 mm male (32343B23).

*[1] 22 mm Male x 22 mm Male Adapter-32001B1, [2] 22 mm Female x 19 mm Male Adapter-32343B23, [3] Rodent Mask-V7105, [4] Breathing Circuit, Silicone Tubing-32372B1 (not pictured), Supply Line Adapter-32343B15 (not pictured)



V7115

Mayo Tray

PRODUCT CODE	DESCRIPTION	FEATURES
V7115	Mayo Tray	<ul style="list-style-type: none"> Constructed of a powder coated aluminum frame with stainless steel tray Tray will accommodate a 1-2 inch pole Dimensions: 15" L x 11" W



Sodasorb® LF CO₂ Absorbent

Sodasorb® LF CO₂ Absorbent

PRODUCT CODE	DESCRIPTION	FEATURES
008855	Sodasorb® LF CO ₂ Absorbent Plastic Bottle (9.3 lbs.) (4.2 kg)	<ul style="list-style-type: none"> Permanent color change when absorbent is exhausted
008865	Sodasorb® LF CO ₂ Absorbent Refill Bag (2.2 lbs.) (1 kg)	<ul style="list-style-type: none"> Specifically formulated to prevent the degradation of anesthesia gases that can lead to problematic byproducts Exhibits high contrast color change, should dessication occur
008875	Sodasorb® LF CO ₂ Absorbent Refill Bag (2.2 lbs.) (1 kg)	<ul style="list-style-type: none"> Long-lasting formulation More efficient use of anesthetic gases Uniform pellet design results in less dust; better absorption.



Sodasorb® CO₂ Absorbent

Sodasorb® CO₂ Absorbent

PRODUCT CODE	DESCRIPTION	FEATURES
008870	Sodasorb® CO ₂ Absorbent Pre Pak (2.2 lbs.) (1 kg)	<ul style="list-style-type: none"> • Long-lasting formulation removes CO₂ for hours while maximizing the use of anesthetic agents • Fewer change-outs and low dust content makes cleaning easier • Specially formulated to help reduce the risk of degradation of anesthetic agents • Proven safety record ensures the comfort and well being of both patients and staff • High contrast color change lets you know when it's time to replace • May be used in conjunction with all the most commonly utilized anesthetic agents
008860	Sodasorb® CO ₂ Absorbent Refill Bag (2.8 lbs.) (1.3 kg)	
008880	Sodasorb® CO ₂ Absorbent 5 Gallon Pail (37 lbs.) (16.8 kg)	



V7306

SurgiVet® V7306 Pop-Off Valve

PRODUCT CODE	DESCRIPTION	FEATURES
V7306	Pop-Off Valve	The Pop-Off Valve maintains a constant passive volume within the anesthesia system allowing the patient to take a breath or the anesthetist to administer a breath. The rebreathing bag will stay passively full while maintaining less than 2 cm H ₂ O pressure. The valve will interface with a properly adjusted active waste gas evacuation system.



V7309

SurgiVet® V7309 Pop-Off Occlusion Valve

PRODUCT CODE	DESCRIPTION	FEATURES
V7309	Pop-Off Occlusion Valve	The Pop-Off Occlusion Valve helps to prevent an accidental pneumothorax due to leaving the Pop-Off Valve closed. By simply depressing a button on top, the valve temporarily blocks the gas flow so that a breath may be given to the patient. Upon release, the valve opens and no adjustment of the Pop-Off Valve is required. A 19 mm to 22 mm adapter (32343B23) is included.



V7114

Absorbers

PRODUCT CODE	DESCRIPTION	FEATURES
V7114	Stand Alone Absorber	The Stand Alone Absorber is an add-on absorber canister for small animals weighing 15-200 lbs that includes inhalation/exhalation valves, manometer, and a Pop-Off Valve. This item is designed for use with the V7025 SurgiVet® LDS 3000 Large Animal Anesthesia Machine.
V711401	Stand Alone Absorber Pole Mount	The Stand Alone Absorber Pole Mount is an add-on absorber canister for small animals weighing 15-200 lbs that includes inhalation/exhalation valves, manometer, and a Pop-Off Valve. The Stand Alone Absorber Pole Mount can fit a 1.5-2" pole.



V7128

Non-Rebreathing Adapter

PRODUCT CODE	DESCRIPTION	NOTE
V7128	Non-Rebreathing Adapter	By placing the adapter in-line, moving from rebreathing to non-rebreathing is quick and easy.

Anesthesia Equipment



V7149

Anesthesia Masks



PRODUCT CODE	DESCRIPTION	DIMENSIONS	REPLACEMENT DIAPHRAGM
32393B5	Small Feline Mask	2" x 2" x 1" diaphragm opening	32368B5
32393B3	Medium Feline Mask	2.25" x 3 5/8" x 1 1/8" diaphragm opening	32368B4
32393B2	Large Feline Mask	3.75" x 3 5/8" x 1 1/8" diaphragm opening	32368B2
32393B4	Small Canine Mask	3.5" x 3 5/8" x 1.25" diaphragm opening	32368B2
32393B1	Large Canine Mask	6" x 5.5" x 2.25" diaphragm opening	32368B1
V7149	5 Mask Set	Includes one of each of the above	
V7105	Rodent Mask	1.5" x 2.25"	No diaphragm included



Rodent Mask Diaphragms

Rodent Mask Diaphragms

PRODUCT CODE	SIZE
V710503	3/4" diaphragm opening
V710502	9/16" diaphragm opening
V710501	7/16" diaphragm opening
32247B1	Without opening (cut to desired size)



V7137

Dual Anesthesia Research Circuit

PRODUCT CODE	DESCRIPTION	FEATURES
V7137	Dual Anesthesia Research Circuit	<ul style="list-style-type: none"> Enables gas to be directed from one station to another or both



V7138

Foal Circuit Set

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED	NOTE
V7138	Foal Circuit Set	32426B5 1-"Y" Piece with Luer Fitting 32364B1 3-Large Animal Endotracheal Tube Funnels 32364B3 3-Funnel Valve Adapters 32346B6 1-3 Liter Bag 32359B5 2-60" 22 mm Hoses 32001B1 3-22 mm x 22 mm Adapters 32343B10 F-Circuit Elbow	Commonly used with the LDS 3000 Large Animal Anesthesia Machine. The hoses included in the set can also be used with SurgiVet® brand small animal anesthesia machines.



V7142

“Y” Piece with Luer Fitting

PRODUCT CODE	DESCRIPTION	FEATURES
V7142	“Y” Piece with Luer Fitting	<ul style="list-style-type: none"> • Stainless Steel “Y” Piece for large animal machine • Fits a 2” diameter circuit • Luer Fitting allows connection for ETCO₂ monitoring



32364B1

Large Animal Endotracheal Tube Funnel Adapter (Silicone)

32364B1	Large Animal Endotracheal Tube Funnel Adapter (Silicone)
---------	--

Miscellaneous Adapters

PRODUCT CODE	DIMENSIONS	MATERIAL	DESCRIPTION/USAGE
32001B1	22 mm Male x 22 mm Male	Clear Acrylic	SAV 2500 Ventilator
32036B1	15 mm Male x 1/4” hose barb	Silver Metal	Fresh Gas Outlet Elbow
32343B15	15 mm Male x 1/4” hose barb	Clear Acrylic	Adapter Supply Line/Modified Jackson Rees Circuit
32343B18	30 mm Female x 19 mm Male	Purple Plastic	SAV 2500 Ventilator/Waste Gas Port
32343B20	22 mm Female x 17 mm Female	Clear Silicone	SAV 2500 Ventilator
32343B23	22 mm Female x 19 mm Male	Translucent White Plastic	Adapts 19 mm tubing to F/Air Canister/V7309
32343B27	13 mm Male x 19 mm Female	Translucent White Plastic	Blue waste gas evacuation hose for MJR, Avian circuit, and AES
32343B6	22 mm Male x 22 mm Female (ribbed cuff)	Translucent White Plastic	Universal F-Circuit (connection at inhalation valve)/V7301 Lab Animal Evacuation Unit inlets
32345B3	22 mm Male x 22 mm Female (ribbed cuff) with 1/8” hose barb	Translucent White Plastic	High Patient Pressure Alarm Line Adapter
V7311	19 mm Male x 1/2” Female	White Plastic	AES (1/2” PVC adapter)
V7313	19 mm Male x 19 mm Male	Translucent White Plastic	Connection of one 19 mm tubing section to another
32426B5	22 mm Male x 22 mm Male x 22 mm Male Wye	Translucent White Plastic	15 mm ID/breathing circuits

Common Features of SurgiVet® CDS 9000 MR Conditional Small Animal Anesthesia Machines



- All SurgiVet® MR Conditional products have been independently tested by an outside agency and have been certified MR Conditional in a static magnetic field of 3-Tesla or less. Each MR conditional product is appropriately labeled according to ASTM standards for MR Conditional products. This allows users to easily identify products that are safe for this use. Contact Smiths Medical for specific technical and performance characteristics for each MR Conditional product.



WWV9000MR

SurgiVet® CDS 9000 MR Conditional Small Animal Anesthesia Machines



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
WWV9000MR	Single Flowmeter (0.2-4 LPM)	<ul style="list-style-type: none"> • 32359B3 F-Circuit with 2 Liter Bag 	<ul style="list-style-type: none"> • Pop-off occlusion valve included
WWV9001MR	Dual Flowmeter (0-1,000 cc, 0.2-4 LPM)	<ul style="list-style-type: none"> • 32359B6 Modified Jackson Rees with 1/2 Liter Bag 	<ul style="list-style-type: none"> • 5 year warranty
WWV9010MR	Single Flowmeter (0.2-4 LPM)	<ul style="list-style-type: none"> • CO₂ Absorbent Refill Bag (2.2 lbs) 	<ul style="list-style-type: none"> • Dimensions: (WWV9000MR, WWV9001MR) 14.5" H x 14.5" W x 10.5" Weight: 23 lbs.
WWV9011MR	Dual Flowmeter (0-1,000 cc, 0.2-4 LPM)	<ul style="list-style-type: none"> • V7309 Pop-Off Occlusion Valve • WWV1932 Operations Manual 	<ul style="list-style-type: none"> • Dimensions: (WWV9010MR, WWV9011MR) 14.5" H x 14.5" W x 10.5" Weight: 15 lbs.



VCT302MR

MR Conditional Vaporizers



PRODUCT CODE	DESCRIPTION	CLASSIC T ^{3™} FEATURES	100 SERIES FEATURES
VCT302MR	SurgiVet® Classic T3™ Isoflurane Funnel Fill Vaporizer	<ul style="list-style-type: none"> • Precision style, out of circuit, single agent use vaporizer • Temperature and flow compensated 	<ul style="list-style-type: none"> • Available in funnel fill style
VCT3K2MR	SurgiVet® Classic T3™ Isoflurane Key Fill Vaporizer	<ul style="list-style-type: none"> • Clearly labeled with the name of agent to be used • Site glass provides easy viewing of liquid volume 	<ul style="list-style-type: none"> • Built from all new parts based on the original Ohio style vaporizer
VCT305MR	SurgiVet® Classic T3™ Sevoflurane Funnel Fill Vaporizer	<ul style="list-style-type: none"> • Available as funnel fill or key fill style 	<ul style="list-style-type: none"> • Three year parts and labor warranty is included
VCT3K5MR	SurgiVet® Classic T3™ Sevoflurane Key Fill Vaporizer	<ul style="list-style-type: none"> • One-piece design fill assembly to eliminate leak points • Three year parts, labor and calibration warranty is included 	<ul style="list-style-type: none"> • One-year warranty on calibration • Annual calibration service is required by an authorized SurgiVet® vaporizer service center to maintain the warranty
V720100MR	SurgiVet® 100 Series Funnel Fill Isoflurane Vaporizer	<ul style="list-style-type: none"> • Funnel Fill Dimensions: 7.25" H x 5" W x 5.5" D (18.42 cm H x 12.7 cm W x 13.97 cm D) 	<ul style="list-style-type: none"> • Funnel Fill Dimensions: 6" L x 5" W x 8" H
V720104MR	SurgiVet® 100 Series Funnel Fill Sevoflurane	<ul style="list-style-type: none"> • Funnel Fill Weight: 16 lbs. (7.27 kg) • Key Fill Dimensions: 7.25" H x 5" W x 6" D (18.42 cm H x 12.7 cm W x 15.24 cm D) • Key Fill Weight: 16 lbs (7.27 kg) 	<ul style="list-style-type: none"> • Weight: 13 lbs. (5.91 kg)



32370B10

Waste Gas Collection



PRODUCT CODE	DESCRIPTION
32370B10	Enviro-Pure® Charcoal Filter Canister
V7305MR	Enviro-Pure® Charcoal Filter Bracket
32086B1	Waste Gas Evacuation Tubing 19 mm
32091B1	Waste Gas Evacuation Tubing 22 mm



V727000

MRI 2550 Small Animal Ventilator



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V727000	MRI 2550 Ventilator	<ul style="list-style-type: none"> • 32359B7 48" 22 mm Cuffed Hose • 32359B8 12" 22 mm Cuffed Hose • 24" O₂ Hose with DISS Female Connectors on each end • 32426B3 DISS Wye Connector (F x M x M Wye) • 32391B1 Test Lung • 124001 Patient Valve • Operations Manual 	<ul style="list-style-type: none"> • Pneumatically powered and controlled • Can be used as an anesthesia ventilator or a long-term respirator • Adjustable controls for inspiratory and expiratory time as well as inspiratory flow • Transportable • MRI compatible • Dimensions: 11.75" H x 9" W x 6.5" D • Weight: 16 lbs.
			OPTIONAL ACCESSORIES
			<ul style="list-style-type: none"> • V7265 MRI 2550 Mounting Bracket • V7266 Roll Stand

*MR Conditional ventilators will work with non MR conditional products. Optional bellows available-pg. 30.



V727001

Pneupac® Respiratory Ventilator



PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V727001	Pneupac® Respiratory Ventilator	<ul style="list-style-type: none"> • 32391B1 Test Lung • 124001 Patient Valve • 32359B8 12" 22 mm Cuffed Hose • 24" O₂ Hose with DISS Female Connectors on Each End • Operations Manual 	<ul style="list-style-type: none"> • Pneumatically powered and controlled • Designed for short or long term ventilation • Confidently control inspiratory and expiratory times, and inspiratory flow • Compressed gas source can be combined with room air • Compact design • Portable • MRI Compatible • Dimensions: 3.5" H x 9" W x 6.5" D • Weight: 9 lbs.

*MR Conditional ventilators will work with non MR conditional products. Optional bellows available-pg. 29.

Ventilators



V725000

SAV 2500 Small Animal Ventilator

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V725000	SAV 2500 Small Animal Ventilator	<ul style="list-style-type: none"> • 32411B16 Bellows Assembly • 32359B8 12" 22 mm Cuffed Hose • 32359B7 48" 22 mm Cuffed Hose • 32001B1 22 mm x 22 mm Adapter • 32343B20 22 mm x 17 mm Silicone Adapter • 32426B3 DISS Wye Connector (F x M x M Wye) • 32391B1 Test Lung • 32343B18 Adapter 30 mm x 19 mm • 32411B10 Power Pack with Cord, 100-220 V, 50/60 Hz • 24" O₂ Hose with DISS Female connectors on each end • Operations Manual 	<ul style="list-style-type: none"> • Electronically controlled, pneumatically powered, time cycled • Adjustable controls for inspiratory flow, inspiratory time and breaths-per-minute • Manual pneumatic breath button in the event of power outage • Portable • Dimensions: 11.5" H x 9.5" W x 12" D • Weight: 12 lbs.
OPTIONAL ACCESSORIES			
		<ul style="list-style-type: none"> • V7310 300 ml Bellows • V7312 3 Liter Bellows (Foal) • V7260 SAV 2500 Mounting Bracket 	



V7275

DHV 1000 Large Animal Ventilator

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V7275	DHV 1000	<ul style="list-style-type: none"> • 32387B5 46" Hose (Attach Bag Port to Ventilator, 2" Diameter) • 2 ft O₂ Hose with DISS Female Connectors on each end • 32411B10 Power Pack with Cord, 100-220 V, 50/60 Hz • 32423C7 Bellows • Operations Manual 	<ul style="list-style-type: none"> • Electronically controlled, pneumatically powered with an adjustable timed cycle • MRI model is pneumatically controlled, pneumatically powered and time cycled • Provides controlled ventilation to anesthetized patients over 250 lbs • Designed with descending bellows • Mount on a cart or shelf • Dimensions: 30" H x 16" W x 16" D • Weight: 45 lbs.
V7275MR	MR Conditional DHV 1000		
OPTIONAL ACCESSORIES			
		<ul style="list-style-type: none"> • V7276 Large Animal Anesthesia and Ventilator Cart • V7276MRI Large Animal Anesthesia and Ventilator Cart MR Conditional • V7025 Large Animal Anesthesia Machine 	



V7310

Pediatric Bellows

PRODUCT CODE	DESCRIPTION	FEATURES
V7310	Pediatric Bellows	<ul style="list-style-type: none"> • 300 ml bellows • Can be used with the MRI 2550 Ventilator and SAV 2500 Small Animal Ventilator • MR Conditional



V7312

Foal Bellows

PRODUCT CODE	DESCRIPTION	FEATURES
V7312	Foal Bellows	<ul style="list-style-type: none"> • 3 liter bellows • Can be used with the MRI 2550 Ventilator and SAV 2500 Small Animal Ventilator • MR Conditional

Vaporizers



VCT302

SurgiVet® Classic T3™ Vaporizer

PRODUCT CODE	DESCRIPTION	FEATURES
VCT302	Classic T3™ Isoflurane Funnel Fill Vaporizer	<ul style="list-style-type: none"> Precision style, out of circuit, single agent use vaporizer
VCT3K2	Classic T3™ Isoflurane Key Fill Vaporizer	<ul style="list-style-type: none"> Temperature and flow compensated Clearly labeled with the name of agent to be used
VCT305	Classic T3™ Sevoflurane Funnel Fill Vaporizer	<ul style="list-style-type: none"> Site glass provides easy viewing of liquid volume
VCT3K5	Classic T3™ Sevoflurane Key Fill Vaporizer	<ul style="list-style-type: none"> Available as funnel fill or key fill style One-piece design fill assembly with no solder joints eliminating leak points Three year parts, labor and calibration warranty is included Funnel Fill Dimensions: 7.25" H x 5" W x 5.5" D (18.42 cm H x 12.7 cm W x 13.97 cm D) Funnel Fill Weight: 16 lbs. (7.27 kg) Key Fill Dimensions: 7.25" H x 5" W x 6" D (18.42 cm H x 12.7 cm W x 15.24 cm D) Key Fill Weight: 16 lbs (7.27 kg)



V720100

SurgiVet® 100 Series Vaporizer

PRODUCT CODE	DESCRIPTION	FEATURES
V720100	Vaporizer 100 Series Funnel Fill Isoflurane	<ul style="list-style-type: none"> Available in funnel fill style Built from all new parts based on the original Ohio style vaporizer
V720104	Vaporizer 100 Series Funnel Fill Sevoflurane	<ul style="list-style-type: none"> Three year parts and labor warranty is included One-year warranty on calibration Annual calibration service is required by an authorized SurgiVet® vaporizer service center to maintain the warranty Funnel Fill Dimensions: 6" L x 5" W x 8" H Weight: 13 lbs. (5.91 kg)



Tech 4 Vaporizer*

Cleaning and Calibration

- To ensure proper calibration of a vaporizer, service is recommended once a year in a controlled environment
- Each vaporizer is carefully taken apart and every component inside is cleaned, gaskets replaced, and wick replaced if necessary
- Servicing Isoflurane or Sevoflurane Tec 3, Tec 4 and Ohio style vaporizers only
- One-year warranty on calibration
- Contact Service & Repair Department to set up this service

*Clean and Cal service only

Vaporizer Conversion

- Cleaning and calibration service with a conversion to a different agent
- Servicing Tec 3, Tec 4 and Ohio style vaporizers only
- Conversion to Isoflurane or Sevoflurane only
- One-year warranty on calibration
- Contact Customer Service to set up this service

Vaporizer Exchange Program

- Eliminates downtime of servicing vaporizer with agent change request
- Place an order for the same style of vaporizer and specify the agent type of current vaporizer and the requested agent type (Isoflurane or Sevoflurane) to be received
- A cleaned and calibrated, remanufactured vaporizer of the same style will be sent with the requested agent
- Shipping label provided for return of the original vaporizer
- Original vaporizer will not be returned
- One-year warranty on calibration
- This service available for an additional fee.
- Contact Customer Service to set up this service.



Ohio Vaporizer



Tech 3 Vaporizer

PRODUCT CODE	DESCRIPTION
V741000	Exchange Ohio Isoflurane Funnel Fill
V741004	Exchange Ohio Sevoflurane Funnel Fill
V740000	Exchange Tech 3 Isoflurane Funnel Fill
V740001	Exchange Tech 3 Isoflurane Key Fill
V740004	Exchange Tech 3 Sevoflurane Funnel Fill
V740005	Exchange Tech 3 Sevoflurane Key Fill
V740100	Exchange 100 Series Isoflurane Funnel Fill
V740104	Exchange 100 Series Sevoflurane Funnel Fill
WWVCT302X	Exchange Classic T3™ Isoflurane Funnel Fill
WWVCT305X	Exchange Classic T3™ Sevoflurane Funnel Fill
WVCT3K2X	Exchange Classic T3™ Isoflurane Key Fill
WWVCT3K5X	Exchange Classic T3™ Sevoflurane Key Fill

Vaporizer Accessories

PRODUCT CODE	DESCRIPTION
V7121	"Z" Bracket allows for vaporizer to be mounted to a 1.5-2" pole
V7122	Universal Vaporizer Mounting Bracket with End Caps allows for vaporizer to be mounted to a 1.5-2" pole
V712201	Universal Vaporizer Mounting Bracket without End Caps allows for vaporizer to be mounted to a 1.5-2" pole
V7123	Tec 3 Vaporizer Mounting Cup
V7132	Dual Vaporizer Mounting Bracket for Vaporstick Plus
32056B2	Female End Cap
32056B1	Male End Cap
32372B1	Silicone Tubing (sold by the foot)



V7121




V7122




V712201



32056B2, 32056B1



V7123



Manifolds

Manifolds

PRODUCT CODE	DESCRIPTION
V7214	Ohio Manifold
V7216	Tec 4 Manifold
32294B7	Tec 3 Manifold



32343B12

Anesthetic Key Filler

PRODUCT CODE	DESCRIPTION	FEATURES
32343B12	Isoflurane	<ul style="list-style-type: none"> Designed to be used exclusively with key fill vaporizers
32343B17	Sevoflurane	<ul style="list-style-type: none"> Agent specific adapter that prevents the inadvertent mixing of agents Only fits the bottle of anesthetic for which it is intended



32343B8


Anesthetic Anti-Spill Adapter

PRODUCT CODE	DESCRIPTION	FEATURES
32343B8	Isoflurane	<ul style="list-style-type: none"> Helps to minimize the loss of agent when filling a funnel style vaporizer
32343B9	Sevoflurane	<ul style="list-style-type: none"> Only fits the bottle of anesthetic for which it is intended

Hard Line Assemblies (Wall Outlets)

Recessed Style Wall Outlet

PRODUCT CODE	DESCRIPTION	FEATURES
32367B4	O ₂ Recessed Rough in Assembly	<ul style="list-style-type: none"> • Best suited for buildings being constructed (copper tubing inside of wall) • Two parts needed to have complete assembly: <ul style="list-style-type: none"> • Recessed Rough in Assembly • Finishing Plate • Two styles of Finishing Plates: DISS Male or Ohio Style Female Quick Connect • Unit is similar to an electrical outlet in appearance when installed
32367B3	Ohio O ₂ Recessed Finishing Plate	
32370B2	DISS O ₂ Recessed Finishing Plate	
32367B2	N ₂ O Recessed Rough in Assembly	
32367B1	Ohio N ₂ O Recessed Finishing Plate	
32370B1	DISS N ₂ O Recessed Finishing Plate	
32367B12	Recessed and Surface Mount Repair Kit	



32367B4



32367B3



32370B2



V7218

Oxygen Manifolds

PRODUCT CODE	DESCRIPTION	FEATURES
V7218	O ₂ Manifold, 2 DISS Male Check Inlets with a 1/4" NPT Male Outlet and a Pressure Relief Valve	<ul style="list-style-type: none"> • Allows for two H-tanks to supply O₂ to facility

Soft Line Assemblies (Hose Assemblies)

Pre-Assembled Hoses



V7162

PRODUCT CODE	DESCRIPTION
V7162	3 ft O ₂ Hose DISS Female x DISS Female
V7156	6 ft O ₂ Hose DISS Female x DISS Female
V7157	8 ft O ₂ Hose DISS Female x DISS Female
V7158	10 ft O ₂ Hose DISS Female x DISS Female
V7163	15 ft O ₂ Hose DISS Female x DISS Female
V7164	20 ft O ₂ Hose DISS Female x DISS Female
V7159	6 ft O ₂ Hose Ohio Male x DISS Female
V7160	8 ft O ₂ Hose Ohio Male x DISS Female
V7161	10 ft O ₂ Hose Ohio Male x DISS Female

Gas Delivery Management

Hose Assemblies

- Custom hose assemblies, including ceiling hose drops
- Variety of quick connects available
- Custom hose orders require the following:
Part Number of one End by Part Number of Hose with length indicated by Part Number of second end plus hose assembly fee (VHASS)
- Ceiling hose drop orders may need to add Ceiling Plate (V716601) to conceal opening in ceiling
- Common pre-assembled hoses available (see pg. 33)


Gas Hoses

- The gas hoses are nonconductive hose and sold by the foot.


PRODUCT CODE	DESCRIPTION
32387B2	Green Oxygen Hose
32387B4	Blue Nitrous Hose
32387B11	Yellow Air Hose
32387B13	White Vacuum Hose
32387B12	Purple Evacuation Hose
32387B14	Gray Carbon Dioxide Hose

Air Hose Connectors (Yellow)

PRODUCT CODE	DESCRIPTION
32370B10	DISS Female x Hose Barb
32406B3	Ohio Male Quick Connect x Hose Barb
32164B2	DISS Male x 1/8" NPT Male x Hose Barb




32370B10



32406B3

Carbon Dioxide Hose Connectors

PRODUCT CODE	DESCRIPTION
32370B7	DISS Female x Hose Barb
32164B3	DISS Male x 1/8" NPT Male



32370B7

Evacuation Hose Connectors (Purple)


PRODUCT CODE	DESCRIPTION
32407B14	Chemtron Male Quick Connect x Hose Barb



32407B14

Miscellaneous Hose Connectors


PRODUCT CODE	DESCRIPTION
32386B1	1/4" NPT Female x Hose Barb
32386B2	1/4" NPT Male x Hose Barb
32386B8	1/4" Hose Coupler, 1/4" Hose Barb x 1/4" Hose Barb
32397B1	1/4" NPT Male x 1/4" NPT Male 90°
32397B7	1/4" NPT Male x Female (Elbow)
32343B31	1/4" NPT Male x Two 1/4" NPT Female Branch Tee
32417B2	1/4" Street T
32501B55	1/4" NPT Female Coupler
32032B1	1/8" NPT Male x Hose Barb
32035B1	1/8" NPT Male x Hose Barb 90°
32371B1	1/8" NPT Male x 1/8" NPT Female 45° elbow
32379B1	1/8" NPT Female x 1/8" NPT Female
32386B6	1/8" NPT Female x Hose Barb
32397B2	1/8" NPT Female Coupler
32397B3	1/8" NPT Male x Female (Elbow)
32397B4	1/8" NPT Male x 1/8" NPT Male 90°
32397B5	1/8" NPT Male x Hex Socket Plug
32417B1	1/8" Street T
32427B1	Wye Three 1/8" NPT Female Inputs
32209B1	Colder Male x Hose Barb
32366B1	Cross Connector Four 1/8" NPT Female



32386B1



32386B2



32035B1



32397B3




32417B1


Gas Delivery Management

Nitrous Hose Connectors (Blue)


PRODUCT CODE	DESCRIPTION
32370B11	DISS Female x Hose Barb
32386B5	DISS Male x Hose Barb
32395B1	DISS Male x 1/8" NPT Male
32406B1	Ohio Female x DISS Female
32406B2	Ohio Female Quick Connect x Hose Barb
32406B4	Ohio Male Quick Connect x Hose Barb
32406B7	Ohio Female x DISS Male
32406B8	Ohio Male x DISS Male




32370B11



32386B5






32406B2



32406B4

PRODUCT CODE	DESCRIPTION
32370B12	DISS Female x Hose Barb
32228B2	Ohio Female Quick Connect x Hose Barb
32228B1	Ohio Male Quick Connect x Hose Barb
32142B1	Chemtron Male Quick Connect x Hose Barb




32228B1




32142B1

Oxygen Hose Connectors With Barbs (Green)


PRODUCT CODE	DESCRIPTION
32370B5	DISS Female x Hose Barb
32370B8	DISS Male x Hose Barb
32407B5	Ohio Female Quick Connect x Hose Barb
32407B6	Ohio Male Quick Connect x Hose Barb
32407B7	Puritan Female Quick Connect x Hose Barb
32407B8	Puritan Male Quick Connect x Hose Barb
32363B8	Schrader Female Quick Connect x Hose Barb
32407B10	Schrader Male Quick Connect x Hose Barb
32407B1	Chemtron Female Quick Connect x Hose Barb
32407B2	Chemtron Male Quick Connect x Hose Barb
32383B2	DISS Handi-Tight x Hose Barb
V716601	Ceiling Plate




32370B5



32370B8




32407B5




32407B6




32407B7



32407B8



32363B8




32407B10




32407B1



32407B2



32383B2




V716601


Gas Delivery Management

Oxygen Hose Connectors (Green)


PRODUCT CODE	DESCRIPTION
32034B1	DISS Male x 1/8" NPT Male
32164B4	DISS Male x Male
32370B9	DISS Male x 1/4" NPT Male
32370B6	DISS Male Demand Valve x 1/8" NPT Male
32370B14	DISS Male Check Valve x 1/4" NPT Male
32421B5	DISS Male Check Valve x 1/8" NPT Male
32343B14	DISS Female Swivel Adapter x 1/8" NPT Male
32386B9	DISS Cap
32407B3	Ohio Female x DISS Female
32407B4	Ohio Male x DISS Female
32407B11	Ohio Female x DISS Male
32407B12	Ohio Male x DISS Male
32407B13	Ohio Male x 1/8" NPT Male
32363B3	Chemtron Male x DISS Male
32363B4	Chemtron Male x DISS Female
32363B6	Chemtron Female x DISS Female
32363B5	Chemtron Female x DISS Male
32363B7	Schrader Male x DISS Male
32407B9	Schrader Female x DISS Male
32406B5	Puritan Male x 1/4 NPT Male
32242B1	Puritan Male x DISS Male




32034B1



32164B4



32370B6



32421B5




32343B14




32386B9



32407B3



32407B4



32407B11



32407B12



32407B13




32363B3



32363B4



32363B6




32363B7




32242B1

E-Tank Accessories

PRODUCT CODE	DESCRIPTION
32426B3	Wye DISS Female x Two DISS Male
32426B4	Wye Ohio Male x Two Ohio Female
V7338	Wye DISS Female x Two DISS Male with Checks
V733803	Wye DISS Male x Two DISS Male with Checks
V711301	Supply T with Demands
V7113	Supply T Dual O ₂ with Checks (For use with two oxygen sources)
32426B2	Wye Chemtron Male x Two Chemtron Female
32426B6	Wye Puritan Male x Two Puritan Female



32426B3



32426B4

Gas Delivery Management

H-Tank Accessories

- All single stage regulators are preset to 50-55 psi.

PRODUCT CODE	DESCRIPTION
V7327	O ₂ H-Tank Regulator
V7334	N ₂ O H-Tank Regulator
V7329	O ₂ H-Tank Regulator and 0-15 lpm Flowmeter
32408B17	O ₂ H-Tank Dual Stage Adjustable Regulator 0-100 psi



V7327



V7329



V7331

Universal Double E-Tank Manifold

PRODUCT CODE	DESCRIPTION	FEATURES
V7331	Universal Double E-Tank Manifold	<ul style="list-style-type: none"> • Internal check valves to allow one or both tanks to be turned on • Ensures anesthesia machine maintains mobility throughout clinic or hospital • Recommend use of V7261 with this manifold for extra stability of tanks



V7261

E-Tank Stabilizer

PRODUCT CODE	DESCRIPTION	FEATURES
V7261	E-Tank Stabilizer	<ul style="list-style-type: none"> • For use with V7331 to ensure stability of E-tanks • Secures tanks in place and eliminates movement which may help reduce leaks from developing and extend the life of the V7331



V7330

E-Tank Holders and Regulator

PRODUCT CODE	DESCRIPTION	FEATURES
V7330	O ₂ E-Tank Regulator, Double Tank Holder, 2 ft O ₂ Hose, E-Tank Wrench	<ul style="list-style-type: none"> • Two sets of rings that slide over the anesthesia machine pole and secures with set screws
V7333	N ₂ O E-Tank Regulator, Double Tank Holder, 2 ft N ₂ O Hose, E-Tank Wrench	<ul style="list-style-type: none"> • One E-tank regulator, wrench and 2 ft hose included • Tight grouping eliminates lateral swing movements • Tanks positioned to rest on legs of base
V7337	Air E-Tank Regulator, Double Tank Holder, 2 ft Air Hose, E-Tank Wrench	
V7335	E-Tank Rings Only	




V7326



32425B1



32066B1



E-Tank Accessories

- All single stage regulators are preset to 50-55 psi.

PRODUCT CODE	DESCRIPTION
V7326	O ₂ E-Tank Regulator
V7332	N ₂ O E-Tank Regulator
V7336	Air E-Tank Regulator
V7328	E-Tank 0-15 lpm Flowmeter
32425B1	E-Tank Wrench
32408B7	E-Tank Regulator Yoke
32408B10	E-Tank Yoke T-Handle
32257B2	E-Tank Yoke Washer

Respiratory Flowmeters

PRODUCT CODE	DESCRIPTION	FEATURES
32066B1	Respiratory Flowmeter 0-3 lpm	• Each has 1/8" NPT Female input fitting on back and DISS Male output fitting below
32375B6	Respiratory Flowmeter 0-5 lpm	
32375B5	Respiratory Flowmeter 0-15 lpm	

Flowmeter Replacement Parts

PRODUCT CODE	DESCRIPTION
32413B2	(A) O ₂ Needle Valve and Knob with Torque Guard
32413B5	O ₂ Needle Valve and Knob (SurgiVet® LDS 3000)
32413B6	N ₂ O Needle Valve and Knob with Torque Guard
32413B7	(A) Flowmeter O ₂ Knob
32413B4	(B) Single Flowmeter Lens Cover
32413B1	(C) O ₂ Thorpe Tube 0.2-4 lpm
32413B3	Flowmeter Adapter, Dual Threaded, Rear Connector
32375B14	Respiratory O ₂ Flowtube Housing

Gas Delivery Management



V7144

Humidifier System

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED	FEATURES
V7144	Humidifier System	<ul style="list-style-type: none"> • 32375B5 0-15 LPM Respiratory Flowmeter • 32388B1 Humidifier Chamber • V7131 Mounting Apparatus • Installation Instructions 	<ul style="list-style-type: none"> • Mounts to a wall or cage • Standard Male O₂ DISS fitting to connect 50 psi oxygen supply • Tapered hose connection on the humidifier chamber to connect cannula line or standard silicone tubing • Oxygen hose not included

Demand Valves

PRODUCT CODE	DESCRIPTION	FEATURES
32421B12	Small Animal Demand Valve	<ul style="list-style-type: none"> • Designed to provide 100% oxygen to breathing or nonbreathing patients
32421B11	Large Animal Demand Valve	<ul style="list-style-type: none"> • Can be used with mask, ET tube, or Tracheotomy tube • Delivery pressure is 60 +/- 5 cm H₂O • Small Animal Demand Valve delivers 40 LPM • Large Animal Demand Valve delivers 160 LPM



V7225

Anesthesia Detection Badge

PRODUCT CODE	DESCRIPTION	FEATURES
V7225	Single Badge	<ul style="list-style-type: none"> Test individuals or a facility for exposure to anesthetic gas
V7226	Five Pack	<ul style="list-style-type: none"> The badge can be used for fifteen minutes or up to a maximum of eight hours in one day The prepaid test includes confidential results that are mailed back within 14 business days



V7300

SurgiVet® AES Gas Evacuation System

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V7300	110 Volt System	<ul style="list-style-type: none"> 32292B1 3 ft and 10 ft Precut Hose (13 mm Diameter tubing with 13 mm x 19 mm adapters) V7311 1/2" PVC Adapter Operation Manual 	<ul style="list-style-type: none"> Product safely removes waste gas associated with everyday inhalation anesthesia Waste gas is collected by passive flow to the AES, then is actively pushed up to 150 feet away Two 19 mm inlets and one exhaust port allows for easy installation Can be mounted onto a wall or anesthesia machine pole Dimensions: 14" H x 5" D Weight: 8 lbs.
V7300220	220 Volt System		
OPTIONAL ACCESSORIES			
		<ul style="list-style-type: none"> V7107 SurgiVet® AES Wall Mount Bracket 	



32343B18

Waste Gas Evacuation Accessories

PRODUCT CODE	DESCRIPTION
32086B1	19 mm Tubing (per foot) (blue corrugated)
32091B1	22 mm Tubing (per foot)
32137B1	Evacuation Tubing Cap (For use with V7166) (white corrugated)
32343B18	Adapter 30 mm x 19 mm (Purple)
32343B23	Adapter 22 mm Female x 19 mm Male
V7166	Ceiling Plate Cover for Evacuation Tubing
V7313	19 mm x 19 mm Male Adapters



32343B23

Waste Gas Management



32373B10

SurgiVet® Enviro-Pure® Filter Canister

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH WWV7318	FEATURES
WWV7318	Enviro-Pure® Charcoal Filter Kit	<ul style="list-style-type: none"> • 1 Enviro-Pure® Charcoal Filter Canister 	<ul style="list-style-type: none"> • Designed to absorb waste anesthetic gas • Larger canister contains 1/3 more charcoal than comparable canisters • Enviro-Pure® Filter canister absorbs twice the agent of typical canisters • Gram for gram, Enviro-Pure® Filter is more cost effective than comparable canisters of charcoal
32373B10	SurgiVet® 6 Pack Enviro-Pure® Charcoal Filter Canisters	<ul style="list-style-type: none"> • V7305 Enviro-Pure® Charcoal Filter Bracket • 32086B1 4 ft Scavenger Tubing 	
			OPTIONAL ACCESSORIES <ul style="list-style-type: none"> • V7305 Enviro-Pure® Charcoal Filter Bracket (For use with Enviro-Pure® Charcoal Filter Canisters Only)

Scavenger Interface Valve For Active Systems

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V7308	Scavenger Interface Valve for Active Systems	<ul style="list-style-type: none"> • 32346B5 2 Liter Rebreathing Bag • 1 U-bolt for mounting 	<ul style="list-style-type: none"> • Improves the control of active draw when using an active suction system • This ensures adequate volume is available for the patient breathing circuit



V7301

SurgiVet® Lab Animal Evacuation System

PRODUCT CODE	DESCRIPTION	ITEMS INCLUDED WITH EACH	FEATURES
V7301	SurgiVet® Lab Animal Evacuation System (110 V Only)	<ul style="list-style-type: none"> • Startup quantity of Activated Charcoal • Operation Manual 	<ul style="list-style-type: none"> • Ideal for facilities without a built-in evacuation system • The system may be vented to the outside via a two inch outlet or activated charcoal may be used to scavenge waste gas • Can be used with the Multi-Station Anesthesia Machine that accommodates up to 6 stations • Dimensions: 9.75" H x 18" W x 9.75" D • Weight 27 lbs.



V7109

Evacuation Charcoal

PRODUCT CODE	DESCRIPTION
V7109	1 Gallon Charcoal



Resuscitation Bags

Resuscitation Bags



PRODUCT CODE	DESCRIPTION	OPTIONAL ACCESSORIES
32346B12	Small Resuscitation (240 ml Volume, Under 7 kg)	WWV7267 Peep Valve 32027B5 Manometer
32346B11	Medium Resuscitation (500 ml Volume, 7-30 kg)	
32346B10	Large Resuscitation (1600 ml Volume, Over 30 kg)	



V7148

Recovery Oxygen Masks



PRODUCT CODE	DESCRIPTION	DIMENSIONS	FEATURES
32490B15	Small Canine Recovery Mask	3" H x 3.5" W	<ul style="list-style-type: none"> Dual vents enable unrestricted inhalation and exhalation Durable polycarbonate construction Supplied with 7 ft oxygen tubing 22 mm oxygen adapter Ideal for rescue operations Latex free
32490B10	Large Canine Recovery Mask	4.5" H x 5.5" W	
32490B16	Feline Recovery Mask	2.5"H x 3.75" W	
V7148	Recovery Mask Kit (Includes One of Each Mask)		

Replacement Mask Diaphragms

PRODUCT CODE	DESCRIPTION
32368B1	Diaphragm Replacement for 32490B10
32368B2	Diaphragm Replacement for 32490B15, 32490B16

Product Index

A			
Absorbents.....	22	CDS 9000 Small Animal	
Adapters		Anesthesia Machines.....	17-26
Airway.....	11	Ceiling Plate.....	37
Large Animal Endotracheal		Chamber, Induction.....	20
Tube Funnel Adapter.....	25	Charcoal, Evacuation.....	44
Miscellaneous Plastic.....	25	Circuits	
Advisor® Vital Signs Monitor.....	6	Non-rebreathing.....	21-22
AES Gas Evacuation System.....	43	Rebreathing.....	21-22
Air		Classic T3™ Vaporizer.....	26, 30
Connectors.....	34	Carbon Dioxide	
Hose.....	34	Absorbent.....	22-23
Regulator.....	40	Connectors.....	34
Airway Adapters.....	11	Hose.....	34
Alarms		Monitor Exhaust Kit.....	11
Low Oxygen Pressure.....	20	Regulator.....	40
High Patient Pressure.....	20	Convective Warmer Accessories....	15
Anesthesia		Convective Warmer.....	15
Absorbent.....	22-23	C-Pram Breathing Circuit.....	21
Bags.....	20	“C” Sensor.....	7
Breathing Circuits.....	21-22		
Resuscitation.....	45-46	D	
Cart.....	18	Data Logger, Veterinary.....	7
Chambers, Induction.....	20	Demand Valves, Oxygen.....	42
Detection Badge.....	43	Detection Badge, Anesthesia.....	43
Machines		DHV 100 Large Animal Ventilator...28	
Research/Pre-Clinical.....	19	Diaphragm, Mask.....	24
Large Animal.....	18	Docking Station.....	12
Small Animal.....	17-19	Dual Anesthesia	
Masks		Research Circuit.....	24
Anesthesia.....	24		
Diaphragm.....	46	E	
Recovery Oxygen.....	46	ECG Leads and Cables.....	8
Rodent.....	24	Equator® Convective	
Non-Rebreathing Circuits.....	21-22	Warming System.....	15
Ventilators.....	28-29	Enviro-Pure® Charcoal	
Anesthetic Anti-Spill Adapter.....	32	Filter Canister.....	27, 44
Anesthetic Key Filler.....	32	Enviro-Pure® Charcoal	
Avian Circuit.....	22	Filter Bracket.....	27, 44
		Enviro-Pure® Charcoal Filter Kit....44	
		E-Tank	
		Accessories.....	41
		Evacuation	
		Charcoal.....	44
		Connectors.....	35
		Hose.....	34
		Lab System.....	44
		Exhaust Kit, CO ₂ Monitor.....	11
		F	
		F-Circuit, Universal.....	21
		Flowmeters	
		Replacement Parts.....	41
		Respiratory.....	41
		Foal	
		Circuit Set.....	24
		G	
		Gas Hoses.....	34
		Graphics Printer.....	12
		H	
		High Patient Pressure Alarm.....	20
		Hoses, gas.....	34
		HOTLINE® Fluid Warmer.....	15
		H-Tank	
		Accessories.....	40
		Humidifier System.....	42
		I	
		Induction Chamber.....	20
		Infrared Printer.....	12
		Invasive Pressure Accessories.....	9
		L	
		Lab Evacuation System.....	44
		Large Animal Endotracheal	
		Tube Funnel Adapter.....	25
		LDS 3000 Large Animal	
		Anesthesia Machine.....	18
		Low Oxygen Pressure Alarm.....	20
		M	
		Manifolds	
		Oxygen.....	33
		Vaporizer.....	32
		Masks	
		Anesthesia.....	24
		Diaphragm.....	45
		Recovery Oxygen.....	45
		Rodent.....	24
		Mayo Tray.....	22
		Mini Clip Sensor.....	7
		Modified Jackson Rees.....	21
		Monitor	
		Accessories.....	7-14
		Capnography.....	4-6
		Multi-Vital Signs.....	6-7
		Non-invasive Blood Pressure.....	3
		Pulse Oximetry.....	2-3, 6-7
		Temperature.....	2-3, 6
		MR Conditional.....	26-27
		MRI 2550 Ventilator.....	27
		Multi-Station Lab Animal	
		Anesthesia Machine.....	19
		N	
		Nitrous	
		Connectors.....	36
		Hose.....	34
		Regulator.....	40-41
		NIBP Cuffs.....	9
		Non-Rebreathing Adapter.....	23

O

Oximetry Extension Cable..... 8
 Oxygen
 Connectors..... 37-39
 Demand Valves..... 42
 Hoses 34
 Manifolds..... 33
 Masks, Recovery 45
 Regulators..... 40-41
 Transfiller..... 41

P

Pneupac® Respiratory Ventilator 27
 Pop-Off Valve..... 23
 Pop-Off Occlusion Valve..... 23
 Power Cords..... 13
 Power Supplies 13
 Pre-Assembled Oxygen Hoses 33
 Printer
 Infrared 12
 Paper..... 12
 Thermal..... 12
 Probes, Temperature 9
 Pulse Oximetry
 Monitors..... 2-7
 Sensors 7

R

Reflectance Sensor..... 7
 Respiration Cable and Sensor 9
 Regulators
 Dual Stage..... 40
 E-Tank 40-41
 H-Tank..... 40
 Single Stage 41-42
 Resuscitation Bags 46
 Rodent
 Circuit..... 22
 Mask..... 24
 Roll Stands
 Advisor 14
 Tabletop Monitor 14

S

Sample Lines 11
 SAV 2500 Small Animal Ventilator .. 28
 Scavenger Interface Valve..... 44
 Sensors
 Pulse Oximetry..... 7
 Respiration..... 9
 Silicone Tubing..... 20, 32
 Snuggle Warm® Warming Blankets.. 16
 SodaSorb® CO₂ Absorbant 22-23
 Stand Alone Absorbers 23
 Standard “Y” Circuit..... 21

T

Tabletop Research
 Anesthesia Machine..... 19
 Tabletop/Wall Mount Small
 Animal Anesthesia Machine 17
 Temperature
 Probes 9
 Cables 9
 Thermal Printer 12
 Transducer
 Cable 9
 Invasive Pressure..... 9
 Tubing
 19 mm 43
 22 mm 43
 Silicone..... 20, 32

U

Universal Double
 E-Tank Manifold 40
 Universal F-Circuit..... 21
 Universal “Y” Sensor 7
 University Anesthesia Machine 18

V

Vacuum
 Connectors..... 36
 Hose 34
 Vaporizer
 Accessories..... 32
 100 Series 26, 30
 Classic T3™ 26, 30
 Cleaning & Calibration 30
 Conversion 30
 Conversion Exchange 31
 Exchange..... 31
 Vaporstick Anesthesia Machine..... 19
 Ventilator
 Cart 18
 Small Animal 27-28
 Large Animal 28
 Veterinary Data Logger..... 7

W

Wall Mount Bracket, Advisor®
 vital signs monitor 14
 Wall Outlets
 Recessed..... 33
 Wye, Oxygen 39

Y

“Y” Circuit, Standard..... 21
 “Y” Piece with Luer Fitting..... 25

Z

Z-Bracket 32

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.

St. Paul, MN 55112, USA
Phone: 1-214-618-0218
Toll-Free USA: 1-800-258-5361

www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

Customer and Clinical Services

USA: 1-214-618-0218
Toll-Free USA: 1-800-258-5361

Customer Service

Canada 1-905-477-2000
Toll-Free Canada: 1-800-387-4346

www.surgivet.com
www.smiths-medical.com

Product(s) described may not be licensed or available for sale outside of the United States.

Advisor, Capnocheck, Classic T3, Enviro-Pure, Pneupac, HOTLINE, Level 1, Snuggle Warm, SurgiVet, Equator and the SurgiVet and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2011 Smiths Medical.
All rights reserved. 02.2011 Equipment Catalog VE195804EN

smiths medical