



KAWASUMI
KAWASUMI LABORATORIES AMERICA, INC.

U.S. Product Catalog





K-Shield Advantage® Winged Blood Collection Sets

CE 0123



- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an “audible” click
- Needle is covered during removal, ensuring confidence in safety
- Superior needle sharpness*
- Latex free - for patients, practitioners and end-users



*NAMS needle sharpness report available upon request.


For product codes, box and case sizes please see page 12

K-Shield Advantage®
Winged Blood Collection Sets
with pre-attached Tube Holder

CE 0123



- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an “audible” click
- Needle is covered during removal, ensuring confidence in safety
- Superior needle sharpness*
- Convenience of a pre-attached vacuum tube holder
- Latex free - for patients, practitioners and end-users




*NAMSA needle sharpness report available upon request.

For product codes, box and case sizes please see page 12




Standard Winged Blood Collection Sets with Multiple Sample Luer Adapter CE 0123

- Low cost winged collection set with multiple sample luer adapter
- Self-sealing synthetic needle sheath
- Flexible wings with textured surface for secure grip
- Easy-grip hub fits tube holders with anti-stick devices
- Superior needle sharpness*
- Latex free - for patients, practitioners and end-users



*NAMS needle sharpness report available upon request.

For product codes, box and case sizes please see page 12



K-Shield Advantage® Small Vein Infusion Sets

- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an "audible" click
- Needle is covered during removal, ensuring confidence in safety
- Superior needle sharpness*
- Latex free - for patients, practitioners and end-users

Standard Small Vein Infusion Sets CE 0123

- Flexible wings with textured surface for secure grip which are color coded for easy identification
- Siliconized needle for easy insertion
- Kink resistant tubing
- Superior needle sharpness*
- Latex free - for patients, practitioners and end-users



Intermittent Injection Site

- Low cost injection site
- Male luer connector
- Allows clinician to perform multiple injections when attached to a winged infusion set
- Latex free - for patients, practitioners and end-users



Intermittent injection Site


*NAMSA needle sharpness report available upon request.



K-Shield Advantage® Blood Drawing Kits

- Low-cost therapeutic blood drawing
- Fast and simple one-step activation with an “audible” click
- Needle-free sampling port
- Does not contain unnecessary anticoagulant solution
- 600 mL blood bag
- Pre-attached 16G venipuncture needle
- Pinch clamp for easy prevention of backflow
- 48-inch tubing length
- Volume graduation marks
- Latex free - for patients, practitioners and end-users

BK-64A




For product codes, box and case sizes please see page 13




K-Shield Advantage® Phlebotomy Sets

- Low-cost design incorporates pre-attached needles for accessing all vacuum bottles
- Fast and simple one-step activation with an “audible” click
- Pre-attached 16G venipuncture needle
- 16G needle with backeye optimizes blood flow
- Needle with wing for improved needle control
- Roller clamp to control blood flow
- 36-inch tubing length
- Latex free - for patients, practitioners and end-users



Standard Phlebotomy Sets

- Sets are available with and without pre-attached 16G venipuncture needle
- No safety device
- Latex free - for patients, practitioners and end-users



For product codes, box and case sizes please see page 13



K-Shield Advantage® Port Access Infusion Sets

CE 0123

- **One step activation** – simple, safe and easy to use for the clinicians
- **Tactile and audible click upon safety activation** – assures you that the needle is safely secured
- **Low profile and small body size** – makes dressing simple and device less visible under patient's clothing
- **Small transparent skin plate** – allows a clear view of the infusion site
- **Non-coring siliconized needle** – minimizes drag on the skin during access and de-access
- **Non-DEHP fluid path** – offers a safe, DEHP-free fluid path
- **Latex free** – for patients, practitioners and end-users

**Non-DEHP
Fluid Path**




Before Activation



After Activation

For product codes, box and case sizes please see page 14



Standard Port Access Infusion Sets CE 0123

- Low-cost port access sets available in a variety of models
- Non-DEHP fluid path
- Non-coring needle
- Needle hub and wing design offers secure and accurate needle placement
- Latex free - for patients, practitioners and end-users

For product codes, box and case sizes please see page 14

**Non-DEHP
Fluid Path**




Gravity IV Administration Sets

- **Flexible** - extremely flexible tubing, does not kink
- **Long Length** - long length sets for accommodating a wide variety of surgical procedures where the patient is distant from the IV pole
- **Needle Safe** - Available with NAC® (Needleless Access Connectors), providing a needle-free environment
- **Available** with .22 micron, low protein binding in-line filter
- **Wide selection** of sets and components
- **Non-DEHP fluid path** - all sets feature a Non-DEHP fluid path including tubing and drip chamber
- **Latex Free** - for patients, practitioners and end-users
- **Economical** - high-quality, low-cost gravity IV administration sets for clinical needs as well as meeting materials management budget requirements


Code #	Description (50 sets per case)	Length	Drip Rate	Priming Volume
IV-009N	Standard set with NAC®	84 inches	10 drops/mL	18 mL
IV-009Y	Set with needle-pierceable Y site	84 inches	10 drops/mL	18 mL
IV-013Y	Micro-drip set with needle-pierceable Y site	84 inches	60 drops/mL	18 mL
IV-016Y	Set with two needle-pierceable Y sites	100 inches	10 drops/mL	23 mL
IV-019	Vented secondary set	40 inches	20 drops/mL	10 mL
IV-020Y	Vented standard set with needle-pierceable Y site	84 inches	20 drops/mL	18 mL
IV-025Y	Vented set with .22 micron filter and needle-pierceable Y site	92 inches	20 drops/mL	22 mL
IV-F02Y	Add-on .22 micron filter set and needle-pierceable Y site	21.5 inches	N/A	6.4 mL

For product codes, box and case sizes please see page 15


One Injection Site Sets




IV-009N Standard set with NAC®



IV-009Y Set with needle-pierceable Y-site




IV-013Y Micro-drip set with needle-pierceable Y-site



IV-020Y Vented standard set with needle-pierceable Y-site

Two Injection Site Sets




IV-016Y Set with two needle-pierceable Y-sites

Secondary Set



IV-019 Vented secondary set

.22 Micron Filter Sets



IV-025Y Vented set with .22 micron filter and needle-pierceable Y-site
 € 0123





IV-F02Y Add-on .22 micron filter set and needle-pierceable Y-site
 € 0123


NAC® (Needleless Access Connector) is a registered trademark of CareFusion


For product codes, box and case sizes please see page 15

Ordering Information

Winged Blood Collection Sets (50 per box / 4 boxes per case)							
 <small>Page 2</small>		 <small>Page 4</small>		Needle Gauge	Wing Color	Needle Length	Tubing Length
Code#							
DBMA-21G	DBM1-21G			21g	Green	3/4"	12"
DBMA-23G	DBM1-23G			23g	Light Blue	3/4"	12"
DBMA-25G	DBM1-25G			25g	Dark Blue	3/4"	12"

With Pre-attached Tube Holder <small>Page 3</small>							
(50 per box / 2 boxes per case)							
				Needle Gauge	Wing Color	Needle Length	Tubing Length
				21g	Green	3/4"	12"
				23g	Light Blue	3/4"	12"
				25g	Dark Blue	3/4"	12"

 <small>Page 5</small>						
K-Shield Advantage® Small Vein Infusion Sets (50 per box / 10 boxes per case)						
Code#	Needle Gauge	Wing Color	Needle Length	Tubing Length	Priming Volume	
D3KA-19G	19g	Brown	3/4"	12"	.34mL	
D3KA-21G	21g	Green	3/4"	12"	.34mL	
D3KA-23G	23g	Light Blue	3/4"	12"	.34mL	
D3KA-25G	25g	Dark Blue	3/4"	12"	.34mL	

 <small>Page 5</small>						
Standard Small Vein Infusion Sets (50 per box / 10 boxes per case)						
Code#	Needle Gauge	Wing Color	Needle Length	Tubing Length	Priming Volume	
D3K2-19G	19g	Brown	7/8"	12"	.34mL	
D3K2-21G	21g	Green	3/4"	12"	.34mL	
D3K2-23G	23g	Light Blue	3/4"	12"	.34mL	
D3K2-25G	25g	Dark Blue	3/4"	12"	.34mL	
D3K1-21G	21g	Green	3/4"	3.5"	.11mL	
D3K1-23G	23g	Light Blue	3/4"	3.5"	.11mL	
D3K1-25G	25g	Dark Blue	3/8"	3.5"	.11mL	

Ordering Information

Intermittent Injection Site

Page 5

(50 per box / 4 boxes per case)

Code #	Description	Length
INT-01	Injection site with male luer connector	

K-Shield Advantage® Blood Drawing Kits

Page 6

(25 per box / 2 boxes per case)

Code #	Description	Length
BK-64A	K-Shield Advantage® 600 mL blood bag with 48 inches 48" tubing, pinch clamp and pre-attached 16G venipuncture needle with safety device, and needle-free sampling port	48 inches

Standard Phlebotomy Sets

Page 7

(25 per box / 2 boxes per case)

Code #	Description	Length
BO-1000	Phlebotomy set with 16G needle vacuum bottle	36 inches
BO-16S1	Phlebotomy set with 16G vacuum bottle needle and pre-attached 16G venipuncture needle with backeye	36 inches



K-Shield Advantage® Phlebotomy Sets

Page 7

(25 per box / 2 boxes per case)

Code #	Description	Length
BA-16S1	K-Shield Advantage® Phlebotomy set with 16G vacuum bottle needle and pre-attached 16G venipuncture needle with backeye and safety device	36 inches

Ordering Information

Port Access Infusion Sets			Page 8		Page 9	
			 (20 per box / 5 boxes per case)		 (25 per box / 4 boxes per case)	
Needle Gauge	Needle Length	Pinch Clamp Color	Microbore Tubing		Microbore Tubing	
			Without NAC®	With NAC®	Without Y-SITE	With Y-SITE
19	3/4"	Brown	PA03090	PA01N90	PI03090	PI01Y90
19	1"	Brown	PA03091	PA01N91	PI03091	PI01Y91
20	3/4"	Yellow	PA03000	PA01N00	PI03000	PI01Y00
20	1"	Yellow	PA03001	PA01N01	PI03001	PI01Y01
22	3/4"	Black	PA03020	PA01N20	PI03020	PI01Y20
22	1"	Black	PA03021	PA01N21	PI03021	PI01Y21
20	1 1/2"	Yellow	-	-	-	PI01Y03
			Standard Tubing		Standard Tubing	
19	1"	Brown	PA04091	-	PI04091	-
20	3/4"	Yellow	PA04000	-	PI04000	-
20	1"	Yellow	PA04001	-	PI04001	-

NAC® (Needleless Access Connector) is a registered trademark of CareFusion

Ordering Information



Gravity IV Administration Sets Page 10, 11				
(50 sets per case)				
Code #	Description	Length	Drip Rate	Priming Volume
IV-009N	Standard set with NAC®	84 inches	10 drops/mL	18 mL
IV-009Y	Set with needle-pierceable Y-site	84 inches	10 drops/mL	18 mL
IV-013Y	Micro-drip set with needle-pierceable Y-site	84 inches	60 drops/mL	18mL
IV-016Y	Set with two needle-pierceable Y-sites	100 inches	10 drops/mL	23 mL
IV-019	Vented secondary set	40 inches	20 drops/mL	10 mL
IV-020Y	Vented standard set with needle-pierceable Y-site	84 inches	20 drops/mL	18 mL
IV-025Y	Vented set with .22 micron filter and needle-pierceable Y-site	92 inches	20 drops/mL	22 mL
IV-F02Y	Add-on .22 micron filter set and needle-pierceable Y site	21.5 inches	N/A	6.4 mL

NAC® (Needleless Access Connector) is a registered trademark of CareFusion



KAWASUMI

KAWASUMI LABORATORIES AMERICA, INC.

www.kawasumiamerica.com

1-800-KAWASUMI (529-2786)

P.O. Box 24355 • Tampa, FL 33623-4355

Phone: 813-630-5554 • Fax: 813-630-5033

e-mail: orders@kawasumiamerica.com



KAWASUMI LABORATORIES, INC.

Shinagawa Intercity Tower B, 2-15-2,

Konan, Minato-ku, Tokyo 108-6109, Japan

Kawasumi continuously improves and develops new products. Visit our website for more details!

PC-03