



ESSENTIAL. EFFECTIVE. EFFICIENT.

[Home](#) [Products](#) [Wound Management](#) [Wound Drains](#) [Gravity Drains](#)

Prescriptive Information

[read more](#)

Literature

[Wound Drainage Brochure](#)

[Closed Wound Drainage Brochure](#)

Training and Education

For additional training and education on a number of topics, visit our Training Center

[read more](#)

CONTACT US

[read more](#)

Gravity Drains



Product Images:

2 of 2

Gravity Drains include:

- Deaver T-Tubes and Cattell T-Tubes which are used following open gallbladder procedure to drain bile. Often attached to BARD® bile bag.
- Whelan-Moss T-Tubes - Stem ID allows stone basket to retrieve stones missed during initial procedure.
- Penrose Drains - Used to promote drainage in an open surgical wound.

Features and Specifications

Description	Part Number	Size	Quantity	Natural Latex?
Drains, Penrose	0912010	Length 12", Width 1/4"	200/case	Yes
Drains, Penrose	0912020	Length 12", Width 3/8"	200/case	Yes
Drains, Penrose	0912030	Length 12", Width 1/2"	200/case	Yes
Drains, Penrose	0912040	Length 12", Width 5/8"	200/case	Yes
Drains, Penrose	0912050	Length 12", Width 3/4"	200/case	Yes
Drains, Penrose	0912060	Length 12", Width 7/8"	200/case	Yes
Drains, Penrose	0912070	Length 12", Width 1"	200/case	Yes
Drains, Penrose	0918010	Length 18", Width 1/4"	200/case	Yes
Drains, Penrose	0918020	Length 18", Width 3/8"	200/case	Yes
Drains, Penrose	0918030	Length 18", Width 1/2"	200/case	Yes
Drains, Penrose	0918040	Length 18", Width 5/8"	200/case	Yes
Drains, Penrose	0918050	Length 18", Width 3/4"	200/case	Yes
Drains, Penrose	0918060	Length 18", Width 7/8"	200/case	Yes
Drains, Penrose	0918070	Length 18", Width 1"	200/case	Yes
Drains, Penrose	0097870	Length 36", Width 1/4"	144/case	Yes
Drains, Penrose	0097880	Length 36", Width 5/16"	144/case	Yes
Drains, Penrose	0097890	Length 36", Width 3/8"	144/case	Yes
Drains, Penrose	0097900	Length 36", Width 1/2"	144/case	Yes
Drains, Penrose	0097910	Length 36", Width 5/8"	144/case	Yes
Drains, Penrose	0097920	Length 36", Width 3/4"	144/case	Yes
Drains, Penrose	0097930	Length 36", Width 7/8"	144/case	Yes
Drains, Penrose	0097940	Length 36", Width 1"	144/case	Yes
Drains, Cattell T-Tube	0100140	12" Crossbar, 11" Stem, 8FR	10/case	Yes
Drains, Cattell T-Tube	0100150	12" Crossbar, 11" Stem, 10FR	10/case	Yes
Drains, Cattell T-Tube	0100160	12" Crossbar, 11" Stem, 12FR	10/case	Yes
Drains, Cattell T-Tube	0100170	12" Crossbar, 11" Stem, 14FR	10/case	Yes
Drains, Cattell T-Tube	0100180	12" Crossbar, 11" Stem, 16FR	10/case	Yes
Drains, Cattell T-Tube	0100190	12" Crossbar, 11" Stem, 18FR	10/case	Yes
Drains, Cattell T-Tube	0100200	12" Crossbar, 11" Stem, 20FR	10/case	Yes
Drains, Cattell T-Tube	0100210	12" Crossbar, 11" Stem, 22FR	10/case	Yes
Drains, Cattell T-Tube	0100220	12" Crossbar, 11" Stem, 24FR	10/case	Yes
Drains, Deaver T-Tube	0098770	5" Crossbar, 12" Stem, 5FR	10/case	Yes

Drains, Deaver T-Tube	0098790	5" Crossbar, 12" Stem, 8FR	10/case	Yes
Drains, Deaver T-Tube	0099800	5" Crossbar, 12" Stem, 10FR	10/case	Yes
Drains, Deaver T-Tube	0099810	5" Crossbar, 12" Stem, 12FR	10/case	Yes
Drains, Deaver T-Tube	0099820	5" Crossbar, 12" Stem, 14FR	10/case	Yes
Drains, Deaver T-Tube	0099830	5" Crossbar, 12" Stem, 16FR	10/case	Yes
Drains, Deaver T-Tube	0099840	5" Crossbar, 12" Stem, 18FR	10/case	Yes
Drains, Deaver T-Tube	0099850	5" Crossbar, 12" Stem, 20FR	10/case	Yes
Drains, Deaver T-Tube	0099860	5" Crossbar, 12" Stem, 22FR	10/case	Yes
Drains, Deaver T-Tube	0099870	5" Crossbar, 12" Stem, 24FR	10/case	Yes
Drains, Hysterectomy T-Tube	0070450	18FR Crossbar	10/case	Yes
Drains, Whelan-Moss T-Tube	0103220	12FR Crossbar	10/case	Yes
Drains, Whelan-Moss T-Tube	0103230	14FR Crossbar	10/case	Yes
Drains, Whelan-Moss T-Tube	0103240	16FR Crossbar	10/case	Yes
Drains, Whelan-Moss T-Tube	0103250	18FR Crossbar	10/case	Yes

Specifications subject to change without notice.



Bard Medical is excited to announce the acquisition of Rochester Medical

**Bard Medical and Rochester Medical are coming together,
combining quality, innovative products and services
to better serve Clinician and Patient needs.**



For information about Bard Medical Home Care products. [Click here](#)



For information about MAGIC 3® Hydrophilic Catheters. [Click here](#)

For Rochester Medical product orders and all other inquiries, please continue to contact Rochester Medical at 800.243.3315 or visit www.rocm.com

Copyright © 2013 C. R. Bard, Inc. All rights reserved. All products on this website may not be available for sale in all geographies. Please contact your local Bard Medical representative for information about products available in your area.

This site is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider.