

Hinged Wraptor

Similar to the Wraptor Ankle Stabilizer with rigid thermoplastic hinge for additional support and stabilization.



Common Examples of Use

- Chronic and acute ankle sprains
- Abnormal eversion and inversion control
- Protective sportswear

Features

- Detachable rigid hinge
- Breathable Spacermesh lining wicks moisture away from the skin
- Sharkskin material at the base of the foot reduces brace slippage
- Latex-free Neoprene at the heel reduces irritation on the Achilles tendon

Coming Soon

Wraptor Ankle Stabilizer

Figure 8 strapping design mimics traditional taping without gapping. High ankle strap provides additional support and secures medial and lateral stirrup straps.



Common Examples of Use

- Chronic and acute ankle sprains
- Abnormal eversion and inversion control
- Protective sportswear

Features

- Breathable Spacermesh lining wicks moisture away from the skin
- Sharkskin material at the base of the foot reduces brace slippage
- Latex-free Neoprene at the heel reduces irritation on the Achilles tendon

Coming Soon!

Ankle Brace Lock (ABL)

The Ankle Brace Lock (ABL) is ideal for providing support for patients with mild to moderate inversion and eversion injuries and chronic ankle instabilities. The body of the ABL is made of a breathable, stretch nylon / lycra weave, allowing it to fit in a compact-style shoe. Two no-stretch stirrup straps provide extra support.



Common Examples of Use

- Mild to moderate inversion and eversion injuries
- Chronic ankle instabilities

Features

- Low profile
- Breathable, stretch nylon / lycra material
- Two no-stretch stirrup straps
- Universal right or left
- XS – XXL sizes

Part # Description

| | |
|----------|--|
| SA7023XX | Hinged Wraptor, Regular Laces w/ Hinge, XS-XXL |
| SA7024XX | Hinged Wraptor, Speed Laces w/ Hinge, XS-XXL |

XX= See size chart above

Wraptor Ankle / Hinged Wraptor

| Size # | Size | Men's shoe size | Women's shoe size |
|---------|------|-----------------|-------------------|
| XX = 01 | XS | 6 – 7 | 7 – 8 |
| XX = 03 | S | 8 – 9 | 9 – 10 |
| XX = 05 | M | 9 - 11 | 10 – 12 |
| XX = 07 | L | 11 – 13 | 12 – 14 |
| XX = 09 | XL | 14+ | 15+ |
| XX = 11 | XXL | 14 W+ | 15W+ |

For wider feet, move up a size.
For narrower feet, move down a size.

Part # Description

| | |
|-------------|---|
| SA7025XX--M | Wraptor Ankle Stabilizer w/ Standard Laces, Black, XS-XXL |
| SA7020XX | Wraptor Ankle Stabilizer w/ Speed Laces, Black, XS-XXL |
| SA7020XX--M | Wraptor Ankle Stabilizer w/ Speed Laces, White, XS-XXL |

XX= See size chart above

Ankle Brace Lock Sizing Chart

| Size # | Size | Men | Women |
|--------|------|-----------|-----------|
| X = 1 | XS | < 5 | < 6 |
| X = 2 | S | 5.5 – 7 | 6.5 – 8 |
| X = 3 | M | 7.5 – 9 | 8.5 – 10 |
| X = 4 | L | 9.5 – 11 | 10.5 – 12 |
| X = 5 | XL | 11.5 – 13 | 12.5 + |
| X = 6 | XXL | 13.5+ | |

Measurements according to shoe size

Part # Description

| | |
|-------|------------------------|
| 1017X | Ankle Brace Lock (ABL) |
|-------|------------------------|

X = See size chart
CE Marking