



Now Available...

4 Fr. Single-Lumen PowerPICC\* Catheter

A Smaller PowerPICC\* Catheter Option



**BAIRD**  
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.\**

**The Power of Purple\***



Experience the efficiency of power injection for contrast enhanced CT scans and PICC therapy with the full line of PowerPICC\* catheters. The wide range of PowerPICC\* catheter options designed to meet diverse patient and clinician needs. The PowerPICC\* catheter family now includes a 4 Fr. single-lumen catheter.

- **Maximum injection rates –**  
allow injection of contrast media for CECT scans at a **maximum rate of 5 ml/sec**
- **Easy identification –**  
catheters are clearly labeled to differentiate them from all non-indicated power injection PICCs
- **Reverse taper design –**  
gently plugs the insertion site and provides superior kink resistance
- **CVP monitoring –**  
all PowerPICC\* catheters are indicated for CVP monitoring
- **Patient satisfaction –**  
improves patient comfort due to decreased needle sticks
- **Multiple kit configurations –**  
available in both microintroducer and safety needle introducer kit configurations. Also available with Sherlock\* stylet.

**Multiple Configurations:**

Single-Lumen	Dual-Lumen	Triple-Lumen
<b>4 Fr.</b>	5 Fr.	6 Fr.
5 Fr.	6 Fr.	

**Ordering Information**

**Nursing:**

Product Code	Tray	Gravity Flow (ml/hr)	Priming Volume (ml)	Gauge
3174108	Full w/MI	1,272	0.67	18
3174118	Basic w/MI	1,272	0.67	18

**Nursing with Sherlock\* Stylet:**

Product Code	Tray	Gravity Flow (ml/hr)	Priming Volume (ml)	Gauge
9174108	Full w/MI	1,272	0.67	18
9174118	Basic w/MI	1,272	0.67	18

**Interventional Radiology:**

Product Code	Guidewire	Gravity Flow (ml/hr)	Priming Volume (ml)	Gauge
3174135	135 cm – stainless steel wire	1,272	0.67	18
3174155	70 cm – stainless steel wire	1,272	0.67	18
3174335	135 cm – nitinol wire	1,272	0.67	18
3174355	70 cm – nitinol wire	1,272	0.67	18



Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

**Learn More** - For more information on the PowerPICC\* Catheter, please contact your Bard Access Systems representative or visit [www.powerpicc.com](http://www.powerpicc.com).  
Call toll-free at 1-800-545-0890, or visit us on the web at [www.bardaccess.com](http://www.bardaccess.com).



ACCESS SYSTEMS

**Bard Access Systems, Inc.**

Salt Lake City, Utah 84116 USA  
801-595-0700  
[www.bardaccess.com](http://www.bardaccess.com)

Customer Service: 800-545-0890  
Clinical Information: 800-443-3385