



pro**glide**4011
ultra procedure chair



The Next Generation Procedure Chair just got better

The **ONLY** Power Procedure Chair with Wheels and ADA Access
Helps you protect patients from infection
Helps give you maximum room flexibility
Helps your staff to be more efficient



CHAIR SPECIFICATIONS

- Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat Max. 37" at seat (93.9cm)
Seat Width: 28" (71.12cm)
Back Width: 26.5" (67.31cm)
Length: 73"-80" with foot section extended (185.42cm – 203.2cm)
Leg Rest: 28" W x 11.5" D (71.12cm W x 29.21cm D)
Stirrups: Four Lateral Positions, Adjustable Length
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60HZ
International Model 4011-650-301: 220/240V, 50/60Hz

ACCESSORIES

- 044 Adjustable/self leveling arms pair - field installable
 050 Removable U-shaped headrest
 051 Removable rectangular headrest
 251 Bierhoff knee crutches
 467 Mounting plate for Welch Allyn Light Bracket #48955 - patient right
 468 Mounting plate for Welch Allyn Light Bracket #48955 - patient left
 830 Patient safety straps
 900 Lifting bracket kit
 920 Easy installation kit
 Fire retardant foam (CAL133) available

KEY FEATURES

model # 4011-650-300

- Chair Designed for Universal Use
 OneTouch Wheelbase® Integrated Central-Locking Caster System
 650 lb. (294.84kg) Patient Weight Capacity
 Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
 Hand Control with OneTouch Patient Positioning®
 2 Programmable Positions & Return-to-Home
 Foot Control with OneTouch Patient Positioning®
 2 Programmable Positions & Return-to-Home
 Removable Seamless Upholstered Cushion
 Double-Articulating Headrest with Adjustable Height
 Interchangeable Headrest Options
 Stirrups
 Stainless Steel Drain Pan
 Trendelenburg, Lithotomy, Semi-Fowler Positions
 Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
 Concealed Paper Roll Holder
 3 Year Limited Warranty

proglide4011

ultra procedure chair



COLORS

Shown with optional Armrests (044)

- River Rock
- Sand Grey
- Mocha
- Adobe
- Steel Blue
- Midnight Blue
- Sage
- Onyx

Why wouldn't you want a procedure chair that moves?

Features that matter - Empowering the Patient Experience

Infection Control

One Touch WheelBase® for cleaning surfaces around and under the Procedure Chair. Seamless upholstery and body surfaces make it fast and easy to wipe clean. No cracks or crevices for mold or bacteria to grow. The upholstered top is removable for deep cleaning and replacement to extend the life of the procedure chair.

OneTouch Patient Positioning®

Hand Control and Foot Control with 2 customizable positions plus a Return-to-Home function designed for uninterrupted Patient-Practitioner Communication. Program each position in less than 2 seconds.

Patient-Centric Design

Accessibility. 650 lb. pt. wt., 18.5" ADA height. Your patient population is changing. You can accommodate change seamlessly with UMF Medical's Power Procedure Chairs. Every power chair has 650 lb. patient weight and 18.5" floor height for wheel-chair transfer designed to help you provide exceptional care to every patient.

Ultimate adjustability, including
Power Seat Tilts up to a 30° angle
Lithotomy Position
Semi-Fowler Position

Hand Control



Articulating, Adjustable Headrest



Foot Control



OneTouch Locks and Unlocks the base.. Glide your procedure chair around with ease.

The OneTouch WheelBase® is a game changer!

Integrated central-locking wheelbase enables staff to easily move the chair **for cleaning** and lock it solidly into place for exams and procedures.



Shown with optional Armrests (044)

“Vulnerable patient populations rely on frequent and intensive use of outpatient care to maintain or improve their health... It is critical that all of this care be provided under conditions that minimize or eliminate risks of healthcare-associated infections (HAI). All healthcare settings, regardless of the level of care provided, must make infection prevention a priority.”

Guide to Infection Prevention For Outpatient Settings: Minimum Expectations for Safe Care, CDC, National Center for Emerging and Zoonotic Diseases. Nov. 2015



The **proglide4011 Procedure Chair** is designed to optimize the patient and caregiver experience. Because your patient population is changing, every detail matters. That's why UMF Medical is setting a **new benchmark** for **infection control, smart efficiencies, optimal accessibility, and lasting value**. By asking what's going to improve your patients' experience, UMF Medical has arrived at a procedure room solution that is the perfect blend of style and function.

THIS IS THE PROCEDURE CHAIR YOU'VE BEEN WAITING FOR!

UMFmedical

1316 Eisenhower Boulevard
Johnstown, PA 15904-3393
Phone: 814-266-8726
Toll Free: 800-638-5322
Fax: 814-266-1870
www.UMFmedical.com

