



Smart.
Connected.
Empowering.
Alere.

Product Catalogue



Welcome to Connected Health.

We believe that better health starts with better information. That's why Alere brings together advanced diagnostic tools and integrated health management solutions that deliver timely, accurate information and drive better outcomes. It's an integrated approach we call connected health — and it's built upon the technology and vision of more than 40 leading diagnostic and health management companies with expertise in over 100 disease categories. It means connecting the widest range of products and services in the widest range of environments — from hospital to home. It means connecting everyone — providers, physicians and the people they treat — to the same reliable, up-to-date data. And it means getting the diagnosis right, selecting the least invasive interventions and allowing health management to continue away from traditional points of care. So whether it's in an emergency room in China, a doctor's office in Spain or a home in New York City, our solutions make it possible to connect with the information needed to make the best health decisions. Because at Alere, we understand that connected health is the smart path to better health. Everywhere.



Smart.

Alere products and solutions are changing how healthcare is delivered—extending diagnosis and health management from the lab, doctor’s office and hospital into people’s homes. The result: improved outcomes and management models that are both more effective and more efficient.

Connected.

Alere actively integrates diagnostic tools with education and health management solutions, making it easier than ever to personalize care. We’re connecting patients to information. Information to diagnostic results. And diagnostic results to healthcare providers, so they have the information needed to make the best decisions.

Empowering.

Alere products and solutions give providers, physicians and patients access to the right information at the right time. That means doctors don’t have to wait for results to provide seamless, continuous care. And individuals have the tools and knowledge to play a more active role in their own health management.



Smart. Connected. Empowered.

Ordering

Alere Inc.
Shared Service Centre
Tel: +1 321 441 7200
Fax: +1 321 441 7400
Email: clientservices@alere.com

Alere International Limited
Parkmore East Business Park
Ballybrit, Galway
Ireland
Tel: + 353 (0) 91 429 900

Please note that individual products may have different specifications in various markets. Also, some products may not be available in all markets worldwide, only in selected markets.

Product Information and General Enquiries

Further information can be obtained from your distributor, or contact
Email: professionaldiagnosics@alere.com
Website: www.alere.com



Ordering and Product Enquiries

Contents

1) Alere	6
2) Alere Actim®	7
3) Alere BinaxNOW®	8
4) Alere NMP22® Bladderchek®	11
5) Alere Cholestech LDX® System	12
6) Alere Determine™	14
7) Alere DxLINK™	16
8) Alere™ Heart Check System	17
9) Alere INRatio® and INRatio®2	18
10) Alere Pima™ CD4	20
11) Alere™ TestPack	21
12) Alere Triage®	22
13) Afinion™ and NycoCard™	28
14) BinaxNOW®	31
15) Clearview®	32
16) DoubleCheckGold™	36
17) epoc®	36
18) ImmunoComb®	37
19) Osteomark®	40
20) Panbio®	41
21) SureStep™	43
22) Techlab®	51
23) Disease/Condition Index	58

MRSA

Alere™ PBP2a

For the rapid qualitative detection of PBP2a in isolates of *Staphylococcus*, as an aid in identifying methicillin-resistant *Staphylococcus aureus* (MRSA).

Product	Product Code	Size
Alere™ PBP2a Culture Colony Test	891-000	25



Women's Health

Alere™ HCG

A qualitative rapid immunoassay for the detection of human Chorionic Gonadotropin (hCG) in urine.

Product	Product Code	Size
Alere™ hCG Cassette - UK ONLY	4531319999	25mIU/mL
	4531309999	

A qualitative rapid immunoassay for the detection of human Chorionic Gonadotropin (hCG) in urine and Serum.

Product	Product Code	Size
Alere™ hCG Combo Cassette	4531299999	25mIU/mL



Preterm Delivery

Alere Actim® Partus

Alere Actim® Partus is a fast and simple bedside test to estimate the risk of preterm delivery. It is an immunochromatographic dipstick test that detects the presence of the phosphorylated form of IGFBP-1 (insulin-like growth factor binding protein-1) in cervical secretions. The level of phosphorylated IGFBP-1 in the cervix increases considerably as the cervix matures.

Product	Product Code	Size
Alere Actim® Partus	31930ETAL	1
Alere Actim® Partus	31931ETAL	10



Alere Actim® PROM

Alere Actim® PROM is a fast immunochromatographic dipstick test that reacts to the presence of amniotic fluid in vaginal secretions. Based on monoclonal antibodies, the Alere Actim® PROM test detects IGFBP-1 (insulin-like growth factor binding protein-1) in a vaginal sample.

Product	Product Code	Size
Alere Actim® PROM	30830ETAL	1
Alere Actim® PROM	30831ETAL	10
Alere Actim® PROM	30832ETAL	20



Infectious Diseases - Respiratory

Influenza A & B

For the qualitative detection and differentiation of influenza A and influenza B viral antigen from nasopharyngeal swab, nasal wash and nasal aspirate specimens.

Product	Product Code	Size
Alere BinaxNOW® Influenza A & B	416-000	22



Legionella

For the qualitative detection of *Legionella pneumophila* serogroup 1 antigen in urine specimens.

Product	Product Code	Size
Alere BinaxNOW® Legionella	852-012	12
	852-000	22
Alere BinaxNOW® Legionella Control Swab Pack	852-010	5 positive & 5 negative
Alere BinaxNOW® Legionella Urinary Antigen EIA	851-000	96 well test kit



Respiratory Syncytial Virus

For the qualitative detection of respiratory syncytial virus (RSV) fusion protein antigen from nasopharyngeal swab and nasal wash specimens in children less than 5 years.

Product	Product Code	Size
Alere BinaxNOW® RSV	430-022	22
	430-000	42
Alere BinaxNOW® Nasopharyngeal Swab Accessory Pack	400-065	20



Streptococcus pneumoniae

For the rapid qualitative detection of *S. pneumoniae* antigen in urine of patients with pneumonia and cerebral spinal fluid (CSF) of patients with meningitis.

Product	Product Code	Size
Alere BinaxNOW® <i>S. pneumoniae</i> Control Swab Pack	710-010	5 positive & 5 negative
Alere BinaxNOW® <i>S. pneumoniae</i> Antigen Card	710-012	12
	710-000	22



Infectious Diseases - Parasitology

Malaria

For the qualitative differential detection of circulating *Plasmodium falciparum* (*P.f.*) antigen and an antigen that is common to all four species of malaria: *Plasmodium falciparum* (*P.f.*), *Plasmodium vivax* (*P.v.*), *Plasmodium ovale* (*P.o.*), and *Plasmodium malariae* (*P.m.*) in whole blood specimens.

Product	Product Code	Size
Alere BinaxNOW® Malaria	66005	5
	660-000	25

Malaria Positive Control

The Alere BinaxNOW® Malaria Positive Control is intended for use as an assayed positive external quality control with the qualitative Alere BinaxNOW® Malaria Test. It is designed for routine use to aid in verifying proper performance of the procedure and the antigen detection reagents of the Alere BinaxNOW® Malaria Test.

Product	Product Code	Size
Alere BinaxNOW® Malaria Positive Control	665-010	1

Infectious Diseases - Healthcare Associated Infections

Staphylococcus aureus

For the presumptive identification of *Staphylococcus aureus*; performed directly from blood culture samples positive for Gram-Positive Cocci in Clusters (GPCC).

Product	Product Code	Size
Alere BinaxNOW® <i>S. aureus</i>	800-000	22

Alere NMP22[®] Bladderchek[®]

Bladder Cancer

The Alere NMP22[®] Bladderchek[®] Test is an in vitro immunoassay intended for the qualitative or quantitative detection of NMP22, a nuclear mitotic apparatus protein (NuMA) in urine of persons with risk factors or symptoms of bladder cancer or with a history of bladder cancer.

The Alere NMP22[®] Bladderchek[®] Test is a painless and non-invasive point-of-care test performed on a single urine sample and is the only in-office test approved by the FDA for the diagnosis of bladder cancer.

Both test formats – POC and lab based ELISA have received FDA and CE approval.



Product	Product Code	Size
Alere NMP22 [®] Bladderchek [®] Test	D1200E	24 Tests
Alere NMP22 [®] Bladderchek [®] Test	D1300E	10 Tests
Alere NMP22 [®] Test Kit EIA	D1100E	96 Well Plate
Urine Stabiliser (only needed for NMP22-EIA)	F4026	100 Bottles

Alere Cholestech LDX[®] System

Lipid Testing

Alere Cholestech LDX[®] System

The Alere Cholestech LDX[®] System delivers the ability to measure a complete lipid profile, with or without glucose, and it does it all in 5 minutes per test cassette.



Product	Quantity	Product Code
Alere Cholestech LDX [®] Analyzer (without Printer) Includes: LDX Analyzer and Universal power plug pack supplied only	1	14-875
LDX System, Int'l		12-282
Includes: Int'l LDX Analyzer, EU transformer, Euro thermal printer, EU Cholestech Transformer, user manual, training video, printer cable & OCC (OUS Kit)		
Alere Universal Printer		55115
Printer Power Supply		14-857 (UK); 14-858 (AUS); 14-873 (North America); 14-856 (EU)
Alere Cholestech LDX [®] Rom Upgrade Kit	1	13-184
MiniPet [®] Pipette (Lipids/Glucose venous whole blood, Lipids/Glucose/controls)	40 μ L	13-014
Pipette Tips	50	11-010
Optics Check Cassette	1	10-228
40 μ L Capillary Tubes (Lipids/Glucose)	50	52193
Capillary Plungers (fits all capillary tubes)	50	10-311
Printer Cable Kit	1	11-787
White Printer Labels	roll of 400	11-785
Carrying Case	1	10-973

Product	Quantity	Product Code
Other Accessories:	1	
RS232 Adapter Kit	1	10-958
Accessory Tray	1	01-401
LDX Power Supply (International Plug Kit)	1	14-698

Alere Cholestech LDX® Lipid & Glucose Tests

Product	Quantity	Product Code
Alere Cholestech LDX® Lipid Profile•GLU	1 box of 10 cassettes	10-991
Alere Cholestech LDX® Lipid Profile	1 box of 10 cassettes	10-989
Alere Cholestech LDX® TC•HDL•GLU	1 box of 10 cassettes	10-990
Alere Cholestech LDX® TC•HDL	1 box of 10 cassettes	10-987
Alere Cholestech LDX® TC•GLU	1 box of 10 cassettes	10-988
Alere Cholestech LDX® TC	1 box of 10 cassettes	10-986
Alere Cholestech LDX® Multianalyte Controls Level 1 & 2 (2mL per vial)	2 ea	88769
Alere Cholestech LDX® Multianalyte Controls Level 1 & 2 (0.25mL per vial)	1 ea	88773
Alere Cholestech LDX® Calibration Verification	1 ea	88770
Alere Cholestech LDX® Lipid Profile (for export only)	1 box of 10 cassettes	99020
Alere Cholestech LDX® Lipid Profile•GLU (for export only)	1 box of 10 cassettes	99021


Infectious Diseases

Hepatitis B

For the qualitative detection of Hepatitis B Surface Antigen (HBsAg) in human serum, plasma and whole blood specimens.

Product	Product Code	Size
Alere Determine™ HBsAg		
Serum & plasma	7D2543	100
Whole blood - chase buffer required	7D2543 + 7D2243	100

Accessories	Product Code
Chase buffer	7D2243
EDTA Capillary Tubes	7D2222



Infectious Diseases - Sexually Transmitted Infections

HIV-1/2

For the qualitative detection of antibodies to HIV-1 and HIV-2 in human serum, plasma and whole blood specimens.

Product	Product Code	Size
Alere Determine™ HIV 1/2		
Serum & plasma	7D2342	20
Serum & plasma	7D2343	100
Whole blood - chase buffer required	7D2352 + 7D2243	20
Whole blood - chase buffer required	7D2353 + 7D2243	100
Serum & Plasma (CE marked)	7D2346	20
Serum & Plasma (CE marked)	7D2347	100
Whole blood - chase buffer required (CE Marked)	7D2346 + 7D2243	20
Whole blood - chase buffer required (CE Marked)	7D2347 + 7D2243	100



Product**Alere Determine™ HIV-1/2 Ag/Ab Combo**

Serum & plasma
 Whole Blood - chase buffer required
 Serum and Plasma (CE Marked)
 Whole Blood - chase buffer required (CE Marked)
 Whole Blood - chase buffer required (CE Marked)

Product Code

7D2643
 7D2643 + 7D2243
 7D2646
 7D2646 + 7D2243
 7D2647 + 7D2243

Size

100
 100
 20
 20
 100

**Accessories**

Chase buffer (CE and Non-CE Marked)
 EDTA Capillary Tubes
 EDTA Capillary Tubes (CE Marked)

Product Code

7D2243
 7D2222
 7D2227

100

Syphilis

For the qualitative detection of antibodies to Treponema pallidum in human serum, plasma and whole blood specimens.

Product**Alere Determine™ Syphilis TP**

Serum & plasma
 Whole blood - chase buffer required

Product Code

7D2443
 7D2453 + 7D2243

Size

100
 100

**Accessories**

Chase buffer (CE + Non-CE Marked)
 EDTA Capillary Tubes
 EDTA Capillary Tubes (CE Marked)

Product Code

7D2243
 7D2222
 7D2227

TB LAM

For the qualitative detection of lipoarabinomannan (LAM) antigen of Mycobacteria in human urine as an aid in the diagnosis of active mycobacterial infection in HIV positive individuals with clinical symptoms of tuberculosis.

Product**Alere Determine™ TB LAM Ag**

Urine

Product Code

7D2740

Size

100



Alere DxLINK™

System Includes Hardware: Scanner, Scanner Stand, Power Cord, AC Adapter and USB Cord Software: Installation CD-ROM – windows based computer with operating systems – Windows XP, Vista or Windows 7 (not provided)
Other: Operator Guide-electronic.



Product

Instrumented Test System with Alere DxLink™ Technology

Product Code

1905190701

Drugs of Abuse

25 Test devices, 25 Droppers, 1 Package insert (available in English, French, German, Spanish and Italian)

Product

Instrumented Test System with Alere DxLink™ Technology

Product Code

4531379012 /
1156050901

Size

1 Kit

Alere™ Heart Check and Alere Influenza A & B

Alere™ Heart Check System

Alere™ BNP Test

The first rapid Point-of-Care test for the quantitative measurement of BNP using fresh capillary whole blood from a fingerstick.

Product	Product Code	Quantity
BNP Test Kit	CAT - 0001	12 Strips
Meters - German	CAT - 0025	
- English	CAT - 0026	
- French	CAT - 0027	
- Spanish	CAT - 0028	
- Italian	CAT - 0029	
- Swedish	CAT - 0030	
- India	CAT - 0031	
- Australia/NZ	CAT - 0032	
Controls	98071EU	



Alere™ Influenza A & B

The Alere™ Influenza A & B Test is an in-vitro immunochromatographic assay for the qualitative detection of influenza A and B nucleoprotein antigens in nasal swab specimens collected from symptomatic patients. It is intended to aid in the rapid differential diagnosis of influenza A and B viral infections.

Product Code	Quantity
Please contact for SKU information	25 tests

Alere INRatio® and INRatio®2

Coagulation Monitoring

Alere INRatio®2 PT/INR Monitoring System

Measures Prothrombin Time (PT)/INR, the critical test performed to maintain the proper “therapeutic range” of oral anticoagulation therapies.



Product	Size	Product Code
Alere INRatio®2 PT/INR Professional Monitoring System	1	0200431
Alere INRatio®2 PT/INR Home Monitoring Kit	1	0200433
Alere INRatio® PT/INR Test Strips	12	0100071
Alere INRatio® PT/INR Monitoring System Test Strips	48	0100139

Alere INRatio®/INRatio®2 Components

Product	Size	Product Code
Alcohol Pads	100	0200086
Data Cable	1	0100263
Gauze Pads	200	0200085
Alere INRatio® Blue Tote Case	1	0200016
Alere INRatio® Instructional DVD, Pro	1	0200264
Alere INRatio® Instructional DVD, PST	1	0200279
Alere INRatio® Printer	1	14-716INR
Alere INRatio® Printer Roll of Labels	400/roll	11-785
Alere INRatio® Pro User Guide	1	0200123
Alere INRatio® Self Test, DVD, Instructional	1	0200279
Alere INRatio® Self Test User Guide, English	1	0200116
Alere INRatio® Quick Reference Guide, English	1	0200047
Alere INRatio® Video, Self-Test, VHS, Instructional	1	0200278
Alere INRatio®2 Carrying Case	1	25275
Alere INRatio®2 Pro User Guide	1	26047EN
Alere INRatio®2 Self Test User Guide	1	26049EN
Alere INRatio®2 Quick Reference Guide, English	1	26052

Alere INRatio®/INRatio®2 Components

Product	Size	Product Code
Alere INRatio®2 Training DVD, Multilingual	1	26050
Lancing Device, Autolet® Impression	1	0200497
Lancet, Unilet® 23G Lancet Needles	12 Pack	18101
Lancet, Unilet® 23G Lancet Needles	100 ea	18102
Lancets, Unistik® 21 Gauge	50 Pack	902009
Log, Patient Result, Generic	1	26057
MicroSafe® Capillary Tubes	50/bag	200235
MicroSafe® Capillary Tubes	10/bag	200455
Power Supply, Boxed 120V USA	1	0100011
Power Supply Plug Set, Worldwide	1	0100260

Alere Pima™ CD4

CD4 Testing

Consisting of the Alere Pima™ Analyser and the Alere Pima™ CD4 cartridge, the Alere Pima™ CD4 test is a unique and breakthrough technology in the enumeration of T-helper cells in whole blood samples. Designed to perform in a laboratory as well as in non-laboratory environments, the Alere Pima™ CD4 test can be performed on external power or using an on-board rechargeable battery and requires no maintenance.



Alere Pima™ CD4 Cartridge Kit

Product

Alere Pima™ CD4 (Cartridge Kit 100x)

Alere Pima™ CD4 (Cartridge Kit 25x)

Product Code

260100100

260100025

Size

100

25

Finger Stick Sample Collection Kit

Product

Finger Stick Sample Collection Kit containing

4 x units of safety lancets

4 x units of gauze swabs

1 x unit of alcoholic swabs

4 x units of plasters

Product Code

260400199

Alere Pima™ CD4 Analyser

Product

Alere Pima™ Analyser

Product Code

260300003

Size

1

The Alere Pima™ Bead Standard is an internal standard for daily quality control (QC) on the Alere Pima Analyser. It comprises two ready-to-use cartridges, Pima Beads “Normal” and Pima Beads “Low”.

Product

Alere Pima™ CD4 Test Bead standard

Alere Pima™ CD4 Test Accessories

Product

Alere Pima™ Analyser Bag

Alere Pima™ Printer

Alere Pima™ Printer Paper I, non adhesive

Alere Pima™ Printer Paper II, coated, adhesive

Alere Pima™ CD4 Bead Standard

Alere™ Connectivity Pack I Samba3G-E-Modem, extension cord

Alere™ Connectivity Pack II CT63 Terminal incl. external FME

antenna, power adaptor incl. cord, USB connection cord

Alere Connectivity Pack III: USB 2.0 Ethernet adapter

Alere Power Pack: external rechargeable battery for the Pima Analyser

Alere Solar Panel: Solar Panel to charge Alere Power Pack

Product Code

260400001

260400007

260400009

260400010

260400011

260400015

260400016

260400046

260400017

260400041

Size

1

10

10

one set

1

1

1

1

1

Infectious Diseases

Strep A

For the qualitative detection of group A Streptococcal antigen from throat swabs.

Product	Product Code	Size
Alere™ TestPack +Plus Strep A	505715J	20
	505796J	40



Women's Health

hCG

A rapid immunoassay for the qualitative detection of human Chorionic Gonadotropin (hCG) in urine or urine and serum specimens.

Product	Product Code	Size
Alere™ TestPack +Plus hCG Combo	505805J	20



A rapid immunoassay for the qualitative detection of human Chorionic Gonadotropin (hCG) in urine.

Product	Product Code	Size
Alere™ TestPack hCG Urine	505798J	20



Alere Triage® System

Alere Triage® Meters are portable fluorescence instruments used to run our expansive Alere Triage® test kit portfolio. The Meters are designed to provide diagnostic results quickly and easily in both laboratory and point of care settings. Quantitative results are displayed and printed in an easy to read format and are conveniently highlighted when the result is above a given cutoff. Qualitative results are displayed as positive (POS) or negative (NEG), eliminating the need for complicated interpretation.



Product	Size	Product Code
Alere Triage® MeterPro	1	55071
Alere Triage® MeterPro User's Manual – English, printed version	1	26213EN
Alere Triage® MeterPro CD – multi-lingual CD	1	26213
Paper Roll for Alere Triage® Meter & Alere Triage® MeterPro	1	52031
Alere Triage® Supervisor Code Chip™ Module	1	80035A
Alere Triage® Meter Power Cord, European for Meter or scanner	1	52059
Alere Triage® Meter Power Cord, UK for Meter or scanner	1	52060
Triage® MeterPro Barcode Scanner Set	1	52142
Y-adapter	1	52098

Alere Triage® Data Management/Connectivity

Product	Size	Product Code
Alere Triage Census® ATM Software Kit (includes the 3 items listed below)	1	55069
Alere Triage® MeterPro Serial Interface Cable (RS-232 DB9M-DB9F)	1	52039
Triage Census® ATM Software Version 3.0x.x User's Guide	1	22643
Triage Census® ATM Software Version 3.0x.x CD	1	52140

Cardiovascular

Alere Triage® Cardiac Panel

Alere Triage® Cardiac Panel is a fluorescence immunoassay, to be used with the Alere Triage® Meters, for quantitative measurements of myoglobin, creatine kinase MB (CK-MB) and troponin I in EDTA whole blood or plasma specimens to be used as an aid in the diagnosis of acute myocardial infarction.

Product	Size	Product Code
Alere Triage® Cardiac Panel (25 tests)		97000HS / 97000HSEU*
Alere Triage® Total 5 Controls 1	5 x 0.25ml	88753
Alere Triage® Total 5 Controls 2	5 x 0.25ml	88754
Alere Triage® Total 5 Calibration Verification		88755



Alere Triage® Natriuretic Peptides

The test is intended to be used as an aid in the diagnosis and assessment of severity of congestive heart failure (also referred to as heart failure). The test also is used for the risk stratification of patients with acute coronary syndromes and for the risk stratification of patients with heart failure.

Product

	Size	Product Code
Alere Triage® BNP Test	25 tests	98000XR / 98000XREU*
Alere Triage® Total 5 Controls 1	5 x 0.25ml	88753
Alere Triage® Total 5 Controls 2	5 x 0.25ml	88754
Alere Triage® Total 5 Calibration Verification		88755
Alere Triage® BNP 15 Test Kit		98010XR / 98010XREU*



The test is to be used as an aid in the diagnosis of individuals suspected of having congestive heart failure (also referred to as heart failure). The test is also to be used as an aid for the risk stratification of patients with heart failure and the risk stratification of patients with acute coronary syndromes (ACS). The test may also serve as an aid in the assessment of increased risk for cardiovascular events and mortality in patients at risk for heart failure who have stable coronary artery disease.

Product

	Product Code
Alere Triage® NT-proBNP Test	98700EU*
Alere Triage® NT-proBNP Control 1	98713EU*
Alere Triage® NT-proBNP Control 2	98714EU*



Alere Triage® Profiler SOB™ Panel

Alere Triage® Profiler SOB™ Panel is a fluorescence immunoassay, to be used with the Alere Triage® Meters, for quantitative measurements of myoglobin, CK-MB, troponin I, B-type natriuretic peptide and D-Dimer in EDTA whole blood or plasma specimens to be used as an aid in the diagnosis of myocardial infarction, aid in the diagnosis and assessment of heart failure disease severity, an aid in risk stratification acute coronary syndromes and as an aid in the assessment and evaluation of thromboembolic events including pulmonary embolism.

Product	Size	Product Code
Alere Triage® Profiler SOB™ Panel (25 tests)		97300 / 97300EU*
Alere Triage® Total 5 Controls 1	5 x 0.25ml	88753
Alere Triage® Total 5 Controls 2	5 x 0.25ml	88754
Alere Triage® Total 5 Calibration Verification		88755



Alere Triage® Troponin I Test

Alere Triage® Troponin I Test is a fluorescence immunoassay, to be used with the Alere Triage® Meters, for the quantitative determination of troponin I in EDTA whole blood or plasma specimens to be used as an aid in the diagnosis of myocardial infarction (injury).

Product	Size	Product Code
Alere Triage® Troponin I Test	25 tests	98600EU*
Alere Triage® Total 3 Control 1		88733*
Alere Triage® Total 3 Control 2		88734*
Alere Triage® Total 3 Calibration Verification		88735*



Alere Triage® Cardio2 Panel

Alere Triage® Cardio2 Panel is a fluorescence immunoassay, to be used with the Alere Triage® Meters, for the rapid quantitative determination of troponin I and B-type natriuretic peptide (BNP) in EDTA whole blood or plasma specimens to be used as an aid in the diagnosis of myocardial infarction, aid in the diagnosis and assessment of heart failure disease severity and as an aid in risk stratification acute coronary syndromes.

Product	Product Code
Alere Triage® Cardio2 Panel	97500EU*
Alere Triage® Total 3 Control 1	88733*
Alere Triage® Total 3 Control 2	88734*
Alere Triage® Total 3 Calibration Verification	88735*



Alere Triage® Cardio3 Panel

Alere Triage® Cardio3 Panel is a fluorescence immunoassay, to be used with the Alere Triage® Meters, for the rapid quantitative determination of troponin I, CK-MB and B-type natriuretic peptide (BNP) in EDTA whole blood or plasma specimens to be used as an aid in the diagnosis of myocardial infarction, aid in the diagnosis and assessment of heart failure disease severity and as an aid in risk stratification acute coronary syndromes.

Product

Alere Triage® Cardio3 Panel
 Alere Triage® Total 3 Control 1
 Alere Triage® Total 3 Control 2
 Alere Triage® Total 3 Calibration Verification

Product Code

97400EU*
 88733*
 88734*
 88735*



Alere Triage® D-Dimer Test

A rapid, quantitative immunoassay to aid in the assessment and evaluation of patients suspected of having disseminated intravascular coagulation (DIC) and thromboembolic events including pulmonary embolism (PE) and deep vein thrombosis (DVT).

Product

Alere Triage® D-Dimer Test
 Alere Triage® Total 5 Control 1
 Alere Triage® Total 5 Control 2
 Alere Triage® Total 5 Calibration Verification

Size

25 tests

Product Code

98100 / 98100EU*
 88753
 88754
 88755



Alere Triage® BNP for Beckman Coulter Immunoassay System

The Alere Triage® BNP test is intended for use with Beckman Coulter Access Family of Immunoassay Systems for the In Vitro quantitative measurement of B-type natriuretic peptide (BNP) in plasma specimens using EDTA as the anticoagulant. The test is intended to be used for the following indications:

- as an aid in the diagnosis of congestive heart failure (also referred to as heart failure)
- as an aid in the assessment of severity of congestive heart failure
- for the risk stratification of patients with acute coronary syndromes
- for the risk stratification of patients with heart failure



Product

Alere Triage® BNP Test for Beckman Coulter® Immunoassay Systems
 Alere Triage® BNP QC Controls for Beckman Coulter® Immunoassay Systems
 Alere Triage® BNP Calibrators for Beckman Coulter® Immunoassay Systems

 Alere Triage® BNP Test on Beckman Coulter® Immunoassay Systems
 Package Insert CD

Product Code

98200
 98201
 98202

 26259

Size

100 tests
 3 2 x 2.5ml
 1 x 4.0ml +
 5 x 0.5ml

Acute Kidney Injury

Alere Triage® NGAL Test and Alere Triage® CardioRenal Panel

The Alere Triage® CardioRenal Panel is a rapid, point of care fluorescence immunoassay to be used with the Alere Triage® Meter for the quantitative determination of neutrophil gelatinase-associated lipocalin (NGAL) and B-type natriuretic peptide (BNP) in EDTA anticoagulated whole blood and plasma specimens. The NGAL test is used as an aid in the diagnosis of acute kidney injury. The BNP Test is used as an aid in the diagnosis and assessment of severity of congestive heart failure (also referred to as heart failure), an aid in the risk stratification of patients with heart failure, and an aid in the risk stratification of patients with acute coronary syndromes.



Product	Product Code
Alere Triage® CardioRenal Panel	98420EU*
Alere Triage® NGAL Test	98430EU*
Alere Triage® CardioRenal Control 1	98423EU*
Alere Triage® CardioRenal Control 2	98424EU*

Placental Growth Factor

Alere Triage® PLGF Test

The Alere Triage® PLGF Test is a rapid immunoassay to be used with the Alere Triage® MeterPro for the quantitative determination of Placental Growth Factor (PLGF) in EDTA anticoagulated plasma specimens.

PLGF is a biomarker of placental function and is found to be at abnormally low levels in pregnant women with preterm pre-eclampsia. Alere Triage® PLGF is used as an aid in the diagnosis of preterm pre-eclampsia and as an aid in the prognosis for preterm delivery in conjunction with other diagnostic and clinical information in women presenting with signs or symptoms of pre-eclampsia before 35 weeks gestation.



Product	Product Code	Size
Alere Triage® PLGF Test	98800EU*	25 tests
Alere Triage® PLGF Control 1	98813EU*	5 x 0.25ml
Alere Triage® PLGF Control 2	98814EU*	5 x 0.25ml

Toxicology

Alere Triage® TOX Drug Screen

The Alere Triage® TOX Drug Screen is a fluorescence immunoassay to be used with the Alere Triage® Meters for the qualitative determination of the presence of drug and/or the major metabolites above the threshold concentrations of up to 10 distinct drug classes, including assays for acetaminophen/paracetamol, amphetamines, methamphetamines, barbiturates, benzodiazepines, cocaine, methadone, opiates, phencyclidine, THC and tricyclic antidepressants in urine.



Product	Product Code	Size
Alere Triage® TOX Drug Screen	94400 / 94400EU*	
Alere Triage® TOX Drug Screen Control 1	94413	5 x 0.25ml
Alere Triage® TOX Drug Screen Control 2	94414	5 x 0.25ml

Alere Triage® Drugs of Abuse Panels

The Alere Triage® Drugs of Abuse Panel is an immunoassay used for the qualitative determination of the presence of the major metabolites of the following drugs: MTD, BZO, COC, AMP/mAMP, OPI, BAR, THC, TCA and gives quality results in about 10 minutes.

Product	Product Code
Alere Triage® 8 Panel for Drugs of Abuse (25 tests)	90008
Alere Triage® Drugs of Abuse Urine Controls, Negative Level 1 (6 x 5mL)	52229
Alere Triage® Drugs of Abuse Urine Controls, Positive Level 2 (6 x 5mL)	52230

* Not available for sale in the US.

** Not CE Marked.

Afinion™

The Afinion™ AS100 Analyzer is designed to enable simple and fast on-the-spot testing, regardless of the sample type, which may be whole blood, plasma or urine, and to deliver accurate results during patient consultation for improved patient management.



Product	Product Code	Country
Afinion™ AS100 Analyzer	1115027	EUR
	1115173	GB
	1115174	JP
	1115176	NO-DK-FI

Afinion™ ACR

Afinion™ ACR is an in vitro diagnostic test for quantitative determination of albumin, creatinine and albumin/creatinine ratio (ACR) in human urine. The measurement of urine albumin, creatinine and ACR aids in the early diagnosis of nephropathy.



Product	Product Code	Country	Size
Afinion™ ACR	1115021	EUR	15
	1115022	NOR	15

Afinion™ CRP

The Afinion™ CRP is a rapid in vitro diagnostic test for the quantitative determination of C-reactive protein (CRP) in human whole blood except from newborns, serum and plasma. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases.



Product	Product Code	Country	Size
Afinion™ CRP	1115013	EUR	15
	1115014	NOR	15
	1115377	JP	
Afinion™ CRP Control	1114993	EUR	

Afinion™ HbA1c

Test for quantitative determination of glycated hemoglobin (HbA1c) in human whole blood, used to monitor metabolic control in patients with diabetes. Aids in the diagnosis of diabetes and in identifying patients who may be at risk for developing diabetes.



Product	Product Code	Country	Size
Afinion™ HbA1c	1115016	EUR	15
	1115017	NOR	15
Afinion™ HbA1c Control	1115052		

Afinion™ Lipid Panel

Test for the quantitative determination of - Total Cholesterol (Chol) - High-Density Lipoprotein cholesterol (HDL) - Low-Density Lipoprotein cholesterol (LDL) - Triglycerides (Trig) - non-HDL - Chol/HDL ratio in whole blood.



Product	Product Code	Country	Size
Afinion™ Lipid Panel	1115779	NOR	15
	1115781	EUR	15
Afinion™ Lipid Panel Control	1115783	EUR	

Afinion™ Data Connectivity Converter

The Afinion™ Data Connectivity Converter (ADCC) is designed to automatically transfer Patient and Control results from the Afinion™ AS100 Analyzer to an electronic journal system (HIS/LIS) at doctor's offices, hospitals or similar. The results are transferred through Ethernet in the standard formats HL7 or AS™ (selectable by configuration).

Product	Product Code	Size
Afinion™ Data Connectivity Converter (ADCC)	1115660	
Afinion™ Microfiber cloth	1115186	
Afinion™ Analyzer Cleaning kit	1115269	Kit 1x5

NycoCard™

The NycoCard™ CRP is an in vitro diagnostic test for the quantitative determination of C-reactive protein (CRP) in whole blood, serum and plasma. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases.

Product	Product Code	Country	Size
NycoCard™ CRP Single Test	1113002		24
	1112999		48
	1114491	JP	48
	1115699	IN	24
	1115698	IN	48
Capillary Holder	1003668		
NycoCard™ Reader II Pen Tip	1047780		1x10 pc

NycoCard™ Reader II Accessories

Product	Product Code	Size
NycoCard™ Reader II Pen Tip	1047780	1
NycoCard™ Reader II Pen Rings	1047781	10 pcs
NycoCard™ Reader II Calibration Device	1047782	1x10 pc
NycoCard™ Reader II, In Field Calibr. Kit	1111562	

NycoCard™

NycoCard™ Test Accessories

Product	Product Code	Size
Capillary Holder	1003668	
Mini-Pet Stand	1113059	
Test Tube Rack	1113060	6 tubes
MINI-PET	1114345	5 µl

NycoCard™ Reader II

Product	Product Code
NycoCard™ Reader II	1113117

NycoCard™ D-Dimer

NycoCard™ D-Dimer is an in vitro test based upon an immunometric flow-through principle for rapid determination of the fibrin degradation product D-dimer in plasma.

Product	Product Code	Country	Size
NycoCard™ D-Dimer Single Test	1114300		24
	1114301		6
	1115697	IN	24

NycoCard™ HbA1c

Test for quantitative determination of glycated hemoglobin (HbA1c) in human whole blood, used to monitor metabolic control in patients with diabetes.

Product	Product Code	Size
NycoCard™ HbA1c Single Test	104218	24
NycoCard™ HbA1c Control	1115271	2x1,5ml

NycoCard™ U-Albumin

Test for the quantitative determination of albumin in human urine, used for early identification of renal disease in patients with diabetes and/or hypertension.

Product	Product Code	Size
NycoCard™ U-Albumin Contro	1114358	2x1x1,0 ml

Infectious Disease

Product	Product Code	Size
MICROSYPH™ TPHA200	FTPHA200	200
MICROSYPH™ TPHA1000	FTPHA1000	1000
Rapid Plasma Reagin (RPR)	FRPR500	500

Women's Health

Product	Product Code	Size
Active-B12 EIA (kit)	FMABT100	96

Infectious Diseases - Respiratory

Strep A

For the qualitative detection of group A Streptococcal antigen from throat swabs.

Product	Product Code	Size
BinaxNOW® Strep A	730-025	25
BinaxNOW® Strep A Control Swabs	730-010	5 positive & 5 negative

Parasitology

Filariasis

For the qualitative detection of *Wuchereria bancrofti* antigen in whole blood, serum and plasma specimens.

Product	Product Code	Size
BinaxNOW® Filariasis	620-000	25

G6PD Deficiency

G6PD Rapid

The BinaxNOW® G6PD (Glucose-6-Phosphate Dehydrogenase) Test is an *in vitro* enzyme chromatographic test for the qualitative detection of G6PD enzyme activity in human venous whole blood, collected in heparin or ethylenediaminetetraacetic acid (EDTA). The BinaxNOW® G6PD Test is a visual screening test used for differentiating normal from deficient G6PD activity levels in whole blood and is intended to aid in the identification of people with G6PD deficiency.

Product	Product Code	Size
BinaxNOW® G6PD Rapid	780-000	25

Cardiovascular

D-dimer

For the qualitative detection of D-dimer in human whole blood, capillary blood and plasma specimens.

Product

Clearview® Simplify D-dimer Test Kit
 Clearview® Simplify D-dimer External Controls

Product Code

505285T
 6101KCVAT

Size

10
 1 Positive and
 1 Negative Control



Infectious Diseases - Gastrointestinal Disorders

Helicobacter pylori

For the qualitative detection of IgG antibodies specific to Helicobacter pylori in human serum, plasma and whole blood specimens.

Product

Clearview® H. pylori

Product Code

CV506757

Size

20



Infectious Diseases - Respiratory

Infectious Mononucleosis

For the qualitative detection of infectious mononucleosis IgM heterophile antibodies in human whole blood, serum and plasma specimens.

Product

Clearview® IM

Product Code

CV506815

Size

20



Influenza A & B

For the qualitative detection and differentiation of influenza A and influenza B virus antigen from nasopharyngeal swabs.

Product

Clearview® Exact Influenza A & B

Product Code

CV506761
 CV506834

Size

20
 5



Respiratory Syncytial Virus

For the qualitative detection of respiratory syncytial virus fusion protein directly from nasopharyngeal swab or nasal aspirate specimens in children less than 6 years and adults over the age of 60.

Product	Product Code	Size
Clearview® RSV	504957	20



Strep A

For the qualitative detection of group A Streptococcal antigen from throat swabs.

Product	Product Code	Size
Clearview® Exact Strep A Cassette	CV506758	25
Clearview® Exact Strep A Dipstick	CV506759	25



Tuberculosis

Immunosorbent assay (ELISA) based system that qualitatively detects the presence of lipoarabinomannan (LAM) antigen of mycobacteria in human urine as an aid in the diagnosis of mycobacterial infection in patients with HIV.

Product	Product Code	Size
Clearview® TB ELISA	M913	96 Wells



Infectious Diseases - Sexually Transmitted Infections

Chlamydia

For the qualitative detection of *Chlamydia trachomatis* antigen from female endocervical swab specimens and male urine specimens.

Product	Product Code	Size
Clearview® <i>Chlamydia</i> MF	CV504524	20
Clearview® <i>Chlamydia</i> Female Specimen Kit	CV506816	20 swabs
Clearview® Workstation	UK 500189 Euro 500185 ROW 500200	



Infectious Diseases - Malaria

Malaria

A rapid, qualitative immunoassay for the detection of HRP-2 (Histidine-rich protein 2) specific to Plasmodium falciparum (p.f.) in capillary or venous whole blood.

Product	Product Code	Size
Clearview® Malaria P.f.	VB01	25
Clearview® Malaria Dual	VB20	25

Malaria Combo

A simple, rapid, qualitative and cost effective method for the differential detection of HRP-2 and Paldolase specific to Plasmodium falciparum (P.f.), Plasmodium vivax (P.v.), Plasmodium malariae (P.m.) and Plasmodium ovale (P.o.) antigens in whole blood.

Product	Product Code	Size
Clearview® Malaria Combo	VB11	25

Tumour Markers

Faecal Occult Blood (Immunochemical)

For the qualitative detection of Human Occult Blood in faecal specimens.

Product	Product Code	Size
Clearview® iFOBT	CV506149	25

Women's Health

hCG

For the qualitative detection of human Chorionic Gonadotrophin (hCG) in urine or urine and serum specimens.

Product	Product Code	Size
Clearview® HCG	CV506788	20
Clearview® EASY HCG	E.S.P.G. UPK009	20
	F.I.G.D. UPK010	20
	N.S.F.D. UPK011	20
Clearview® hCG COMBO	CV506760	40



Vaginal pH

For the measurement of vaginal pH as part of a self-examination programme to help reduce the risk of pre-term birth.

Product	Product Code	Size
Careplan® VpH Test Gloves	505006	20
	504102	50



DoubleCheckGold™

HIV 1&2

HIV 1&2

For the qualitative detection of antibodies to HIV-1 and HIV-2 in human serum, plasma and whole blood specimens.



Product	Product Code	Size
DoubleCheckGold™ HIV 1&2 Serum/Plasma	70632020	20
DoubleCheckGold™ HIV 1&2 Serum/Plasma	70632100	100
DoubleCheckGold™ HIV 1&2 Whole Blood/Serum/Plasma	70633020	20
DoubleCheckGold™ HIV 1&2 Whole Blood/Serum/Plasma	70633100	100
DoubleCheckGold™ Ultra HIV 1&2	70635020	20
DoubleCheckGold™ Ultra HIV 1&2	70635100	100

epoc® Blood analysis system – The epoc® point of care blood analysis system is a wireless bedside testing solution that uses SmartCard technology. This breakthrough platform provides state-of-the-art lab results directly to a hand-held PDA. It's easy to use and interfaces with existing wireless systems. The epoc® BGEM test card is a blood gas, electrolyte and metabolite panel consisting of the following:



- measured analytes: pH, PCO₂, PO₂, Na, K, iCa, Glu, Lac, Crea, Cl, Hct
- calculated values: HCO₃, TCO₂, Base Excess, sO₂, Hb, eGFR, AGap

Product	Item No.	Pack Size
epoc® System Hardware		
epoc® Reader (Incl. Power Supply)	HR-1002-00-00	1
epoc® Host ² (Motorola MC55A0-HC)	HH-1009-00-00	1
epoc® Data Management		
epoc® Enterprise Data Manager		
epoc® System Data Manager Software	SD-0001-00-00	1
epoc® Accessories		
epoc® Care-Film Capillary Tubes (90uL)		
Mylar-wrapped, glass, capillary tubes	AT-0006-00-00	50
Paper, Thermal Printer		
Portable printer paper rolls	AP-0002-00-00	50
Power Supply, epoc Reader	AR-0007-00-00(EU)	1
Power Supply, epoc Reader	AR-0008-00-00 (AU)	1
epoc® Host Accessories (Motorola)		
Stylus, epoc Mobile Computer Stylus		
for epoc Host	AH-1007-00-00	3
Battery, 3600 MAH, replacement	AH-1008-00-00	1
Protector, Screen, replacements	AH-1009-00-00	3
Charger, 4 Slot, Battery	AH-1010-00-00	1
Power Supply, Charger (US)	AH-1012-00-00	1
Cord, Line, 3 wire for Charger	AH-1011-00-00	1

Product	Item No.	Pack Size
epoc[®] Test Cards		
epoc BGEM Test Card (pH, pCO ₂ , pO ₂ , Na ⁺ , K, Ca ⁺⁺ , Hct, Glu, Lac)	CT-1004-00-00	50
epoc BGEM2 Test Card (pH, pCO ₂ , pO ₂ , Na ⁺ , K, Ca ⁺⁺ , Hct, Glu, Lac, Crea, Cl ⁻)	CT-1006-00-00	50
epoc[®] Calibration and Control Fluids		
Eurotrol [®] GAS-ISE Metabolites Level 1	179.001.010	10 ampoules x 2.5mL
Eurotrol [®] GAS-ISE Metabolites Level 2	179.002.010	10 ampoules x 2.5mL
Eurotrol [®] GAS-ISE Metabolites Level 3	179.003.010	10 ampoules x 2.5mL
Eurotrol [®] epoc [®] Calibration Verification Fluids Levels 1–5	183.000.005	4 vials each level, 2.5 mL
Eurotrol [®] epoc [®] Hct Control, Level A	195.002.010	10 vials x 2.5 mL
Eurotrol [®] epoc [®] Hct Control, Level B	195.004.010	10 vials x 2.5 mL
Eurotrol [®] epoc [®] Hct Control, Level C	195.003.010	10 vials x 2.5 mL
Eurotrol [®] epoc [®] Hct Linearity Verification Material Levels 1–5	190.000.005	4 vials each level, 2.5 mL
epoc[®] User Documentation		
epoc [®] System Manual System Manual, English	UM-0001-US-EN	1
epoc[®] Service Agreements		
epoc [®] Blood Analysis System SA Host Mobile Computer & Reader, 12 Months	WB-0001-00-00	
epoc [®] Reader SA epoc Reader, 12 Months	WB-0002-00-00	
epoc [®] Host Mobile Computer SA Host Mobile Computer, 12 Months	WB-0003-00-00	
epoc [®] Enterprise Data Manager SA EDM Software, 12 Months	WB-0004-00-00	
epoc [®] Mobile Printer SA Thermal printer, 12 Months	WP-0006-00-00	
Wireless Printers		
Epson TM-P60, IEEE820.11b (Wi-Fi) with WPA2 security, Portable Thermal Printer Local Purchase Recommended		
Epson TM-P60, Bluetooth, Portable Thermal Printer Local Purchase Recommended		
Zebra MZ 320, Portable Thermal Printer Local Purchase Recommended		
Sanei BLM-80BT, Mobile Printer Local Purchase Recommended		

Please contact your local representative for full ordering details.

Infectious Diseases

A rapid EIA for the detection of antibodies or antigens in human serum or plasma



Product	CE-Mark	Product Code	Size
ImmunoComb® Chlamydia trachomatis IgA	√	50411002	36
ImmunoComb® Chlamydia trachomatis IgG	√	50410002	36
ImmunoComb® II Chlamydia trach. Monovalent IgA	√	60412002	36
ImmunoComb® Chlamydia Bivalent IgG	√	50416002	36
ImmunoComb® II HAV IgM	√	60457002	36
ImmunoComb® II HAV Ab	√	60456002	36
ImmunoComb® II Anti HBs		60453002	36
ImmunoComb® II HBs Ag	√	60451002	36
ImmunoComb® II HBs Ag 90'		60454002	36
ImmunoComb® II HCV	√	60455002	36
ImmunoComb® II CMV IgG	√	60460002	36
ImmunoComb® II CMV IgM	√	60461002	36
ImmunoComb® II Rubella IgG	√	60400002	36
ImmunoComb® II Rubella IgM	√	60401002	36
ImmunoComb® Toxo IgG		50440002	36
ImmunoComb® Toxo IgM		50441002	36
ImmunoComb® II HIV 1& 2 BiSpot	√	60432002	36
ImmunoComb® II HIV 1& 2 TriSpot Ag-Ab		60433002	36
ImmunoComb® II HIV 1& 2 CombFirm		60434002	18
ImmunoComb® II Helicobacter pylori	√	60425002	36
ImmunoComb® II Chagas Ab		60481002	36
ImmunoComb® II HTLV 1 & 2		60482002	36
ImmunoComb® II Dengue IgG & IgM BiSpot		60483002	36

Infectious Diseases

For early detection of HIV and p24 antigen: ImmunoComb® II HIV 1&2 TriSpot Ag-Ab

A 4th generation HIV assay for the qualitative and differential detection of antibodies to Human Immunodeficiency Viruses types 1 and 2 (HIV-1 and HIV-2) and for the simultaneous detection of HIV p24 antigen in human serum or plasma.

Product	Product Code	Size	Detectable HIV Types
ImmunoComb® II HIV 1&2 TriSpot Ag-Ab	60433002	36	HIV 1&2, Sub-type O Ab and p24 Ag

For HIV Blood Screening: ImmunoComb® II HIV 1&2 BiSpot

A rapid test intended for the qualitative and differential detection of IgG antibodies to HIV 1&2 in human serum or plasma.

Product	Product Code	Size	Detectable HIV Types
ImmunoComb® II HIV 1&2 BiSpot	60432002	36	HIV 1&2 and Subtype O

For confirmation of HIV positive result: ImmunoComb® II HIV 1&2 CombFirm

A quick confirmatory test for the presence of HIV antibodies in samples found positive by a screening test. The test employs five distinct recombinant and synthetic HIV antigens for the rapid serological determination of the HIV antibody profile.

Product	Product Code	Size	Detectable HIV Types
ImmunoComb® II HIV 1&2 CombFirm	60434002	18	HIV 1&2 and Subtype O

Equipment

CombScan™

An optical reader for Origenics ImmunoComb® kits

The CombScan™ provides objective quantitative or qualitative results of ImmunoComb® tests

- High sensitivity
- High reliability of results
- Clinical results, Absolute and Relative absorbance values can be viewed, as well as the ImmunoComb® image
- Results may be saved and easily retrieved, printed and viewed on a computer's screen

The CombScan™ package includes:

- Scanner - reads the comb by using optical pattern recognition algorithm on the scanned image.
- CombScan™ Software - interprets and displays the results on the computer's screen
- Validation accessories

Minimum requirements for installation and setup of CombScan™:

- Pentium PC with a USB connection
- Microsoft Windows operation systems: '98, Millennium, 2000, XP

Product	Product Code
CombScan™	60590001

Bone Health

Cross-linked N-telopeptides of bone type I collagen (NTx)

For the quantitative detection of NTx in urine and serum specimens as an indicator of human bone resorption.

Product

Osteomark® NTx Serum ELISA

Product Code

504836

Size

96 well plate with
12 x 8 microwell strips



Osteomark® NTx Urine ELISA

504837

96 well plate with
12 x 8 microwell strips



ELISA Kits



Product	Product Code	Size
Barmah Forest Virus IgG	E-BFV01G	96 tests
Barmah Forest Virus IgM	E-BFV01M	96 tests
<i>Bordetella pertussis</i> IgA	E-BPB01A	96 tests
<i>Bordetella pertussis</i> IgG	E-BPB01G	96 tests
<i>Coxiella burnetii</i> (Q Fever) IgG	E-QFB01G	96 tests
<i>Coxiella burnetii</i> (Q Fever) IgM	E-QFB01M	96 tests
Dengue IgG Indirect	E-DEN01G	96 tests
Dengue IgM Capture	E-DEN01M	96 tests
Dengue Early (2nd Generation)	E-DEN02P	96 tests
Dengue IgG Capture	E-DEN02G	96 tests
Dengue Duo IgM Capture and IgG Capture	E-DEN01D	96 tests
Epstein-Barr Virus VCA IgG	E-EBV01G	96 tests
Epstein-Barr Virus VCA-p18 IgG	E-EBV04G	96 tests
Epstein-Barr Virus VCA IgM	E-EBV01M	96 tests
Epstein-Barr Virus EBNA IgG	E-EBV02G	96 tests
EBV VCA M-One	E-EBV03M	96 tests
Human Herpesvirus-6 IgG	E-HV601G	96 tests
Human Herpesvirus-6 IgM	E-HV601M	96 tests
Japanese Encephalitis - Dengue IgM Combo	E-JED01C	96 tests
<i>Leptospira</i> IgM	E-LEP01M	96 tests
Ross River Virus IgG	E-RRV01G	96 tests
Ross River Virus IgM	E-RRV01M	96 tests
Varicella-zoster IgG	E-VZV01G	96 tests
West Nile Virus IgM Capture	E-WNV02M	96 tests

Rapid Tests

Dengue Early Rapid Test

The Panbio® Dengue Early Rapid is a dengue NS1 antigen immunochromatographic assay. It is for the qualitative detection of NS1 antigen in serum, used as an aid in the clinical laboratory diagnosis of patients with clinical symptoms consistent with dengue fever.

Product	Product Code	Size
Dengue Early Rapid test	R-DEN01P	25 tests

Dengue Duo Cassette

The Panbio® Dengue Duo Cassette is for the qualitative presumptive detection of IgM and elevated high titre IgG antibodies to dengue virus in human serum, plasma and whole blood. The assay can be used for the presumptive differentiation between primary and secondary infection.

Product	Product Code	Size
Dengue Duo Cassette	R-DEN03D	25 tests

Miscellaneous Reagents

Product	Product Code	Size
Buffered Avidity Reagent	O-AVR02G	25 mL

SureStep™ Drug of Abuse Tests

Drug Names - Abreviation & Cutoffs

Product	Abbreviation	Urine		Oral Fluid
		Standard Cut-Off (ng/mL)	Alternative Cut-Off (ng/mL)	Standard Cut-Off (ng/mL)
Amphetamine	AMP	1000	500, 300	50
Barbiturates	BAR	300	-	-
Benzodiazepines	BZO	300	200	10
Buprenorphine	BUP	10	-	-
Cocaine	COC	300	150	20
Cotinine*	COT	100	-	30
Ecstasy	MDMA	500	-	-
Fentanyl	FTY	20	-	-
Ketamine	KET	1000	-	-
Marijuana	THC	50	150, 20	100
Methadone	MTD	300	-	30
Methadone Metabolite	EDDP	300	100	-
Methamphetamine	MET	1000	500, 300	50
Morphine	MOP	300	-	-
Opiate	OPI	2000	-	40
Oxycodone	OXY	100	-	20
Phencyclidine	PCP	25	-	10
Propoxyphene	PPX	300	-	-
Tramadol	TRA	100	-	-
Tricyclic Antidepressant	TCA	1000	-	-

Single Drug Format: Strip Test In Sealed Pouch



Product Format: Cut-Off Level	Test Per Kit
SureStep™ AMP One Step Amphetamine Test Strip (AMP - 1000 ng/mL)	50
SureStep™ AMP 500 One Step Amphetamine Test Strip (AMP - 500 ng/mL)	50
SureStep™ AMP 300 One Step Amphetamine Test Strip (AMP - 300 ng/mL)	50
SuteStep™ BAR One Step Barbiturates Test Strip (BAR - 300 ng/mL)	50
SureStep™ BZO One Step Benzodiazepines Test Strip (BZO - 300 ng/mL)	50
SureStep™ BZO 200 One Step Benzodiazepines 200 Test Strip (BZO - 200 ng/mL)	50
SureStep™ BUP One Step Buprenorphine Test Strip (BUP - 10 ng/mL)	50
SureStep™ COC One Step Cocaine Test Strip (COC - 300 ng/mL)	50
SureStep™ COC 150 One Step Cocaine Test Strip (COC - 150 ng/mL)	50
SureStep™ COT One Step Cotinine Test Strip (COT - 100 ng/mL)	50
SureStep™ MDMA One Step Ecstasy Test Strip (MDMA - 500 ng/mL)	50
SureStep™ FTY One Step Fentanyl Test Strip (FTY - 20 ng/mL)	50
SureStep™ KET One Step Ketamine Test Strip (KET - 1000 ng/mL)	50
SureStep™ THC One Step Marijuana Test Strip (THC - 50 ng/mL)	50
SureStep™ THC 150 One Step Marijuana Test Strip (THC - 150 ng/mL)	50
SureStep™ THC 20 One Step Marijuana Test Strip (THC - 20 ng/mL)	50
SureStep™ MTD One Step Methadone Test Strip (MTD - 300 ng/mL)	50
SureStep™ EDDP One Step (Methadone Metabolite) Test Strip (EDDP - 300 ng/mL)	50
SureStep™ EDDP 100 One Step (Methadone Metabolite) Test Strip (EDDP - 100 ng/mL)	50
SureStep™ MET One Step Methamphetamine Test Strip (MET - 1000 ng/mL)	50
SureStep™ MET 500 One Step Methamphetamine Test Strip (MET - 500 ng/mL)	50
SureStep™ MET 300 One Step Methamphetamine Test Strip (MET - 300 ng/mL)	50
SureStep™ MOP One Step Morphine Test Strip (MOP - 300 ng/mL)	50
SureStep™ OPI One Step Opiate Test Strip (OPI - 2000 ng/mL)	50
SureStep™ OXY One Step Oxycodone Test Strip (OXY - 100 ng/mL)	50
SureStep™ PCP One Step Phencyclidine Test Strip (PCP - 25 ng/mL)	50
SureStep™ PPX One Step Propoxyphene Test Strip (PPX - 300 ng/mL)	50
SureStep™ TRA One Step Tramadol Test strip (TRA - 100 ng/mL)	50
SureStep™ TCA One Step Tricyclic Antidepressants Test Strip (TCA - 1000 ng/mL)	50

Single Drug Format: Device Format



Product Format: Cut-Off Level	Test Per Kit
SureStep™ AMP One Step Amphetamine Test Device (AMP - 1000 ng/mL)	25
SureStep™ AMP 500 One Step Amphetamine Test Device (AMP - 500 ng/mL)	25
SureStep™ AMP 300 One Step Amphetamine Test Device (AMP - 300 ng/mL)	25
SureStep™ BAR One Step Barbiturates Test Device (BAR - 300 ng/mL)	25
SureStep™ BZO One Step Benzodiazepines Test Device (BZO - 300 ng/mL)	25
SureStep™ BZO 200 One Step Benzodiazepines Test Device (BZO - 200 ng/mL)	25
SureStep™ BUP One Step Buprenorphine Test Device (BUP - 10 ng/mL)	25
SureStep™ COC One Step Cocaine Test Device (COC - 300 ng/mL)	25
SureStep™ COC 150 One Step Cocaine Test Device (COC - 150 ng/mL)	25
SureStep™ COT One Step Cotinine Test Device (COT - 100 ng/mL)*	25
SureStep™ MDMA One Step Ecstasy Test Device (MDMA - 500 ng/mL)	25
SureStep™ FTY One Step Fentanyl Test Device (FTY - 20 ng/mL)	25
SureStep™ KET One Step Ketamine Test Device (KET - 1000 ng/mL) *	25
SureStep™ THC One Step Marijuana Test Device (THC - 50 ng/mL)	25
SureStep™ THC 150 One Step Marijuana Test Device (THC - 150 ng/mL)	25
SureStep™ THC 20 One Step Marijuana Test Device (THC - 20 ng/mL)	25
SureStep™ MTD One Step Methadone Test Device (MTD - 300 ng/mL)	25
SureStep™ EDDP One Step (Methadone Metabolite) Test Device (EDDP - 300 ng/mL)	25
SureStep™ EDDP 100 One Step (Methadone Metabolite) Test Device (EDDP - 100 ng/mL)	25
SureStep™ MET One Step Methamphetamine Test Device (MET - 1000 ng/mL)	25
SureStep™ MET 500 One Step Methamphetamine Test Device (MET - 500 ng/mL)	25
SureStep™ MET 300 One Step Methamphetamine Test Device (MET - 300 ng/mL)	25
SureStep™ MOP One Step Morphine Test Device (MOP - 300 ng/mL)	25
SureStep™ OPI One Step Opiate Test Device (OPI - 2000 ng/mL)	25
SureStep™ OXY One Step Oxycodone Test Device (OXY - 100 ng/mL)	25
SureStep™ PCP One Step Phencyclidine Test Device (PCP - 25 ng/mL)	25
SureStep™ PPX One Step Propoxyphene Test Device (PPX - 300 ng/mL)	25
SureStep™ TRA One Step One Step Tramadol Test Device (TRA - 100 ng/mL)	25
SureStep™ TCA One Step Tricyclic Antidepressants Test Device (TCA - 1000 ng/mL)	25

Single Drug Format: Strip Test In Closed Canister*



Product Format: Cut-Off Level	Test Per Kit
SureStep™ AMP One Step Amphetamine Test Strip (AMP - 1000 ng/mL)	25
SureStep™ AMP 500 One Step Amphetamine Test Strip (AMP - 500 ng/mL)	25
SureStep™ AMP 300 One Step Amphetamine Test Strip (AMP - 300 ng/mL)	25
SuteStep™ BAR One Step Barbiturates Test Strip (BAR - 300 ng/mL)	25
SureStep™ BZO One Step Benzodiazepines Test Strip (BZO - 300 ng/mL)	25
SureStep™ BZO 200 One Step Benzodiazepines 200 Test Strip (BZO - 200 ng/mL)	25
SureStep™ BUP One Step Buprenorphine Test Strip (BUP - 10 ng/mL)	25
SureStep™ COC One Step Cocaine Test Strip (COC - 300 ng/mL)	25
SureStep™ COC 150 One Step Cocaine Test Strip (COC - 150 ng/mL)	25
SureStep™ MDMA One Step Ecstasy Test Strip (MDMA - 500 ng/mL)	25
SureStep™ THC One Step Marijuana Test Strip (THC - 50 ng/mL)	25
SureStep™ THC 150 One Step Marijuana Test Strip (THC - 150 ng/mL)	25
SureStep™ THC 20 One Step Marijuana Test Strip (THC - 20 ng/mL)	25
SureStep™ MTD One Step Methadone Test Strip (MTD - 300 ng/mL)	25
SureStep™ EDDP One Step (Methadone Metabolite) Test Strip (EDDP - 300 ng/mL)	25
SureStep™ EDDP 100 One Step (Methadone Metabolite) Test Strip (EDDP - 100 ng/mL)	25
SureStep™ MET One Step Methamphetamine Test Strip (MET - 1000 ng/mL)	25
SureStep™ MET 500 One Step Methamphetamine Test Strip (MET - 500 ng/mL)	25
SureStep™ MET 300 One Step Methamphetamine Test Strip (MET - 300 ng/mL)	25
SureStep™ MOP One Step Morphine Test Strip (MOP - 300 ng/mL)	25
SureStep™ OPI One Step Opiate Test Strip (OPI - 2000 ng/mL)	25
SureStep™ OXY One Step Oxycodone Test Strip (OXY - 100 ng/mL)	25
SureStep™ PCP One Step Phencyclidine Test Strip (PCP - 25 ng/mL)	25
SureStep™ PPX One Step Propoxyphene Test Strip (PPX - 300 ng/mL)	25
SureStep™ TRA One Step Tramadol Test Strip (TRA - 100 ng/mL)	25
SureStep™ TCA One Step Tricyclic Antidepressants Test Strip (TCA - 1000 ng/mL)	25
Urines Adulteration Test Strip	25

Single Drug Format: Dip Card Format



Product Format: Cut-Off Level	Test Per Kit
SureStep™ AMP One Step 1 Drug Screen Test Panel (AMP - 1000 ng/mL)	25
SureStep™ AMP 500 One Step 1 Drug Screen Test Panel (AMP - 500 ng/mL)	25
SureStep™ AMP 300 One Step 1 Drug Screen Test Panel (AMP - 300 ng/mL)	25
SureStep™ BAR One Step 1 Drug Screen Test Panel (BAR - 300 ng/mL)	25
SureStep™ BZO One Step 1 Drug Screen Test Panel (BZO - 300 ng/mL)	25
SureStep™ BZO 200 One Step 1 Drug Screen Test Panel (BZO - 200 ng/mL)	25
SureStep™ BUP One Step 1 Drug Screen Test Panel (BUP - 10 ng/mL)	25
SureStep™ COC One Step 1 Drug Screen Test Panel (COC - 300 ng/mL)	25
SureStep™ COC 150 One Step 1 Drug Screen Test Panel (COC - 150 ng/mL)	25
SureStep™ MDMA One Step 1 Drug Screen Test Panel (MDMA - 500 ng/mL)	25
SureStep™ KET One Step 1 Drug Screen Test Panel (KET - 1000 ng/mL)	25
SureStep™ THC One Step 1 Drug Screen Test Panel (THC - 50 ng/mL)	25
SureStep™ THC 150 One Step 1 Drug Screen Test Panel (THC - 150 ng/mL)	25
SureStep™ THC 20 One Step 1 Drug Screen Test Panel (THC - 20 ng/mL)	25
SureStep™ MTD One Step 1 Drug Screen Test Panel (MTD - 300 ng/mL)	25
SureStep™ EDDP One Step 1 Drug Screen Test Panel (EDDP - 300 ng/mL)	25
SureStep™ EDDP 100 One Step 1 Drug Screen Test Panel (EDDP - 100 ng/mL)	25
SureStep™ MET One Step 1 Drug Screen Test Panel (MET - 1000 ng/mL)	25
SureStep™ MET 500 One Step 1 Drug Screen Test Panel (MET - 500 ng/mL)	25
SureStep™ MET 300 One Step 1 Drug Screen Test Panel (MET - 300 ng/mL)	25
SureStep™ MOP One Step 1 Drug Screen Test Panel (MOP - 300 ng/mL)	25
SureStep™ OPI One Step 1 Drug Screen Test Panel (OPI - 2000 ng/mL)	25
SureStep™ OXY One Step 1 Drug Screen Test Panel (OXY - 100 ng/mL)	25
SureStep™ PCP One Step 1 Drug Screen Test Panel (PCP - 25 ng/mL)	25
SureStep™ PPX One Step 1 Drug Screen Test Panel (PPX - 300 ng/mL)	25
SureStep™ TRA One Step 1 Drug Screen Test Panel (TRA - 100 ng/mL)	25
SureStep™ TCA One Step 1 Drug Screen Test Panel (TCA - 1000 ng/mL)	25

Multi-Drug Test: Multi-Line Test Format: Multi-Line Test Device



Product Name	Code Extension	Test Per Kit	AMP 1000*	BAR 300	BZO 300*	BUP 10	COC 300*	MDMA 500	THC 50	MTD 300	EDDP 300*	MET 1000*	MOP 300	TCA 1000
SureStep™ Multi-Drug One Step Multi-Line 2 Drug Screen Test Device	011	25					■		■					
SureStep™ Multi-Drug One Step Multi-Line 5 Drug Screen Test Device	151	25	■		■		■		■				■	
SureStep™ Multi-Drug One Step Multi-Line 5 Drug Screen Test Device	171	25	■				■	■	■				■	
SureStep™ Multi-Drug One Step Multi-Line 5 Drug Screen Test Device	111	25	■				■		■			■	■	
SureStep™ Multi-Drug One Step Multi-Line 6 Drug Screen Test Device	071	25	■		■		■		■			■	■	
SureStep™ Multi-Drug One Step Multi-Line 6 Drug Screen Test Device	tbt	25	■			■	■		■			■	■	
SureStep™ Multi-Drug One Step Multi-Line 7 Drug Screen Test Device	011	25	■		■		■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 8 Drug Screen Test Device	061-xa	25	■ AMP 300	■	■		■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 8 Drug Screen Test Device	091 **	25	■ AMP 500			■	■	■	■	■		■ MET 500	■	
SureStep™ Multi-Drug One Step Multi-Line 8 Drug Screen Test Device	081-a **	25	■		■	■	■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 9 Drug Screen Test Device	061	25	■ AMP 300	■	■		■		■	■		■ MET 500	■	■
SureStep™ Multi-Drug One Step Multi-Line 9 Drug Screen Test Device	***	25	■		■	■	■		■	■	■	■	■	
SureStep™ Multi-Drug One Step Multi-Line 10 Drug Screen Test Device	041-b	25	■	■	■		■	■	■	■		■	■	■
SureStep™ Multi-Drug One Step Multi-Line 10 Drug Screen Test Device	081-a	25	■ AMP 500		■	■	■	■	■	■		■ MET 500	■	■
SureStep™ Multi-Drug One Step Multi-Line 11 Drug Screen Test Device	031-a	25	■	■	■	■	■	■	■	■		■	■	■
SureStep™ Multi-Drug One Step Multi-Line 11 Drug Screen Test Device	021-xa	25	■ AMP 300	■	■	■	■	■	■	■		■ MET 300	■	■

Other drug combinations on demand (also including combination with other drugs as OXY, PCP, PPX, TRA
 * Alternative cut offs (ng/mL): AMP 500, 300 - BZO 200 - COC 150 - THC 150, 20 - MET 500, 300 - EDDP 100
 ** Combination validated but needs SureStep silk-printing / label development
 *** Combination which needs to be validated - requires annual forecast

Multi-Drug Test Panel Format: Test Panel & Test Panel With Specimen Validity Test (Adulteration) Multi-Drug Screen Test Panel

Product Name	Code Extension	Test Per Kit	pH SG*** Oxi***	Cre*** Glut*** Nit***	AMP 1000*	BAR 300	BZO 300*	BUP 10	COC 300*	COT 100	FTY 20	KET 1000	MDMA 500	THC 50*	MTD 300	EDDP 300*	MET 1000*	MOP 300	OPI 2000	OXY 100	PCP 25	PPX 300	TRA 100	TCA 1000
SureStep™ Multi-Drug One Step 4 Drug Screen Test Panel with Adulteration Strips	See below **	25	■	■	4 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 5 Drug Screen Test Panel with Adulteration Strips	**	25	■	■	5 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 6 Drug Screen Test Panel with Adulteration Strips	**	25	■	■	6 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 7 Drug Screen Test Panel with Adulteration Strips	**	25	■	■	7 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 8 Drug Screen Test Panel with Adulteration Strips	**	25	■	■	8 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 9 Drug Screen Test Panel with Adulteration Strips	**	25	■	■	9 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 10 Drug Screen Test Panel with Adulteration Strips	**	25	■	■	10 drugs among 20 to select - All combinations are allowed																			

Multi-Drug Test Panel



Product Name	Code Extension	Test Per Kit	AMP 1000*	BAR 300	BZO 300*	BUP 10	COC 300*	COT 100	FTY 20	KET 1000	MDMA 500	THC 50*	MTD 300	EDDP 300*	MET 1000*	MOP 300	OPI 2000	OXY 100	PCP 25	PPX 300	TRA 100	TCA 1000
SureStep™ Multi-Drug One Step 2 Drug Screen Test Panel	See below **	25	2 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 3 Drug Screen Test Panel	**	25	3 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 4 Drug Screen Test Panel	**	25	4 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 5 Drug Screen Test Panel	**	25	5 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 6 Drug Screen Test Panel	**	25	6 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 7 Drug Screen Test Panel	**	25	7 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 8 Drug Screen Test Panel	**	25	8 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 9 Drug Screen Test Panel	**	25	9 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 10 Drug Screen Test Panel	**	25	10 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 11 Drug Screen Test Panel	**	25	11 drugs among 20 to select - All combinations are allowed																			
SureStep™ Multi-Drug One Step 12 Drug Screen Test Panel	**	25	12 drugs among 20 to select - All combinations are allowed																			

* Alternative cut offs (ng/mL): AMP 500, 300 - BZO 200 - COC 150 - MET 500, 300 - THC 150, 20 - EDDP 100

** Please order the Multi-Drug Test Panel by using the Product Code followed by the drugs combination

*** SG: Specific Gravity - Oxi: Oxidant - Cre: Creatinine - Glut: Glutaraldehyde - Nit: Nitrite

Multi-Drug Test With Integrated Cup II (A) With Or Without Specimen Validity Test (Adulteration): Drug Test & Cup II System With Adulteration Strips



Product Name	Code Extension	Test Per Kit	pH SG**	Cre**	SG**	AMP 1000*	BAR 300	BZO 300*	BUP 10	COC 300*	MDMA 500	THC 50	MTD 300	EDDP 300*	MET 1000*	MOP 300	TCA 1000
SureStep™ Multi-Drug One Step Multi-Line 5 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II with Adulteration Strips	See Below**	25	■	■		5 drugs among 20 to select - All combinations are allowed											
SureStep™ Multi-Drug One Step Multi-Line 6 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II with Adulteration Strips	651-a	25	■	■	■					■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 7 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II with Adulteration Strips	131	25	■	■	■	AMP 500				■	■	■	■		■	AMP 500	■
SureStep™ Multi-Drug One Step Multi-Line 8 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II with Adulteration Strips	121	25	■	■	■	AMP 500			■	■	■	■	■		■	AMP 500	■
SureStep™ Multi-Drug One Step Multi-Line 9 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II with Adulteration Strips	****	25	■	■	■			■	■	■		■	■	■	■	■	■
SureStep™ Multi-Drug One Step Multi-Line 10 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II with Adulteration Strips	042-a***	25	■	■	■	■	■			■	■	■	■		■	■	■

To order Surestep products please contact your local customer service representative

Multi-Drug Test With Cup II System



Product Name	Code Extension	Test Per Kit	AMP 1000*	BAR 300	BZO 300*	BU 10	COC 300*	MDMA 500	THC 50	MTD 300	EDDP 300*	MET 1000*	MOP 300	TCA 1000
SureStep™ Multi-Drug One Step Multi-Line 2 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	See Below**	25												
SureStep™ Multi-Drug One Step Multi-Line 3 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	**	25												
SureStep™ Multi-Drug One Step Multi-Line 4 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	**	25												
SureStep™ Multi-Drug One Step Multi-Line 5 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	**	25												
SureStep™ Multi-Drug One Step Multi-Line 6 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	**	25												
SureStep™ Multi-Drug One Step Multi-Line 7 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	021	25	■		■		■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 7 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	121	25	■ AMP 500				■	■	■	■		■ AMP 500	■	
SureStep™ Multi-Drug One Step Multi-Line 7 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	171	25	■		■	■	■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 8 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	101-a***	25	■		■	■	■		■	■		■	■	
SureStep™ Multi-Drug One Step Multi-Line 8 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	123	25	■ AMP 500				■	■	■	■		■ AMP 500	■	
SureStep™ Multi-Drug One Step Multi-Line 9 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	151-a	25	■		■	■	■		■	■	■	■	■	
SureStep™ Multi-Drug One Step Multi-Line 10 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	041	25	■	■	■		■	■	■	■		■	■	■
SureStep™ Multi-Drug One Step Multi-Line 10 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	191	25	■		■	■	■	■	■	■		■	■	■
SureStep™ Multi-Drug One Step Multi-Line 11 Drug Screen Test Panel with Integrated E-Z Split Key™ Cup II	051-a	25	■	■	■	■	■	■	■	■		■	■	■

Other drug combinations on demand (also including combination with other drugs as COT, FTY, KET, OXY, PCP, PPX, TRA)

* Alternative cut offs (ng/mL): AMP 500, 300 - BZO 200 - COC 150 - MET 500, 300 - THC 150, 20 - EDDP 100

** Please order the Multi-Drug Test Panel by using the Product Code followed by the drugs combination

*** Combination validated but needs SureStep silk-printing / label development

**** Combination which needs to be validated - requires annual forecast

(A) Multi-Drug Test with Cup 009 is still available

Multi-Drug Oral Fluid Test Multi-Drug Multi-Line Twist Test Device



Product Name	Code Extension	Test Per Kit	AMP 50	BZO 10	COC 30	THC 100	MTD 30	MET 50	MOP 40
SureStep™ Oral fluid Multi-Drug Multi-Line 6 Twist Screen Test Device	021	25	■		■	■	■	■	■
SureStep™ Oral fluid Multi-Drug Multi-Line 6 Twist Screen Test Device	031	25	■	■	■	■		■	■
SureStep™ Oral fluid Multi-Drug Multi-Line 6 Twist Screen Test Device	041	25	■	■	■		■	■	■

Gastrointestinal Disorders - Rapid Tests Lateral Flow Immunochromatography

Clostridium difficile Complete

For the simultaneous qualitative detection of *Clostridium difficile* antigen (Glutamate Dehydrogenase) and Toxins A & B in faecal specimens.

Product	Product Code	Size
<i>C. DIFF QUIK CHEK COMPLETE</i> ®	T30525C	25
	T30550C	50



Clostridium difficile

A rapid enzyme immunoassay for the qualitative detection of *Clostridium difficile* antigen (Glutamate Dehydrogenase) in faecal specimens.

Product	Product Code	Size
<i>C. DIFF QUIK CHEK</i> ™	30390	25



Clostridium difficile toxins A & B

For the qualitative detection of *Clostridium difficile* toxins A and B in faecal specimens.

Product	Product Code	Size
<i>TOX A/B QUIK CHEK</i> ®	T30394	25

For the simultaneous qualitative detection and differentiation of Shiga toxin 1 and Shiga toxin 2 in human faecal specimens and cultures derived from faecal specimens.

Product	Product Code	Size
<i>SHIGA TOXIN QUIK CHEK</i>	T30625	25




For the rapid simultaneous qualitative detection of shiga toxins 1 & 2 in human faecal specimens and broth cultures.

Product	Product Code	Size
<i>SHIGA TOXIN CHEK</i>	T30696	96

Lactoferrin

For the qualitative detection of elevated levels of faecal lactoferrin, a marker for faecal leukocytes and an indicator of intestinal inflammation.

Product	Product Code	Size
<i>IBD-EZ VUE</i> ®	30353	25
<i>LEUKO EZ VUE</i> ®	30355	25



Gastrointestinal Disorders - ELISA

Clostridium difficile

For the qualitative detection of *Clostridium difficile* antigen (Glutamate Dehydrogenase) in faecal specimens.

Product	Product Code	Size
<i>C. DIFFICILE CHEK™</i> - 60	30392	96 well plate

Clostridium difficile toxins A & B

For the qualitative detection of *Clostridium difficile* toxins A and B in faecal specimens.

Product	Product Code	Size
<i>C. DIFFICILE TOX A/B II™</i>	30397	96 well plate

Clostridium perfringens

For the qualitative detection of *Clostridium perfringens* enterotoxin in faecal specimens.

Product	Product Code	Size
<i>C. PERFRINGENS ENTEROTOXIN TEST</i>	RUO	96 well plate

Cryptosporidium

For the qualitative detection of *Cryptosporidium* oocyst antigen in faecal specimens.

Product	Product Code	Size
<i>CRYPTOSPORIDIUM II</i>	30406	96 well plate

Entamoeba Histolytica

For the qualitative rapid detection of *Entamoeba Histolytica* adhesin in human faecal specimens.

Indicates the presence of the organism in faeces as a diagnostic aid.

Product	Product Code	Size
<i>E. Histolytica II</i>	30404	96 well plate

Giardia lamblia

For the qualitative detection of *Giardia lamblia* cyst antigen in faecal specimens.

Product	Product Code	Size
<i>GIARDIA II</i>	30405	96 well plate

Lactoferrin

For the qualitative detection of elevated levels of lactoferrin as a marker of intestinal inflammation in faecal specimens.

Product	Product Code	Size
<i>IBD-CHEK®</i>	30352	96 well plate

For the quantitative detection of elevated levels of lactoferrin as a marker of intestinal inflammation in faecal specimens.

Product	Product Code	Size
<i>IBD-SCAN®</i>	30351	96 well plate

Disease/Condition Index

ACR	Afinion	28
Acute Coronary Syndromes	Alere Triage CardioProfilER	23
Acute Kidney Injury	Alere Triage NGAL	26
Barmah Forest virus	Panbio	41
Bladder cancer	Alere BladderChek	11
Blood gases	epoc	36
Bordetella Pertussis	Panbio	41
CD4	Alere PIMA CD4	20
Chagas	ImmunoComb	37
Chlamydia	Clearview	33
	ImmunoComb	37
Cholesterol/Lipids	Alere Cholestech LDX	12
Cholesterol/Lipids/Glucose	Alere Cholestech LDX	12
Clostridium Difficile	Techlab	52
Colon cancer	Clearview	34
Coxiella burnetti - Q Fever	Panbio	41
CRP	Afinion	28
Cytomegalovirus	Techlab	53
Dengue	ImmunoComb	37
	Panbio	41
Diabetes	Afinion	28
	NycoCard	30
Drugs of Abuse	DxLink	16
	SureStep	43
	Alere Triage	27
Entamoebae Histolytica	Techlab	51
Epstein-Barr virus	Panbio	41
Filariasis	Alere BinaxNOW	9
G6PD	Alere BinaxNOW	10
Giardiasis	Techlab	51
Giardiasis/Cryptosporosis	Techlab	51

Disease/Condition Index

HbA1c	Afinion	28
	NycoCard	31
Heart Failure	Alere Heart Check	17
	Alere Triage	22
Helicobacter pylori	Clearview	32
	ImmunoComb	37
Hepatitis	Alere Determine	14
	ImmunoComb	37
Herpes Simplex Virus	Wampole	41
HTLV	ImmunoComb	37
Human Herpesvirus-6	Panbio	41
Human Immunodeficiency Virus (HIV) 1/2	Alere Determine	14
	Alere Determine Combo	14
	DoubleCheckGold	36
	ImmunoComb	37
Infectious Mononucleosis (IM)	Clearview	32
Influenza	Alere BinaxNOW	8
	Clearview	32
INR	Alere INRatio	18
Irritable Bowel Disorder (IBD)	Techlab	51
Japanese Encephalitis	Panbio	41
Lactoferrin	Techlab	51
Legionella	Alere BinaxNOW	8
Leptospira	Panbio	41
Liver enzymes	Alere Cholestech LDX	12
Luteinising Hormone	Clearview	35
Malaria	Alere BinaxNOW	9
	Clearview	34
Myocardial Infarction & Heart Failure	Alere Triage CardioProfiler	23
	Alere Triage Cardiac Panels	22

Osteoporosis	Osteomark NTx	40
Pneumonia	Alere BinaxNOW	8
Pregnancy	Clearview	35
	Alere TestPack	21
PBP2a	Alere	6
Premature Birth	Careplan Vp18H Test Gloves	35
	Actim Partus	7
	Actim PROM	7
Pulmonary Embolism / Deep Vein Thrombosis	Clearview	32
	Alere Triage D-Dimer	24
	NycoCard	30
Respiratory Syncytial Virus (RSV)	Alere BinaxNOW	8
	Clearview	33
Ross River virus	Panbio	41
Shortness of breath	Alere Triage Profiler SOB	24
Staphylococcus aureus	Alere BinaxNOW	10
Streptococcus A (Strep A)	BinaxNOW	31
	Clearview	33
	Alere TestPack	21
Syphilis	Alere Determine	15
Toxoplasma gondii	ImmunoComb	37
Tuberculosis	Clearview	15
Varicella-zoster	Panbio	41
West Nile virus	Panbio	41



Alere Worldwide Locations

North America

Corporate Headquarters

Alere Inc.
51 Sawyer Road, Suite 200
Waltham, MA 02453, USA
Tel: (781) 647 3900

Alere North America, Inc.
30 S. Keller Rd
Orlando, FL 32810, USA
Tel: (877) 441 7440

Alere San Diego, Inc.
9975 Summers Ridge Road
San Diego, CA 92121, USA
Tel: (858) 805 4808

Alere North America, Inc.
2 Research Way
Princeton, NJ 08540, USA
Tel: (609) 627 8000

Canada

Alere Inc.
1-57 Iber Road, Ottawa,
Ontario, K2S 1E7, Canada
Tel: 1-800-818-8335

Latin America & Africa

Argentina

Biolinker S.A.
14 de julio 618, Buenos Aires
Tel: (54 11) 4554-4007
Fax: (54 11) 4553-2141
www.biolinker.com.ar

Brazil

Alere S.A.
Rua does Pinheiros 498
7th Floor 05422-000 - São Paulo
SP Brazil
South America
Tel: +55 11 2131 5100
Fax: +55 11 2131 5181
www.alerebrasil.com.br

Bioeasy Diagnóstica Ltda.
Avenida Barão Homem de Melo, 567 - Nova Granada
30431-285 - Belo Horizonte - MG
Tel: +55 31 3048-0008
http://www.bioeasy.com.br

Prodímol Biotecnologia S.A.
Praça Carlos Chagas, 49 - 3º andar - Santo Agostinho
30170-020 - Belo Horizonte - MG
Tel: +55 31 2122 2900
http://www.prodímol.com.br

Colombia

Biosystems S.A.
Av. 9 Norte No. 17A -24
Cali, Colombia
PBX: +57 (2) 661 8797-661 8916
Fax: +57 (2) 661 8917
e-mail: info@biosystems.com.co

South Africa

Alere Healthcare (Pty) Ltd
11B Riley Road,
Eastwood Office Park,
Bedfordview 2007, South Africa
Tel: +27 11 4504411
Fax: +27 11 4502267
E-mail: yeswecan@alere.com
www.alerehealthcare.co.za

Asia Pacific

Australia

Alere
532 Seventeen Mile Rocks Rd,
Sinnamon Park, QLD, 4073, Australia
Tel: +61 7 3363 7100
1800 622 642 (Freecall within Australia)
Fax: +61 7 3363 7199
Email: au.enquiries@alere.com

China

Alere Shanghai
Suite 4306, Wheelock Square, No. 1717,
West Nanjing Road,
Shanghai, 200040,
China
Tel: +86-21-5131-5356
Fax: +86-21-3383-0280

India

Alere Medical Pvt. Ltd.
No.113, BPTP Park CentraNH-8,
Opposite 32nd Milestone
Gurgaon - 122001
(Haryana) India
Tel: +91 124 456 9000
Fax: +91 124 456 9011/022

Indonesia

PT. Alere Health
MT. Haryono Square No. OF 03/06 Jl. Otto Iskandarinata
No. 390
Jakarta Timur - 13330
Indonesia
Tel: +62.21.8590 6006 (ext. 106)
Fax: +62.21.8590 9321

Japan

Alere Medical Co Limited
18F STEC Jyoho Building,
1-24-1 Nishi-Shinjuku,
Shinjuku-ku, Tokyo160-0023, Japan
Tel: +81 (0)3- 5326- 7747
www.alere.co.jp

Malaysia

Alere Health Sdn Bhd (986219-V)
B-11-01, Level 11, Block B, PJ8, Jalan Barat,
46050 Petaling Jaya, Selangor, Malaysia
Tel: +603 7960 1500
Fax: +603 7932 2500

Singapore

Alere Pte Ltd
1 Maritime Square, Harbour Front Centre #09-64
Singapore 099253
Tel: +65 6643 5002

South Korea

Alere Healthcare Inc.
6F, Yumin Technotown, 437-4 Bora-dong,
Giheung-gu, Yongin-si, Gyeonggi-do, Korea
Tel: +82-31-899-2957
Fax: +82-32-899-2840

Taiwan

Alere Health Corp
3F, No 64, Tacheng Street
Taipei City
Taiwan
Tel: +886-2-2552-2605
Email: tw.info@alere.com

New Zealand

Alere Limited
PO Box 105653, Auckland City,
Auckland 1143, New Zealand
Tel: 0800 466 951
Fax: 0800 466 958
Email: nz.enquiries@alere.com

International

Ireland

Alere International Limited
Parkmore East Business Park
Ballybrit, Galway, Ireland
Tel: +353 (0) 91 429 900
e-mail: professionaldiagnostics@alere.com

Europe & Middle East

Austria

Alere GmbH
Estermannstrasse 17, A 4020
Linz, Austria
Tel: +43 (0)732 77 10 77 - 0

Belgium

Alere Health bvba/sprl
Derbystraat 359, 9051 Gent,
St. Denijs Westrem, Belgium
Tel: +32 9 242 88 44
Fax: +32 9 242 88 45
Email: belux@alere.com

France

Alere France
21, rue Albert Calmette-Bâtiment B4
78 350 Jouy-en-Josas, France
Tel: 01 39 46 83 18
Fax: 01 39 46 64 20

Germany

Alere GmbH
Am Wassermann 28,
50829 Köln, Germany
Tel: +49 221 271 43-0
+49 221 569 730 (Gabmed GmbH)
Fax: +49 221 271 43-400

Ireland

Alere Health Limited
c/o Barr Pomeroy,
21 Herbert Place,
Dublin 2, Ireland
Tel: 048 90 469992
Fax: 048 90 462154

Italy

Alere srl
Via Drizzagno, 11
30037 Scorzè (VE), Italy
Tel: +39 041 5841546
Fax: +39 041 5841872

Via Eustachi, 36

20129 Milano, Italy
Tel: +39 02 2774131
Fax: +39 02 27741323

Israel

Orgenics Ltd.
PO Box 360, 70650 Yavne, Israel
Tel: +972 (8) 9429201
Fax: +972 (8) 9438758
Email: info@orgenics.com

Netherlands

Alere Health B.V.
Charles Stulemeijerweg 15
5026 RS TILBURG, The Netherlands
Tel: (+31) (0)13 467 27 00
Fax: (+31) (0)13 467 80 16

Spain

Alere Healthcare, S.A.U.
Botanica 146, 08908
L'Hospitalet de Llobregat,
Barcelona, Spain
Tel: +34 936008000
Fax: +34 936008005

Switzerland

Alere GmbH
Grossmattstrasse 9
CH-8902 Urdorf, Switzerland
Tel: +41 (0) 41 761 72 72
Fax: +41 (0) 41 761 72 73

United Kingdom

Alere Limited
Pepper Road,
Hazel Grove, Stockport,
SK7 5BW, United Kingdom
Tel: +44 (0)161 483 5884
Fax: +44 (0)161 483 5778