

# Test Tubes and Closures

## Plug-Type Closures/Screw Caps/Centrifuge Tubes - Conical

### Polyethylene Tube Caps

STOCKWELL - Use polyethylene thumb caps for resealing evacuated glass and glass culture tubes. Choose from a variety of colors.



Cat. No.	Mfr. No.	Description	Qty.
<b>For 12mm glass culture tubes and 13mm glass evacuated tubes</b>			
8569B	8569B	Blue	1,000/bg
8569G	8569G	Green	1,000/bg
8569L	8569L	Lavender	1,000/bg
8569R	8569R	Red	1,000/bg
8569	8569	White	1,000/bg
8569Y	8569Y	Yellow	1,000/bg

### For 13mm glass culture tubes

8573B	8573B	Blue	1,000/bg
8573SI	8573SI	Gray	1,000/bg
8573G	8573G	Green	1,000/bg
8573L	8573L	Lavender	1,000/bg
8573R	8573R	Red	1,000/bg
8573	8573	White	1,000/bg

### For 16mm glass tubes

8579B	8579B	Blue	1,000/bg
8579G	8579G	Green	1,000/bg
8579R	8579R	Red	1,000/bg
8579	8579	White	1,000/bg

### Corning® Disposable Screw Caps

CORNING - These caps are made of a phenolic compound. Each cap includes a white rubber liner.



Cat. No.	Mfr. No.	Description	Qty.
T1355-750	99999-13	Screw cap with liner 13-415	1,000/cs
T1355-800	99999-15	Screw cap with liner 15-415	1,000/cs

### Corning® Reusable Phenolic Screw Caps

CORNING - These caps are fabricated from a resin specially formulated to resist the effects of temperature and steam autoclaving.



Cat. No.	Mfr. No.	Description	Qty.
<b>With Teflon®-faced rubber liner</b>			
T1358-1	9998-13	Screw cap with liner 13-415	288/cs
T1358-2	9998-15	Screw cap with liner 15-415	288/cs

### KIM-KAP™ Polypropylene Closures

KIMBLE CHASE - Use these closures with a variety of disposable culture tubes to keep contents sterile without the need for cotton plugging. They are autoclavable and can be removed and replaced without contaminating tube contents.



Cat. No.	Mfr. No.	Description	Qty.
T1291-13A	73660 13	13mm O.D., natural	1,000/cs

### Nonsterile Screw Caps for Sample Tubes

SIMPORT PLASTICS - These nonsterile screw caps are made of polypropylene. Caps feature a long skirt and fast thread design, enabling you to easily remove or seal the caps with a single turn. Caps do not have silicone washer.



Cat. No.	Mfr. No.	Description	Qty.
T502N	T502N	Natural color	1,000/cs



### Sterile Disposable Centrifuge Tubes

CARDINAL HEALTH - Our disposable centrifuge tubes are designed to fit standard centrifuge apparatuses for general laboratory use. Conical tubes feature plug seal cap with leak-resistant seal.

Translucent polypropylene tubes have black printed graduations and greater strength for high-speed centrifugation. Available in poly-foam rack.

Cat. No.	Description	Qty.
SP3920-80	15mL, plug cap w/rack	50/pk, 10 pk/cs
SP3920-15A	15mL, plug seal bulk	50/pk, 10 pk/cs
SP3920-15B	15mL, plug seal w/rack	50/pk, 10 pk/cs
SP3920-50A	50mL, plug cap bulk	25/pk, 20 pk/cs
SP3920-50	50mL, plug cap w/rack	25/pk, 20 pk/cs
SP3921-50A	50mL, top cap w/rack	25/pk, 20 pk/cs

### Corning® Self-Standing Centrifuge Tubes

CORNING - These disposable sterile polypropylene tubes are ideal for sample storage and for working in laminar flow hoods. The tubes have printed graduations and can withstand a maximum RCF of 3,000 x g. Tubes feature conical bottom with attached skirt, so they can stand independently. Tubes are available with plug seal or flat top cap.



Cat. No.	Mfr. No.	Description	Qty.
C3974-10	430897	Plug seal	500/cs
C3974-11	430921	Flat top	500/cs