

417 Power Podiatry
Treatment Chair
Designer Series

MIDMARK®



Power to th

Midmark presents the 417 Power Podiatry Chair—made just for podiatrists.

Before we designed the 417 we listened to what podiatrists like you wanted. Then we built the Midmark 417 from the ground up. You'll find it's "just what the podiatrist ordered." It'll put your patients' feet exactly where you want them.

Convenient low height

Your patients will appreciate how easy it is to get on and off the 417, with its low minimum height of 22 1/2". They'll like the sculptured top, built-in pillow and smooth, rounded chair arms that move as the back section moves. Makes your patients feel secure.

The 417's seat section goes as high as 41"—easy for you to work standing or seated at whatever height is most comfortable for you.



It's beautiful

The first things you'll notice about the 417 are its good looks, contemporary styling and clean lines. And it's easy to *keep* your 417 looking good. Dirt and stain-resistant vinyl back, seat and foot sections are fastened with Velcro, so they're easy to remove. Makes care, cleaning and changing colors a breeze! And you can choose from thirteen upholstery colors to coordinate with Midmark Modular Casework.

User friendly, simultaneous movement

As easy to program as a car radio, the 417 is the fastest power podiatric treatment chair available. All you do is touch a button on the foot control or optional Hand-held Commander to achieve one of your pre-programmed positions—height, tilt and back move *simultaneously*.

The standard foot control helps you maintain the sterile field while positioning patients during surgical procedures. Use it yourself... or have your assistant do it.

The auto return function brings the 417 back to the chair position. It'll save you precious seconds on each procedure.

Foot control

A touch of your toe is all it takes! You control height, tilt and back with three power positioning choices and two programmed positions. Place the lightweight foot control anywhere around the table—it won't slip, slide or scuff your floors. The plastic molding is durable and easy to clean. The control's standard telephone cord connection plugs into the front of the chair's base or left side hand control outlet.



Hand-held Commander

The Hand-held Commander is an accessory and offers you three power positioning choices—height, tilt and back—and four programmed positions, two more than with the foot control. But it's just as easy to use.

When not using the Commander, just place it in its holder mounted on either side of the foot section.



e Podiatrist!



Debris tray

It's conveniently located underneath the foot section. Just pull it out when you want to use it. It tilts 180° so you can dispose of waste into your waste receptacle without having to remove the tray. The tray tilts at various angles, depending on the tilt of the chair and can extend 8 3/4" beyond the foot section so it's not in your way. All done? Slide tray underneath foot section. Rounded corners make cleaning easy. Just remove by rotating tray knobs.



Electric receptacles, on/off switch

Two duplex hospital grade electrical receptacles—one on each side of the chair—let you use electrical instruments and lights without tangling or tripping over cords.

Flip the 417's on/off switch to "off" during delicate procedures to prevent any accidental movement. Power is still supplied to the receptacles.

Built to last

The 417 will help increase your income for years to come. Heavy-duty construction and state-of-the-art electronics ensure a long-lasting relationship between you and your 417. The durable electro-mechanical drive offers smooth, quiet and reliable movement throughout positioning.

Foot extension

Do a lot of orthotics? The 417's adjustable foot extension makes orthotic procedures easier by positioning patients' feet exactly where you need them. It extends 8 7/8" and comfortably accommodates patients of all sizes.

Prefer to work standing? With the 417 at maximum tilt, the foot extension can be as high as 58". Locks automatically to avoid accidental movement.



417 Power Podiatry Treatment Chair Designer Series

Specifications

Minimum height: 22.5" (57.2 cm)
Maximum height: 41" (106.7 cm)
Length w/back reclined: 70" (177.8 cm)
Length w/foot section fully extended:
78.69" (199.9 cm)
Foot section: 21.25" (53.9 cm) wide tapering
to 17.75" (45.1 cm) x 14" long (35.6 cm)
Back section width: 22.75" (57.8 cm)
tapering to 17.5" (44.4 cm)
Seat section width: 22.75" (57.8 cm)
tapering to 21.5" (54.6 cm)
Minimum height at foot section w/o tilt:
19.5" (49.5 cm)
Maximum height at foot section w/o tilt:
38.5" (97.8 cm)
Maximum height at foot section
w/maximum tilt: 38.5" (97.8 cm)
Maximum height at foot section
w/maximum tilt and maximum lift:
55.25" (140.3 cm)
Maximum height at foot section
w/maximum tilt, maximum lift and foot
section fully extended: 58" (147.3 cm)

Total width including arms: 28.12" (71.4 cm)
Arm rest width: 3.12" (7.9 cm)
Debris tray: Polyethylene, extends 8.75"
beyond foot section, 8.75" W x
12.5" L x 1.75" D (22.2 x 31.8 x 4.4 cm)



Electrical

115 VAC, 60 Hz, 12 AMPS
Two duplex hospital grade electrical
receptacles. Rated 115 VAC 60 Hz,
6 AMPS maximum.

Options

Stainless steel debris tray, 9.5" W x 11.5" L
x 2.5" D (24.1 x 29.2 x 6.4 cm)

Accessories

Hand-held Commander control, vision block
screen w/clamp, instrument tray, utility tray,
foot section cover

Colors

Molded housings and painted steel surfaces
are Pearl Grey. Base is black. Upholstery
colors: Blueberry, Silver Sage, Tea Green,
Fossil Grey, Cashmere Blue, Terra Cotta,
Black, Belagio, Siena, Tuscany, Firenze,
Capri and Treviso.

Pearl Grey modular units available with
door and drawer pulls to match 417
upholstery colors.

Midmark is an ISO 9001 Certified Company.

For more information
or a demonstration, contact
your Midmark dealer or
Midmark Corporation,
Versailles, Ohio.
Call: 1-800-MIDMARK
Fax: 1-800-365-8631
midmark.com



Because we care.

© Midmark Corporation 1994
Products subject to improvement changes without notice.
Litho in U.S.A. 003-0493-00 Rev. A1 (03/07)