

# Health Design



POST OP SHOES • ABDOMINAL AND BACK SUPPORT • KNEE BRACES • PROTECTION GEAR • WALKERS • WRIST BRACES





## Health Design

A division of DARCO International

Health Design is committed to providing healthcare professionals with high-quality, affordable medical products.

While maintaining an inventory of orthopedic soft goods and footwear products, we also provide competitive OEM and private label services.

Please visit our website at [www.healthdesign.net](http://www.healthdesign.net) for additional information.

<b><i>PRODUCTS</i></b>	<b><i>PAGE</i></b>
<i>ABDOMINAL AND RIB.....</i>	<i>.4</i>
<i>CLAVICLE AND BACK SUPPORT.....</i>	<i>.5-6</i>
<i>WRIST BRACES.....</i>	<i>.7-9</i>
<i>ARM SLINGS.....</i>	<i>.10-12</i>
<i>SHOULDER IMMOBILIZER.....</i>	<i>.13-14</i>
<i>KNEE BRACES.....</i>	<i>.15-16</i>
<i>CERVICAL SUPPORTS.....</i>	<i>.17</i>
<i>TRACTION AND PROTECTION GEAR.....</i>	<i>.18-19</i>
<i>CAST SHOES.....</i>	<i>.20-21</i>
<i>POST OP SHOES.....</i>	<i>.22-24</i>
<i>DIMENSION GUIDE.....</i>	<i>.25-26</i>

## **ABDOMINAL BINDER**

The Abdominal Binder provides circular compression and support to the abdomen. It may also be used to hold postoperative dressings in place. The binder is constructed of cotton lined elastic with a tri laminate abdominal pad. This pad receives the 1.5" hook allowing for ease of application. This product is latex free. Available in white.

### **9" ABDOMINAL BINDER**

<i>SIZE</i>	<i>HIP CIRCUM.</i>	<i>CAT. NO.</i>
small/med	30"- 45"	202-0050-35
med/large	46"- 62"	202-0050-45
X-large	62"+	202-0050-55

### **12" ABDOMINAL BINDER**

<i>SIZE</i>	<i>HIP CIRCUM.</i>	<i>CAT. NO.</i>
small/med	30"- 45"	202-0051-35
med/large	46"- 62"	202-0051-45
X-large	62"+	202-0051-55



ABDOMINAL BINDER

## 9" DUAL ADJUSTMENT BACK SUPPORT

The Dual Adjustment Back Support provides maximum support to both the abdomen and lumbosacral areas. Dual adjustment panels of elastic increase or decrease support and tension to the desired level. Two posterior stays provide increased support and stability against the lower back. Two lateral anti-roll stays maintain the shape of the belt and prevent it from rolling up. Available in white.

SIZE	HIP CIRCUM.	CAT. NO.
small	26"- 40"	232-0010-30
med	36"- 49"	232-0010-40
large	44"- 59"	232-0010-50
X-large	52"- 66"	232-0010-60



9" DUAL ADJUSTMENT BACK SUPPORT

## 9" ELASTIC CRISS CROSS BACK SUPPORT

Designed for the treatment of lower back pain from muscle tension, spasm or overuse, the Criss-Cross Brace features a solid elastic panel, criss-crossing elastic straps, and lateral supports. The brace uses these features to provide intra-abdominal circular compression and support. Criss-cross elastic straps provide additional support to the lumbosacral area. Available in white.

SIZE	HIP CIRCUM.	CAT. NO.
X-small	25"- 30"	232-0050-20
small	30"- 33"	232-0050-30
med	36"- 38"	232-0050-40
large	37"- 41"	232-0050-50
X-large	41"- 44"	232-0050-60
XX-large	44"- 48"	232-0050-70



9" ELASTIC CRISS CROSS BACK SUPPORT

## CLAVICLE HEAVY DUTY SUPPORT WITH BUCKLE

The Heavy-Duty Clavicle is comfortably padded throughout with medium density foam. It features extra-wide underarm pads and fully adjustable dual safety-lock tooth buckle adjustment closures. Available in beige.

SIZE	HIP CIRCUM.	CAT. NO.
Pediatric	13"- 20"	432-0030-20
small	16"- 23"	432-0030-30
med	22"- 30"	432-0030-40
large	28"- 39"	432-0030-50
X-large	37"- 45"	432-0030-60



CLAVICLE HEAVY DUTY SUPPORT

**Clavicle strap sizing:** Measure from center of the back around the arm to the same point on the back.

## **INDUSTRIAL UTILITY BELT**

This heavy-duty belt is designed to provide maximum support of the abdomen and lumbosacral area during repetitive lifting tasks. The belt is constructed of sturdy, 8" black, knit elastic. A panel consisting of heavy, black, nylon woven loop is bonded to foam padding & tricot. Additional features include adjustable suspenders and a full-front pressure sensitive closure. Lateral elastic double-pull controls allow the patient to adjust additional pressure against the lower back as needed. Available in black.

<i>SIZE</i>	<i>HIP CIRCUM.</i>	<i>CAT. NO.</i>
small	34"- 38"	231-0080-30
med	38"- 42"	231-0080-40
large	42"- 46"	231-0080-50
X-large	46"- 50"	231-0080-60
XX-large	50"- 54"	231-0080-70



INDUSTRIAL UTILITY BELT

## WRIST SPLINT DELUXE

This splint's three-layer bonded construction eases maintenance and assures durability. The interior is fully padded with foam and lined with tricot for comfort. The exterior is covered with a sturdy nylon oxford. A dorsal aluminum support assures rigidity while a contoured aluminum palmar support maintains the palmar arch and permits use of the hand. A circumferential Velcro® strap firmly secures the splint to the wrist. Available in both sized and universal versions. Sized versions are one-piece construction. Universal splints feature a unique two-piece construction which allow them to be worn on either the left or right hand. Available in black.

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
Pediatric	2.0" - 2.5"	452-1010-20L	20R
small	2.5" - 3.0"	452-1010-30L	30R
med	3.0" - 3.5"	452-1010-40L	40R
large	3.5" - 4.0"	452-1010-50L	50R
X-large	4.0" - 4.5"	452-1010-60L	60R
<b>SIZE CAT. NO.</b>			
Universal		452-1010-10	



WRIST SPLINT DELUXE

## BLACK FOAM WRIST SPLINT

This premium wrist splint features a tri laminate material with a white cotton lining, and a brushed nylon exterior. A dorsal aluminum support and a split contoured aluminum palmar stay help immobilize and support the girth changes in the forearm. One adjustable Velcro® strap, and one loop lock strap firmly secure the splint to the wrist. A foam strip provides increased comfort over the palmar stay. Available in black.

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
X-small	2.0" - 2.5"	451-4010-20L	20R
small	2.5" - 3.0"	451-4010-30L	30R
med	3.0" - 3.5"	451-4010-40L	40R
large	3.5" - 4.0"	451-4010-50L	50R
X-large	4.0" - 4.5"	451-4010-60L	60R



BLACK FOAM WRIST SPLINT

## WRIST AND FOREARM SPLINT PLUSH BEIGE FOAM

This Plush Beige Foam Wrist and Forearm Splint is the same as the Plush Beige Foam Wrist Splint, except it provides an additional 2" of length for more immobilization and support of the wrist and lower arm. A third circumferential Velcro® strap ensures a proper. Please specify left or right when ordering.

SIZE	CAT. NO. (left)	(right)
Universal	452-1120-10L	10R



WRIST AND FOREARM SPLINT PLUSH FOAM

## WRIST BRACES

### PERFORATED VINYL WRIST BRACE

This heavy-duty wrist support is constructed of durable perforated vinyl, lined with a soft cotton backing. A dorsal aluminum support and contoured aluminum palmar support will maintain rigidity and immobilize the wrist. One adjustable Velcro® strap, and one loop lock firmly secure the splint to the wrist. An inner elastic liner protects the arm from coming in contact with the Velcro® straps. Please specify left or right when ordering. Available in tan.



PERFORATED VINYL WRIST BRACE

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
X-small	2.0" - 2.5"	451-0351-20L	20R
small	2.5" - 3.0"	451-0351-30L	30R
med	3.0" - 3.5"	451-0351-40L	40R
large	3.5" - 4.0"	451-0351-50L	50R
X-large	4.0" - 4.5"	451-0351-60L	60R

### 8" BILATERAL WRIST BRACE TAN

Designed for mild wrist support. This brace has a wide palmar stay that is reversible to accommodate a right or left hand. The brace has hook and loop closures on the dorsal side.



8" BILATERAL WRIST BRACE TAN

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
small	2.75" - 3.00"	252-1030-30	
med	3.00" - 3.50"	252-1030-40	
large	3.50" - 4.00"	252-1030-50	
X-large	4.00" - 4.25"	252-1030-60	

### BUNGEE WRIST SPLINT

The Bungee Wrist Splint has an easy to use bungee strap allowing the wearer to make quick adjustments without assistance. Designed to help provide compression. The hook & loop closure allows for easy adjustment. Removable aluminum palmar spoon with an aluminum dorsal stay. Black exterior with blue interior.



BUNGEE WRIST SPLINT

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
X-small	2.0" - 2.5"	9BW0L	9BW0R
small	2.5" - 3.0"	9BW1L	9BW1R
med	3.0" - 3.5"	9BW2L	9BW2R
large	3.5" - 4.0"	9BW3L	9BW3R
X-large	4.0" - 4.5"	9BW4L	9BW4R



## WRIST BRACES

### PLUSH WRIST SPLINT WITH THUMB SPICA

Designed to help provide compression and mild wrist support and immobilize the thumb. The interior is fully padded for comfort. The three hook & loop straps provide uniform support. Includes a removable aluminum palmar spoon with two aluminum dorsal stays with a sturdy aluminum stay to immobilize the thumb. Black exterior with blue interior.

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
X-small	2.0" - 2.5"	9PT0L	9PT0R
small	2.5" - 3.0"	9PT1L	9PT1R
med	3.0" - 3.5"	9PT2L	9PT2R
large	3.5" - 4.0"	9PT3L	9PT3R
X-large	4.0" - 4.5"	9PT4L	9PT4R



PLUSH WRIST SPLINT WITH THUMB SPICA

### PLUSH WRIST SPLINT

Designed to help provide compression and mild wrist support. The interior is fully padded for comfort. The three hook & loop straps provide uniform support. Includes a removable aluminum palmar spoon with two aluminum dorsal stays. Black exterior with blue interior.

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
X-small	2.0" - 2.5"	9PW0L	9PW0R
small	2.5" - 3.0"	9PW1L	9PW1R
med	3.0" - 3.5"	9PW2L	9PW2R
large	3.5" - 4.0"	9PW3L	9PW3R
X-large	4.0" - 4.5"	9PW4L	9PW4R



PLUSH WRIST SPLINT

### LACE UP WRIST SPLINT

Beige canvas wrist splint with a removable palmar spoon to immobilize the wrist. The cotton laces provide great compression.

SIZE	PALM WIDTH	CAT. NO. (left)	(right)
X-small	2.0" - 2.5"	E59185	E59180
small	2.5" - 3.0"	E59186	E59181
med	3.0" - 3.5"	E59187	E59182
large	3.5" - 4.0"	E59188	E59183
X-large	4.0" - 4.5"	E59189	E59184



LACE UP WRIST SPLINT

## CRADLE SLING DELUXE

Designed to help prevent elbow pressure, the Cradle Sling is constructed of one-piece, navy blue duck for maximum durability. The sling features an easy-release buckle adjustment strap along with a thumb loop to help eliminate hand drop and to maintain positioning in the sling.

SIZE	ENVELOPE SIZE	CAT. NO.
X-small	4.50" x 12"	422-0030-20
small	5.50" x 14"	422-0030-30
med	6.50" x 16"	422-0030-40
large	7.00" x 18"	422-0030-50
X-large	7.00" x 18"	422-0030-60*

\*The X-large has an elongated shoulder strap of 5"



CRADLE SLING DELUXE

## UNIVERSAL SLING DELUXE

Constructed of heavy-duty 10 oz. black duck for maximum durability, the Universal Sling features a 1.5" easy-release Velcro® closure shoulder strap and a slide buckle adjustment.

SIZE	ENVELOPE SIZE	CAT. NO.
Adult Univ.	7" - 17.5"	422-0040-10



UNIVERSAL SLING DELUXE

## CRADLE SLING ECONOMY

Constructed from one-piece heavy muslin material the Economy Cradle Sling features a slide-buckle shoulder strap adjustment. Available in Beige.

SIZE	ENVELOPE SIZE	CAT. NO.
X-small	4.50" x 12"	422-0025-20
small	5.50" x 14"	422-0025-30
med	6.50" x 16"	422-0025-40
large	7.00" x 18"	422-0025-50
X-large	7.00" x 18"	422-0025-60*

\*The X-large has an elongated shoulder strap of 5"



CRADLE SLING ECONOMY

## UNIVERSAL SLING ECONOMY

Constructed for lightweight 7 oz., navy blue duck, this Universal Sling features a one inch webbing shoulder strap and a slide buckle adjustment.

SIZE	ENVELOPE SIZE	CAT. NO.
Universal	7"x 17.5"	422-0041-10



UNIVERSAL SLING ECONOMY

## FOAM SLING AND SWATHE

Constructed of heavy, foam-lined, pressure sensitive material, the Foam Sling and Swathe features a durable 1.5" strap which evenly distributes weight across the patients shoulders, and thereby maximizes comfort. Other features include pressure sensitive closures on the swathe and shoulder straps for easy adjustment and application. Available in black.

SIZE	ENVELOPE SIZE	CAT. NO.
S/M	9.5" x 12"	422-0060-35
L/XL	9.5" x 12"	422-0060-55



FOAM SLING AND SWATHE

## TENNIS ELBOW WRAP

This Tennis Elbow Wrap is recommended for extensor muscle pain associated with tennis elbow. It features heavy-duty twin pressure sensitive loop for universal adjustment and a soft, brushed nylon interior for comfort. Available in beige.

SIZE	CAT. NO.
Universal	422-0110-10



TENNIS ELBOW WRAP

## **ELEVATED ICE BAG SLING**

Constructed of blue cotton material, the sling is a dual purpose sling, it provides two ice pockets that can be used for cold therapy. The sling may be used as a regular arm sling or as an elevated ice bag sleeve.

<i>SIZE</i>	<i>ENVELOPE SIZE</i>	<i>CAT. NO.</i>
small	5.50" - 14"	422-0045-30
large	7.00" - 18"	422-0045-50



ELEVATED ICE BAG SLING

## **SHOULDER IMMOBILIZER ECONOMY**

Constructed from one piece navy blue muslin. Features a slide buckle shoulder strap for easy adjustments, neck pad and thumb loop.

<i>SIZE</i>	<i>ENVELOPE SIZE</i>	<i>CAT. NO.</i>
large	7.00" - 18"	422-0072-50



SHOULDER IMMOBILIZER ECONOMY

# SHOULDER IMMOBILIZER

## SLING AND SWATHE DELUXE CLOSED ELBOW

Manufactured from black duck, The Sling and Swathe (Closed Elbow) features an easy release buckle closure for simple adjustment, dual thumb loops, and a neck comfort pad. The sling features an angled, hand-heart position for patient comfort which places minimal pressure on the elbow. The sling is designed for universal (left or right) fit.

SIZE	CHEST CIRCUM.	CAT. NO.
small/med	30"- 46"	422-0055-35
large/X-large	32"- 52"	422-0055-55



SLING AND SWATHE DELUXE CLOSED ELBOW

## UNIVERSAL SHOULDER IMMOBILIZER

Manufactured from plush 3/8" foam laminated to a soft, white, brushed nylon pressure sensitive material. The wrist and upper arm cuffs may be positioned as desired on either arm with Velcro® attachments. Universally sized to fit male and female patients.

SIZE	ENVELOPE SIZE	CAT. NO.
Universal	up to 60"	422-0075-10



UNIVERSAL SHOULDER IMMOBILIZER

## SHOULDER IMMOBILIZER ECONOMY

Manufactured from one-piece, heavy-duty black duck, the Economy Shoulder Immobilizer features a body strap with a slide buckle adjustment, an easy-release buckle shoulder strap, and a thumb loop.

SIZE	PALM WIDTH	CAT. NO. (left)
Pediatric	4.50"- 12"	422-0070-20
small	5.50"- 14"	422-0070-30
med	6.50"- 16"	422-0070-40
large	7.00"- 18"	422-0070-50
X-large	7.00"- 18"	422-0070-60*
XX-large	7.00"- 18"	422-0070-70**

\*The X-large has an elongated shoulder strap of 5"

\*\*The XX-large has an elongated shoulder strap of 10"



SHOULDER IMMOBILIZER ECONOMY

## SHOULDER IMMOBILIZER

### **SHOULDER IMMOBILIZER ELASTIC**

Constructed of 6" firm-ribbed elastic, the Elastic Shoulder Immobilizer features heavy, sewn-in padding on the humerus and wrist straps. Pressure sensitive closures and a support strap to prevent hand drop, this product will immobilize the shoulder by anchoring the elbow to the body. Available in white.

<i>SIZE</i>	<i>PALM WIDTH</i>	<i>CAT. NO. (left)</i>
Pediatric	18" - 22"	222-0065-20
small	27" - 31"	222-0065-30
med	31" - 35"	222-0065-40
large	35" - 39"	222-0065-50
X-large	39" - 43"	222-0065-60
XX-large	43" - 47"	222-0065-70
XXX-large	47" - 51"	222-0065-90



SHOULDER IMMOBILIZER ELASTIC

## **DELUXE THREE-PANEL KNEE IMMOBILIZER**

The Deluxe Three-Panel Knee Immobilizer features 3 aluminum stays located medially, laterally, and posteriorly to provide support and stability. Three-panel universal construction offers custom fitting regardless of size (greatly reducing costly inventory and shelf space requirements). Additional features include Velcro® loop-lock straps with removable adjustment tabs and a popliteal pad.

### **12" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0021-10

### **14" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0022-10

### **16" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0023-10

### **18" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0024-10

### **22" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0026-10

### **26" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0029-10

### **20" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0025-10

### **24" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0028-10



DELUXE THREE-PANEL KNEE IMMOBILIZER

## **DELUXE POST-OP KNEE IMMOBILIZER**

Designed for maximum postoperative support and stability, the Deluxe Post-Op Knee Immobilizer is constructed of high density foam and pliable duck with 3 aluminum stays located medially, laterally, and posteriorly. The Immobilizer is x-ray lucent. Pressure sensitive loop-lock closures and removable locking tabs assure proper fit. A bottom fitting foam strip eliminates gap areas.

### **16" DELUXE POST-OP KNEE IMMOBILIZER**

SIZE	CIRCUM.	CAT. NO.
small	14" - 16"	462-0012-30
med	16" - 18"	462-0012-40
large	18" - 20"	462-0012-50
X-large	20" - 22"	462-0012-60
XX-large	22" - 24"	462-0012-70

### **20" DELUXE POST-OP KNEE IMMOBILIZER**

SIZE	CIRCUM.	CAT. NO.
small	14" - 16"	462-0014-30
med	16" - 18"	462-0014-40
large	18" - 20"	462-0014-50
X-large	20" - 22"	462-0014-60
XX-large	22" - 24"	462-0014-70



DELUXE POST-OP KNEE IMMOBILIZER

## KNEE BRACES

### **SUPER-LITE KNEE IMMOBILIZER**

Provides removable contoured aluminum stays located medially, laterally and posteriorly, three wide hook & loop straps to help secure and provide immobilization. Provides comfort and maintains proper positioning of the leg.

#### **12" SUPPER-LITE KNEE IMMOBILIZER**

SIZE	CAT. NO.
small	E59410
med	E59411
large	E59412
X-large	E59413

#### **16" SUPPER-LITE KNEE IMMOBILIZER**

SIZE	CAT. NO.
small	E59414
med	E59415
large	E59416
X-large	E59417
XX-large	E59418

#### **20" SUPPER-LITE KNEE IMMOBILIZER**

SIZE	CAT. NO.
small	E59419
X-large	E59420
XX-large	E59421

#### **24" SUPPER-LITE KNEE IMMOBILIZER**

SIZE	CAT. NO.
small	E59422
med	E59423
large	E59424
X-large	E59425



SUPER-LITE KNEE IMMOBILIZER

### **12" DELUXE THREE-PANEL KNEE IMMOBILIZER**

The Deluxe Three-Panel Knee Immobilizer features 3 unique air-foiled aluminum stays located medially, laterally, and posteriorly. Together they provide support and stability. Additional features include pressure sensitive loop-lock closures with removable adjustment tabs and a popliteal pad. Available in black.

#### **12" DELUXE THREE-PANEL KNEE IMMOBILIZER**

SIZE	CAT. NO.
Universal	462-0021-10



12" DELUXE THREE-PANEL KNEE IMMOBILIZER



## **CERVICAL COLLAR MEDIUM DENSITY**

The Medium Density Standard Cervical Collar is manufactured from medium density poly-foam. All collars are manufactured from medical grade poly-foam covered in an orthopedic stockinette. They feature a contoured design, and an adjustable rear Velcro® closure. The Firm Density Standard Cervical Collar is manufactured with two layers of orthopedic stockinette, and features an adjustable comfort plug for protruding Adam's apple. Available in beige.

### **3.00" CERVICAL COLLAR MEDIUM DENSITY**

SIZE	NECK CIRCUM.	CAT. NO.
small	12" - 14"	482-0021-30
med	14" - 16"	482-0021-40
large	16" - 18"	482-0021-50
X-large	18" - 20"	482-0021-60

### **3.50" CERVICAL COLLAR MEDIUM DENSITY UNIVERSAL**

SIZE	NECK CIRCUM.	CAT. NO.
Universal	12" - 20"	482-0025-10



CERVICAL COLLAR MEDIUM DENSITY

## **ECONOMY CERVICAL COLLAR**

This 2.5" collar features a serpentine design covered in an orthopedic stockinette, with a pressure sensitive closure. The Cervical Collar is manufactured from low density poly-foam. Available in beige.

### **2.5" ECONOMY CERVICAL COLLAR**

SIZE	NECK CIRCUM.	CAT. NO.
Universal	12" - 20"	482-0048-10



2.5" ECONOMY CERVICAL COLLAR

## **ECONOMY CERVICAL COLLAR**

This 3" collar features a serpentine design covered in an orthopedic stockinette, with a pressure sensitive closure. The Cervical Collar is manufactured from low density poly-foam. Available in beige.

### **3" ECONOMY CERVICAL COLLAR**

SIZE	NECK CIRCUM.	CAT. NO.
Universal	12" - 20"	482-0050-10




3" ECONOMY CERVICAL COLLAR

## **BUCKS TRACTION WITH LINER**

Bucks Traction is used to provide traction to fractures of the shaft of the femur. Flexible lateral and medial stays are structured to fit leg contours. The elastic crossover foot strap eliminates pressure on the heel and prevents foot drop.

<b>SIZE</b>	<b>CAT. NO.</b>
Universal	872-0012-10



BUCKS TRACTION WITH LINER

## **TUMMYGRIP 33 FOOT ROLL**

Tummy Grip is a wide, elastic enhanced, stockinette that is used clinically to support the abdomen after surgery, to retain dressing placement, and as a way to keep electrodes in place during child birth without the need for adhesives. Tummy Grip stockinette is reinforced with a circumferential elastic thread knitted into the product which provides even consistent compression.

### **12" TUMMYGRIP**

<b>SIZE</b>	<b>WIDTH</b>	<b>CAT. NO.</b>
33 FOOT	12"	CT120



TUMMYGRIP 33 FOOT ROLL

## **HED-LOC**

Head Lock immobilizer for spine board. Long securing straps, re-useable, easy to clean

<b>SIZE</b>	<b>CAT. NO.</b>
Universal	HL-112



HED-LOC

## ADD-A-SPLINT

Fast and effective first responder splint applications. The splints are re-useable and economically priced to fit the tightest budgets. Several of the splints are designed for multi-use applications. Easy to apply and remove. Easy to use hook & loop closure no buckles or straps, universal left and right design, X-ray and MRI compatible, suitable for extreme hot or cold temperatures.

SIZE	DESCRIPTION	CAT. NO.
Universal	Complete Kit	AD-100
Universal	Wrist/Lower Arm	AD-101
Universal	Half Arm	AD-102
Universal	Full Arm	AD-103
Universal	Sing/Swathe Immobilizer	AD-104
Universal	Full Leg	AD-105
Universal	Ankle/Elbow	AD-106
Universal	Carrying Case	AD-107



ADD-A-SPLINT KIT

## POLYCAST CAST BOOT

Ideal for use with thinner casting materials such as polyester or fiberglass. Conforms to any casting material.

### Features:

- One-buckle closure on each strap
- Non-skid outsole
- Light weight high top upper

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
Pediatric	13 - 3	-	-	CB-PC1
X-small	-	3 - 5	-	CB-PC2
small	-	5.5 - 8	-	CB-PC3
med	-	8.5 - 10	6 - 8	CB-PC4
large	-	-	8.5 - 11	CB-PC5
X-large	-	-	11.5 - 14	CB-PC6



POLYCAST CAST BOOT

## LOW PROFILE CAST SHOE

The Low Profile boot provides a durable canvas upper, with an open toe and heel. Padded EVA insole protection. The lower outsole compensates for casting thickness. Available in six unisex sizes.

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
small	-	3 - 4.5	-	CB-TP1
small/med	-	5 - 7.5	4.5 - 6	CB-TP2
med	-	8 - 9	6.5 - 8	CB-TP3
med/large	-	9.5 - 12	8.5 - 11	CB-TP4
large	-	-	11.5 - 13	CB-TP5
X-large	-	-	13.5 - 16	CB-TP6



LOW PROFILE CAST SHOE

## CANVAS / OPEN TOE CAST SHOE

This lightweight cast boot features a sturdy canvas upper. The three layer EVA rocker outsole provides excellent shock absorption.

### Features:

- Adjustable heel strap
- Open toe and heel
- Five sizes: XS-XL

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
X-small	-	2 - 3	-	CB-CV1
small	-	3.5 - 5	-	CB-CV2
med	-	5.5 - 9.5	7 - 9	CB-CV3
large	-	10 - 12	9.5 - 12	CB-CV4
X-large	-	-	12.5 - 15	CB-CV5



CANVAS / OPEN TOE CAST SHOE

## VINYL CLOSED TOE CAST BOOT

The choice for inclement weather. This boot provides both closed toe and heel with a moisture resistant vinyl upper.

Features:

- Three layer EVA rocker outsole
- Excellent shock absorption

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
X-small	-	2 - 3	-	CB-VN1
small	-	3.5 - 5	-	CB-VN2
med	-	5.5 - 9.5	7 - 9	CB-VN3
large	-	10 - 12	9.5 - 12	CB-VN4
X-large	-	-	12.5 - 15	CB-VN5



VINYL CLOSED TOE CAST BOOT

## CLASSIC POST OP SHOE

This popular post op shoe provides quality care at an outstanding value. The semi-rigid outsole provides stability and minimizes flexure.

Features:

- Semi-rigid outsole provides medial/lateral stability
- Breathable nylon mesh upper
- Simple, one buckle closure on each strap

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
Pediatric	10 - 1	-	-	PO-CL1
WS	-	4 - 6	-	PO-CL2
WM	-	6.5 - 8	-	PO-CL3
WL	-	8.5 - 10	-	PO-CL4
MS	-	-	6 - 8	PO-CL5
MM	-	-	8.5 - 10	PO-CL6
ML	-	-	10.5 - 12	PO-CL7
MXL	-	-	12.5 - 14	PO-CL8
MXXL	-	-	14.5 - 16.5	PO-CL9



CLASSIC POST OP SHOE

## EZY CLOSE

The EZY Close post op shoe's flap-on-flap closure provides a secure fit and can accommodate even the bulkiest of dressing. Perfect for older patients with limited hand dexterity.

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
Pediatric	10 - 1	-	-	PO-EZ1
WS	-	4 - 6	-	PO-EZ2
WM	-	6.5 - 8	-	PO-EZ3
WL	-	8.5 - 10	-	PO-EZ4
MS	-	-	6 - 8	PO-EZ5
MM	-	-	8.5 - 10	PO-EZ6
ML	-	-	10.5 - 12	PO-EZ7
MXL	-	-	12.5 - 14	PO-EZ8



EZY CLOSE

## HI-TOP POST OP SHOE

The design of the Hi-Top ankle strap allows the foot to be held back in the heel cup of the shoe counter. The front is cut back which relieves pressure and compression on the toes. Features breathable nylon upper, adjustable ankle strap and two adjustable top straps.

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
Pediatric	10 - 1	-	-	PO-HR1
WS	-	4 - 6	-	PO-HR2
WM	-	6.5 - 8	-	PO-HR3
WL	-	8.5 - 10	-	PO-HR4
MS	-	-	6 - 8	PO-HR5
MM	-	-	8.5 - 10	PO-HR6
ML	-	-	10.5 - 12	PO-HR7
MXL	-	-	12.5 - 14	PO-HR8



HI-TOP POST OP SHOE

## SOFT TOUCH POST OP SHOE

This economically priced semi-rigid post op shoe is ideal for the sensitive foot. A soft cotton-poly blend breathable upper will keep your patient both cool and comfortable. The hook and loop closure provides even tension and compression. The reinforced heel counter provides additional comfort and support for your patient. The black TPR outsole is flared in the heel area to provide better medial/lateral stability. The outsole has the famous Health Design notch to give it a more shoe-like profile. Ideal for use following surgery, as well as for post trauma, dislocations, burns and the foot sensitive patient.



SOFT TOUCH POST OP SHOE

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
Pediatric	10 - 1	-	-	PO-ST1
WS	-	4 - 6	-	PO-ST2
WM	-	6.5 - 8	-	PO-ST3
WL	-	8.5 - 10	-	PO-ST4
MS	-	-	6 - 8	PO-ST5
MM	-	-	8.5 - 10	PO-ST6
ML	-	-	10.5 - 12	PO-ST7
MXL	-	-	12.5 - 14	PO-ST8

## FOREFOOT RELIEVER

Designed for forefoot relief and healing. This shoe has a unique wedge outsole and a soft nylon upper. We also added the hi-top upper to keep the foot held back and reduce slippage. Available in five sizes XS - XL.



FOREFOOT RELIEVER

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
X-small	-	3 - 5.5	-	PO-FR1
small	-	6 - 8.5	-	PO-FR2
med	-	9 - 11.5	7 - 9.5	PO-FR3
large	-	-	10 - 12.5	PO-FR4
X-large	-	-	13 - 15	PO-FR5

## POSI WEDGE

The Posi Wedge is a half shoe designed to give forefoot relief. This shoe uses three hook and loop straps for closure. Available in four sizes S - XL.



POSI WEDGE

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
small	-	5 - 7.5	-	PO-PW1
med	-	8 - 11.5	6 - 9.5	PO-PW2
large	-	-	10 - 12.5	PO-PW3
X-large	-	-	13 - 15	PO-PW4

## HD LOW WALKER

This low cost walker features a non-skid rocker outsole with a cushioned EVA innersole. The low profile design minimizes leg length discrepancy. Three hook & loop closures with D-rings makes strapping quick and easy. Steel uprights are wider between the struts to accommodate swollen ankle and bulky bandaging. Designed for easy application and patient comfort. Replaceable pads.

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
X-small	-	4 - 5.5	-	LWS-MM0
small	-	6 - 8.0	4.5 - 7	LWS-MM1
med	-	8.5 - 11.5	7.5 - 10.5	LWS-MM2
large	-	12 - 13.5	11 - 12.5	LWS-MM3
X-large	-	-	13 - 15	LWS-MM4



HD LOW WALKER

## HD TALL WALKER

This low cost walker features a non-skid rocker outsole with a cushioned EVA innersole. The low profile design minimizes leg length discrepancies. Five hook & loop closures with D-rings makes strapping quick and easy. Steel uprights are wider between the struts to accommodate swollen ankles and extra bandaging. This extra width insures easy application and better patient comfort. Replaceable pads.

SIZE	PEDIATRIC	WOMEN	MEN	CAT. NO.
X-small	-	4 - 5.5	-	LWL-MM0
small	-	6 - 8.0	4.5 - 7	LWL-MM1
med	-	8.5 - 11.5	7.5 - 10.5	LWL-MM2
large	-	12 - 13.5	11 - 12.5	LWL-MM3
X-large	-	-	13 - 15	LWL-MM4



HD TALL WALKER



## 1. CERVICAL NECK

SIZE	NECK CIRCUM.
X-small	10" - 13"
small	12" - 14"
medium	14" - 16"
large	16" - 18"
X-large	18" - 20"

## 2. Rib Belt & Shoulder Immobilizer

SIZE	CHEST CIRCUM.
X-small	> 15.5
small	25.5" - 35.5"
medium	35.5" - 43.5"
large	43" - 51"

## 3. WRIST SUPPORT

SIZE	WRIST CIRCUM.
small	< 6"
medium	6" - 7.5"
large	7.5" - 8.5"
X-large	> 8.5"

## 4. ARM SLING


SIZE	LENGTH	DIMENSION
small	7"	13.5" - 15"
medium	7"	14.5" - 17"
large	8"	16.5" - 19"
X-large	8"	> 19"

## 5. KNEE BRACES

SIZE	THIGH CIRCUM.
small	14" - 16"
medium	16" - 18"
large	18" - 20"
X-large	20" - 22"


## 5. NEOPRENE KNEE SUPPORT

SIZE	KNEE CIRCUM.
X-small	12" - 13"
small	13" - 14"
medium	14" - 15"
large	15" - 16"
X-large	16" - 17"
XX-large	17" - 18"
XXX-large	18" - 19"



**1. CLAVICLE BACK SUPPORT**

<i>SIZE</i>	<i>CHEST CIRCUM.</i>
Pediatric	13" - 20"
small	16" - 23"
medium	22" - 30"
large	28" - 39"
X-large	37" - 45"



**2. BACK SUPPORTS  
ABDOMINAL BINDER**

<i>SIZE</i>	<i>WAIST CIRCUM.</i>
small/med.	16" - 40"
large	40" - 62"

**3. ANKLE SUPPORT**

<i>SIZE</i>	<i>MALE SHOE SIZE</i>	<i>FEMALE SHOE SIZE</i>
X-small	3 - 5	4 - 6
small	6 - 7	7 - 9
medium	8 - 10	10 - 11
large	11 - 12	12 - 13
X-large	13 - 14	14+





**Health Design**

810 Memorial Blvd.

Huntington, WV 25701 | USA

Toll Free voice 888-847-8090

Toll Free Fax 888-847-8089

Voice 304-781-5234

[www.healthdesign.net](http://www.healthdesign.net)

