

Henry Schein PRO Digital Thermometer

From A Name You Can Trust.



Features & Benefits

- This unit is intended to measure body temperature in oral, axillary, and rectal
- 60 second digital thermometer
- 5 probe covers

(570-2055)Each

Physical Characteristics:

Battery:	One 1.5V button size battery (LR 41)
Power consumption:	0.15 MW in measurement mode
Battery life:	More than 200 hours of continuous operation
Weight:	Approx. 10~12 grams including battery
Dimensions:	130mm x 20mm x 14mm

Description of Performance and Characteristics:

Range:	32.0 ~ 43.0 °C (89.6°F ~ 109.4°F)
Accuracy:	±0.2 °C (±0.3 °F) Less than 35.5 °C (95.9 °F)
(tested in Laboratory)	±0.1 °C (±0.2 °F) 35.5 °C to 42.0 °C (95.9 °F to 107.6 °F) ±0.2 °C (±0.3 °F) 42.0 °C to 43.0 °C (107.6 °F to 109.4 °F)
Display:	Liquid crystal display 3 1/2 digits (19 mm x 6 mm)
Operation Condition:	15°C ~ 40°C (59°F ~ 104°F) / 15% – 90% RH
Storage Condition:	-25°C ~ 55°C (-13°F ~ 131°F) / 15% - 90% RH



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!