

One Step+ H.Pylori Test Kit

From A Name You Can Trust.

Features & Benefits

Rapid immunoassay for the qualitative detection of anti-Helicobacter pylori antibodies in human whole blood, serum or plasma

- CLIA Waived for whole blood
- Sensitivity of 95.9% with 91.6% overall agreement to biopsy histology
- Easy-to-use single reagent procedure
- Results in 10 minutes
- Requires only 25µl of whole blood or 10µl serum/plasma



Item Number: 570-2778

Contents:

- 20 Test devices containing a membrane strip coated with H. pylori antigen and a pad with indicator conjugates in a protein matrix)
- 20 Capillary tubes
- 1 Developer solution (containing 0.1% sodium azide)
- 1 Positive control
- 1 Negative control
- 1 Instructional insert



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!