

Connectors, Adapters and Tees

Page

42-48 Connectors, Adapters and Tees

49-51 Adapters

52-54 Couplers

55..... Block Assemblies

56..... DISS Fittings

57-58 Hose Assemblies

Connectors, Adapters and Tees Tees



AirLife™ trach tee adapter 001500

22 mm OD x 22 mm OD
15 mm ID base
50/case



AirLife™ tee adapter for small nebulizers 001503

22 mm OD x 22 mm OD
18 mm base
50/case



AirLife™ tee adapter 001505 shown with mouthpiece

22 mm ID x 22 mm OD
18 mm ID base and mouthpiece
50/case



AirLife™ tee adapter with mouthpiece 001506

22 mm OD x 22 mm OD
18 mm base with mouthpiece
50/case



AirLife™ baffled tee adapter 001508

22 mm OD x 22 mm OD
18 mm ID base
50/case



AirLife™ tee 004000

22 mm OD x 22 mm ID
15 mm base
50/case



AirLife™ tee 004001

22 mm ID x 22 mm ID
15 mm base
50/case



AirLife™ tee 004011

22 mm OD x 22 mm OD
15 mm base
50/case



AirLife™ tee with base 004225

22 mm ID x 22 mm OD
18 mm ID base
50/case



AirLife™ tee w/one way valves 004050

22 mm OD x 22 mm OD and 15 mm base
with oxygen stem and thermometer port
50/case



AirLife™ tee w/one way valves 004051

22 mm OD x 22 mm OD and 15 mm base
with oxygen stem
50/case



AirLife™ tee 004052

22 mm OD x 22 mm OD and 15 mm base
with thermometer port
50/case

Connectors, Adapters and Tees *Valved Tees and Straight Connectors*



AirLife™ Valved Tee adapter 002058

15 mm ID x 15 mm OD
30/case



AirLife™ Valved Tee adapter 002060

22 mm OD x 22 mm OD
30/case



AirLife™ Valved Tee adapter 002061

22 mm ID x 22 mm OD
30/case



AirLife™ straight connector 004078

22 mm ID x 22 mm OD
with oxygen stem and thermometer port
50/case



AirLife™ straight connector 004080

22 mm OD x 22 mm OD
with oxygen stem and thermometer port
50/case



AirLife™ straight connector 004081

22 mm ID x 22 mm OD
with oxygen stem
50/case



AirLife™ straight connector 004083

22 mm OD x 22 mm OD
with oxygen stem
50/case



AirLife™ straight connector 004084

22 mm ID x 22 mm OD
with thermometer port
50/case



AirLife™ straight connector 004086

22 mm OD x 22 mm OD
with thermometer port
50/case



AirLife™ straight connector 004154

22 mm ID x 22 mm OD
with oxygen stem and thermometer port
50/case



AirLife™ straight connector 004204

22 mm ID x 22 mm OD
with oxygen stem
50/case



AirLife™ temperature probe adapter 5979-504

22 mm ID x 22 mm OD
with temperature port
50/case

Connectors, Adapters and Tees *Connectors and Adapters*



AirLife™ elbow 004301

15 mm ID base x 22 mm ID
50/case



AirLife™ elbow 004303

15 mm ID base x 22 mm ID
with temperature port and oxygen stem
50/case



AirLife™ ventilator elbow 5930-504

15 mm OD x 15 mm OD
50/case



AirLife™ ventilator elbow 5999-504

22 mm ID x 22 mm OD
50/case



AirLife™ elbow 001550

22 mm ID x 22 mm OD inlet arm
50/case



AirLife™ elbow 001551

22 mm ID x 22 mm ID inlet arm
50/case



AirLife™ Y connector 001830

22 mm OD x 22 mm OD x 15 mm ID
50/case



AirLife™ Y connector 001831

22 mm OD x 22 mm OD x 15 mm ID
with integral thermometer and
pressure-sensing ports
50/case



AirLife™ Y connector 5924-504

Standard infant Y connector with
suction port and strap plug
30/case



AirLife™ Y connector 5925-504

Standard infant Y connector
(no ports)
30/case



AirLife™ pediatric connector 3215

Omni-Flex™ connector, pediatric
50/case



AirLife™ adult connector 3222

Omni-Flex™ connector, adult
50/case

Connectors, Adapters and Tees *Connectors and Adapters*



AirLife™ connector 5908-504

0.165" OD x 0.258" ID
50/case



AirLife™ connector 5911-504

22 mm OD x 22 mm ID
150/case



AirLife™ connector 5912-504

7.5 mm OD x 0.258" OD
150/case



AirLife™ connector 5913-504

22 mm ID x 22 mm ID
150/case



AirLife™ connector 5915-504

7.5 mm OD x 0.165" OD
150/case



AirLife™ connector 5916-504

15 mm OD x 22 mm OD
150/case



AirLife™ oxygen plug 5917-504

7.5 mm OD plug
150/case



AirLife™ connector 5918-504

3-way
150/case



AirLife™ oxygen plug 5920-504

4.5 mm plug
150/case



AirLife™ connector 5923-504

15 mm ID x 22 mm ID
150/case



AirLife™ connector 5939-504

15 mm OD x 15 mm OD
150/case



AirLife™ oxygen tubing adapter 5906-504

15 mm OD x 0.258" OD/22 mm OD
50/case

Connectors, Adapters and Tees *Connectors and Adapters*



AirLife™ drive line adapter 001810

50/case




AirLife™ circuit hanger 001832

50/case



AirLife™ pressure adapter 001801

50/case



AirLife™ large one-way valve 001802

22 mm OD x 22 mm OD
50/case



AirLife™ straight connector 001800

22 mm OD x 22 mm OD
50/case



AirLife™ trach tee oxygenator 001507

with flextube and oxygen stem connector
50/case



AirLife™ weaning parameter kit 001504

50/case

Connectors, Adapters and Tees *Connectors and Adapters*



AirLife™ intubation adapter 001820
22 mm OD/15 mm ID both ends
50/case



AirLife™ cuff connector 001821
22 mm ID x 22 mm ID
50/case



AirLife™ connector 5940-504
22 mm OD x 22 mm OD
30/case



AirLife™ tubing connector 001823
22 mm OD x 22 mm OD
One end tapered
50/case




AirLife™ tubing connector 001822
22 mm OD x 22 mm OD
Both ends tapered
50/case



AirLife™ spirometer adapter 001812
For use with Puritan-Bennett MA-1™
and MA-2™ spirometers
50/case



AirLife™ oxygen tubing connector 001811
50/case



AirLife™ 001842
Oxygen tubing swivel connector
Male/female
25/case



AirLife™ male/male connector 001841
Oxygen tubing swivel connector
Male/male
25/case



AirLife™ two-step tubing adapter 001824
22 mm OD x 15 mm OD
2-step tubing adapter
100/case




AirLife™ oxygen stem with plug 001803
50/case



AirLife™ oxygen nut/nipple connector 001840
50/box 250/case

Connectors, Adapters and Tees *Connectors, Adapters and Tees*



AirLife™ connector 5991-504
15 mm OD x 22 mm ID
150/case



AirLife™ tee-type infant temperature probe adapter 5941-504
30/case



AirLife™ dog-leg-type infant temperature probe adapter 5921-504
30/case



AirLife™ infant neb tee with cap 5903-504
30/case



AirLife™ infant neb tee shown with nebulizer (not included) 5900-504
30/case



AirLife™ oxygen sensor adapter kit 5928-504
30/case



AirLife™ universal mouthpiece 001000
500/case
















AirLife™ Brand UIAdapt-it™ exhalation manifold body 004381
22 mm OD x 22 mm OD



AirLife™ Brand UIAdapt-it™ exhalation manifold body 004379
22 mm ID x 22 mm OD

Ohmeda® Adapters





<p>Please contact us for information about MRI options</p>	 Oxygen	 Medical Air	 Vacuum	 N ₂ O	 WAGD	 CO ₂
 1/8" NPT Male	ADOHOM2	ADOHAM2	ADOHVM2	ADOH2M2	ADOHWM2	ADOHCM2
 1/4" Hose Barb 5/16" Hose Barb	ADOHOHB4	ADOHAHB4	ADOHVHB4 ADOHVHB5	ADOH2HB4	ADOHWHB4 ADOHWHB5	ADOHCHB4
 DISS Male	ADOHODM	ADOHADM	ADOHVDM	ADOH2DM	ADOHWDM	ADOHCDM
 1/4" NPT Male	ADOHOM4	ADOHAM4	ADOHVM4	ADOH2M4	ADOHWM4	ADOHCM4
 1/8" NPT Female	ADOHOF2	ADOHAF2	ADOHVF2	ADOH2F2	ADOHWF2	ADOHCF2
 DISS Handtight	ADOHODH	ADOHADH	ADOHVDH	ADOH2DH	ADOHWDH	ADOHCDH
 DISS Nut and Nipple	ADOHODN	ADOHADN	ADOHVDN	ADOH2DN	ADOHWDN	ADOHCDN

Connectors, Adapters and Tees *Adapters*





Chemetron® Adapters

<p>Please contact us for information about MRI options</p>	 Oxygen	 Medical Air	 Vacuum	 N ₂ O	 WAGD	 CO ₂
 1/8" NPT Male	ADCHOM2	ADCHAM2	ADCHVM2	ADCH2M2	ADCHWM2	ADCHCM2
 1/4" Hose Barb 3/16" Hose Barb	ADCHOHB4	ADCHAHB4	ADCHVHB4 ADCHVHB5	ADCH2HB4	ADCHWHB4 ADCHWHB5	ADCHCHB4
 DISS Male	ADCHODM	ADCHADM	ADCHVDM	ADCH2DM	ADCHWDM	ADCHCDM
 1/4" NPT Male	ADCHOM4	ADCHAM4	ADCHVM4	ADCH2M4	ADCHWM4	ADCHCM4
 1/8" NPT Female	ADCHOF2	ADCHAF2	ADCHVF2	ADCH2F2	ADCHWF2	ADCHCF2
 DISS Handtight	ADCHODH	ADCHADH	ADCHVDH	ADCH2DH	ADCHWDH	ADCHUCDH
 DISS Nut and Nipple	ADCHODN	ADCHADN	ADCHVDN	ADCH2DN	ADCHWDN	ADCHCDN

Oxequip® and Medstar Adapters





<i>Please contact us for information about MRI options</i>	Oxygen	Medical Air	Vacuum	N ₂ O
 1/8" NPT Male, Oxequip	ADXQOM2	ADXQAM2	ADXQVM2	ADXQ2M2
 1/4" Hose Barb, Oxequip 5/16" Hose Barb, Oxequip	ADXQOHB4	ADXQAHB4	ADXQVHB4 ADXQVHB5	ADXQ2HB4
 1/8" NPT male, Medstar	ADMSOM2	ADMSAM2	ADMSVM2	ADMS2M2
 1/4" Hose Barb, Medstar 5/16" Hose Barb Medstar	ADMSOHB4	ADMSAHB4	ADMSVHB4 ADMSVHB5	ADMS2HB4

Puritan Bennett® Adapters and Couplers

<i>Please contact us for information about MRI options</i>	Oxygen	Medical Air	Vacuum	N ₂ O	WAGD
 1/8" NPT Male	ADPBOM2 COPBOM2	ADPBAM2 COPBAM2	ADPBVM2 COPBVM2	ADPB2M2 COPB2M2	ADPBWM2 COPBWM2
 1/4" Hose Barb 5/16" Hose Barb	ADPBOHB4 COPBOHB4	ADPBAHB4 COPBAHB4	ADPBVHB4 ADPBVHB5 COPBVHB4 COPBVHB5	ADPB2HB4 COPB2HB4	ADPBWHB4 ADPBWHB5 COPBWHB4 COPBWHB5
 1/4" NPT Male	ADPBOM4 COPBOM4	ADPBAM4 COPBAM4	ADPBVM4 COPBVM4	ADPB2M4 COPB2M4	ADPBWM4 COPBWM4
 1/8" NPT Female	ADPBOF2 COPBOF2	ADPBAF2 COPBAF2	ADPBVF2 COPBVF2	ADPB2F2 COPB2F2	ADPBWF2 COPBWF2

Connectors, Adapters and Tees *Couplers*

Schrader® Adapters and Couplers














<i>Please contact us for information about MRI options</i>	Oxygen	Medical Air	Vacuum	N ₂ O	Nitrogen
 1/8" NPT Male	ADSHOM2 COSHOM2	ADSHAM2 COSHAM2	ADSHVM2 COSHVM2	ADSH2M2 COSH2M2	ADSHNM2 COSHNM2
 1/4" Hose Barb 3/16" Hose Barb	ADSHOHB4 COSHOB4	ADSHAHB4 COSHAHB4	ADSHVHB4 COSHVHB4	ADSH2HB4 COSH2HB4	ADSHNHB4 COSHNB4
 1/4" NPT Male	ADSHOM4 COSHOM4	ADSHAM4 COSHAM4	ADSHVM4 COSHVM4	ADSH2M4 COSH2M4	ADSHNM4 COSHNM4
 1/8" NPT Female	ADSHOF2 COSHOF2	ADSHAF2 COSHAF2	ADSHVF2 COSHVF2	ADSH2F2 COSH2F2	ADSHNF2 COSHNF2

Important Safety Reminders

- Gas-specific adapters were developed to help prevent delivering the wrong gas to a patient. Before using an adapter, verify the source to ensure it is the appropriate gas for the application. Never use adapters to change the connectors or fittings on cryogenic vessels.
- To reduce the risk of fire, ensure that substances such as dirt, oil, grease or hand creams do not contaminate connections between medical devices and medical gas cylinders.

Visit the Compressed Gas Association website at <http://cganet.com> for free safety posters and a variety of safety and educational materials.

Ohmeda® Couplers

<p>Please contact us for information about MRI options</p>	 Oxygen	 Medical Air	 Vacuum	 N ₂ O	 WAGD	 CO ₂
 1/8" NPT Male	COOHOM2	COOHAM2	COOHVM2	COOH2M2	COOHWM2	COOHCM2
 1/4" Hose Barb 5/16" Hose Barb	COOHOHB4	COOHAHB4	COOHVHB4 COOHVHB5	COOH2HB4	COOHWHB4 COOHWHB5	COOHCHB4
 DISS Male	COOHODM	COOHADM	COOHVDM	COOH2DM	COOHWDM	COOHCDM
 1/4" NPT Male	COOHOM4	COOHAM4	COOHVM4	COOH2M4	COOHWM4	COOHCM4
 1/8" NPT Female	COOHOF2	COOHAF2	COOHVF2	COOH2F2	COOHWF2	COOHCF2
 DISS Handtight	COOHODH	COOHADH	COOHVDH	COOH2DH	COOHWDH	COOHCDH
 DISS Nut and Nipple	COOHODN	COOHADN	COOHVDN	COOH2DN	COOHWDN	COOHCDN

Connectors, Adapters and Tees *Couplers*

Chemetron® Couplers

<p>Please contact us for information about MRI options</p>	 Oxygen	 Medical Air	 Vacuum	 N ₂ O	 WAGD	 CO ₂
 1/8" NPT Male	COCHOM2	COCHAM2	COCHVM2	COCH2M2	COCHWM2	COCHCM2
 1/4" Hose Barb 5/16" Hose Barb	COCHOHB4	COCHAHB4	COCHVHB4 COCHVHB5	COCH2HB4	COCHWHB4 COCHWHB5	COCHCHB4
 DISS Male	COCHODM	COCHADM	COCHVDM	COCH2DM	COCHWDM	COCHCDM
 1/4" NPT Male	COCHOM4	COCHAM4	COCHVM4	COCH2M4	COCHWM4	COCHCM4
 1/8" NPT Female	COCHOF2	COCHAF2	COCHVF2	COCH2F2	COCHWF2	COCHCF2
 DISS Handtight	COCHODH	COCHADH	COCHVDH	COCH2DH	COCHWDH	COCHCDH
 DISS Nut and Nipple	COCHODN	COCHADN	COCHVDN	COCH2DN	COCHWDN	COCHCDN

Block Assemblies

<p><i>Please contact us for information about MRI options</i></p>			
Oxygen	YBC00H00H	YBCOCHOCH	YBDMDCODH
Air	YBC00HA0H	YBCOCHACH	YBDMDCADH
VAC	YBC00HVOH	YBCOCHVCH	YBDMDCVDH
N ₂ O	YBC00H20H	YBCOCH2CH	YBDMDC2DH
CO ₂	YBC00HCOH	YBCOCHCCH	YBDMDCCDH
NIT	N/A	N/A	YBDMDCNDH













Who is CareFusion?

Our proven product family of brands include Pyxis® for medication and supply dispensing, Alaris® for infusion, AVEA® ventilators, Jaeger® for respiratory diagnostic instruments, AirLife™ for respiratory consumables, AVAmax® and PleurX® for interventional procedures, V. Mueller® and Snowden-Pencer® surgical instruments, Chloraprep® skin antiseptic and Nicolet™ neurodiagnostic devices.

To make global healthcare better, we partner with our customers to help them improve medication management, lower costs in procedural areas, reduce risk of infection, advance the care of ventilated patients and turn the endless amount of data generated in healthcare into actionable information.











Connectors, Adapters and Tees *DISS Fittings*

DISS Fittings

<i>Please contact us for information about MRI options</i>		Oxygen	Medical Air	Vacuum	N ₂ O	Nitrogen	WAGD	CO ₂	Instrument Air
1/8" NPT Male		BAOM2	BAAM2	BAVM2	BA2M2	BANM2	BAWM2	BACM2	BAIM2
		DA150M2DN	DA15AM2DN	DA15VM2DN	DA152M2DN	DA15NM2DN	DA15WM2DN	DA15CM2DN	DA15IM2DN
		DA150M2DH	DA15AM2DH	DA15VM2DH	DA152M2DH	DA15NM2DH	DA15WM2DH	DA15CM2DH	DA15IM2DH
1/4" & 5/16" Hose Barb		BAOHB4	BAAHB4	BAVHB4 BAVHB5	BA2HB4	BANHB4	BAWHB4 BAWHB5	BACHB4	BAIHB4
		DAH40DN	DAH4ADN	DAH4VDN DAH5VDN	DAH42DN	DAH4NDN	DAH4WDN DAH5WDN	DAH4CDN	DAH4IDN
		DAH40DH	DAH4ADH	DAH4VDH DAH5VDH	DAH42DH	DAH4NDH	DAH4WDH DAH5WDH	DAH4CDH	DAH4IDH
1/4" NPT Male		BAOM4	BAAM4	BAVM4	BA2M4	BANM4	BAWM4	BACM4	BAIM4
		DA150M4DN	DA15AM4DN	DA15VM4DN	DA152M4DN	DA15NM4DN	DA15WM4DN	DA15CM4DN	DA15IM4DN
		DA150M4DH	DA15AM4DH	DA15VM4DH	DA152M4DH	DA15NM4DH	DA15WM4DH	DA15CM4DH	DA15IM4DH
1/8" NPT Female		BAOF2	BAAF2	BAVF2	BA2F2	BANF2	BAWF2	BACF2	BAIF2
		DA150F2DN	DA15AF2DN	DA15VF2DN	DA152F2DN	DA15NF2DN	DA15WF2DN	DA15CF2DN	DA15IF2DN
		DA150F2DH	DA15AF2DH	DA15VF2DH	DA152F2DH	DA15NF2DH	DA15WF2DH	DA15CF2DH	DA15IF2DH

Connectors, Adapters and Tees *Hose Assemblies*

Hose Assemblies

Please contact us for information about MRI options		 Ohmeda Adapter	 Chemetron Adapter	 DISS Handtight	 DISS Hex Nut	 DISS Male with Demand Check
 Ohmeda Coupler	HSXXYZOFOMC4	HSXXYZOFCMC4	HSXXYZDHOFC4	HSXXYZDFOFC4	HSXXYZOFDCC4	
 Chemetron Coupler	HSXXYZCFOMC4	HSXXYZCFCMC4	HSXXYZCFDHC4	HSXXYZCFDFC4	HSXXYZCFDCC4	
 DISS Handtight	HSXXYZDHOMC4	HSXXYZDHCMC4	HSXXYZDHDHC4	HSXXYZDFDHC4	HSXXYZDHDCC4	
 DISS Hex Nut	HSXXYZDFOMC4	HSXXYZDFCMC4	HSXXYZDFDHC4	HSXXYZDFDFC4	HSXXYZDFDCC4	
 DISS Male with Demand Check	HSXXYZDCOMC4	HSXXYZCMDCC4	HSXXYZDHDCC4	HSXXYZDFDCC4	HSXXYZDCDCC4	

Hose Configurator:

Gas Type:

Oxygen:	O
Medical air:	A
Medical vacuum:	V
N ₂ O:	Z
Nitrogen:	N
WAGD:	W
CO ₂ :	C
Instrument Air:	I
O ₂ He:	J

For correct hose matrix, follow these rules when selecting inlet and outlet connections:

1. Female fittings first
2. Numbers before letters
3. In alphabetical order

Color Code:

USA:	U
ISO:	I
Black non-conductive:	K
Pin stripe (ISO air):	P
Stripe:	S

Length in feet

HS-XYZ-XXYYXX

Hose Type:

Conductive 1/4":	C4
Conductive 5/16":	C5
Emission hose 1/4":	E4
General purpose 1/4":	G4
Nonconductive 1/4":	N4
Nonconductive 5/16":	N5

Inlet (XX) and Outlet (YY) Connections:

<p>Female:</p> <table border="0"> <tr><td>Air M16 adapter:</td><td>16</td><td>Medstar female:</td><td>MF</td></tr> <tr><td>Australian handtight:</td><td>AH</td><td>NIST female:</td><td>NF</td></tr> <tr><td>British female:</td><td>BF</td><td>Ohmeda female:</td><td>OF</td></tr> <tr><td>Chemetron female:</td><td>CF</td><td>Puritan female:</td><td>PF</td></tr> <tr><td>DISS female</td><td></td><td>Schrader female:</td><td>SF</td></tr> <tr><td>with backbody:</td><td>DB</td><td>Italian female:</td><td>UM</td></tr> <tr><td>DISS hex nut:</td><td>DF</td><td>AGA female:</td><td>ZF</td></tr> <tr><td>DISS handtight:</td><td>DH</td><td></td><td></td></tr> <tr><td>1/8" FNPT:</td><td>F2</td><td></td><td></td></tr> <tr><td>1/4" FNPT:</td><td>F4</td><td></td><td></td></tr> <tr><td>German female:</td><td>GF</td><td></td><td></td></tr> <tr><td>Japanese female:</td><td>JF</td><td></td><td></td></tr> </table>	Air M16 adapter:	16	Medstar female:	MF	Australian handtight:	AH	NIST female:	NF	British female:	BF	Ohmeda female:	OF	Chemetron female:	CF	Puritan female:	PF	DISS female		Schrader female:	SF	with backbody:	DB	Italian female:	UM	DISS hex nut:	DF	AGA female:	ZF	DISS handtight:	DH			1/8" FNPT:	F2			1/4" FNPT:	F4			German female:	GF			Japanese female:	JF			<p>Male:</p> <table border="0"> <tr><td>British male:</td><td>BM</td><td>German male straight:</td><td>GM</td></tr> <tr><td>Chemetron male:</td><td>CM</td><td>Japanese male:</td><td>JM</td></tr> <tr><td>DISS male with Demand check:</td><td>DC</td><td>1/8" MNPT:</td><td>M2</td></tr> <tr><td>DISS male with no Demand check:</td><td>DM</td><td>1/4" MNPT:</td><td>M4</td></tr> <tr><td>French male angled:</td><td>FE</td><td>Medstar male:</td><td>MM</td></tr> <tr><td>French male:</td><td>FM</td><td>NIST male:</td><td>NM</td></tr> <tr><td>German male</td><td></td><td>Ohmeda male:</td><td>OM</td></tr> <tr><td>no Check angled:</td><td>GE</td><td>Puritan male:</td><td>PM</td></tr> <tr><td>German male check</td><td></td><td>Schrader male:</td><td>SF</td></tr> <tr><td>Angled:</td><td>GC</td><td>Oxequip male:</td><td>SM</td></tr> <tr><td></td><td></td><td>Italian male:</td><td>UM</td></tr> <tr><td></td><td></td><td>AGA male:</td><td>XM</td></tr> <tr><td></td><td></td><td>Blank:</td><td>XX</td></tr> </table>	British male:	BM	German male straight:	GM	Chemetron male:	CM	Japanese male:	JM	DISS male with Demand check:	DC	1/8" MNPT:	M2	DISS male with no Demand check:	DM	1/4" MNPT:	M4	French male angled:	FE	Medstar male:	MM	French male:	FM	NIST male:	NM	German male		Ohmeda male:	OM	no Check angled:	GE	Puritan male:	PM	German male check		Schrader male:	SF	Angled:	GC	Oxequip male:	SM			Italian male:	UM			AGA male:	XM			Blank:	XX
Air M16 adapter:	16	Medstar female:	MF																																																																																																		
Australian handtight:	AH	NIST female:	NF																																																																																																		
British female:	BF	Ohmeda female:	OF																																																																																																		
Chemetron female:	CF	Puritan female:	PF																																																																																																		
DISS female		Schrader female:	SF																																																																																																		
with backbody:	DB	Italian female:	UM																																																																																																		
DISS hex nut:	DF	AGA female:	ZF																																																																																																		
DISS handtight:	DH																																																																																																				
1/8" FNPT:	F2																																																																																																				
1/4" FNPT:	F4																																																																																																				
German female:	GF																																																																																																				
Japanese female:	JF																																																																																																				
British male:	BM	German male straight:	GM																																																																																																		
Chemetron male:	CM	Japanese male:	JM																																																																																																		
DISS male with Demand check:	DC	1/8" MNPT:	M2																																																																																																		
DISS male with no Demand check:	DM	1/4" MNPT:	M4																																																																																																		
French male angled:	FE	Medstar male:	MM																																																																																																		
French male:	FM	NIST male:	NM																																																																																																		
German male		Ohmeda male:	OM																																																																																																		
no Check angled:	GE	Puritan male:	PM																																																																																																		
German male check		Schrader male:	SF																																																																																																		
Angled:	GC	Oxequip male:	SM																																																																																																		
		Italian male:	UM																																																																																																		
		AGA male:	XM																																																																																																		
		Blank:	XX																																																																																																		

* Not all configurations are possible



*Connecting and Empowering clinicians with technology,
products and services to help improve patient care.*