

Ferno® Model 12 Emergency Stretcher

- › Exclusive slimline center hinges enable stretcher to fold for compact storage
- › Heavy-duty vinyl-coated nylon cover
- › Lightweight aluminum frame
- › Can be used with Model 385-A Storage/Carrying Case
- › Includes two 5' (1524 mm) quick-release patient restraints

- Model 12, Burgundy** #0107742
- Model 12, Orange** #0101200
- Storage Case, Burgundy Unlettered** #0311650
- Storage Case, Burgundy Lettered** #0311640
- Storage Case, Orange Unlettered** #0311670
- Storage Case, Orange Lettered** #0311630



Specifications

	Unfolded Length	Folded Length	Width	Unfolded Height	Folded Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
English	74 in.	37 in.	18 in.	1 in.	3 in.	8 lb	350 lb	12 lb	38x21x4 in.
Metric	1873 mm	940 mm	464 mm	29 mm	60 mm	4 kg	159 kg	5 kg	965x533x102 mm

Ferno® Model 12-C Emergency Breakaway Stretcher

- › Low-profile frame for easier logrolling and lifting
- › Stretcher halves separate for easy removal without moving the patient
- › Frame adjusts to align with ambulance-cot backrest, enabling patient to have head elevated during transport
- › Includes two 5' (1524 mm) quick-release patient restraints
- › Heavy-duty vinyl-coated nylon cover

- Model 12-C, Burgundy** #0107743
- Model 12-C, Orange** #0107744



Specifications

	Unfolded Length	Folded Length	Width	Unfolded Height	Folded Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
English	74 in.	47 in.	9 in.	1 in.	3 in.	10 lb	350 lb	15 lb	75x19x2 in.
Metric	1870 mm	1190 mm	220 mm	30 mm	60 mm	5 kg	159 kg	7 kg	1905x483x51 mm

Ferno® 107 Series Combination Stretcher Chairs

- › Vinyl-coated nylon cover
- › Folds in half for compact storage
- › High-back chair design
- › Folding foot-end handles provide excellent grip
- › Stationary or hinged 4" wheels
- › Folding head-end posts double as lift handles and support legs
- › Legs fold underneath the stretcher for compact storage
- › Foot-end wheels fold flat for storage and snap back into place
- › Fits squad bench post-and-wheel cups for use as a second stretcher
- › Includes two 5' (1524 mm) quick-release patient restraints

Model 107, Burgundy #0107714
Model 107, Orange #0101070



Specifications

	Stretcher Length	Depth Chair	Depth Folded	Width	Height Stretcher	Height Chair	Height Folded	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
English	73 in.	30 in.	5 in.	20 in.	9 in.	53 in.	39 in.	24 lb	350 lb	27 lb	41x21x6 in.
Metric	1850 mm	760 mm	130 mm	510 mm	230 mm	1350 mm	980 mm	11 kg	159 kg	11 kg	1041x533x152 mm

Ferno® Model 107-B4

- › Fixed head-end posts are support legs
- › Fixed posts are in position when chair is unfolded
- › Fixed foot-end wheels allow easy maneuvering
- › Folding foot-end handles provide excellent grip
- › Can be used with Model 385-F carrying case
- › Includes two 5' (1524 mm) quick-release patient restraints

Model 107-B4, Burgundy #0107715
Model 107-B4, Orange #0101074
Storage Case, Lettered, Burgundy 385-F #0311762
Storage Case, Lettered, Orange 385-F #0311760



Specifications

	Stretcher Length	Depth Chair	Depth Folded	Width	Height Stretcher	Height Chair	Height Folded	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
English	82 in.	26 in.	6 in.	20 in.	6 in.	54 in.	40 in.	20 lb	350 lb	22 lb	41x21x6 in.
Metric	2080 mm	650 mm	140 mm	500 mm	150 mm	1370 mm	1010 mm	9 kg	159 kg	10 kg	1041x533x152 mm

Ferno® Model 107-C

- › Fixed foot-end posts provide stability when the stretcher is used flat
- › Fixed wheels at both ends permit it to be rolled like a chair or loaded into a specialty ambulance
- › Folding handles at both ends enable it to be used as a standard pole stretcher
- › Includes two 5' (1524 mm) quick-release patient restraints

Model 107-C, Burgundy #0107716
Model 107-C, Orange #0101072
Storage Case, Lettered, Burgundy 385-F #0311762
Storage Case, Lettered, Orange 385-F #0311760



Specifications

	Stretcher Length	Depth Chair	Depth Folded	Width	Height Stretcher	Height Chair	Height Folded	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
English	82 in.	26 in.	6 in.	20 in.	6 in.	54 in.	40 in.	25 lb	350 lb	30 lb	41x21x8 in.
Metric	2080 mm	650 mm	160 mm	500 mm	140 mm	1370 mm	1010 mm	11 kg	159 kg	14 kg	1041x533x203 mm