

BARD®

# WOUND DRAINAGE

PRODUCTS TO MEET YOUR DIFFERENT PROCEDURAL NEEDS



Closed Wound Suction  
Gravity Drains  
Sump Drains

# Channel Drains

The BARD® Channel Drain is a radiopaque closed wound suction drain distinguished by four free-flowing channels that offers optimized drainage for any procedure.

Latex-Free, Single-use, 10/case.

1. Improves drainage with less clotting
2. Minimizes tissue invagination
3. Reduces pain and tissue trauma for the patient upon removal
4. Features sizes for every surgery



Part Number	Size			Drain		
	Size Fr.	In.	mm	Shape	Fluted	Trocar
072186	10	1/8	3.2	Round	Full	
072187	10	1/8	3.2	Round	Full	1/8
072221	10	1/8	3.2	Round	3/4	1/8
072226	10	1/8	3.2	Round Hubless	3/4	
072227	10	1/8	3.2	Round Hubless	3/4	1/8
072188	15	3/16	4.7	Round	Full	
072189	15	3/16	4.7	Round	Full	3/16
072223	15	3/16	4.7	Round	3/4	3/16
072228	15	3/16	4.7	Round Hubless	3/4	
072229	15	3/16	4.7	Round Hubless	3/4	3/16
072190	19	1/4	6.3	Round	Full	
072191	19	1/4	6.3	Round	Full	1/4
072192	19	1/4	6.3	Round	Full	1/4
072225	19	1/4	6.3	Round	3/4	1/4
072230	19	1/4	6.3	Round Hubless	3/4	
072231	19	1/4	6.3	Round Hubless	3/4	1/4
072234	24	5/16	8.0	Round Hubless	Full	
072290	28	3/8	9.3	Round Hubless	Full	
072292	32	7/16	10.7	Round Hubless	Full	
072328	28	3/4	9.3	Round Hubless	Full, Ext. Length	
072332	32	7/16	10.7	Round Hubless	Full, Ext. Length	
072210	n/a	1/4	7	Flat	3/4	
072211	n/a	1/4	7	Flat	Full	
072212	n/a	1/4	7	Flat	Full	1/4
072216	n/a	1/4	7	Flat	3/4	1/4
072213	n/a	3/8	10	Flat	3/4	
072214	n/a	3/8	10	Flat	Full	
072215	n/a	3/8	10	Flat	Full	3/8
072217	n/a	3/8	10	Flat	3/4	3/8

# Wound Drains

## Silicone Drains - Round

Perforated silicone round drains are available with and without trocars in an assortment of sizes. Drains feature a radiopaque stripe for x-ray detection. Latex Free, Single-use, Sterile. 10/case.

Part Number	Shape	Size (in)	Hole Pattern	Perforation	Trocar
0070200	Round	3/32	10 in	End	
0070210	Round	1/8	10 in	End	
0070220	Round	3/16	10 in	End	
0070230	Round	1/4	10 in	End	
0070310	Round	1/8	10 in	End	Yes
0070320	Round	3/16	10 in	End	Yes
0070330	Round	1/4	10 in	End	Yes
0070710	Double Round	1/8	12.25 in	Mid	Yes
0070720	Double Round	3/16	12.25 in	Mid	Yes
0070730	Double Round	1/4	12.25 in	Mid	Yes



## Silicone Drains - Flat

Hubless silicone flat drains are available in 7mm and 10mm sizes with 3/4 and full perforations. Drains are radiopaque for x-ray detection. Latex Free, Single-use, sterile. 10/case.

Part Number	Shape	Size (mm)	Perforation	Trocar
0070370	Flat- Hubless	10	3/4	No
0070430	Flat- Hubless	7	Full	No
0070440	Flat- Hubless	10	Full	No



## PVC Drains

Siliconized PVC round drains with 12 in. hole pattern and radiopaque stripe. Available with or without trocar in an assortment of sizes. Latex Free, Single-use, sterile. 10/case.

Part Number	Shape	Size (in)	Perforation	Hole Pattern	Trocar
0034660	Double Round	1/8	Mid	12	No
0034670	Double Round	1/4	Mid	12	No
0034680	Double Round	3/16	Mid	12	No
0080700	Double Round	3/32	Mid	12	Yes
0080710	Double Round	1/8	Mid	12	Yes
0080720	Double Round	3/16	Mid	12	Yes
0080730	Double Round	1/4	Mid	12	Yes



# Evacuators

## Silicone Bulb Evacuator Kits

100cc evacuators offer lightweight, low level suction. Available with round or flat drain, with or without trocar. Latex Free, Single-use, sterile. 20/case.

Part Number	Evacuator	Drain				Trocar
		Drain Shape	Length	Width	Perforation	
0071360	100cc	Flat Hubless	20 cm	7 mm	3/4	No
0071370	100cc	Flat Hubless	20 cm	10 mm	3/4	No
0071430	100cc	Flat Hubless	20 cm	7 mm	Full	No
0071440	100cc	Flat Hubless	20 cm	10 mm	Full	No
0072360	100cc	Flat Hubless	20 cm	7 mm	3/4	Yes
0072370	100cc	Flat Hubless	20 cm	10 mm	3/4	Yes
0072430	100cc	Flat Hubless	20 cm	7 mm	Full	Yes
0072440	100cc	Flat Hubless	20 cm	10 mm	Full	Yes
0073310	100cc	Round		1/8 in		Yes
0073320	100cc	Round		3/16 in		Yes



## RELIAVAC® Evacuator Kits (100cc)

100cc RELIAVAC® PVC evacuators offer a more constant suction profile than silicone evacuators. Bottom emptying port reduces chance of contamination. Available in an assortment of sizes, round or flat drain, with or without trocar. Latex Free, Single-use, sterile. 20/case.

Part Number	Evacuator			Drain			Y Connector	Trocar
	Size	Suction	Material	Shape	Size	Perforation		
0085310	100cc	Medium (120-140mm Hg)	Silicone	Round	1/8 in		Yes	Yes
0085320	100cc	Medium (120-140mm Hg)	Silicone	Round	3/16 in		Yes	Yes
0085360	100cc	Medium (120-140mm Hg)	Silicone	Flat	7 mm	3/4	Yes	No
0085370	100cc	Medium (120-140mm Hg)	Silicone	Flat	10 mm	3/4	Yes	No
0085380	100cc	Medium (120-140mm Hg)	Silicone	Hubless Flat	7 mm	Full	Yes	No
0085390	100cc	Medium (120-140mm Hg)	Silicone	Hubless Flat	10 mm	Full	Yes	No
0090360	100cc	Low (40-80 mm Hg)	Silicone	Flat	10 mm	Full	Yes	No
0090370	100cc	Low (40-80 mm Hg)	Silicone	Flat	10 mm	Full	Yes	No
0090380	100cc	Low (40-80 mm Hg)	Silicone	Flat	10 mm	Full	Yes	No
0090390	100cc	Low (40-80 mm Hg)	Silicone	Flat	10 mm	Full	Yes	No



## Three-Spring Evacuator Kits

400cc Three-spring evacuators provide reliable suction for surgical cases with larger amounts of drainage. Latex-Free, Single-use, sterile. 10/case.

Part Number	Evacuator	Drain Material	Drain Size (in)	Hole Pattern	Y Connector	Trocar
0043600	400cc	PVC	3/32	12 in	Yes	Yes
0043610	400cc	PVC	1/8	12 in	Yes	Yes
0043620	400cc	PVC	1/4	12 in	Yes	Yes
0043630	400cc	PVC	3/16	12 in	Yes	Yes
0043660	400cc	Silicone	1/8	10 in	Yes	Yes
0043670	400cc	Silicone	1/4	10 in	Yes	Yes



## RELIAVAC® Evacuator Kits

Offers a more constant suction profile than three-spring evacuators, resulting in the potential for improved tissue approximation and reduced risk of hematoma. 400cc RELIAVAC® kits are available in an assortment of sizes, with or without trocar. Single-use, sterile. 10/case.

Part Number	Evacuator	Drain Material	Shape	Size (in)	Hole Pattern	Y Connector	Trocar
0034600	400cc	PVC	Round	3/32		Yes	Yes
0034610	400cc	PVC	Round	1/8		Yes	Yes
0034620	400cc	PVC	Round	1/4		Yes	Yes
0034630	400cc	PVC	Round	3/16		Yes	Yes
0070480	400cc	Silicone	Round	1/8	10 in	Yes	Yes
0070490	400cc	Silicone	Round	3/16	10 in	Yes	Yes
0070510	400cc	Silicone	Round	1/4	10 in	Yes	Yes
0070520	400cc	Silicone	Flat	1/4	5.4 in	Yes	No
0070610	400cc	Silicone	Flat	1/4	5.4 in	Yes	No




## RELIAVAC® Evacuators

Part Number	Description	Pkg
0034830	Evacuator only (400cc)	10/cs
0070800	Evacuator with Y-Connector (400cc)	10/cs
0070850	Bulb Evacuator with Y-Connector (100cc)- Medium	20/cs
0070900	Bulb Evacuator with Y-Connector (100cc)- Low	20/cs

## Evacuators Only

Offers lightweight low level suction and is easy for a patient to carry.

Part Number	Evacuator	Material	Shape	Pkg
0070740	100cc	Silicone	Bulb	20/cs
0071200	200cc	Silicone	Bulb	10/cs
0070400	400cc	Silicone	Bulb	10/cs
0043650	400cc	PVC	3 Spring	10/cs



# Gravity and Sump Drains

## Sump Drains

Abramson Triple-Lumen Sump Drains are used in the peritoneal cavity to remove large volumes of thick exudates. The triple-lumen sheath allows an air vent, suction lumen and irrigation/medication lumen. The sheath that can be left in place as a gravity drain. It has an anti-bacterial filter that allows air to enter without contamination. Latex Free, Single-use, sterile. 5/case.

Part Number	Material	Shape	Length (in)	Width (in)
0080400	Silicone	3 Lumen	15	5/8
0080420	Silicone	3 Lumen	15	3/8
0080450	Silicone	3 Lumen with Sheath (.67" x 15")	15	3/8



## Sump Drain Accessories

Part Number	Component	Pkg
0080500	Filter	10/cs

## Gravity Drains

- Penrose Drains – Used to promote drainage in an open surgical wound.

Part Number	Shape	Length (in)	Width (in)	Pkg/case	
0912010	Flat	12	1/4	200/cs	Sterile
0912020	Flat	12	3/8	200/cs	Sterile
0912030	Flat	12	1/2	200/cs	Sterile
0912040	Flat	12	5/8	200/cs	Sterile
0912050	Flat	12	3/4	200/cs	Sterile
0912060	Flat	12	7/8	200/cs	Sterile
0912070	Flat	12	1	200/cs	Sterile
0918010	Flat	18	1/4	200/cs	Sterile
0918020	Flat	18	3/8	200/cs	Sterile
0918030	Flat	18	1/2	200/cs	Sterile
0918040	Flat	18	5/8	200/cs	Sterile
0918050	Flat	18	3/4	200/cs	Sterile
0918060	Flat	18	7/8	200/cs	Sterile
0918070	Flat	18	1	200/cs	Sterile
0097870	Flat	36	1/4	144/cs	Sterile
0097880	Flat	36	5/16	144/cs	Non-Sterile
0097890	Flat	36	3/8	144/cs	Non-Sterile
0097900	Flat	36	1/2	144/cs	Non-Sterile
0097910	Flat	36	5/8	144/cs	Non-Sterile
0097920	Flat	36	3/4	144/cs	Non-Sterile
0097930	Flat	36	7/8	144/cs	Non-Sterile
0097940	Flat	36	1	144/cs	Non-Sterile



For other types of specialty gravity drains refer to the Bard Medical Division at [www.bardmedical.com](http://www.bardmedical.com).

## Accessories

### Adapters

Part Number	Fits	Fits	Pkg
0034850	Reliavac Suction		10/cs
9019500	Silicone	1/8 Round	10/cs
9019510	PVC/Silicone Round/Flat	3/16, 1/4	10/cs
072251	Channel Drains	24/28 Fr	20/cs
072252	Channel Drains	32 Fr	20/cs

### Connectors

Part Number	Type	Connects	Size (in)	Pkg
0034800	Y	PVC	1/8	10/cs
0034810	Y	PVC	1/4	10/cs
0034820	Y	PVC	3/16	10/cs
0070780	Y	Silicone	1/8	10/cs
0070790	Y	Silicone Round	3/16, 1/4	10/cs
0070790	Y	Silicone Flat	10mm, 12.7mm	10/cs
0070770	Universal Claw	PVC and Silicone	Universal	10/cs

### Trocars

Part Number	Size (in)	Pkg
0034760	1/8	10/cs
0034770	1/4	10/cs
0034780	3/16	10/cs



0034800



0070790



0070770



0034760  
0034780  
0034770



# WOUND DRAINAGE

PRODUCTS TO MEET YOUR DIFFERENT PROCEDURAL NEEDS

Closed Wound Suction  
Gravity Drains  
Sump Drains

For more information about the **Wound Drainage Products** and other Bard products, please contact your Bard sales representative, call 800.526.4455, or visit [www.bardmedical.com](http://www.bardmedical.com).

**BARD** | MEDICAL

**C. R. Bard, Inc.** Covington, GA 30014 [www.bardmedical.com](http://www.bardmedical.com)

Bard, Davol and Reliavac are trademarks and/or is a registered trademarks of C. R. Bard, Inc.

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use.

©2011 C. R. Bard, Inc. All Rights Reserved.  
1009-62 R03/11 THP P04/11 1M