



Stax Finger Splints

Sturdy, lightweight splints that protect and support fingertip and nail bed injuries.

Number
69990000

HARTMANN

Stax Finger Splints

Support & Protection

Starter Kit

4-size 1	4-size 5
4-size 2	4-size 5 1/2
4-size 3	4-size 6
4-size 4	2-size 7

30 Splints Per Kit

HARTMANN USA, Inc. Rock Hill, SC 29730
1-800-243-2294 www.hartmanninfo.com






helps healing.

HARTMANN USA, Inc.
Rock Hill, SC 29730
1-800-243-2294
www.hartmanninfo.com

Stax Finger Splints

Latex-Free Support and Protection



Features

- Composed of sturdy plastic for support and protection
- Lightweight, perforated and shaped
- Adhesive strip elevates distal joint
- Made to withstand wet or soiled conditions

Benefits

- Strong support for distal joint, while permitting unrestricted movement of proximal interphalangeal joint
- Designed for patient comfort
- Effective alignment of injured finger to proper extension
- Quick and easy to apply, and can be changed easily by patient

Stax Finger Splints

Order No.	Size	Splints Per Pack
69100000	Size 1	6
69200000	Size 2	6
69300000	Size 3	6
69400000	Size 4	6
69500000	Size 5	6
69550000	Size 5.5	6
69600000	Size 6	6
69700000	Size 7	6
69990000	Starter Kit With Assorted Sizes (4 each size 1-6, 2 each size 7)	30

Applications

- Normally used without padding, and is applied to middle phalanx
- For use in treating nail bed injuries and baseball finger
- In cases where the finger is too tender, a larger splint with padding on the inside should be used until it can be replaced with an appropriate sized Stax Finger Splint