

BLICKMAN

SPEC SHEET

PP Bouffant Cover



Product Description

Bouffant covers should be worn to cover all of the hair, scalp, and ears to minimize skin and hair shedding and to protect patients from bacteria. Bouffant covers are not to be reused. Available in three sizes.

PRODUCT SPECIFICATIONS

Intended Use	Cover all of the hair, scalp, and ears in medical settings.
Material	Spunbound polypropylene (PP), single layer of nonwoven fabric, elastic band.
Liquid Resistance	Meant for minimal fluid exposure
Sizes Available	Regular, 64x13cm (R); Large, 64x14cm(L); and Extra-Large, 64x15cm (XL)
Single-use Only	Yes
Warranty	Blickman's limited lifetime warranty
QTY/Box	100 bouffant covers per box