



PASSION. PARTNERSHIP. POSSIBILITIES.

# Composite Island Dressings



## Description

A convenient multiple layer dressing which provides essential functions such as absorption, adhesion, and bacterial barrier. It allows the exchange of moisture vapor and can be used as a primary or secondary dressing.

## Features & Benefits

- Non-adherent contact layer that covers the wound site
- Soft absorbent pad that collects exudate and protects the wound from trauma
- Non-woven backing with adhesive border that is gentle and holds pad in place
- Waterproof backing to keep wounds dry and keep external contaminants out

## Common Applications:

- Light to moderately draining wounds
- Pressure ulcers
- Partial and full thickness wounds
- Incision sites
- Venous Ulcers
- Lacerations and abrasions

Item#	Description	Packaging
4080	4" x 10"	25/bx, 4 bx/cs
4081	4" x 14"	25/bx, 2 bx/cs
4082	4" x 4"	25/bx, 8 bx/cs
4083	6" x 6"	25/bx, 4 bx/cs
4085	4" x 6"	25/bx, 4 bx/cs
4086	6" x 8"	25/bx, 4 bx/cs

## Connect with us!

Find our full line of products online at:  
[www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)

To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779

