

## *High Frequency Desiccator*



**Item Number: 570-0583**

40 watt high frequency desiccator

Less to do and less to buy makes this the most cost effective and user friendly unit on the market

Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar

More accurate power with Universal Power Control 100-240 VAC

The Henry Schein Desiccator is a ground reference unit which means no grounding pad is required.



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**



**STANDARD EQUIPMENT**

- 570-0583 (1) Henry Schein High Frequency Generator
- 895-0124 (1) 3-Button Autoclavable Handpiece
- 112-7171 (5) non-sterile sharp dermal tips
- 112-7172 (2) sterile sharp dermal tips
- 112-7173 (5) non-sterile blunt dermal tips
- 112-7175 (2) sterile blunt dermal tips
- 895-0127 (2) non-sterile DERM-Elite™ Premium blunt dermal tips
- 895-0126 (2) sterile DERM-Elite Premium blunt dermal tips
- 112-7176 (2) sterile disposable handpiece sheath
- 570-0331 (2) non-sterile disposable handpiece sheath
- 926-3686 (1) wall mount kit
- 110 VAC Hospital-grade power cord\*
- User's guide\*



895-0124



112-7171, 112-7172, 112-7173,  
112-7175, 895-0127, 895-0126



926-3686



112-7176



570-0331

**COMPLEMENTARY PRODUCTS †**


† Sold Separately



102-7058  
Table Top Stand †



600-5943  
Footswitch †



5700888, 5700890  
Sharp H-Type Tips †



5700891, 5700889  
Blunt H-Type Tips †



118-3875  
Force Cable †




147-2819  
Reusable Bipolar Forceps †



682-5062  
Ball Electrode †



650-3112  
Needle Electrode †



329-9379  
Top Tray for Mobile Stand †



330-8972  
Mobile Stand †

**PRODUCT SPECIFICATIONS**

**OUTPUT POWER**  
TO 40 WATTS (0-10 WATTS IN  
1/10TH WATT INCREMENTS  
ON BOTH MONOPOLAR  
AND BIPOLAR)

**DIMENSIONS**  
HEIGHT: 188MM (7.4 IN.)  
WIDTH: 228MM (8.98 IN.)  
DEPTH: 105MM (4.13 IN.)  
WEIGHT: < 2.26 KG (< 5 LBS.)

**OUTPUT FREQUENCY**  
328 – 492KHZ

**LINE VOLTAGE**  
570-0583: 100-240 VAC

**LINE FREQUENCY**  
50 - 60 HZ

**WARRANTY**  
FOUR YEARS



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**