

Instrument Lubricant

Instrument Lubricant contains an all-natural corrosion inhibitor and lubricant specially formulated to lubricate and protect surgical instruments, hinges, box locks and general equipment during autoclaving



RECOMMENDED USES:

- Medical/surgical instruments

ADVANTAGES:

- Contains an anticorrosive agent to help prevent corrosion
- The gentle formula promotes ease of use
- Extends instrument life by leaving a protective coating
- Bacteriostatic

TECHNICAL DATA:

- pH 6.5 to 7.5
- Non-corrosive
- Milk Bath Dilution 1 part solution to 6 parts water
- Non-Silicone
- Automated Equipment 1/8 to 1/2 ounce per gallon of water

PACKAGING:

- 4x1 Gallon-2000980
- 2x2.5 Gallons-2000981
- 5 Gallons-2000990
- 15 Gallons-2000995

ENZIME SOLUTIONS®
INCORPORATED

2105 Forrest Park Drive
Garrett, IN 46738
800-523-1323 Phone • 260.497.0512 Fax
www.enzymesolutions.com