

Medical Enzyme Detergent

Medical Enzyme Detergent is a low foaming multi-enzyme detergent presoak cleaner for medical instruments and apparatuses. It is a blend of enzymes that are effective in the removal of all proteins, blood, fatty tissue, etc. from all medical instruments, equipment and supplies.

Medical Enzyme Detergent will clean surgical instruments, endoscopes, respiratory therapy equipment and glassware. It will also remove stains in surgical scrubs, linens and medical clothing.

Medical Enzyme Detergent is also available in a 32 ounce Gel or Foam Transport Spray.



ENZYME SOLUTIONS®
INCORPORATED

2105 Forrest Park Drive • Garrett, IN 46738
800-523-1323 Phone • 260.553.9104 Fax
www.enzymesolutions.com

RECOMMENDED USES:

- General Cleaning
- Medical/surgical instruments
- Endoscopes
- Respiratory Therapy Equipment
- Evacuation Systems

ADVANTAGES:

- Provides quick and thorough cleaning in manual, ultrasonic and automatic cleaning applications
- pH neutral-instrument and equipment safe
- Multi-enzyme formula effectively works within minutes to disintegrate and carry away organic soils

TECHNICAL DATA:

- pH neutral
- 100 % bio-degradable
- Nonflammable
- Non-toxic
- Non-Caustic
- Dilution 1/2 to 1 oz. per gallon of water

PACKAGING:

- 12x32 oz. Ready To use Gel Spray
- 12x32 oz. Ready To Use Foam Spray
- 4x1 Gallon • 2x2.5 Gallons • 5 Gallons
- 15 Gallons • 30 Gallons

