

# Medical Enzyme Detergent PLUS

## With high level biofilm removal

Medical Enzyme Detergent PLUS is a low foaming multi-enzyme detergent presoak cleaner for medical instruments and apparatuses. It has double the concentration of enzymes than the regular concentrated formula for cleaning heavily soiled and stained instruments and supplies.

Medical Enzyme Detergent PLUS with high level biofilm removal is a proprietary blend of enzymes effective against proteins, lipids, carbohydrates, mucopolysaccharides and biofilm. It will clean surgical instruments, endoscopes, respiratory therapy equipment and glassware and will also remove stains in surgical scrubs, linens and medical clothing.



### RECOMMENDED USES:

- Medical/surgical instruments
- Fiberoptics
- Endoscopes
- Respiratory Therapy Equipment
- Evacuation Systems

### ADVANTAGES:

- Provides quick and thorough cleaning in manual, ultrasonic and automatic cleaning applications
- pH neutral-instrument and equipment safe
- Multi-enzyme formula effectively works within minutes to disintegrate and carry away organic soils

### TECHNICAL DATA:

- pH neutral
- 100 % bio-degradable
- Nonflammable
- Non-toxic
- Non-Caustic
- Dilution 1/2 to 1 oz. per gallon of water

### PACKAGING:

- 4x1 Gallon
- 2x2.5 Gallons
- 5 Gallons
- 15 Gallons
- 30 Gallons

**ENZYME SOLUTIONS®**  
INCORPORATED

2105 Forrest Park Drive  
Garrett, IN 46738  
800-523-1323 Phone • 260.497.0512 Fax  
[www.enzymesolutions.com](http://www.enzymesolutions.com)