



ESSENTIALS 776HS

An economical aneroid sphygmomanometer for today's cost conscious healthcare market.

The Essentials 776HS series features:

- Black enamel 300mmHg CE no-pin stop manometer
- Navy blue nylon cuff with range markings to facilitate selection of correct cuff size; artery label, and gauge holder
- Latex free inflation bladder and bulb
- Chrome-plated deflation valve
- Complete with leatherette zippered carrying case and operating instructions
- Assembled, inspected, and packaged in the U.S.A. from foreign components



PRIVATE LABELING	
COMPONENT	PACKAGING
N/A	1 Pc. Box/Label

Feature	Benefit
Gauge is 100% visually inspected for 0 point, 2% re-inspected for accuracy	To ensure conformance with National Standard ANSI SP 10
Cuff has a properly positioned artery mark at the bladder midpoint	Ensures proper compression of the brachial artery. Most low cost imports fail to meet this requirement from ANSI SP 10
Leak tested Valve	Valves are tested to ensure conformance with ANSI SP 10 for dump rate and leakage
100% inspected, assembled and packaged in the USA	Ensures unmatched consistency, quality, and performance in a value priced BP instrument
5 year calibration, 1 year system warranty	To protect the consumer's investment

ORDERING INFORMATION					
Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
1126076	776CHS	ESSENTIALS 776, Child	Latex Free	ea.	20
1126077	776HS	ESSENTIALS 776, Adult	Latex Free	ea.	20
1126080	776XHS	ESSENTIALS 776, Lrg. Adult	Latex Free	ea.	20