

Fast, Accurate, Efficient Lipid and Glucose Results at the Point of Care

Recognizing the strong connection between heart attack, stroke, and diabetes, the Henry Schein CardioChek® Plus analyzer provides rapid on-site blood chemistry values and measurements to help identify at-risk individuals.

It's one of the fastest, most cost-efficient, and user-friendly analyzers. The Henry Schein CardioChek Plus system accurately tests lipid profile and glucose simultaneously with one fingerstick. In as little as 90 seconds, it provides on-site results using the same technology as clinical laboratories.



Measures:

- Total Cholesterol
- HDL Cholesterol
- Triglycerides
- Glucose

Calculates:

- LDL Cholesterol
- TC/HDL Ratio
- LDL/HDL Ratio
- Non-HDL Cholesterol



Certified Accuracy

- CLIA-waived
- Meets NCEP guidelines for accuracy and precision
- FCC and ISO certified wireless communications
- CRMLN-certified for CHOL and HDL tests
- FDA-cleared

Fast

- On-site lipid and glucose results in as little as 90 seconds
- Simple quality control testing

Portable

- Compact and lightweight, easy to transport and store
- Handheld, battery-operated
- Room temperature test strips

Patient-Friendly

- Fingerstick with small sample

Wireless Communication

- Easy data transfer

High ROI

- Rapid throughput and low waste, with minimal start-up cost
- Fully reimbursable



Simultaneous testing of lipid profile and glucose with Henry Schein PTS Panels® test strips.