

Digital Thermometers

From A Name You Can Trust.



Item Number	Item Description	Size
900-4856	Digital Thermometer	Ea
100-9361	Probe Covers for Digital Thermometers	100/Bx

Features & Benefits

- Records temperature in approximately 60 seconds
- Beeps when reading is completed
- Waterproof
- Memory recall of last reading
- Comes with 5 Probe Covers
- Auto-off function conserves battery life
- Compact size - under 5" long, weighs 11g
- Replaceable 1.55V (LR41) battery provides approximately 1,500 measurements
- Range: 90°F - 109.9°F +/-0.2°F
32°C - 43.9°C +/-0.1°C



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!