

# Ankle Brace Stirrup

## From A Name You Can Trust.

*Our Ankle Brace Stirrup contains all the features necessary to ensure your patients heal with a proper alignment, for a quick and complete recovery, while having the best patient experience available.*



Ankle Brace Stirrup	
Item Number	Size
(570-1016)	Universal

### Features & Benefits:

*Designed and manufactured to provide strength, protection, and value.*

- Provides medial and lateral stability while allowing plantar and dorsiflexion
- Air chambers decrease edema and help break the paincycle
- Anatomically contoured semi-rigid shells conform closely to the ankle with minimal bulk in the shoe
- Adjustable heel strap ensures comfortable and functional fit
- Supports injured and unstable ankles



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**