

# Scales & Stadiometer

## From A Name You Can Trust.

*Henry Schein Brand Scales and Stadiometer are a perfect fit for any practice setting!*



### Digital Waist High Scale

Item Number: 112-7155

#### Features & Benefits:

- Includes integrated mechanical height rod
- Large backlit LCD display with a capacity of up to 500 lbs (225 kg)
- Low profile for easy access, platform is only 2.5" high
- Contains BMI Calculator and a clock
- 2 year Limited Warranty
- 22" x 15" x 36"




### Eye Height Digital Scale

Item Number: 112-7154

#### Features & Benefits:

- Includes integrated mechanical height rod
- Large backlit LCD display with a capacity of up to 500 lbs (225 kg)
- Low profile for easy access, platform is only 2.5" high
- Contains BMI Calculator and a clock
- 2 year Limited Warranty
- 22" x 15" x 53"



### Small Countertop Pediatric Scale

Item Number: 112-7156

#### Features & Benefits:

- Curved design to support and protect the infant
- Sta-bilt stabilizing software ensures accurate measurement
- Large backlit LCD display measuring in lbs and ounces
- Weight capacity of 45 lbs
- 2 year Limited Warranty
- 16.75" x 27.75" x 4"



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**

# Scales & Stadiometer

## From A Name You Can Trust.

*Henry Schein Brand Scales and Stadiometer are a perfect fit for any practice setting!*



### Standing Mechanical Beam Scale

Item Number: 112-7157

#### Features & Benefits:

- Includes built in height rod
- Weight capacity of 450 lbs (200 kg)
- Low profile for easy access, platform is only 2.5" high
- 1 year Limited Warranty
- 22.25" x 18.25" x 58.25"




### Digital Remote Indicator Scale

Item Number: 570-0596

#### Features & Benefits:

- 1.5" Backlit Remote Display w/ 6' Long Cable Connection to Base
- Stabilizing Software ensures accurate weight even when patient moves
- Auto-Hold Feature Shows Weight Even After Patient Leaves Scale
- 15" x 14.5" x 2.3", 500 lb. Max Capacity
- 2 Year Limited Warranty




### Digital Floor Scale

Item Number: 112-7153

#### Features & Benefits:

- Large backlit LCD display
- Weighs accurately in both lbs and kg, up to 400 lbs (200 kg)
- Features anti-slip feet and a low profile for added safety
- Automatically turns on
- 1 year Limited Warranty
- 14" x 11.625" x 1"



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**

# Scales & Stadiometer

## From A Name You Can Trust.

*Henry Schein Brand Scales and Stadiometer are a perfect fit for any practice setting!*



### Digital Wheelchair Scale

Item Number: 570-0597

#### Features & Benefits:

- Single, stable platform allows scale access with Wheelchair
- Easily Folded to Minimize Footprint and moved with casters
- Individual Wheelchair weight can be programmed for easy recall
- 2 Year Limited Warranty
- 32" x 33" x 2.5", 1000 lb. Max Capacity



### Digital Handrail Scale w/ Height Rod

Item Number: 570-0598

#### Features & Benefits:

- Low Profile Design for Easy Access, Platform is Only 2.75" High
- Contains BMI Calculator and Clock
- 28" x 23.5" x 2.75", 1000 lb. Max Capacity
- 2 Year Limited Warranty
- Handrails Allow Patients to Hold Steady During Weighing



### Stadiometer

Item Number: 570-0606

#### Features & Benefits:

- Easy Installation on Wall
- Mounting hardware included
- Height range from 2' to 6' 6.75"



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**