

Hextio Desktop Air Purifier by Radic8

Designed for Reception Desks and Small Offices

NOW
AVAILABLE
through
Henry Schein



Radic8 Viruskiller[™] technology was funded in 2004 by the Korean Government in response to SARS. The Viruskiller[™] technology which combines air purification and photocatalytic oxidation is now in over 80% of all hospitals and schools in South Korea with over 370,000 installations.

The Viruskiller[™] Technology has a **99.9999% kill rate—in a single pass—against viruses** similar to SARS-CoV-2 or harder-to-kill, such as Coronavirus DF2, Adenovirus, Influenza virus, and Poliovirus. The patented Hextio chamber utilizes ultimate reflection and refraction to amplify the power of the UV lamp by 10 times. Even though Hextio is small in size, it is big in capability.

Viruskiller[™] technology has 14 patents granted in South Korea.

Purification Methods

- HEPA filter
- Carbon filter
- UVGI
- Photocatalytic oxidation
- 1 UVC lamp
- 10 chromed TiO₂ hexagon filters
- Covers spaces up to 215 square feet*
- 40–49 dB sound level
- 5" x 4" x 13"
- 2.6 lbs
- Freestanding or wall-hung option
- Portable and can be easily transported

Top Innovation Company Award 2018 for Infection Control as determined by InfectionControl.tips

Available Now from Henry Schein



*Assuming 8' ceiling.

This product is not FDA cleared or approved. Marketed under FDA's Enforcement Policy for Air Purifiers, during the COVID-19 public health emergency. Tested on viruses that are representative of SARS-CoV-2, but not SARS-CoV-2 itself. This product has not been tested on SARS-CoV-2. This product is for medical purposes and should be used as an adjunct to personal protective equipment (PPE) and other infection control procedures. When using Viruskiller[™], basic safety precautions should be observed, and instructions should be read to reduce the risk of fire, electric shock, or injury; the appliance should not be used by children or adults with reduced capabilities. Looking directly at an ultra-violet light may cause permanent damage to the eyes. Please visit <https://henryscheindigital.com/radic8>

© 2020 Henry Schein, Inc. All rights reserved. No copying without permission. Not responsible for typographical errors. 20BS8613