

ESSENTIALS™

Henry Schein Essentials combine lightweight comfort, sensitive acoustics, and handsome styling, with truly affordable pricing.

The ESSENTIALS Series features:


- Choice of diaphragm only or combination chestpiece
- Lightweight aluminum chestpiece construction
- Non-chill diaphragm rim for patient comfort
- Proprietary diaphragm enhances acoustic response by 50% over traditional designs
- Adjustable aluminum binaurals with PVC eartips
- Flexible 22" PVC tubing
- Spare pair of black mushroom eartips included
- Inspected and packaged in the U.S.A.
- One Year Warranty

ESSENTIALS 660 series features:

- Non-chill diaphragm only nurse style chestpiece (1³/₄" dia.)
- Ideally suited for blood pressure measurement
- Weighs 2.5 oz. / Overall length 32"

ESSENTIALS 670 series features:

- Combination diaphragm/bell chest piece for wide frequency response
- Diaphragm diameter 1³/₄", bell diameter 3³/₄"
- Non-chill bell ring and diaphragm for patient comfort
- Weighs 3 oz. / Overall length 32.5"



ESSENTIALS 675, Pediatric Dual-head
ESSENTIALS 676, Infant Dual-head

PRIVATE LABELING	
COMPONENT	PACKAGING
N/A	1 Pc. Box/Label

Feature	Benefit
660, 670, 675, 676 Exclusive component sourcing	Unsurpassed quality and consistency
660, 670, 675, 676 Aluminum chestpiece and binaurals	Outstanding acoustics in an ultra-lightweight design
660, 670 Proprietary diaphragm/rim design	Improves acoustic performance by up to 50% over competing brands with threaded rims
660, 670, 675, 676 Non chill rim	For patient comfort
660, 670, 675, 676 Binaural up to 1/2" longer than competing brands	More practitioner comfort
660, 670, 675, 676 Y tubing maintained at tight specs	For greater consistency
660, 670, 675, 676 Spare pair of mushroom eartips	For greater comfort
660, 670, 675, 676 Inspected and packaged in the USA	For more consistent quality
660, 670, 675, 676 1 year warranty	To protect the consumer's investment

ORDERING INFORMATION

Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
9004818	660BKHS	ESSENTIALS Nursescope, Black	Latex-Free	ea.	50
9004819	660GHS	ESSENTIALS Nursescope, Gray	Latex-Free	ea.	50
9004820	660RBHS	ESSENTIALS Nursescope, Royal Blue	Latex-Free	ea.	50
9004821	660RHS	ESSENTIALS Nursescope, Red	Latex-Free	ea.	50
9004809	670BKHS	ESSENTIALS Dualhead, Black	Latex-Free	ea.	50

ORDERING INFORMATION

Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
9004810	670GHS	ESSENTIALS Dualhead, Gray	Latex-Free	ea.	50
9004811	670RBHS	ESSENTIALS Dualhead, Royal Blue	Latex-Free	ea.	50
9004812	670RHS	ESSENTIALS Dualhead, Red	Latex-Free	ea.	50
9004813	675BKHS	ESSENTIALS Pediatric, Black	Latex-Free	ea.	50
9004814	676GHS	ESSENTIALS Infant, Gray	Latex-Free	ea.	50