

**One less worry ...**



**... when it comes to choosing  
KIMVENT\* Closed Suction Systems.**

## Closed Suction Systems / Accessories

### Paediatric Endotracheal Tubes

#### ***KIMBERLY-CLARK\* KIMVENT\* Multi-Access Port Closed Suction System***


- Rotating Manifold
- Single Use Cassette
- Turbulent Cleaning Chamber
- Sleeve Tether
- MDI Drug Delivery
- Replacement Catheters

<b>Code</b>	<b>Description</b>	<b>Packaging</b>
240128-5	KIMVENT* Multi-Access Port Closed Suction System For Adults, Endotracheal, MDI, 12 F	20/Case
240148-5	KIMVENT* Multi-Access Port Closed Suction System For Adults, Endotracheal, MDI, 14 F	20/Case
2401238-5	KIMVENT* Multi-Access Port Closed Suction System For Adults, Tracheostomy, DSY, MDI, 12 F	20/Case
2401438-5	KIMVENT* Multi-Access Port Closed Suction System For Adults, Tracheostomy, DSY, MDI, 14 F	20/Case
240125-5	KIMVENT* Multi-Access Port Replacement Catheter For Adults, Endotracheal, 12 F	20/Case
240145-5	KIMVENT* Multi-Access Port Replacement Catheter For Adults, Endotracheal, 14 F	20/Case
2401253-5	KIMVENT* Multi-Access Port Replacement Catheter For Adults, Tracheostomy, 12 F	20/Case
2401453-5	KIMVENT* Multi-Access Port Replacement Catheter For Adults, Tracheostomy, 14 F	20/Case
2401286-5	KIMVENT* WET PAK* Multi-Access Port Closed Suction System For Adults, Endotracheal, DSY, MDI, 12 F	20/Case
2401486-5	KIMVENT* WET PAK* Multi-Access Port Closed Suction System For Adults, Endotracheal, DSY, MDI, 14 F	20/Case
24013-5	KIMVENT* Seal Cassette	20/Case

## KIMVENT\* Closed Suction Systems

### **KIMBERLY-CLARK\* Adult TRACH CARE\* 24 Hour Double Swivel Elbows**

- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility Without Connectors
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes




Code	Outer Diameter	Length	Description	Packaging
2160-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow	20/Case
21603-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
21609-5	12 Fr / 4mm	22.2" / 56cm	Double Swivel Elbow Directional Tip	20/Case
2210-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow	20/Case
22103-5	14 Fr / 4.6mm	21.3" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
22107-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Dual Lumen	20/Case
221037-5	14 Fr / 4.6mm	21.3" / 30.5cm	Double Swivel Elbow Tracheostomy Dual Lumen	20/Case
221079-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Dual Lumen Directional Tip	20/Case
22108-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow MDI Built-In	20/Case
221038-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy MDI Built-In	20/Case
22109-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Directional Tip	20/Case
2260-5	16 Fr / 5.3mm	21.3" / 54cm	Double Swivel Elbow	20/Case

### **KIMBERLY-CLARK\* Adult TRACH CARE\* 24 Hour Double Swivel Elbow WET PAK\***

Code	Outer Diameter	Length	Description	Packaging
21606-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline	20/Case
216069-5	12 Fr / 4mm	22.2" / 56cm	Double Swivel Elbow WET PAK* Directional Tip, includes 12 x 15ml Vials Saline	20/Case
216036-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline Tracheostomy	20/Case
22106-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow WET PAK*	20/Case
221069-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow WET PAK* Directional Tip, includes 12 x 15ml Vials Saline	20/Case
221036-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow WET PAK* Tracheostomy	20/Case
22606-5	16 Fr / 5.3mm	21.3" / 54cm	Double Swivel Elbow WET PAK*	20/Case

**KIMBERLY-CLARK® Adult TRACH CARE®  
24 Hour T-Pieces**


- Gamma Sterilised
- Lavage / Irrigation Port
- Locking Thumb Port
- T-Piece
- Multiple French Sizes



Code	Outer Diameter	Length	Description	Packaging
2105-5	10 Fr / 3.3mm	21.3" / 54cm	Clear T-Piece	20/Case
2155-5	12 Fr / 4mm	21.3" / 54cm	Clear T-Piece	20/Case
2205-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece	20/Case
220135-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece Tracheostomy Swivel Adapter	20/Case
2255-5	16 Fr / 5.3mm	21.3" / 54cm	Clear T-Piece	20/Case
25-2210-5	14 Fr / 4.6mm	21.3" / 54cm	24hr Double Swivel Elbow With Safety Lavage Port	20/Case

**KIMBERLY-CLARK® Adult TRACH CARE®  
24 Hour T-Piece WET PAK®**

- Available in T-Piece, Elbow and Double Swivel Elbow
- Lavage / Irrigation Port
- Catheters Packaged Separate from Vials
- 12 Vials Per Catheter
- Vials Designed for Respiratory Use




Code	Outer Diameter	Length	Description	Packaging
22056-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece WET PAK	20/Case
220568-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece WET PAK MDI Adapter	20/Case
22012356-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece WET PAK Tracheostomy Swivel Adapter 15mm x 22mm Adapter	20/Case
220123568-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece WET PAK Tracheostomy Swivel Adapter 15mm x 22mm Adapter MDI Adapter	20/Case
22556-5	16 Fr / 5.3mm	21.3" / 54cm	Clear T-Piece WET PAK	20/Case

## Turbo Adult Closed Suction Systems

### **KIMBERLY-CLARK® Adult TRACH CARE® 72 Hour Double Swivel Elbows**

- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility Without Connectors
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes




2160, 21603, 2210

Code	Outer Diameter	Length	Description	Packaging
227-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Endotracheal	20/Case
22703-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
22716-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow Endotracheal	20/Case
2271603-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case

### **KIMBERLY-CLARK® Adult TRACH CARE® 72 Hour Double Swivel Elbows WET PAK®**

Code	Outer Diameter	Length	Description	Packaging
227166-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow WET PAK®, includes 12 x 15ml Vials Saline	20/Case
22716036-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow WET PAK®, includes 12 x 15ml Vials Saline Tracheostomy	20/Case
2276-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow WET PAK®, includes 12 x 15ml Vials Saline	20/Case
227036-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow WET PAK®, includes 12 x 15ml Vials Saline Tracheostomy	20/Case

### **KIMBERLY-CLARK® Adult TRACH CARE® 72 Hour T-Pieces**



Code	Outer Diameter	Length	Description	Packaging
22705-5	14 Fr / 4.6mm	21.3" / 54cm	T-Piece Endotracheal	20/Case
2270135-5	14 Fr / 4.6mm	12" / 30.5cm	T-Piece Tracheostomy	20/Case
22715-5	12 Fr / 4mm	21.3" / 54cm	T-Piece Endotracheal	20/Case
227135-5	12 Fr / 4mm	12" / 30.5cm	T-Piece Tracheostomy	20/Case

**KIMBERLY-CLARK® Adult TRACH CARE®  
72 Hour T-Piece WET PAK®**

- Available in T-Piece, Elbow and Double Swivel Elbow
- Lavage / Irrigation Port
- Catheters Packaged Separate from Vials
- 12 Vials Per Catheter
- Vials Designed for Respiratory Use




<b>Code</b>	<b>Outer Diameter</b>	<b>Length</b>	<b>Description</b>	<b>Packaging</b>
227156-5	12 Fr / 4.0mm	21.3" / 54cm	Adult WET PAK* T-Piece Standard	20/Case
2271356-5	12 Fr / 4.0mm	12" / 30.5cm	Adult WET PAK* T-Piece Tracheostomy	20/Case
227056-5	14 Fr / 4.6mm	21.3" / 54cm	Adult WET PAK* T-Piece Standard	20/Case
22701356-5	14 Fr / 4.6mm	12" / 30.5cm	Adult WET PAK* T-Piece Tracheostomy	20/Case

## Specialty Closed Suction Systems

### **KIMBERLY-CLARK® Adult TRACH CARE® 24 Hour Dual Lumens**

- Gamma Sterilised
- Dual Lavage / Irrigation Ports
- Available in T-Piece and Double Swivel Elbow
- Available in Directional Tip
- Allows Lavage to Take Place at the Patient End of the Catheter



Code	Outer Diameter	Length	Description	Packaging
22107-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Dual Lumen	20/Case
221037-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Dual Lumen Tracheostomy	20/Case
221079-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Dual Lumen Directional Tip	20/Case

### **KIMBERLY-CLARK® Adult TRACH CARE® 24 Hour Directional Tips**

- Gamma Sterilised
- Lavage / Irrigation Port
- Available in T-Piece, Elbow and Double Swivel Elbow
- Allows for Suctioning of the Right or Left Main Stem Bronchus
- Radiopaque Line Provides for Observation of Catheter Path Under Fluoroscopy
- Colored Center Line Allows Maintaining Manipulation of the Catheter Tip to Access Proper Location of the Right or Left Main Stem Bronchus

Code	Outer Diameter	Length	Description	Packaging
21559-5	12 Fr / 4mm	22.2" / 56cm	Clear T-Piece Directional Tip	20/Case
22059	14 Fr / 4.6mm	22.2" / 56cm	Clear T-Piece Directional Tip	20/Case
220569-5	14 Fr / 4.6mm	22.2" / 56cm	Clear T-Piece WET PAK* Directional Tip	20/Case
22559-5	16 Fr / 5.3mm	22.2" / 56cm	Clear T-Piece Directional Tip	20/Case
22109-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Directional Tip	20/Case
221079-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Directional Tip Dual Lumen	20/Case
225579-5	16 Fr / 5.3mm	22.2" / 56cm	Closed T-Piece Directional Tip	20/Case

## Paediatric / Neonatal Closed Suction Systems

### ***KIMBERLY-CLARK Neonatal and Paediatric TRACH CARE Elbows***



Code	Outer Diameter	Length	Description	Packaging
202-5	5 Fr / 1.66mm	12" / 30.5cm	Elbow Endotracheal Tube Adapters (2mm, 2.5mm)	20/Case
206-5	6 Fr / 2mm	12" / 30.5cm	Elbow Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	20/Case
207-5	7 Fr / 2.3mm	12" / 30.5cm	Elbow Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm, 4mm)	20/Case
208-5	8 Fr / 2.6mm	14" / 35.5cm	Elbow Endotracheal Tube Adapters (3mm, 3.5mm, 4mm)	20/Case
210-5	10 Fr / 3.3mm	16" / 40.5cm	Elbow	20/Case


### ***KIMBERLY-CLARK Neonatal and Paediatric TRACH CARE "Y"***



Code	Outer Diameter	Length	Description	Packaging
195-5	5 Fr / 1.66mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (2.0mm, 2.5mm)	20/Case
196-5	6 Fr / 2mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (2.5mm, 3.0mm, 3.5mm)	20/Case
197-5	7 Fr / 2.3mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (3.0mm, 3.5mm, 4mm)	20/Case
198-5	8 Fr / 2.6mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (3.5mm, 4.0mm, 4.5mm)	20/Case
1910-5	10 Fr / 3.3mm	16" / 40.5cm	"Y" Endotracheal Tube Adapters (4.5mm, 5.0mm, 5.5mm)	20/Case
1912-5	12 Fr / 4mm	16" / 40.5cm	"Y" Endotracheal Tube Adapters (5.5mm, 6.0mm, 6.5mm)	20/Case



## **KIMBERLY-CLARK Neonatal and Paediatric MAC Catheters**



<b>Code</b>	<b>Outer Diameter</b>	<b>"Y" Size / Plug Size</b>	<b>Description</b>	<b>Packaging</b>
1900	5 Fr / 1.66mm	2.5mm / 7.5mm, 3.0mm / 7.5mm, 3.5mm / 7.5mm	MAC Catheter	5/Case
1940	5 Fr / 1.66mm	4.0mm / 7.5mm	MAC Catheter	5/Case

## **KIMBERLY-CLARK TRACH CARE Paediatric / Neonatal Component Kits**

- Gamma Sterilised
- Neonatal Sizes
- Available in T-Pieces, Elbows and Double Swivel Elbows
- Multiple Configurations
- Customer-Oriented Presentation Designed to Meet Unique Needs

<b>Code</b>	<b>Outer Diameter</b>	<b>Length</b>	<b>Description</b>	<b>Packaging</b>
1968127	6 Fr / 2mm	12" / 30.5cm	Neonatal "Y" WET PAK* (6 Vials) READY CARE* (5 Ea, 8 Fr Oral Catheters) 9', 274cm Suction Tubing	20/Case
1988128	8 Fr / 2.6mm	12" / 30.5cm	Neonatal "Y" WET PAK* (6 Vials) READY CARE* (5 Ea, 10 Fr Oral Catheters) 9', 274cm Suction Tubing	20/Case

## The KIMBERLY-CLARK ADVANTAGE\*

### KNOWLEDGE NETWORK<sup>®</sup> Clinical Education

Online at [HAWatch.com](http://HAWatch.com)  
On-site rep-facilitated programs

### Knowledgeable Customer Support

In-service training  
Product technical support  
Unsurpassed customer service

### Expert Sales Force

Healthcare Industry Representative Credentialed  
On-site trained in hospitals

Infection prevention website:

[www.HAWatch.com](http://www.HAWatch.com)



### Tools & Best Practices

Infection Prevention Education Toolkits  
Product Use and Selection Tools  
Best Practices/ Industry Guidelines  
Cost-in-use Calculators  
Inventory Control Tools

### Clinical Research

Staff medical professionals to advise and direct clinical research  
Clinical Scientific research to measure clinical outcomes in patients  
Peer-to-peer consultation

### Commitment to Excellence

*If, for any reason, our products do not meet your expectations, please let us now your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.*

*At Kimberly-Clark, we deliver innovative healthcare solutions that you can depend on to meet the demands of your fast-paced world, supported by in-service training, clinical research and accredited education. Whenever your needs involve infection prevention, surgical solutions, digestive health, or Respiratory Health Solutions, with Kimberly-Clark solutions you'll have one less worry.*



Protection & Infection Prevention



Surgical Solutions



Respiratory Health Solutions



Digestive Health

Visit our website at [www.kchealthcare.com](http://www.kchealthcare.com)

#### Australia:

Kimberly-Clark Australia Pty. Limited ABN 65 000 032 333  
52 Alfred Street, Milsons Point, NSW 2061  
To Order: 1800 810 077  
Advisory: 1800 101 021  
Fax: 1800 012 255

#### New Zealand:

Kimberly-Clark New Zealand  
86 Plunket Avenue, Manukau, 2104 Auckland  
Tel: 0800 482 211  
Fax: +649 250 2511

Kimberly-Clark makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

\* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc.

© 2012 KCWW

KCHC2012WFKC369



**Kimberly-Clark**

*Trusted Clinical Solutions\**