



Nitrile

## Nitrile exam glove with HYDRASOFT™ Skin Moisturizing Technology

- Designed with our patented HYDRASOFT™ Skin Moisturizing Technology which is proven to allow skin to retain 2X more moisture\* than gloves without HYDRASOFT™
- Soft nitrile formulation provides a comfortable fit and feel
- Thin and lightweight to ensure excellent dexterity



NOW FEATURING  
 **HYDRASOFT™**  
Skin Moisturizing Technology

### Industries

- Healthcare
- Life Sciences
- Food Processing

### Recommended For

- Sample taking and processing
- Protection from Type I latex allergy in HCW's or patients
- Standard patient examinations
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Oral examinations
- Food contact
- Laboratory Analysis
- Testing
- Helping protect skin from drying or cracking
- Pharmaceutical manufacturing
- General dentistry



# MICROFLEX®

ADVANCED DISPOSABLE HAND PROTECTION

## Soft White Nitrile TQ-601

Comfortably protects hands while replenishing moisture and preventing dryness

### TECHNICAL DATA SHEET

#### PRODUCT INFORMATION

Soft White Nitrile TQ-601	
Material	Nitrile
Color	White
Glove Design	Chlorinated, Coated with HYDRASOFT™ Technology, Powder-Free, Textured Fingers
Cuff	Beaded
Regulatory/Standards Compliance	ASTM D6319, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 13485, FDA21 CFR 177-2600, US QSR/FDA 510(k) Medical Examination Grade
Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Country of Origin	Indonesia
User Needs Segment	High Touch
Available sizes	XS, S, M, L, XL

#### PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	240 / 9.5		ASTM D3767,EN420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 455-1
Palm Thickness (mm/inches)	0.07 / 2.8		ASTM D3767,EN420
Finger Thickness (mm/inches)	0.10 / 3.9		ASTM D3767,EN420
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	24	26	ASTM D412
Elongation at Break (%)	500	400	ASTM D412
Force at break (N)	7	7	

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

#### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
Riverside Business Park  
Blvd International, 55  
1070 Brussels, Belgium  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

#### Asia Pacific Region

Ansell Global Trading Center  
(Malaysia) Sdn Bhd  
Prima 6, Prima Avenue  
Block 3512, Jalan Teknokrat 6  
63000 Cyberjaya  
Selangor, Malaysia  
T: +603 8310 6688  
F: +603 8310 6699

#### North America Region

Ansell Healthcare Products LLC  
111 Wood Avenue South,  
Suite 210  
Iselin, NJ 08830, USA  
T: +1 800 800 0444  
F: +1 800 800 0445

#### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
Blvd. Bernardo Quintana No. 7001-C,  
Q7001 Torre II,  
Suites 1304, 1305 y 1306,  
Col. Centro Sur, c.p. 76079  
Queretaro, Qro. Mexico  
T: +52 442 248 1544 / 248 3133

#### Australia

Ansell Limited  
Level 3, 678 Victoria Street,  
Richmond, Vic, 3121  
Australia  
T: +61 1800 337 041  
F: +61 1800 803 578

#### Russia

Ahcenn PYC  
123610 Poccnr, Mockba  
Kpachonpechehckar  
Ha6. 12, n.3, op. 1304A  
Ten. +7 495 258 13 16

Technology  
  
Skin Moisturizing Technology

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2019 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

