


NEW



## CaviWipes Bleach Disinfecting Towelettes for hard, non-porous surfaces.

With patent pending Enhanced Chlorine Technology™ (ECT), these new multi-purpose cleaning wipes are gentle on surfaces, yet tough on pathogens. They kill *C. diff* spores in 3 minutes—and do not require PPE\*. They also meet the CDC recommended 1:10 bleach dilution and are compatible with most medical device equipment\*\*. So when it comes to disinfecting hard non-porous surfaces, leave the gloves behind and reach for CaviWipes Bleach.



- Kills *C. diff* in 3 minutes - no other surface wipe is faster



- No PPE required\*



- Compatible with most medical device equipment\*\*



- One-step cleaner and disinfectant



- Low Odor

Trust the surface disinfectant brand used by thousands of healthcare facilities every day.

\*Special instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores, HIV-1, HBV, HCV.  
\*\*See medical device equipment manual for compatibility information.

CaviWipes Bleach have demonstrated effectiveness against the following 41 relevant microorganisms in 3 minutes:

**3** MINUTE KILL TIME

## Bacteria

*Clostridium difficile* spores  
*Staphylococcus aureus*  
*Pseudomonas aeruginosa*  
*Salmonella enterica*  
*Burkholderia cepacia*  
*Escherichia coli*  
*Klebsiella pneumoniae*  
*Legionella pneumophila*  
*Listeria monocytogenes*  
*Serratia marcescens*  
*Streptococcus pyogenes*  
*Neisseria gonorrhoeae*

## Mold, Mildew, Fungi

*Candida albicans*  
*Aspergillus brasiliensis*  
*Trichophyton interdigitale*

## Non-enveloped Viruses

Adenovirus  
 Hepatitis A Virus  
 Norovirus (as Feline Calicivirus)  
 Poliovirus  
 Rhinovirus  
 Enterovirus (EV-D68)  
 Rotavirus

## Drug Resistant Bacteria

*Enterobacter cloacae*  
 Multidrug Resistant *Acinetobacter baumannii* (MRAB)  
 Carbapenem-resistant *Klebsiella pneumoniae* (CRKP) (CRE)  
 Methicillin Resistant *Staphylococcus aureus* (MRSA)  
 Penicillin Resistant *Streptococcus pneumoniae*  
 Vancomycin Resistant *Staphylococcus aureus* (VRSA)  
 Vancomycin Resistant *Enterococcus faecalis* (VRE)  
 ESBL *Escherichia coli*

*Mycobacterium tuberculosis* var: *bovis* (BCG) (TB)

## Enveloped Viruses

Herpes simplex virus type 1  
 Herpes simplex virus type 2  
 Human Coronavirus  
 Influenza A virus Strain A/Hong Kong 8/68 (H3N2)  
 Influenza B virus Strain B/Hong Kong/5/72  
 HBV  
 HCV  
 HIV-1  
 Respiratory syncytial virus (RSV)

## Parvovirus

Canine Parvovirus



**NO PPE REQUIRED**

Metrex offers complimentary inservicing and wall charts for your education needs.



## Part No. Product Description

13-9100 CaviWipes Bleach 6" x 10.5", 90 wipes per canister  
 13-1175 CaviWipes Bracket

## Quantity

12 canisters/case  
 12 brackets/case

Visit [Metrex.com](http://Metrex.com) to learn more about CaviWipes Bleach.