

# MaxiCide® OPA 28

**From A Name You Can Trust.**

*The Fastest and Longest Lasting  
OPA on the Market Today!*

## WORKS FASTER

Take quicker action against microorganisms. High-level disinfectant disinfects devices manually in just 10 minutes at room temperature (20°C minimum), to achieve bactericidal, fungicidal, and tuberculocidal kills. Studies indicate MaxiCide® OPA 28 effectively inactivates:

- MRSA
- VRE
- TB
- Hepatitis Viruses
- HIV

OPA's faster manual disinfection allows you to reprocess one extra scope per hour as compared to other OPA products.

## LASTS LONGER

*Use for up to 28 days*

It's important for health care providers to manage budgets carefully while still delivering high-quality patient care. MaxiCide® OPA 28 has been tested to provide 84 reuse cycles\* along with a reuse period of up to 28 days—twice the reuse period of other OPA brands.

*\*Minimum recommended concentration*

**(112-6976)** MaxiCide® OPA 28 1-Gallon Btl.  
**(112-6977)** MaxiCide® Test Strips 50 Strips/Btl.



## The MaxiCide® OPA 28 Advantage by the numbers

**5**  
minutes

AER disinfection  
time at 25°C

**10**  
minutes

MANUAL disinfection  
time at 20°C minimum

**84**  
reuse cycles

Up to 10–24 additional  
AER reuse cycles†

**28**  
days

reuse period, twice as long  
as other OPA brands

**2**  
final rinses

needed in  
MEDIVATORS® AER

**3**  
final rinses

needed in manual  
applications



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**