

# DIGITAL ECG SOLUTION BY MIDMARK



IQ  
ecg





The IQecg® from Midmark is an easy-to-use, digital, 12-lead resting ECG solution that helps you quickly administer testing and access patient information you need to make informed cardiac assessments and provide more efficient patient care.

# IQecg®: An ECG that grows with your practice

Functions as a standalone product with our IQmanager® software or seamlessly integrates with today's top EHR/EMR solutions.

## INNOVATIVE FEATURES HELP IMPROVE PATIENT CARE

### Efficiency

- Lightweight, portable design, weighing less than 11 oz.
- View and monitor live, 12-lead ECG in selectable formats on-screen
- Capture abnormal ECG as it appears at the touch of a button
- Midmark 12-lead Resting Electrocardiogram Analysis Program automatically analyzes and interprets ECG data for neonate, pediatric and adult patients
- Review and edit ECG results on-screen with free text or customizable quick edit statements
- Detects possible arm lead reversal during the test so the technician can take appropriate action while the patient is still connected

### Technology

- Save time and reduce data entry errors by populating patient information directly from your EHR/EMR software
- Facilitate a paperless office strategy and eliminate the costs of thermal paper and the hassles of scanning and filing patient charts
- Stores binary ECG waveforms and displays in digital data instead of static images, resulting in smaller file sizes and providing accurate viewing and measurements
- Quickly assess ECG changes with direct overlay serial comparison of up to four ECG tests
- Enhance your diagnosis with electronic tools like zoom and on-screen calipers

### Adaptability

- Functions as a standalone product with our IQmanager® software or seamlessly integrates with today's top EHR/EMR solutions
- Integration allows you to perform tests within the EHR/EMR and store the information in the patient's electronic chart
- IQpath™ software solution offers uninterrupted testing, reduced latency issues and smooth data flow for facilities operating in a thin-client environment
- Part of a family of products created to grow with your practice



Lightweight, portable design, weighing less than 11 oz.




An easy-to-use, 12-lead, digital ECG that helps to quickly administer testing.



Lightweight, portable design makes moving from room to room quick and easy.



Review and edit ECG results on-screen.



Detects and warns possible arm lead reversal while in live ECG screen.

Make your practice more efficient with an IQ Digital Diagnostic System  
**Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to grow as your practice grows.**

Diagnostic Products				
	EHR/EMR Connectivity	Single Platform (IQmanager®)	Standalone	Thin - Client (IQpath™)
IQecg®	■	■		■
IQholter®	■	■		■
IQholter® EX	■	■		■
IQholter® EP	■	■		■
IQstress®	■	■		
IQspiro®	■	■		■
IQvitals®	■	■	■	■
IQvitals® PC	■	■		■



## RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working with you to redefine the future of the clinical space.

[midmarkclinicalsolutions.com](http://midmarkclinicalsolutions.com)

### Minimum System Requirements

- Windows® 8, Professional and Enterprise; Windows® 7; Windows® Vista; Windows® XP Pro; 32-bit and 64-bit
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- USB port
- Mouse, keyboard and high resolution display monitor

Contact us or visit [midmark.com](http://midmark.com) for more detailed system requirements.

Windows is a registered trademark of Microsoft Corporation.  
 Intel Core is a trademark of Intel Corporation in the U.S. and/or other countries.

### Thin-Client (IQpath™) System Requirements

Please contact your local Midmark Sales Representative for detailed system requirements.

### Ordering Information

4-000-0061 IQecg® Digital ECG USB Version

For information, pricing or to place an order, contact us at:  
 1-800-624-8950 or  
 1-310-516-6050 or  
[midmark.com](http://midmark.com)

Midmark is an ISO 13485 Certified Company  
 Midmark Corporation  
 690 Knox St, Suite 100, Torrance, CA 90502  
[midmark.com](http://midmark.com)  
 1-800-624-8950

© 2011 Midmark Corporation  
 60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286  
 Products subject to improvement changes without notice  
 Litho in U.S.A. 47-80-0001 Rev. A (3/13)

### Accessories



Universal ECG Clips (3mm, 4mm & Snap)  
 3-047-0001



Clear ECG Clips (4mm & Snap)  
 3-047-0005



IQcart®  
 3-004-1000



Disposable EZ-Trode ECG Electrodes  
 2-100-0205 (Box of 500)  
 2-100-0206 (10 boxes of 500)



IQmanager® Software  
 4-100-1200



Because we care.