

Fecal Immunochemical Testing

From a Name you can Trust.

OneStepPro+ FIT tests are a simple strip procedure that provides qualitative results in 5 minutes. Requires no dietary or medicinal restrictions, while ensuring quick and accurate results.



OneStepPro+ FIT Tests

Personal Use Kit (20/Bx)

570-0607

20 sampling bottles
20 mailing envelopes

Manual Kit (50/Bx)

570-0609

50 test strips
50 sampling bottles
50 collection papers

Strips Kit (50/Bt)

570-0610

50 test strips

Sample Bottles

570-1172

50 bottles

Control Kit

570-1171

1 positive control
1 negative control

Intended Use: Complete GI

Dietary Restrictions: None

Patient Sample Collection Method: Sampling bottle and collection paper

Number of Specimens Needed: 1 sample

Sample Stability: 15 days room temperature and 30 days when store
2-8 degrees Celsius

Sensitivity: 50 ng Hb/ml

Time for Results: 5 minutes

Reimbursement Codes: G0328QW: Colorectal cancer screening with immunochemical FOBT

82274QW: Blood, occult, by fecal hemoglobin determination with immunochemical FOBT



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!